

CytoScan XON data analysis in Chromosome Analysis Suite (ChAS) software

A quick reference guide for CytoScan XON data analysis using ChAS software

Software version

For analysis of Applied Biosystems™

CytoScan™ XON microarray files, you must install:

- Applied Biosystems™ ChAS software v3.3 or higher
 - [Link to download ChAS v3.3 or higher](#)
- CytoScan XON library files
- Browser/NetAffx™ Genomic Annotations for **hg19** and/or **hg38** (you can only use one at a time)

Supported sample types

Default QC metrics are for postnatal blood samples and could vary for prenatal and cancer specimens.

Suggested files to use

CytoScan XON Suite (.xnchp) sample file(s): Open

(ensure that NetAffx Genomic Annotations version is correct)

XON Region Levels: **On** or **Off** per user preference (can use **Filters** to select gene levels 1–4)

Default Histogram/XON DGV Histogram: **On** (CytoScan XON_aDGV is derived from 1,855 peripheral blood samples from phenotypically normal individuals)

Genes or OMIM® Genes: **On**


AED/BED files: **On** if relevant to phenotype or NGS results. In addition to disease-specific AED/BED files, the following Database of Genomic Variants AED files may be useful:

- CytoScan_XON_DGV_Gains_hg19_20180521.aed
CytoScan_XON_DGV_Losses_hg19_20180521.aed


Or

CytoScan_XON_DGV_Gains_hg38_20180521.aed
CytoScan_XON_DGV_Losses_hg38_20180521.aed

- If desired, you can load one or more additional AED/BED files. Right-click on each file individually and select “**Include in CytoRegions**”. Multiple AED/BED files will be treated as one large CytoRegion, with vertical gray bars running through the data, indicating clinically significant gene regions.
- To perform truly targeted analysis, you can restrict the analysis to only those CytoRegions, and all other regions will be hidden from the view:

– **View** menu > **Restrict to CytoRegions**
(or  button)

VCF files: **On** if NGS findings are relevant. (VCF files must contain only genotype or indel data.)



Database: These aDGV XON files can be useful for comparison to segment calls in your samples:

- ChASDB_XON_aDGV_20180626.backup
ChASDB_XON_aDGV_20180628_hg38.backup

Or

You may prefer to use your own database of abnormal segments in conjunction with the XON_DGV_Gains.aed and XON_DGV_Losses.aed files mentioned above.

- If you use the XON_aDGV as the database, you can hover your mouse over the histogram (**Default Histogram**), and it will show the number and percentage of samples in the database with that same segment.
- If you use the XON_DGV_Gains/Losses.aed files, the frequency of positive samples will appear above the segment in the AED file, but only if the frequency is greater than 1%. Thus, if you use both the database and the AED files, there may be a discrepancy in frequency.

Named Setting section

Select **XON-Level 1** or user-defined **Named Setting** as shown here.

Data Types section

Gain (XON Region) and Loss (XON Region): **On**

Log2 Ratio: **On**



Summarized Log2Ratio: **On**

- If you have **XON Merging** turned **On** in the **User Configuration**, you won't see the **Summarized Log2Ratio** values for merged segments, because **Summarized Log2Ratio** is specific to each exon region.

LOH: **Off**, unless looking for large regions (>5 Mb)

SmoothSignal: **Off**, unless looking for larger regions, because there's no distinct difference for exon-level changes. If you choose to have **SmoothSignal On**, adjust the graph display to **Bar** (rather than **Line**) for easier visualization and interpretation.


Allele Difference and **BAF** (B Allele Frequency): **Off**, unless looking for large regions



Segments Filters

Gene levels:

- Level 1: Genes of high clinical relevance
- Level 2: ClinVar genes not in Level 1
- Level 3: Other OMIM genes
- Level 4: Other regions from RefSeq, UCSC, Ensembl, and/or Leiden Open Variation Database (LOVD)
- **Start with Level 1 genes to minimize number of calls for analysis**
 - In ChAS, **Level 1** genes are displayed in **yellow blocks**
 - In ChAS, **Levels 2–4** are displayed in **brown blocks**
- **Gain:**
 - Marker Count: **Off**
 - Size (kbp): **Off**
- **Loss:**
 - Marker Count: **Off**
 - Size (kbp): **Off**
- **LOH:**
 - Size (kbp): **On**, minimum of 5 Mb




User Configuration

Segment Data tab:

- Use default segment data rules configuration: **Off**
- XON Merging: **Off**

DB Query tab:

Additional information on coverage and overlap functions can be found in the ChAS User Guide.



User Configuration

Segment Data | QC Thresholds | Color Rules | Misc | Vocabularies | DB Query | Filtered DB Query | Exports


CytoScan XON

For CHP files of each array type, edit the configuration of Smoothing, Joining or Merging to create Copy Number Segments. This processing happens at file loading and at configuration saving (always before any Segment Filtering).

☒ Use default segment data rules configuration

XON Merging [merges XON segments of the same type (Gain, or Loss) which not separated by any normal state data]

☒ Merge adjacent XON Gain or XON Loss segments.



User Configuration

Segment Data | QC Thresholds | Color Rules | Misc | Vocabularies | DB Query | Filtered DB Query | Exports

Copy Number Query Parameters

Minimum % CN Coverage: 50

Minimum % CN Overlap: 50

☒ Match only same gain/loss type ☐ Include LOH ☒ Include XON Regions

LOH Query Parameters

Minimum % LOH Coverage: 50

Minimum % LOH Overlap: 50

XON Regions Query Parameters

Minimum % Coverage: 50

Minimum % Overlap: 50


☒ Match only same gain/loss type ☐ Include Copy Number Segments

Minimum % Coverage: The minimum percentage that the current segment is covered by other similar segments from the database.

Minimum % Overlap: The minimum percentage that the current segment overlaps similar segments in the database.

Segments Table column selection

- Ensure that **Summarized Log2Ratio** and **DB Count Both** are in the **Chosen** section (left in the screenshot at right) for visualization in the table view.
- Set up **Table Sorting** as shown at right (or in a user-defined way).



Select Columns

Choose columns to display by dragging between 'Available' and 'Chosen'. Re-order entries in 'Chosen' by dragging.

Chosen

- File
- CN State
- Type
- Chromosome
- Cytoband Start
- Size (kbp)
- Marker Count
- Summarized Log2Ratio
- DB Count Both
- Genes
- Gene Count
- CytoRegions
- DB Coverage Count
- DB Overlap Count
- Call

Available

- Allele Origin
- BACs
- Call Approval
- Condition Category
- Condition ID Type
- Condition ID Value
- Confidence (ReproSeq)
- Curation Time
- Ensembl Gene Count
- Ensembl Genes
- Filtered DB Count Both
- Filtered DB Coverage Count
- Filtered DB Overlap Count
- FISH Clones
- Fluorescing Marker

Filter 'Available':

Table Sorting

Optionally select one or more columns for sorting data in the table.

Sort by: DB Count Both ASCENDING Clear



Then by: Summarized Log2Ratio ASCENDING Clear


Then by: <No Selection> ASCENDING Clear

OK Cancel


CytoScan XON sample analysis

Open a CytoScan XON sample file:

1. Look at the **QC Metrics**; thresholds are shown below.
 - MAPD ≤ 0.20
 - SNPQC ≥ 10
 - Waviness SD ≤ 0.08
2. Go to **Filter Settings** and turn **On** “**Level 1 Genes**” only (check box).
3. Go to **Whole Genome View** to see if there’s anything obvious in the tracks (large **Loss/Gain**). 
4. Go to **Segments Table**.
 - There are two new columns in this table:
 - **Summarized Log2Ratio**
 - **XON Region Level**
 - Note that the software doesn’t give **Copy Number State**. It only calls **Gain** or **Loss**.
 - Note the number of calls in the upper right corner of the **Segments Table**:  is shown as an example
 - **Sort on DB Count Both**: Sort on “0”s first (large arrow up); focus on losses/gains that are <10% of what’s in the DB. In the aDGV database of 1,855 phenotypically normal samples, that would be ~18.
 - **Sort on Summarized Log2Ratio**: Sort on most negative values (small arrow up) for **Loss**; values below 0.5 might be more “real”. Sort on most positive values (further down the list) for **Gain**. The larger the number, the more likely it is to be a true call.
5. Apply **Gene Lists/AED/BED** files (DECIPHER, disease-specific files, etc.), especially when there’s a distinct phenotype. Right-click on the name of the file and select **Set a CytoRegion for Targeted XON analysis**.
6. If you’ve set a CytoRegion for targeted analysis (step 5), you can also restrict the view to these areas using the **Restricted Mode** button in the toolbar, or go to the **View** menu > **Restrict to CytoRegion**. This will hide all data with the exception of what’s in the targeted region file.

7. Using **Edit** mode , you can manually edit the calls to merge exon segments into one larger call, or adjust the breakpoints. You can assign a copy number if you feel confident about it, but it's not necessary. As mentioned previously, merging segments will eliminate the **Summarized Log2Ratio** for an exon call.
8. Some probes (SNPs) are used for **Allele Difference** or **B Allele Frequency**, but not for XON calls. These probes are shown as gray dots in the tracks, not as colored dots.
9. **Smoothing** and **Joining** are optional and should be turned **Off** in **User Configuration**.

General comments

- For a **targeted approach**, as an orthogonal method to confirm NGS or data from other microarrays:
 - Search for a gene or region, and/or load targeted **AED/BED** file(s), and right-click on file name to create a CytoRegion for targeted XON analysis.
 - Restrict the region of analysis as described above, if desired.
 - Scan through the **Loss** and **Gain** segments in the vertical gray bars, examining the size, number of markers, number of genes, gene significance, etc., and confirm whether the calls are significant.
 - Examine any large (≥ 5 Mb) regions of LOH to determine relevance.
- For a **nontargeted approach**, for samples where you may be blinded to other findings, or there's no accompanying case history or phenotype:
- Ultimately, you should sort data based on **DB Count Both (ascending)** and then on **Summarized Log2Ratio (ascending)**. You can do this in two different ways:
 - Click on the **Table** icon at the top right corner of the **Segments Table** tab:  The **Select Columns** table appears, in which you can drag attributes from the right (**Available**) column to the left (**Chosen**) column. You can then drag the attributes up or down to place them in the desired order in the **Segments Table**.
 - Alternatively, you can click on the column header and double-sort by dragging the columns to the desired location.
 - Ultimately, you should focus on those segments that have a small **DB Count Both**, but a large **Summarized Log2Ratio** (most positive and most negative numbers) to help gauge whether it's likely to be a true call vs. a false-positive call, or a normal copy number variant. Opening other samples (.xnchp files) at the same time can provide additional insight as to veracity.
 - **Hint:** For exploratory studies, use an arbitrary cutoff value of 15–18 counts for **DB Count Both** when using the aDGV as the database. (This would represent less than 1% of the samples in the database.)

Notes

[illegible]

Notes

[illegible]