## EN (English)

### Instructions for Safe Use of IVD/CE marked F1-ClipTip pipettes and ClipTip tips

The CE mark covers the combination of CE marked F1-ClipTip pipette and CE marked ClipTip tip. (F1-ClipTip - ClipTip Compatibility Table, see Instructions for Use). Any other pipette - tip combination has to be validated by the user. In IVD use, the F1-ClipTip is to be used by **trained personnel**, who has adequate laboratory skills. The pipette shall be **regularly serviced** according to the instructions for use.

**Intended Use:** Thermo Scientific F1-ClipTip pipettes used in conjunction with the corresponding compatible IVD/CE-marked ClipTip tips form a pipetting system which is designed and constructed to transfer liquids originated from human body or reagents related to an IVD application.

Therefore, the ClipTip pipetting system serves as an IVD accessory under the terms of the Directive 98/79/EC. These products shall be used in clinical laboratory environment by trained personnel, who have adequate clinical laboratory skills.

As this pipetting system is part of the analyzing system of the end user, the user is responsible for validation of the whole system to enable production of reliable and safe results. If the assay performance is essential to the analysis, the test result has to be assured with additional means e.g. using internal quality controls or an alternative test.

**Performance:** When leaving the factory, the pipette is tested against the **factory specifications**. The performance in normal use may vary (due to: temperature, pressure, humidity, operator, liquid, type of tip etc.). The user may also check the performance of the pipette with water against **limits of the standard ISO 8655** (See Instructions for Use). The user shall determine the **performance and applicability** of the pipette and tip combination **for the IVD application** and determine the required interval for performance checking. If the pipetting performance is critical to the medical diagnosis the diagnostic test result has to be ensured with an alternative test.

 **Prior to use:** Read the Instructions for Use. The performance of each new pipette needs to be checked and documented prior to use. The performance of each pipette after service needs to be checked and documented prior to use. Check that the package of sterile tips is unbroken (to ensure sterility). Check that the tip is complete and unbroken. Check that the filter tips have filters in place and that the filters are unbroken or not tilted. Do not reuse the tips. Before use, check that there is a tip attached to the pipette.

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