

## smart H<sub>2</sub>O

for you, and your science






# Thermo Scientific Barnstead LabTower EDI water purification systems

Pure water is elemental to the success of your experiments. But finding the right water system for your research goes deeper than water quality alone. You need a smart choice that supports both your science and your budget – on day 1, and with every cartridge and filter change. One that reflects 130 years of innovations, like effortless cartridge changeouts. With a Thermo Scientific™ Barnstead™ lab water system, the only thing you'll have in your water is confidence.

	LabTower EDI
Type 1 Water	
Resistivity at 25°C, MΩ•cm	18.2
Conductivity, µS/cm	0.055
TOC, ppb	1 - 5
Bacterial content, CFU/ml	<1
Particles, 0.22 µm/ml	<1
Flow Rate at dispenser, L/min	1.5
Type 2 Water (from storage tank)	
Pure water production at 15°C, L/hr	15 or 30
Resistivity at 25°C, MΩ•cm	10 - 1
Conductivity, µS/cm	0.1 - 1.0



### Two systems in one

- Unique system with pre-treatment plus polisher to produce both Type 1 and Type 2 water
- Type 2 water is stored in the integrated 100 L high purity water reservoir





### Barnstead LabTower EDI System Features

### Smart consumables

- Validated 0.2 μM final filter with folded membrane can be autoclaved up to 5 times
- Standard UV photo-oxidation at 185/254 nm reduces organic compounds in the water to ultra-low levels as well as microorganisms and their metabolites

### Easy-to-operate

- Dispensing is easy and features variable speed controlled flow
- The display can be tilted for optimal reading

### Compact, stylish, mobile

- Free-standing unit takes up no bench space
- Easily relocated with bottom-mounted rollers
- Continuous monitoring of all critical parameters
- Recirculation pump protects purified water in tank from bacterial growth during standstill


### **GLP-compliant**


- Real-time clock and code-protected operating system prevents unauthorized changes to system settings
- RS-232 interface with adjustable send-interval for safe data transfer of all measured data, faults, date and time to a PC computer or log printer
- Digital microprocessor control automatically monitors and stores faults from the last four weeks
- USP-compliant conductivity measurement with temperature compensation that can be switched on or off

### 100 L Integrated polyethylene reservoir safely stores Type 2 water

- Integrated reservoir provides 100 L of purified water storage
- Automatic water recirculation in the reservoir moves water across a special polisher module ensuring fresh Type 2 water on demand
- High-purity water reservoir outlet for convenient drainage
- Efficient cleaning and disinfection enhanced by conical bottom design for completely emptying reservoir
- Sterile vent filter and reservoir overflow prevents internal contamination by microorganisms from the surrounding air
- Optional CO<sub>2</sub> adsorber prevents CO<sub>2</sub> adsorption into the water, to help save cartridge life
- Highly visible reservoir volume display on LabTower system controller
- Customize volume stored in the reservoir via the controller

### **EDI Technology**





Electrodeionization, EDI, unites two proven technologies for producing ultrapure water: Electrodialysis and ion exchange. In contrast to conventional ion exchange in which resins must be either chemically regenerated or the cartridge discarded, EDI utilizes an electric current for continual resin regeneration.

### Technologies that keep you one step ahead

- Constant high efficiency of the ion exchange layer bed produces consistently high quality water
- and no disposal of cartridges, making this technology friendly for the environment
- it with no wait time

### How EDI works


Several layers of ion selective membranes are situated between an anode and a cathode. Layered ion exchange resin beds and concentrate chambers are alternately positioned between them.

On applying an electric voltage, water (H<sub>2</sub>0) is split into H+ and OH- in the cell.


The H<sup>+</sup> and Na<sup>+</sup> cations can migrate through the cation permeable membranes, anions through the anion permeable membranes.

The ions migrate in the direction of the applied voltage, i.e. anions to the positive pole (anode), cations to the negative pole (cathode). The water ions that migrate through an ion exchange chamber displace salt ions retained by the ion exchange resins and so continually regenerate the resins.

The salt ions migrate through the appropriate ion selective membranes into the concentrate chambers and are flushed out by water. As all concentration chambers are flushed, excess H<sup>+</sup> and OH<sup>-</sup> ions can again combine to form H<sub>2</sub>0.



### The following diagram illustrates the flow of water through the components of a Barnstead LabTower EDI water system.



#### **LabTower EDI**

Product Dimensions H x W x D, mm (in.)	1500 x 450 x 580 (59 x 18 x 23)
Product Weight kg (lbs.) w/o water	66 (146)
<b>Shipping Dimension</b> s H x W x D, mm (in.)	1210 x 1200 x 800 (48 X 47 X 32)
Shipping Weight kg (lbs.)	77 (170)

#### Feed Water Requirements\*

Source	Potable tap water softened or hardness stabilized.
Conductivity, µS/cm	< 1000
Silt Density Index (SDI)	< 3
pH Range	4 - 11
Temperature, °C	2 - 35
Pressure, psi (bar)	29 - 87 (2 - 6)

<sup>\*</sup> Please see user manual for complete list of feed water requirements.

### ordering information

LabTower EDI System Options	LabTower EDI 15	LabTower EDI 30
All systems include an ultrapure polisher cartridge, UV lamp, RO membrane, RO pretreatment filters (5 µm filter and hardness stabilizer), EDI module, sterile 0.2 µm filter, 10" 1 µm filter (tank outlet) and pressure regulator.	50132395	50132396

Required Accessories			
Sterile overflow for reservoir Prevents reservoir contamination of bacteria and other micro-organisms.		50132714	
CO <sub>2</sub> adsorber + sterile filter, 0.2 μm Combination sterile filter with CO <sub>2</sub> adsorber is designed to prevent CO <sub>2</sub> from entering the tank, to help save cartridge life.		06.5002	
Mix Multi Mini water softener Required for hard feed water and/or if silt density index (SDI) is greater than 3. Also required for purchase with a softener is the softener salt, hardness detection kit, and the 5 µm filter with carbon.	120V, 50/60Hz	50129892	
	240V, 50 Hz	50130297	
Softener salt Required for use with the water softener.	Europe	50129893	
	North America/Latin America Asia Pacific	06.2000	
5 μm filter with carbon cartridge Required when Mixed Multi Mini is purchased.		06.5201	
Hardness detection kit – required with purchase of softener Alerts user when water is no longer softened.	Europe	06.1000	
	North America/Latin America Asia Pacific	50134335	

Optional Accessories			
Hand dispenser kit Hand dispenser with 3 meter cord that connects to tank (to dispense Type 2 water). Ships with a 0.2 µm final filter.		50138221	
Sterile 0.2 µm filter for reservoir outlet		06.5555	
<b>Disinfection cartridge</b> Reusable cartridge needed to clean the water purification system.		09.2201	
Qualification documents (IQOQ)		IOQD0CE50133916	
<b>Printer</b> Utilizes RS-232 interface for safe documentation of all measured values and faults with date and time in compliance with GLP-Guidelines.	120V, 50/60Hz	AY1137X1	
	230V, 50Hz	09.2207	

Replacement Consumables			
Ultrapure polisher cartridge		09.2005	
CO <sub>2</sub> adsorber + sterile filter, 0.2 μm, for tank		06.5002	
Hardness stabilizer, 5" with 5 µm filter		06.5204	
Sterile 0.2 µm filter for system		09.1003	
UV lamp for the system		09.2002	
Oleaning colution	Europe/Asia Pacific	09.2202	
Cleaning solution	North and South America	CMX25	
Reverse osmosis membrane		22.0046 (requires 2)	22.0087 (requires 2)



#### thermoscientific.com/purewater

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All (other) trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Australia +61 39757 4300 Austria +43 1 801 40 0 Belgium +32 53 73 42 41 China +800 810 5118 or +400 650 5118 France +33 2 2803 2180

France +33 2 2803 2180
Germany national toll free 0800 1 536 376
Germany international +49 6184 90 6000

India toll free 1800 22 8374 India +91 22 6716 2200 Italy +32 02 95059 552 Japan +81 3 5826 1616 Netherlands +31 76 579 55 55 New Zealand +64 9 980 6700 Nordic/Baltic/CIS countries +358 9 329 10200 Russia +7 812 703 42 15 Spain/Portugal +34 93 223 09 18 Switzerland +41 44 454 12 22 UK/Ireland +44 870 609 9203 USA/Canada +1 866 984 3766

Other Asian countries +852 2885 4613 Countries not listed +49 6184 90 6000

