# **Application Note 020**

# Measuring pH in Soil

Water Analysis Instruments, Thermo Fisher Scientific

#### **Key Words**

pH, soil, alkalinity, acidity strength, agriculture, environmental testing.

#### Goal

This application note explains how to measure the pH (alkalinity) of soil using a Thermo Scientific<sup>™</sup> Orion<sup>™</sup> ROSS<sup>™</sup> Sure-Flow<sup>™</sup> combination pH electrode.



#### Meter Setup

Connect the pH electrode, stirrer and ATC probe to the Orion Star A211 pH meter. Set measurement mode to pH. In the meter setup mode, set resolution to 0.01, Buffer Set to USA and read type to continuous. If all steps were followed correctly the meter display will show units as "pH" and a value with two decimal places. The temperature will also be displayed on the screen.

### **Electrode Setup**

See the electrode manual for preparation of the electrode.

#### **Electrode Performance Check**

Check slope at least daily according to the electrode manual. Drift may be checked by comparing a 1-minute to 2-minute reading. Results should agree with desired criteria. See troubleshooting section of manual if slope and/or drift are not acceptable.



# Introduction

The pH of soil indicates its alkalinity or acidity strength. This measurement is important for determining which plants or crops would thrive the most in the soil. This measurement can be made by suspending soil in water and measuring the pH with an epoxy body pH electrode. In this application note, we used a Thermo Scientific Orion ROSS Sure-Flow combination pH electrode with an epoxy-body, which is designed to measure pH in soil and other difficult samples.

#### Equipment

- Thermo Scientific<sup>™</sup> Orion Star<sup>™</sup> A211 pH benchtop meter kit (Cat. No. STARA2119), which includes an Orion 9165BNWP ROSS Sure-Flow pH electrode, an Orion 927007MD ATC temperature probe, a meter-attachable electrode stand, pH 4/7/10 buffers and standard storage solution
- Orion automatic stirrer probe (Cat. No. 096019)

#### Solutions

- Orion pH electrode fill solution (Cat. No. 900011)
- Deionized water (DI)

**Solutions Preparation** 

None required.

# Electrode Storage, Soaking, and Rinsing

See electrode manual for storage:

- 1. Between measurements
- 2. Overnight, and
- 3. For long periods of time.

Between measurements, rinse the electrode with DI water and blot dry with lint free tissue before measuring the next sample.

#### **Sample Preservation**

None required.

# **Sample Preparation**

Weigh 20g of soil into a beaker. Using a graduated cylinder, add 20mL of deionized water to the soil, mix for a 30 minute period and let stand for 1 hour. For precise measurements, allow all the standards and the samples to reach the same temperature before analysis.

### Calibration

Perform a two point calibration using pH 7.00 and 10.01 buffers as these bracket the expected sample pH. The electrode slope should be between 92 and 102%. Read a fresh portion of buffers to verify calibration. If readings are not acceptable and/or slope is not within range, see the troubleshooting section of your manual.

### Analysis

Rinse electrode, ATC probe and stirrer with DI water and blot dry. Place all probes in sample and measure. The pH value will be displayed. When a stable reading is achieved, the screen will stop flashing and the "ready" message will appear.

# Comments

The pH of soil samples can vary greatly depending on the source of the soil. If the pH of the first sample measured is not in the 7.00-10.01 range, it is recommended to recalibrate the electrode using buffers which will bracket the sample pH.

# Quality Control (QC)

Recommended QC procedures include: calibration and calibration verification, sample duplicates, slope, drift and lab control sample.

# Results

Soil Sample	рН
Sample 1	7.96
Sample 2	7.94
Sample 3	7.91
Sample 4	7.98
Sample 5	7.92
Mean	7.94
Standard Deviation	0.03

### **Result Statistics**

# Trials	Average	SD
5	7.94	$\pm$ 0.03 pH units

To purchase Thermo Scientific laboratory products, please contact your local equipment distributor and reference the part numbers listed below:

Product	Description	Part Number
Meters	Thermo Scientific Orion Star A211 pH benchtop meter kit, which includes an Orion 9165BNWP Sure-Flow epoxy-body pH electrode, an Orion 927007MD stainless steel ATC probe, an Orion 910199 pH buffer and solution kit and electrode stand	STARA2119
Electrodee	Thermo Scientific Orion Sure-Flow combination pH electrode	
Electrodes	Thermo Scientific Orion ATC probe	927007MD
Solutions	Thermo Scientific Orion pH electrode filling solution	900011
	Thermo Scientific Orion ROSS pH electrode storage solution	810001
	Thermo Scientific Orion pH 4.01, 7.00, and 10.01 buffers	910104, 910107 and 910110
	Thermo Scientific Orion automatic stirrer probe	096019
Accessories	Thermo Scientific Orion swing arm electrode stand	090043
	Thermo Scientific Nalgene Griffin polypropylene low-form beaker, 30mL	1201-0030
	Thermo Scientific Nalgene polypropylene graduated cylinder, 25mL	3662-0025
Deionized Water	Thermo Scientific Barnstead water purification system	Multiple

#### thermoscientific.com/water

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Water and Lab Products

#### North America

Toll Free: 1-800-225-1480 Tel: 1-978-232-6000 info.water@thermofisher.com Germany Tel: (49) 6184-90-6000 info.water.uk@thermofisher.com

China Tel: (86) 21-68654588 wai.asia@thermofisher.com India Tel: (91) 22-4157-8800 wai.asia@thermofisher.com

Singapore Tel: (65) 6778-6876 wai.asia@thermofisher.com Japan Tel: (81) 045-453-9175 wai.asia@thermofisher.com

Australia Tel: (613) 9757-4300 In Australia: (1300) 735-295 InfoWaterAU@thermofisher.com



A Thermo Fisher Scientific Brand