

**Overview** The concentration of acetic acid in canned products was determined by the Orion 960 Titrator PLUS using first derivative technique and using 0.1 N sodium hydroxide as titrant.

---

<b>Industry</b>	Food and Beverage
<b>Species Measured</b>	Acetic Acid
<b>Sample</b>	Canned Products
<b>Sample Size</b>	0.9g
<b>Typical Concentration</b>	2 % w/w
<b>Technique</b>	# 6 First Derivative
<b>Electrode</b>	Ross Sure-Flow Combination pH Electrode 8172BN
<b>Solutions</b>	Electrode Fill Solution 810007. 0.1N NaOH.
<b>Sample Prep</b>	Mix sample well and weigh 50g. Blend the sample with 500g of deionized water for 1 min. Weigh 10g of mixture to a 150 mL beaker. Add 50 mL DI water to the same beaker. Titrate.

#### Statistics

<b># of Trials</b>	3	<b>Mean</b>	98% Acetic Acid w	<b>%CV</b>	0.19%
<b>Analysis Time</b>	13.7				
<b>Comments</b>	Rinse the electrodes, stirrer, and dispenser probe between measurements with deionized water.				