2007 Chemiluminescent Product Guide





1 Chemiluminescent Substrates and Chemiluminescent Enhancers

1	Chemiluminescent Substrates and
	Chemiluminescent Enhancers1
2	Reporter Gene Assays and Reagents9
3	Immunodetection Products25
4	Nucleic Acid Membrane-Based Detection
	Products
5	Reagents and Accessories for
	Chemiluminescence

Introduction1
CDP-Star [®] Substrate and CSPD [®] Substrates
Galacton® / Galacton-Plus® / Galacton-Star® Substrates
Glucuron [®] Substrate5
Glucon [™] Substrate5
NA-Star [®] Substrate6
Solution-based Luminescence Enhancers
Sapphire [™] , Sapphire-II [™] , Emerald [™] , Emerald-II [™] , and Ruby [™] 7
Membrane-based Luminescence Enhancers
Nitro-Block™ and Nitro-Block-II™8

Introduction

Chemiluminescence

Chemiluminescence is the conversion of chemical energy to light energy. Several different chemical reactions, including some enzyme-catalyzed reactions, result in the production of visible light. Chemiluminescence reactions occur naturally (bioluminescence) in a wide variety of organisms, including beetles, jellyfish, bacteria, and many marine organisms. In addition, there are several classes of synthetic chemical structures that upon chemical or enzymatic cleavage produce light emission. Chemiluminescent reactions are employed in a wide variety of applications, including but not limited to biological assays, clinical diagnostic assays, biosensors, hygiene monitoring, and commercial low-level lighting.

Principles of Enzyme-activated Chemiluminescence

1,2-Dioxetane substrates emit visible light upon enzyme-catalyzed decomposition. Chemiluminescent detection of biomolecules with 1,2-dioxetane enzyme substrates is extremely sensitive as a result of low background luminescence coupled with high-intensity light output (due to enzyme cleavage and turnover of a large number of substrate molecules). The energy for light emission is generated internally upon dioxetane decomposition. In comparison, fluorescence requires an external light source for excitation energy, which must be filtered to discriminate the fluorescent signal emission. This limits the sensitivity and introduces complexity into the instrumentation and data analysis.

Glow Kinetics of 1,2-Dioxetanes

The decomposition of CDP-*Star** substrate is shown in Figure 1. Upon dephosphorylation of the substrate by alkaline phosphatase, a metastable phenolate anion intermediate is formed that decomposes and emits light with a maximum intensity at a wavelength of 475 nm. A delay in reaching maximum light intensity results, the length of which depends upon anion structure and the surrounding environment. Film or simple instrumentation may be used to quantitate the chemiluminescent signal that is produced as a steady glow arising from the reaction kinetics of the system.

1,2-Dioxetane Chemiluminescent Substrates

Tropix[®] 1,2-dioxetane chemiluminescent substrates enable extremely sensitive detection of biomolecules by producing visible light that is detected with film or instrumentation. Applied Biosystems offers several different enzyme-activated substrates, including CDP-*Star*[®] and CSPD[®] substrates for alkaline phosphatase; Galacton,[®] Galacton-Plus,[®] and Galacton-*Star*[®] substrates for β -galactosidase; Glucuron[®] substrate for β -glucuronidase; Glucon[®] substrate for β -glucosidase; and NA-*Star*[®] substrate for neuraminidase. Use of these reagents in immunoassays, enzyme assays, reporter gene assays, membrane-based protein detection, and nucleic acid detection on membranes or in tube or microplate assay formats offers substratial benefits compared to colorimetric, fluorescent, or isotopic detection.

Advantages of 1,2-Dioxetane Substrates

Applied Biosystems offers a wide selection of 1,2-dioxetane enzyme substrates that meet the challenging demands of a broad variety of applications. 1,2-Dioxetane substrates are non-isotopic and provide high intensity signal, low background, high sensitivity, wide dynamic range, rapid results, and are compatible with multiple assay formats under physiologically relevant conditions. The high quality, purity, and lot-to-lot consistency of Tropix substrates enable excellent reproducibility. In membrane-based assays, multiple film exposures can be acquired with standard X-ray film for over 24 hours following substrate addition.



Figure 1. Light Emission Mechanism of CDP-Star® Substrate.

Instrumentation for Chemiluminescent Signal Detection

Luminescent signal quantitation with microplate or tube solution-based assays is performed with a luminometer, which measures light being emitted from a sample with either a PMT-based or photodiode detector, or with a CCD camera detector. Commercial luminometers include dedicated (single-mode) luminescence detection platform, or multi-mode platforms that are equipped for luminescence detection. Luminescence can also be quantitated with a scintillation counter, with the appropriate counting mode selected, or with a fluorometer, without the presence of filters that would absorb light emission. Detection with non-luminometer instrumentation does not provide optimal sensitivity, as these instruments are frequently not light-tight, and are not optimized for luminescence detection. For some substrates and assays, automatic injection capabilities are required when the kinetics of light emission are very rapid.

Signal detection with membrane-based assays, including Western blotting and Southern blotting applications, can be performed with X-ray film, photographic film, or with CCD camera-based imaging platforms. Phosphorimaging platforms can be used if the appropriate screens, requiring sensitivity to visible light emission, are available.

Microplates

Luminescence assays in a microplate format are ideally performed in opaque white microplates, including 96-, 384- or 1536-well format, depending on assay design and capability, and on instrumentation capability. Black microplates can be used, but a significant absorption of the emitted light can reduce sensitivity. For cell-based assays performed in the culture plate, white plates with a clear bottom can be utilized, permitting microscopic observation of cells prior to start of assay. With these plates, a white backing sheet is applied to the bottom of the plate prior to reading, to prevent absorption of light by the black plate platform. The continuous clear plate bottom can introduce well-well signal crosstalk, so wells with anticipated large signal differences should not be adjacent. Clear microplates are not appropriate, as well-to-well crosstalk will obscure results.

Applications

Alkaline Phosphatase Substrates

CDP-Star[®] and CSPD[®] substrates are used in both solution and membranebased detection of alkaline phosphatase or alkaline phosphatase conjugates as superior alternatives to radioisotopic, colorimetric, or fluorimetric methods. For immunoassays and DNA probe assays, 1,2-dioxetane substrates improve the sensitivity [1,2,5,8,9,10,11] and generate results faster than radioimmunoassays and colorimetric detection methods. Chemiluminescent

Introduction, continued

detection of biomolecules labeled directly or indirectly with alkaline phosphatase (AP) is performed on membranes for Southern, Northern, and Western blotting, and DNA sequencing [3,4,9]. These substrates are also ideally suited for the detection of placental alkaline phosphatase (PLAP) and secreted placental alkaline phosphatase (SEAP) in reporter gene assays [6,12]. 1,2-Dioxetane substrates for alkaline phosphatase are used widely in immunoassay platforms for analytes, in identification of small molecules and viral nucleic acid detection. Additional potential applications include pasteurization monitoring, alkaline or acid phosphatase biomarker detection (stem cells, tumor markers) and sensitive quantitation of purified protein phosphatase activity [13].

β-Galactosidase Substrates

Galacton," Galacton-Plus," and Galacton-Star" substrates are widely utilized and have become the gold standard for sensitive quantitation of β -galactosidase in reporter gene assays in both mammalian and yeast cells. The chemiluminescent assay for β -galactosidase exhibits over three orders of magnitude greater sensitivity than colorimetric assays. In addition to reporter gene assays, 1,2-dioxetane substrates for β -galactosidase are used for β -galactosidase enzyme complementation assays. β -galactosidase enzyme complementation has a wide variety of applications, including both in vitro detection of biomolecules and cell-based assay systems that monitor cellular functions, such as intracellular protein-protein interactions, receptor activation, protein translocation, and cell fusion. Additional potential applications for chemiluminescent Galacton substrates include measurement of β -galactosidase activity in bacteria, for either gene expression monitoring or direct coliform detection, and immunoassay detection with β -galactosidase-labeled detection reagents.

β -Glucuronidase Substrate

Glucuron* substrate is a highly sensitive substrate for quantitating β -glucuronidase in reporter gene assays in plants or mammalian cells using the bacterial β -glucuronidase (GUS) gene. Chemiluminescent detection with Glucuron substrate exhibits greater sensitivity compared to fluorescence detection. Reporter gene assays employing Glucuron substrate are simple and convenient to perform. Additional potential applications include specific microbe detection and quantitation of mammalian β -glucuronidase. Please inquire regarding availability.

β -Glucosidase Substrate

Glucon[™] substrate is utilized for highly sensitive detection of β -glucosidase. The use of Glucon substrate provides researchers with another quantitative tool for the emerging widespread use of glycosidic enzymes in environmental and biomedical testing, clinical evaluation, toxicology, and pharmaceutical screening. Please inquire regarding availability.

Neuraminidase Substrate

NA-Star[®] substrate was designed for highly sensitive chemiluminescent detection of viral neuraminidase. NA-Star substrate is a sensitive chemiluminescent replacement for fluorescent neuraminidase substrates, and has been used widely for global monitoring of influenza strains for resistance to neuraminidase inhibitors.

Comparison of Chemiluminescent Systems and Substrates

An important feature of Tropix[®] alkaline phosphatase substrates is the longlived signal, especially on membranes. The chemiluminescent signal from CSPD[®] and CDP-*Star[®]* substrate may persist for up to several days on nylon membrane[7]. Since film exposure times range from seconds to several hours, multiple images may be acquired. Varying the film exposure time enables the user to optimize signal-to-noise. Other chemiluminescent systems, such as enhanced luminol, generate shorter-lived signals, making multiple film exposures difficult.

In solution assays (such as an immunoassay), the kinetics of CSPD and CDP-*Star* substrates are similar. Both substrates exhibit peak light emission within 30 minutes after adding substrate solution to a reaction well. Once the maximum signal is reached, it will be maintained as long as substrate is available (at least 60-90 minutes, depending on the amount of alkaline phosphatase present).

Membrane-based Detection Assays

CDP-Star* chemiluminescent substrate combines high intensity with rapid kinetics of light emission. This feature, coupled with up to five-fold higher signal intensity compared to CSPD substrate, makes CDP-Star substrate the ideal choice when rapid imaging exposures are needed. CDP-Star Substrate is used with Nitro-Block-II^{**} enhancer for nitrocellulose membranes. (see Chemiluminescence Enhancers, page 7).

Solution-based Detection Assays

CDP-Star[®] substrate also produces a higher signal in solution assays. CDP-Star substrate with either Sapphire-II[™] or Emerald-II[™] enhancer produces a signal that is nearly five-fold higher than the signal produced by CSPD[®] substrate with enhancer. For five-fold maximum sensitivity, CDP-Star Substrate solutions will provide the highest signal intensity and the greatest sensitivity.

REFERENCES (For complete reference lists, please see system-specific reference lists.)

- Albrecht S, Ehle H, Schollberg K, Bublitz R, Horn A (1991) Chemiluminescent enzyme immunoassay of human growth hormone based on adamantly dioxetane phenyl phosphate substrate, *Bioluminescence and Chemiluminescence: Current Status*, p 115-118.
- Bronstein I, Edwards B, Voyta JC (1989) 1,2-dioxetanes: novel chemiluminescent enzyme substrates: Applications to immunoassays, J. Biolum. Chemilum. 4:99-111.
- Bronstein I, Voyta JC, Lazzari KG, Murphy OJ, Edwards B, Kricka LJ (1990) Rapid and sensitive detection of DNA in Southern blots with chemiluminescence, *Bio Techniques* 8:310-313.
- Bronstein I, Voyta JC, Murphy OJ, Bresnick L, Kricka LJ (1992) Improved chemiluminescent western blotting procedure, *BioTechniques* 12:748-753.
- Bronstein I, Voyta JC, Thorpe GH, Kricka LJ, Armstrong G (1989) Chemiluminescent assay of alkaline phosphatase applied in an ultrasensitive enzyme immunoassay of thyrotropin, *Clin. Chem.* 35:1 441-1446.
- Bronstein I, Fortin JJ, Voyta JC, Juo RR, Edwards B, Olesen CE, Lijam N, Kricka LJ (1994) Chemiluminescent reporter gene assays: Sensitive detection of the GUS and SEAP gene products, *BioTechniques* 17:172-178.
- Bronstein I, Olesen CEM, Martin C, Schneider G, Edwards B, Sparks A, and Voyta JC (1994b) Chemiluminescent detection of DNA and protein with CDP and CDP-Star 1,2-dioxetane enzyme substrates, Bioluminescence and Chemiluminescence: Fundamentals and Applied Aspects. p 269-272.
- Edwards B, Sparks A, Voyta A, Bronstein I (1994) New chemiluminescent dioxetane enzyme substrates, *Bioluminescence and Chemiluminescence: Fundamentals and Applied Aspects*. p 56-59.
- Martin C, Bresnick L, Juo RR, Voyta JC, Bronstein I. (1991) Improved chemiluminescent DNA sequencing, *BioTechniques* 11:110-113.
- Martin CS, Butler L, Bronstein I (1995). Quantitation of PCR products with chemiluminescence, BioTechniques 18:908-913.
- Nishizono I, lida S, Suzuki N, Kawada H, Murakami H, Ashihara Y, Okada M. (1991) Rapid and sensitive chemiluminescent enzyme immunoassay for measuring tumor markers, *Clin. Chem.* 37:1639-1644.
- O'Connor KL, Culp LA. (1994) Quantitation of two histochemical markers in the same extract using chemiluminescent substrates, *BioTechniques* 17:502-509.
- Zhao S, Zhu Q, Somerville RL (2000) The σ(70) transcription factor TyrR has zinc-stimulated phosphatase activity that is inhibited by ATP and tyrosine. J Bacteriol. 182(4):1053-1061.

CDP-Star[®] and CSPD[®] Substrates

Alkaline phosphatase chemiluminescent substrates

Description

CDP-*Star*[®] and CSPD[®] chemiluminescent substrates for alkaline phosphatase enable light-based detection of alkaline phosphatase and alkaline phosphatase-labeled molecules with unparalleled sensitivity, speed and ease. Chemiluminescent detection with 1,2-dioxetane substrates for alkaline phosphatase offers a highly sensitive alternative to fluorescent, isotopic, and colorimetric detection methods.

Membrane-based Applications

Tropix[®] chemiluminescent substrates exhibit high sensitivity detection of biomolecules labeled either directly or indirectly with alkaline phosphatase (AP) in membrane-based applications such as Southern, Northern, and Western blotting. On membranes, the faster emission kinetics and higher light intensity of CDP-*Star* substrate enable exposure times up to 10-fold shorter than CSPD substrate. Exposure times range from 1 second to 15 minutes with CDP-*Star* substrate and from 2 minutes to 2 hours with CSPD substrate with standard X-ray film. CDP-*Star* is ideal when quick film exposures are desired and in applications that require long film exposures, such as single copy gene detection by Southern analysis. Note: With nitro-cellulose membrane, we recommend the use of CDP-*Star* substrate with Nitro-Block-II[™] (see page 8 for substrate formulations with Nitro-Block and Nitro-Block-II).

Solution-based Applications

CDP-Star and CSPD substrates (Figure 2) are also used in solution-based assays such as immunoassays, nucleic acid probe assays, and direct enzyme assays. They are ideally suited for the detection of placental alkaline phosphatase (PLAP) and secreted placental alkaline phosphatase (SEAP) in reporter gene assays [1, 2]. They are used widely in immunoassay platforms for analytes identification, small molecules and viral nucleic acid detection. Maximum light levels are reached at approximately 30 minutes and glow emission persists for several hours. Emission kinetics are similar for CSPD and CDP-Star substrates in solution. Additional potential applications include pasteurization monitoring, alkaline or acid phosphatase biomarker detection (stem cells, tumor markers) and sensitive quantitation of purified protein phosphatase activity [3].

Product Configuration

CDP-Star and CSPD substrates are available as concentrates or as Ready-to-Use solutions, with or without luminescence enhancers. CDP-Star and CSPD Substrate concentrates are supplied at 5 or 25 mM (respectively) in aqueous buffer. The recommended working concentration for these substrates is 0.4 mM for solution-based assays and 0.25 mM for membrane-based assays. Ready-to-Use formulations of 0.25 mM substrate for membrane-based applications are listed below. Ready-to-Use substrate/enhancer formulations for solution-based assays are listed below and described on page 7. The use of enhancers is necessary in solution assays for optimal light output and sensitivity (see page 7 for enhancers and Ready-to-Use substrate/enhancer formulations).

REFERENCES

- Bronstein I, Fortin JJ, Voyta JC, Juo RR, Edwards B, Olesen CE, Lijam N, Kricka LJ (1994) Chemiluminescent reporter gene assays: Sensitive detection of the GUS and SEAP gene products, *BioTechniques* 17:172-178.
- O'Connor KL, Culp LA. (1994) Quantitation of two histochemical markers in the same extract using chemiluminescent substrates, *BioTechniques* 17:502-509.
- Zhao S, Zhu Q, Somerville RL (2000) The σ(70) transcription factor TyrR has zinc-stimulated phosphatase activity that is inhibited by ATP and tyrosine. J Bacteriol. 182(4):1053-1061.

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
CDP- <i>Star</i> ® Chemiluminescent Substrate for Alk		on a
5 mM Concentrate	1.0 mL	T2304
5 mM Concentrate	2.0 mL	T2305
5 mM Concentrate	5.0 mL	T2306
5 mM Concentrate	10 mL	T2307
5 mM Concentrate	20 mL	T2308
5 mM Concentrate	50 mL	T2309
5 mM Concentrate	200 mL	T2310
Ready-to-Use (0.25 mM)	50 mL	T2145
Ready-to-Use (0.25 mM)	100 mL	T2146
Ready-to-Use (0.25 mM)	250 mL	T2147
Ready-to-Use (0.25 mM) with Nitro-Block-II™ Enhancer	100 mL	T2218
Ready-to-Use (0.4 mM) with Sapphire-II™ Enhancer	100 mL	T2214
Ready-to-Use (0.4 mM) with Emerald-II™ Enhancer	100 mL	T2216
Larger sizes of CDP-Star® are available upon request	-	
or Research Use Only. Not for use in diagnostic proced	lures.	

APPLIED BIOSYSTEMS ORDER INFO

CSPD® Chemiluminescent Substrate for Alkaline Phosphatase

	. noopnataoo	
25 mM Concentrate	25 mL	T2044
25 mM Concentrate	0.5 mL	T2040
25 mM Concentrate	1.0 mL	T2098
25 mM Concentrate	2.5 mL	T2041
25 mM Concentrate	5.0 mL	T2042
25 mM Concentrate	10 mL	T2043
Ready-to-Use (0.25 mM)	50 mL	T2141
Ready-to-Use (0.25 mM)	100 mL	T2142
Ready-to-Use (0.25 mM)	250 mL	T2143
Ready-to-Use (0.25 mM)) with Nitro-Block™ Enhancer	100 mL	T2217
Ready-to-Use (0.4 mM) with Sapphire-II™ Enhancer	100 mL	T2210
Ready-to-Use (0.4 mM) with Emerald-II™ Enhancer	100 mL	T2212
Larger sizes of CSPD® are available upon request		

SIZE

CAT#

Larger sizes of CSPD® are available upon request.

For Research Use Only. Not for use in diagnostic procedures.



Panel A.



Panel B.

Figure 2. Chemical Structures of CSPD® (Panel A) and CDP-Star® (Panel B) Substrates.

Substrates and Enhancers

Galacton[®], Galacton-Plus[®] and Galacton-Star[®] Substrates β-Galactosidase chemiluminescent substrates

Description

Galacton[°], Galacton-Plus[°], and Galacton-Star[®] chemiluminescent substrates (Figure 3) for β -galactosidase provide chemiluminescent detection of β -galactosidase and β -galactosidase-labeled molecules. These substrates offer significant sensitivity improvement over colorimetric and fluorescent detection methods.

Applications

Galacton,[®] Galacton-Plus,[®] and Galacton-*Star*[®] substrates have been very widely used for reporter gene assays in many different organisms, including mammals, yeast (including two hybrid assays), bacteria, protozoans, fish cells, and frog oocytes. Galacton-*Star* substrate has been employed in β -galactosidase enzyme complementation technology, which provides a novel cell-based assay capability for intracellular monitoring of protein-protein interactions, cell fusion assays, and intracellular protein translocation. Additional applications include detection of β -galactosidase conjugates in immunoassay and immunodetection formats, and detection of endogenous bacterial β -galactosidase for coliform enumeration.

Substrate Selection

Galacton substrate is a first-generation substrate that enables highly sensitive detection of β -galactosidase in reporter gene assays and direct enzyme assays. Galacton-Plus[®] chemiluminescent substrate for β -galactosidase provides prolonged emission kinetics and superior sensitivity compared to Galacton substrate. Light emission is maintained for several minutes with Galacton-Plus substrate, which is necessary when light signal is measured in a microplate luminometer without automatic injection capabilities. Figure 4 shows the kinetics of detection of purified β -galactosidase with Galacton substrate and Galacton-Plus substrate. Galacton-Star chemiluminescent substrate enables detection in a single-step reaction format. Light emission with Galacton-Star substrate typically reaches maximum in 60 minutes and exhibits glow kinetics for nearly one hour (Figure 5). Galacton-Star substrate provides the most facile detection protocol, requiring only a single reagent addition to sample, and in addition, provides the highest sensitivity detection of β -galactosidase. For new users, the Galacton-Star substrate, included in the Galacto-Star and Gal-Screen® reporter gene assay systems, is recommended for all applications. Enhancers and/or accelerators are required for optimum performance with each of these substrates, and are provided as kit components with each of the β -galactosidase reporter gene assay systems (see page 12 for Gal-Screen system; page 14 for Galacto-Star system; page 16 for Galacto-Light[™] system and Galacto-Light Plus[™] system; page 21 for Dual-Light® system).

Product Configuration

Galacton and Galacton-Plus substrates are supplied as 100X concentrates. Galacton substrate is also available as a component of the Galacto-Light system. Galacton-Plus substrate is also available in the Dual-Light and Galacto-Light Plus systems. Galacton-Star substrate is supplied as a 10 mM concentrate, and is diluted 50-fold for use. Galacton-Star substrate is also available in the Galacto-Star and Gal-Screen systems. For reporter gene assays, we recommend purchase of the complete systems, so that all necessary reagents and protocol are obtained.



Figure 3. Chemical Structures of β -Galactosidase Substrates.



Figure 4. Detection of β -Galactosidase with Galacton® and Galacton-Plus® Substrates.



Figure 5. Detection of β -Galactosidase with Galacton-Star[®] Substrate.

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Galacton® Chemiluminescent Substrate for β -O	Galactosidase	
100X Concentrate	2 mL	T2094
100X Concentrate	4 mL	T2121
100X Concentrate	10 mL	T2095
Galacton-Plus® Chemiluminescent Substrate	for β -Galactosidase	
100X Concentrate	2 mL	T2118
100X Concentrate	4 mL	T2119
100X Concentrate	10 mL	T2120
Galacton-Star® Chemiluminescent Substrate	for β -Galactosidase	
10 mM Concentrate	2 mL	T2264
10 mM Concentrate	4 mL	T2265
10 mM Concentrate	10 mL	T2266

For Research Use Only. Not for use in diagnostic procedures

Substrates and Enhancers

Glucuron[®] Substrate

 β -Glucuronidase chemiluminescent substrate

Description

Glucuron° chemiluminescent substrate enables chemiluminescent detection of β -glucuronidase (GUS). Glucuron substrate is a highly sensitive alternative to colorimetric and fluorescent detection methods. The superior sensitivity achieved with Glucuron substrate makes it the ultimate choice for GUS reporter gene assays in plant or mammalian cells.

Applications

Quantitation of glycosidic enzymes has wide-spread application, including reporter gene assays, environmental testing, biomedical research, toxicology, and pharmaceutical screening. Glucuron substrate can be employed for sensitive chemiluminescent detection in these varied/various applications, and has primarily been used for GUS reporter gene assays [1]. Measurement of GUS activity has been used in a novel bioassay with *E. coli* to evaluate toxicity of metal ions in environmental samples [4]. Additional potential applications include specific microbe detection and detection of mammalian β -glucuronidase. Mammalian β -glucuronidase activity has been measured in assays for mast cell degranulation [2] and to monitor periodontal disease [3].



Glucuron

Figure 6. Chemical Structure of Glucuron® Substrate.

REFERENCES

- Bronstein I, Fortin JJ, Voyta JC, Juo RR, Edwards B, Olesen CE, Lijam N, Kricka LJ (1994) Chemiluminescent reporter gene assays: Sensitive detection of the GUS and SEAP gene products, *BioTechniques* 17:172-178.
- Dreskin SC, Pribluda VS, Metzger H (1989) IgE receptor-mediated hydrolysis of phosphoinositides by cytoplasts from rat basophilic leukemia cells, J Immunol. 142:4407-4415.
- Lamster IB, Holmes LG, Gross KB, Oshrain RL, Cohen DW, Rose LF, Peters LM, Pope MR (1994) The relationship of beta-glucuronidase activity in crevicular fluid to clinical parameters of periodontal disease. Findings from a multicenter study. J Clin Periodontol. 21:118-127.
- Mariscal A, Garcia A, Carnero M, Gomez E, Fernandez-Crehuet J (1994) New toxicity determination method that uses fluorescent assay of Escherichia coli, *Biotechniques* 16:888-892.

Product Configuration

Glucuron substrate is supplied as a 100X concentrate. Please inquire regarding availability.

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Glucuron [®] Substrate		
Tropixcustom (GL005) 100X Concentrate	0.5 mL	Tropixcustom
Tropixcustom (GL010) 100X Concentrate	1.0 mL	Tropixcustom
Tropixcustom (GL025) 100X Concentrate	2.5 mL	Tropixcustom

For Research Use Only. Not for use in diagnostic procedures.

Glucon[™] Substrate

β -Glucosidase chemiluminescent substrate

Description

Glucon^{∞} chemiluminescent substrate enables chemiluminescent detection of β -glucosidase. Glucon, is a highly sensitive alternative to colorimetric and fluorescent detection methods.

Applications

β-Glucosidase has been used in conjunction with β-glucuronidase in a physiologic test for the rapid identification and differentiation of enterococci and streptococci [2]. Vectors containing the sequence for β-glucosidase have been developed [1]. In addition, thermostable β-glucosidase has been cloned and isolated [1] and could provide an extremely sensitive, robust reporter assay system that will eliminate background from endogenous enzyme activity. Glucon[™] substrate can also be used for detection of mammalian β-glucosidase activity.



Figure 7. Chemical Structure of Glucon™ Substrate.

REFERENCES

- Gabelsberger J, Liebl W, Schleifer K (1993) Purification and properties of recombinant β-glucosidase of the hyperthermophilic bacterium Thermotoga maritime, *Appl Microbiol Technol* 40:44-52.
- Kirby R, Ruoff K (1995) Cost-effective, clinically relevant method for rapid identification of betahemolytic streptococci and enterococci, J Clin Microbiol 33:1154-1157.

Product Configuration

Glucon substrate is supplied as a 10 mM concentrate. The suggested working concentration for solution-based assays is 0.1 mM. Please inquire regarding availability.

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Glucon™ Substrate		
Tropixcustom	10 mM concentrate	Tropixcustom

For Research Use Only. Not for use in diagnostic procedures.

Substrates and Enhancers

NA-Star[®] Substrate

Neuraminidase chemiluminescent substrate

Description

NA-Star[®] chemiluminescent substrate enables sensitive detection of neuraminidase (NA) activity. This substrate is a highly sensitive replacement for the widely used fluorogenic substrate, methylumbelliferyl N-acetylneuraminic acid (MUNANA). NA-Star substrate has been incorporated into the NA-Star Influenza Neuraminidase Inhibitor Resistance Detection Kit. The NA-Star Influenza Neuraminidase Inhibitor Resistance Detection Kit includes everything needed to quantitate neuraminidase activity and neuraminidase inhibitor resistance in avian, equine, human (types A and B), and porcine influenza viruses. The kit's fast and easy protocol and convenient 96-well plate format make it ideal for monitoring influenza virus neuraminidase inhibitor resistance, as well as high-throughput inhibitor compound screening.

Applications

NA-Star chemiluminescent substrate has applications in viral research for the detection and characterization of viral neuraminidases. NA-Star substrate has been applied to the detection of influenza virus neuraminidase activity in clinical isolates providing up to 60-fold higher sensitivity than with the fluorescence assay [1]. The chemiluminescent-based detection technology provides a wide dynamic range—greater than four orders of magnitude of neuraminidase concentration (two orders of magnitude greater than fluorescent MUNANA-based assays), enabling accurate quantitation of neuraminidase inhibitor resistance levels over a broad range of virus concentration and neuraminidase activity without having to test multiple virus dilutions. To date, this substrate has been used primarily for global screening of flu strains for neuraminidase inhibitor sensitivity. Additional applications include high throughput screening for identification of new neuraminidase inhibitor anti-viral therapeutics, detection of neuraminidase in other organisms, including bacteria.



10101

Figure 8. Chemical Structure of NA-Star® Substrate.



Figure 9. Sensitivity Comparison of Chemiluminescent Assay Detection with NA-Star® Substrate to Fluorescent Assay Detection with MUNANA Substrate. Dilutions of Influenza Type B (ATCC VR-1535) virus culture supernatant (cultured on MDCK cells) were assayed at different temperatures and Signal/Noise (S/N) ratio calculated using uninfected MDCK cell supernatant (Noise). The lower limit of detection (S/N = 2) is at least 30-fold lower with the chemiluminescent NA-Star assay, and S/N is approx. 50-fold higher with NA-Star assay. The dynamic range of detection with the NA-Star assay with virus samples is three orders of magnitude.

The fast and easy protocol enables you to perform assays in less than 1.5 hours. Simply incubate your virus samples with dilutions of neuraminidase inhibitor, add NA-*Star* chemiluminescent substrate, incubate, and then inject/add the accelerator solution, which triggers light emission from the reaction product. Light signal is measured with a luminometer, including multi-mode instruments that include a luminometer mode. For best results, use a luminometer with an automatic injector to add the accelerator solution, although a multichannel pipettor can be used if plate is read immediately.

Product Configuration

The NA-Star Influenza Neuraminidase Inhibitor Resistance Detection Kit includes NA-Star chemiluminescent substrate for neuraminidase, all necessary assay reagents, and microplates everything needed for fast, accurate quantitation of neuraminidase inhibitor resistance in influenza virus isolates.

REFERENCES

- Buxton RC, Edwards B, Juo RR, Voyta JC, Tisdale M, Bethell RC (2000) Development of a sensitive chemiluminescent neuraminidase assay for the determination of influenza virus susceptibility to Zanamivir, Anal. Biochem. 280:291-300.
- Boivin G, Goyette N (2002) Susceptibility of recent Canadian influenza A and B virus isolates to different neuraminidase inhibitors, *Antiviral Res.* 54:143-147.
- Aymard M, Ferraris O, Gerentes L, Jolly J, Kessler N (2003) Neuraminidase assays, Developmental Biology 115:75-83.
- McKimm-Breschkin J, Trivedi T, Hampson A, Hay A, Klimov A, Tashiro M, Hayden F, Zambon M (2003) Neuraminidase sequence analysis and susceptibilities of influenza virus clinical isolates to Zanamivir and Oseltamivir, *Antimicrobial Agents and Chemotherapy* 47(7):2264-2272.
- Mungall BA, Xu X, Klimov A (2003) Assaying susceptibility of avian and other influenza A viruses to Zanamivir: Comparisons of fluorescent and chemiluminescent neuraminidase assays, Avian Diseases 47:1141-1144.
- Wetherall NT, Trivedi T, Zeller J, Hodges-Savola C, McKimm-Breschkin JL, Zambon M, Hayden FG (2003) Evaluation of neuraminidase enzyme assays using different substrates to measure susceptibility of influenza virus clinical isolates to neuraminidase inhibitors: Report of the Neuraminidase Inhibitor Susceptibility Network, J. Clin. Microbiol. 41(2):742-750.
- Mungall BA, Xu X, Klimov A (2004) Surveillance of influenza isolates for susceptibility to neuraminidase inhibitors during the 2000-2002 influenza seasons, Virus Research 103:195-197.
- Hurt AC, McKimm-Breschkin JL, McDonald M, Barr IG, Komadina N, Hampson AW (2004) Identification of a human influenza type B strain with reduced sensitivity to neuraminidase inhibitor drugs, Virus Research 103:205-211.
- Jackson D, Barclay W, Zurcher T (2005) Characterization of recombinant influenza B viruses with key neuraminidase inhibitor resistance mutations, J. Antimicrobial Chemotherapy 55(2):162-169.
- Neuraminidase Inhibitor Susceptibility Network (2005) Evolution of HSNI avian influenza viruses in Asia, Emerging Infectious Diseases 11(10):1515-1521.

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
NA- <i>Star</i> ® Influenza Neuraminidase	960 assay wells +	
Inhibitor Resistance Detection Kit	10 x 96-well microplates	4374422
NA- <i>Star</i> ® Influenza Neuraminidase		
Inhibitor Resistance Detection Reagent Set	960 assay wells	4374348
NA-Star® Detection Microplates	10 x 96-well microplates	4374349
NA- <i>Star</i> ® Influenza Neuraminidase		
Inhibitor Resistance Detection Kit Protocol	1 protocol	4375714

Sapphire[™], Sapphire-II[™], Emerald[™], Emerald-II[™] and Ruby[™] Enhancers Luminescence enhancers for solution-based assays

Chemiluminescent Enhancement

Applied Biosystems has developed macromolecular signal-enhancement technology for use with 1,2-dioxetane enzyme substrates in solutionbased assays. Aqueous environments reduce the chemiluminescent signal intensity by water-induced quenching. Addition of Tropix[®] luminescence enhancers increases the emission efficiency of light production by partitioning the water away from the site of signal production. Tropix enhancers such as Sapphire[™], Emerald[™], Ruby[™], Sapphire-II[™], and Emerald-II[™] enhancers are essential components of solution-based assays. Tropix enhancers provide signal enhancement with minimal delay of light-emission kinetics. For example, adding Sapphire enhancer to chemiluminescent substrate in 0.1 M diethanolamine (DEA) (pH 10) increases the chemiluminescent half-time to plateau (T1/2) by only 15-20% [1].

Choosing an Enhancer

Sapphire[™], Emerald[™], Ruby[™], Sapphire-II[™], or Emerald-II[™] Enhancers Tropix enhancers shift the wavelength of light emission. Sapphire and Sapphire-II enhancers slightly shift the light emission maximum from 475 nm observed for the dioxetane alone to 461 nm. Emerald and Emerald-II enhancers shift the light emission maximum to 542 nm, while Ruby[™] enhancer shifts the light emission maximum to 620 nm. Although the signal intensity obtained with Emerald and Emerald-II enhancers is much greater compared to other enhancers, use of Sapphire and Sapphire-II enhancers produces a wider dynamic range since photodetector saturation is less likely to occur. The Emerald or Emerald-II enhancer is the optimum choice for applications requiring maximum signal intensity. Applied Biosystems offers an enhancer variety pack to help match the best enhancer for your application and instrumentation.

Applications

Tropix enhancers are essential components in immunoassays, nucleic acid hybridization assays performed in a microplate or similar format, reporter gene assays, and enzyme assays because of improved signal-to-noise performance. Enhancers are included in the ELISA-Light[™] system (page 32) and reporter gene assay system reagents (pages 10-24).



Enhancer Variety Pack

Sample Sizes of Sapphire[™], Emerald , Sapphire-II[™] and Emerald-II[™] Enhancers

Optimizing signal intensity with the appropriate enhancer can be easily accomplished with a selection of enhancers. The Enhancer Variety Pack contains a sample of Sapphire, Sapphire-II, Emerald, and Emerald-II enhancers for solution-based applications in convenient 5 mL sizes.

Ready-to-Use Substrate/Enhancer Formulations

CSPD[®] Ready-to-Use Formulations (with Sapphire-II[™] or Emerald-II[™] Enhancers), CDP-Star[®] Ready-to-Use Formulations (with Sapphire-II[™] or Emerald-II[™] Enhancers)

Ready-to-Use substrate/enhancer formulations are supplied at 0.4 mM substrate and 1X enhancer in a stable buffer system, eliminating the need for substrate dilution buffer preparation and substrate and enhancer dilution. Ready-to-Use formulations simplify procedures for immunoassays, nucleic acid hybridization assays in microplate formats, enzyme assays, and other solution-based assays.

Enhancer	Maximum Light Emission	S/N		
Sapphire-II [™] Enhancer	461 nm	HIGH		
Sapphire [™] Enhancer	461 nm	MODERATE		
Comments: Provides optimum signal-to-noise ratio (S/N) and best dynamic range for highly sensitive luminometer measurement				
Emerald-II™ Enhancer	542 nm	HIGH		
Emerald [™] Enhancer	542 nm	MODERATE		
Comments: Provides highest signal intensity for assays performed in black microplates; for camera imaging systems requiring high intensity signal				
Ruby™ Enhancer	620 nm	MODERATE		
Comments: Provides wavelength shift that may be compatible with some solid state detectors				

Figure 11. Enhancer Comparison for Solution-based Assays

REFERENCES

 Erve V, Voyta JC, Edwards B, Kricka LJ, Bronstein I (1993) Influence of reaction conditions on the chemiluminescent dephosphorylation of AMPPD, *Bioluminescence and Chemiluminescence:* Status Report p 306-311.

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#		
Sapphire™ 10X Concentrate	25 mL	T2011		
Sapphire-II™ 10X Concentrate	25 mL	T2113		
Emerald™ 10X Concentrate	25 mL	T2002		
Emerald-II™ 10X Concentrate	25 mL	T2115		
Ruby™ 10X Concentrate	25 mL	Tropixcustom		
Enhancer Variety Pack	4 x 5 mL	Tropixcustom		
Ready-to-Use Substrate/Enhancer Formulations				
CSPD [®] Ready-to-Use with Sapphire-II™	100 mL	T2210		
CSPD [®] Ready-to-Use with Emerald-II™	100 mL	T2212		
CDP- <i>Star</i> [®] Ready-to-Use with Sapphire-II™	100 mL	T2214		
CDP- <i>Star</i> ® Ready-to-Use with Emerald-II™	100 mL	T2216		
For Research Use Only. Not for use in diagnostic procedures.				

Figure 10. Chemiluminescence Enhancement.

Nitro-Block[™] and Nitro-Block-II[™]

Luminescence enhancers for membrane-based assays

Chemiluminescent Enhancement on Membranes

Membranes commonly used in nucleic acid and protein detection applications include nylon, nitrocellulose, and PVDF (polyvinylidene fluoride). The signal intensity generated from 1,2-dioxetane substrates varies on each of the membranes.

Nylon Membranes

Nylon membranes are essentially "self-enhancing" since the surface contains hydrophobic micro-domains compatible with the chemiluminescence process. The use of enhancer on these membranes is not recommended.

Nitrocellulose and PVDF Membranes

Nitrocellulose membranes provide an inefficient environment for chemiluminescence, resulting in very low signal intensity. Tropix® Nitro-Block™ and Nitro-Block-II™ membrane enhancers increase signal intensity on both nitrocellulose and PVDF membranes [1]. Nitro-Block and Nitro-Block-II enhancers generate a hydrophobic environment on the membrane surface that increases the intensity of chemiluminescence. The effects of Nitro-Block enhancer treatment on nitrocellulose and PVDF membranes are shown (Figure 12). Without Nitro-Block enhancer, the signal on nitrocellulose membranes is weak, making detection either impossible or requiring extremely long exposure times. With Nitro-Block enhancer, short exposures of 10 to 45 minutes are achieved. While PVDF membranes do not require the use of Nitro-Block enhancer, exposure times are tenfold faster with Nitro-Block enhancer. This permits very short exposure times on PVDF, ranging from 15 seconds to 15 minutes for Western blots.

Nitro-Block[™] Enhancer Compared to Nitro-Block-II[™] Enhancer

Nitro-Block-II membrane enhancer is a modified version of Nitro-Block membrane enhancer that has been developed to further optimize the performance of chemiluminescent 1,2-dioxetane substrates. Nitro-Block or Nitro-Block-II enhancers are required for use with CSPD* substrate on nitrocellulose membranes. Nitro-Block-II Enhancer is required for use with CDP-*Star** substrate on nitrocellulose membranes. The use of Nitro-Block-II Enhancer is not necessary with CDP-*Star* substrate on PVDF membranes.

Nitrocellulose



Figure 12. Chemiluminescent Signal Enhancement on Nitrocellulose and PVDF Membranes with Nitro-Block[™] Enhancer. Xray images of Western blots detected using CSPD[®] substrate.

Ready-to-Use Substrate/Enhancer Formulations

CSPD* Ready-to-Use Substrate with Nitro-Block**, CDP-Star* Ready-to-Use Substrate with Nitro-Block-II**

Ready-to-Use formulations for membrane-based applications are supplied with enhancer pre-mixed with substrate for convenience and ease of use. The shelf-life for all the Ready-to-Use formulations is one year at 4°C. Ready-to-Use formulations with Nitro-Block and Nitro-Block-II Enhancer are intended for use with nitrocellulose membranes to increase light intensity. Substrate/enhancer formulations for membrane-based applications are supplied at 0.25 mM substrate and 1X enhancer concentrations. Note: prior to using Ready-to-Use formulations, a high pH wash of the membrane with assay buffer is recommended.

SIZE	CAT#
20 mL	T2026
20 mL	T2184
100 mL	T2217
100 mL	T2218
	20 mL 100 mL

For Research Use Only. Not for use in diagnostic procedures.

REFERENCES

 Bronstein I, Fortin J, Voyta JC, Kricka LJ (1992) Nitro-Block enhancement of AMPPD chemiluminescent signal in the detection of DNA, *Biotechniques* 12: 500-502.

Substrate	Enhancer	Performance
Nitrocellulose		
CSPD [®] Substrate	Nitro-Block [™] or Nitro-Block-II [™]	
CDP-Star® Substrate	Nitro-Block-II™	Optimal
PVDF		
CSPD [®] Substrate	Nitro-Block [™] or Nitro-Block-II [™]	
CDP-Star® Substrate	—	Optimal
Nylon		
CSPD [®] or CDP-Star [®]	_	Optimal

Figure 13. Enhancer Comparison for Membrane-based Assays.

Reporter Gene Assays and Reagents

Introduction
Gal-Screen® Assay System 12
Galacto-Star™ Assay System14
Galacto-Light™, Galacto-Light Plus™ Assay Systems
Luc-Screen® Assay System
Luciferase Assay System 20
Dual-Light® Assay System
Phospha-Light™ Assay System

Introduction

Reporter Gene Assays

Reporter gene assays are invaluable for studying regulation of gene expression, both by cis-acting factors (gene regulatory elements) or trans-acting factors (transcription factors or exogenous regulators). In these assays, the reporter gene acts as a surrogate for the coding region of the gene under study. The reporter gene construct contains one or more gene regulatory elements being analyzed, the structural sequence of the reporter gene, and the sequences required for the formation of functional mRNA. Upon introduction of the reporter construct into cells, expression levels of the reporter gene are monitored through a direct assay of the reporter protein's enzymatic activity. The sensitivity of each reporter gene assay is a function of several factors including detection method, reporter mRNA and protein turnover, and endogenous (background) levels of the reporter activity. Both protein turnover and levels of endogenous background vary with each reporter protein and the cell line used. Commonly used detection techniques utilize isotopic, colorimetric, fluorometric, or luminescent enzyme substrates and immunoassay-based procedures with isotopic, colorimetric, or chemiluminescent end points.

Common Reporter Genes

Below is a list of the most common reporter genes, detection methods for the reporter protein, and corresponding detection limits. β -Galactosidase and luciferase are among the most widely used reporter genes to date. β -Galactosidase is often used in conjunction with other reporter genes to normalize transfection efficiency. β -Galactosidase is traditionally detected with the colorimetric substrate o-nitrophenyl β -D-galactopyranoside (ONPG) [1]. As indicated in table, this colorimetric assay is less insensitive compared to many other reporter gene assays. With Tropix[®] 1,2-dioxetane chemiluminescent substrates for β -galactosidase, the sensitivity is increased dramatically [3,4]. The assay for chloramphenicol acetyl transferase (CAT) exhibits only moderate sensitivity, suffers from a narrow dynamic range, and usually incorporates radioisotopes. This reporter gene assay method is now only infrequently used. Assays with limited dynamic range, such as CAT, require testing of several sample dilutions to verify that sample values are within the linear range.

Secreted placental alkaline phosphatase (SEAP) is secreted by cells directly into the culture media, and can be assayed simply by taking samples of cell culture media. Secreted reporter proteins enable non-destructive assay of cell culture medium, preserving cells for additional assays and enabling time-course monitoring of gene expression. SEAP is detected with both colorimetric and chemiluminescent substrates. Phospha-Light[™] Assay System, the chemiluminescent SEAP reporter gene assay, exhibits remarkable sensitivity and ease of use. Human growth hormone (hGH), another secreted reporter protein, is detected using a radioimmunoassay (RIA) procedure, has disadvantages associated with the use of radioisotopes, and exhibits only a moderate detection limit.

Luciferase has become increasingly popular as a reporter gene, especially for co-transfection experiments where it is important to normalize transfection efficiency. As with the Tropix chemiluminescent reporter gene assays, the bioluminescent luciferase assay offers high sensitivity and a simple assay procedure. The high level of sensitivity attained with this assay is partly due to the lack of luciferase activity in most cell types (see pages 18-22 for luciferase assay systems). Dual-Light[®] Assay System, a combined reporter gene assay system for the sequential detection of luciferase and β -galactosidase, enables the user to perform both measurements from a single aliquot of cell extract in the same reaction well or tube (see pages 21 and 22 for Dual-Light kit), minimizing experimental error.

Reporter Gene	Detection Method	Detection Limit	Advantages	Disadvantages
Chloramphenicol Acetyl Transferase (CAT)	Isotopic ELISA	5x10 ⁷ molecules 1x10 ⁹ molecules	Widely Used No Radioactivity	Radioactive High Cost Low Dynamic Range Labor Intensive Low Dynamic Range High Cost/Assay Labor Intensive
β -Galactosidase	ONPG (Color) MUG (Fluorescence) Galscreen [®] System Galacto-Star™ System Galacto-Light Plus [™] System and Dual-Light [®] System (Luminescence)	3x10 ⁸ molecules 6x10 ⁵ molecules 3x10 ³ molecules	Widely Used High Sensitivity Wide Dynamic Range Simplicity	Poor Sensitivity Autofluorescence
Human Growth Hormone	Radioimmunoassay	3x10 ⁸ molecules	Secreted into Media	Radioactivity High Cost/Assay Low Sensitivity
Luciferase	Luc-Screen [®] System Tropix Luciferase Assay Kit Dual-Light [®] System (Luminescence)	10 ³ -10 ⁴ molecules	Assay Simplicity High Sensitivity Wide Dynamic Range	Protein Instability
β-Glucuronidase	MUG (Fluorescence) Glucuron® Substrate (Luminescence)	2x10 ^s molecules 5x10 ^s molecules	Autofluorescence High Sensitivity Simplicity Wide Dynamic Range	Protein Quenching
Secreted Placental Alkaline Phosphatase	pNPP (Color) Phospha-Light™ System (Luminescence)	1x10 ⁸ molecules 3x10 ⁴ molecules	Secreted into Media Secreted into Media High Sensitivity Wide Dynamic Range	Poor Sensitivity

Figure 14. Comparison of Reporter Gene Assays

Introduction, continued

Bacterial β -glucuronidase (GUS) is one of the most widely used reporter genes in plant genetic research. It is also used to a lesser extent in mammalian cells. In plant cells, β -glucuronidase activity is absent or present at very low levels. Although GUS is present in mammalian cells, its pH profile is dramatically different compared to the transfected bacterial form, enabling discrimination between the two [10]. GUS is commonly measured in extracts using the fluorescent substrate MUG. Glucuron[®] chemiluminescent substrate (as a component of the discontinued GUS-Light[™] reporter gene assay system) has been used for highly sensitive detection in GUS reporter gene assays in plant [6,7,8], mammalian [4,12], yeast [11] and arthropod cells [9]. For information on running GUS reporter gene assays with Glucuron substrate, please contact Applied Biosystems Technical Support.

GENERAL REFERENCES

- Alam J, Cook JL (1990) Reporter genes: application to the study of mammalian gene transcription, Anal Biochem 188:245-254.
- 2 Bronstein I, Martin CS, Fortin JJ, Olesen CE, Voyta JC (1996) Chemiluminescence: Sensitive detection technology for reporter gene assays, *Clin Chem* 42(9):1542-1546.
- Bronstein I, Fortin J, Stanley PE, Stewart GS, Kricka LJ (1994) Chemiluminescent and bioluminescent reporter gene assays, Anal Biochem 219(2):169-181.

B-Glucoronidase (Glucuron® substrate)

 Alonso S, Sola I, Teifke JP, Reimann I, Izeta A, Balasch M, Plana-Duran J, Moormann RJ, Enjuanes L (2002) In vitro and in vivo expression of foreign genes by transmissible gastroenteritis coronavirus-derived minigenomes, J Gen Virol 83:567-579.

- Bronstein I, Fortin JJ, Voyta JC, Juo RR, Edwards B, Olesen CE, Lijam N, Kricka LJ (1994a) Chemiluminescent reporter gene assays: Sensitive detection of the GUS and SEAP gene products, *BioTechniques* 17:172-178.
- Brouwer C, Bruce W, Maddock S, Avramova Z, Bowen B (2002) Suppression of transgene silencing by matrix attachment regions in maize: a dual role for the maize 5' ADH1 matrix attachment region, The Plant Cell 14:2251-2264.
- Denekamp M, Smeekens SC (2003) Integration of wounding and osmotic stress signals determines the expression of the AtMYB102 transcription factor gene, *Plant Physiology* 132:1415-1423.
- Diaz I, Vicente-Carbajosa J, Abraham Z, Martinez M, Isabel-La Moneda I, Carbonero P(2002) The GAMYB protein from barley interacts with the DOF transcription factor BPBF and activates endosperm-specific genes during seed development, The Plant Journal 29(4):453464.
- Kimmick MW, Afanasiev BN, Beaty BJ, Carlson JO (1998) Gene expression and regulation from the p7 promoter of Aedes densonucleosis virus, J Virol 72(5):4364-4370.
- Martin, T, Wohner, RV, Hummel S, Willmitzer L, Frommer WB (1992b) The GUS reporter system as a tool to study plant gene expression, GUS Protocols: Using GUS Gene as a Reporter of Gene Expression, San Diego: Academic Press, p 23-43.
- Noueiry AO, Chen J, Ahlquist P (2000) A mutant allele of essential, general translation initiation factor DED1 selectively inhibits translation of a viral mRNA, *Proc Natl Acad Sci USA* 97(24):12985-1990.
- Spicher A, Guicherit OM, Duret L, Aslanian A, Sanjines EM, Denko NC, Giaccia AJ, Blau HM (1998) Highly conserved RNA sequences that are sensors of environmental stress, *Mol Cell Biol* 18(12):7371-7382.

Reporter Assay System	Reporter Enzyme(s)	Key Attributes	Substrate	Duration of Light Emission	Reagent Injection	Special Features
Gal-Screen® see page 12	β-Galactosidase	Homogeneous assay does not require removal of culture medium.	Galacton-Star®	60-90 min (+)	Not required	Designed for assays on cells cultured in luminometer plates for high-throughput screening. Choice of reaction buffers for lysis of mammalian or yeast and mammalian cells.
Galacto-Star™ see page 14	β-Galactosidase	Single step addition of sub- strate and enhancer	Galacton–Star®	1 hr (+)	Not required	Useful for high throughput applications. Choice of lysis buffer for yeast and mammalian model systems. Recommended for all new users.
Galacto-Light Plus™ see page 16	β-Galactosidase	Extended light emission and better sensitivity than Galacto-Light™	Galacton-Plus®	30-60 min	Recommended	Established system with high sensitivity.
Galacto-Light™ see page 16	β-Galactosidase	First chemiluminescent system for detection of β-galactosidase	Galacton®	Several minutes	Required	Established system with good sensitivity.
Luc-Screen® see page 18	Firefly Luciferase	Extended-glow homogeneous assay does not require removal of culture medium	Luciferin	90 min (+)	Not required	Designed for assays on cells cultured in luminometer plates for high-throughput screening.
Luciferase Assay Kit see page 20	Firefly Luciferase	Enhanced signal	Luciferin	Several minutes	Required	Lysis buffer compatible with other Tropix [®] reporter systems.
Dual-Light® see page 21	Firefly Luciferase, β-Galactosidase	Detection of luciferase and β-galactosidase from same cell extract	Luciferin Galacton-Plus®	5 min/30 min	Required	Two assays from one cell extract enables better preci- sion for normalizing transfection efficiency.
Phospha-Light™ see page 23	Secreted placental alkaline phosphatase	No cell lysis required	CSPD®	1-2 hr (+)	Not required	Cells remain viable; useful for stable transfectants, time course studies, etc.

Figure 15. Selection Guide for Reporter Gene Assay Systems. Applied Biosystems offers several highly sensitive reporter gene assay systems. This chart describes the differences and relative merits of each system to help you choose the most appropriate system.

Gal-Screen[®] Assay System

Homogeneous β -galactosidase reporter gene assay system for mammalian or yeast cells

Description

The Gal-Screen[®] assay system combines direct cell lysis with rapid ultra-sensitive chemiluminescent detection of β -galactosidase reporter enzyme. This homogeneous assay is ideally suited for screening applications where assay automation is required. Gal-Screen system uses Galacton-*Star*^{*} chemiluminescent substrate and Sapphire-II[®] luminescence enhancer. A single reagent, providing cell lysis and chemiluminescent enzyme substrate, is added to cells in the presence of culture medium with or without phenol red. Light emission reaches maximum in 60-90 minutes and remains constant for 45–90 minutes.

The Gal-Screen system utilizes a simple protocol that can be used with either mammalian or yeast cells (Figure 18). A single assay reagent, which provides cell lysis and contains all reaction components, is prepared by adding Galacton-*Star* substrate to Reaction Buffer A or B. The assay reagent is added to an equal volume of cells in culture medium in either 96-, 384- or 1536-well microplates. After an incubation period, light emission is measured in a luminometer.

Advantages

The Gal-Screen assay provides greater sensitivity than colorimetric or fluorescent assays. With a lower limit of detection of 1 picogram in the presence of culture medium (Figure 16), this assay provides excellent sensitivity compared to other reporter systems. Colorimetric reporter gene assays cannot rival the dynamic range of the Galacton-*Star** substrate in the Gal-Screen system. The wide dynamic range spans five orders of magnitude, from picogram to nanogram levels, enabling detection of a wide range of reporter enzyme concentration in cells. The assay protocol was developed particularly for use with automation. The Gal-Screen assay protocol is adaptable for use in 96-, 384- or 1536-well microplate formats (Figure 19), with either mammalian or yeast cells.

Applications

The Gal-Screen assay system is widely used for traditional reporter gene assays in transiently and stably-transfected mammalian cells [7, 11], including assays for studying viral infectivity and function [8, 12, 15]. It is widely used for reporter gene assays in yeast cells [3, 9, 10], including quantitative yeast two hybrid analysis [13, 14, 16]. In addition, it has been used for reporter gene assays in fish cells [6] and bacterial cells [5]. Gal-Screen assays provide highly sensitive detection for β -galactosidase complementation assays used for intracellular monitoring of protein-protein interactions [4], protein translocation [17], and receptor dimerization/activation [1, 2], including for high throughput compound screening for receptor activation [18].

Product Configuration

The Gal-Screen β -galactosidase reporter gene assay system is formatted with two alternative Reaction Buffers for lysis of mammalian or yeast cells. Reaction Buffer A is for use with mammalian cells, while Reaction Buffer B is for yeast or mammalian cells. In comparison on mammalian cells (Figure 18), Reaction Buffer A provides a faster time to peak light emission with a shorter duration of peak signal, while Reaction Buffer B requires a longer time to peak and provides a longer duration of glow light emission.



Figure 16. Sensitivity of Gal-Screen® Assay Detection. Gal-Screen® assays were performed with the indicated amount of purified β -Galactosidase diluted in PBS. Assays were measured on a microplate luminometer.



Seed cells into tissue culture (TC) - treated luminometer plates and culture or treat as desired.

Add Reaction Buffer.

Incubate 60-90 minutes.



Place plate in luminometer and measure $\beta\mbox{-}galactosidase$ light emission.

Figure 17. $\beta\mbox{-}Galactosidase$ Reporter Gene Assay with Gal-Screen® System.

Gal-Screen® Assay System, continued

Homogeneous β -galactosidase reporter gene assay system for mammalian or yeast cells

The system is available in three sizes:

Gal-Screen® System Standard Size

T1029 – Reaction Buffer A

- T1032 Reaction Buffer B
- 0.8 mL Galacton-Star® Substrate
- \bullet 19.2 mL Reaction Buffer A (for mammalian cells) or 19.2 mL Reaction Buffer B (for yeast or mammalian cells)

Gal-Screen[®] System Large Size

T1027 – Reaction Buffer A

- T1030 Reaction Buffer B
- 4 mL Galacton-Star® Substrate
- 96 mL Reaction Buffer A (for mammalian cells) or 96 mL Reaction Buffer B (for yeast or mammalian cells)
- Gal-Screen® System Screening Size

T1028 – Reaction Buffer A

- T1031 Reaction Buffer B
- 42 mL Galacton-Star® Substrate
- 1 L Reaction Buffer A (for mammalian cells)
- or 1 L Reaction Buffer B (for yeast or mammalian cells)

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Gal-Screen® System for Mammalian Cells	200 assays	T1029
Gal-Screen® System for Mammalian Cells	1,000 assays	T1027
Gal-Screen® System for Mammalian Cells	10,000 assays	T1028
Gal-Screen® System for Yeast or Mammalian Cells	200 assays	T1032
Gal-Screen® System for Yeast or Mammalian Cells	1,000 assays	T1030
Gal-Screen® System for Yeast		
or Mammalian Cells	10,000 assays	T1031

For Research Use Only. Not for use in diagnostic procedures.

β-Galactosidase (Gal-Screen[®] Reporter Gene Assay System; for a complete reference list, please see the Applied Biosystems website.)

- Blakely BT, Rossi FM, Tillotson B, Palmer M, Estelles A, Blau HM (2000) Epidermal growth factor receptor dimerization monitored in live cells, *Nature Biotechnology* 18:218-222.
- Buensuceso C, de Virgilio M, Shattil SJ (2003) Detection of Integrin allbb3 clustering in living cells, J Biol Chem 278(17):15217-15224.
- Carrigan PE, Riggs DL, Chinkers M, Smith DF (2005) Functional comparison of human and Drosophila Hop reveals novel role in steroid receptor maturation, J Biol Chem 280(10):8906-8911.
- Higuchi T, Orita T, Katsuya K, Yamasaki Y, Akiyama K, Li H, Yamamoto T, Saito Y, Nakamura M (2004) MUC20 suppresses the Hepatocyte Growth Factor-induced Grb2-Ras pathway by binding to a multifunctional docking site of Met., *Mol Cell Biol* 24(17):7456-7468.
- Kim CC, Falkow S (2004) Delineation of upstream signaling events in the Salmonella Pathogenicity Island 2 transcriptional activation pathway, J Bacteriol 186(14):4694-4704.
- Lopez A, Fernandez-Alonso M, Rocha A, Estepa A Coll JM (2001) Transfection of epithelioma paulosum cyprini (EPC) carp cells, *Biotechnology Letters* 23:481-487.
- Ludewig AH, Kober-Eisermann C, Weitzel C, Bethke A, Neubert K, Gerisch B, Hutter H, Antebi A (2004) A novel nuclear receptor/coregulator complex controls C. elegans lipid metabolism, larval development, and aging, Genes & Devel 18:2120-2133.
- Mori T, O'Keefe BR, Sowder RC 2nd, Bringans S, Gardella R, Berg S, Cochran P, Turpin JA, Buckheit RW Jr, McMahon JB, Boyd MR (2005) Isolation and characterization of Griffithsin, a novel HV-inactivating protein, from the red alga Griffithsia sp., J Biol Chem 280(10):9345-9353.
- Natarajan L, Witwer NE, Eisenmann DM (2001) The divergent Caenorhabditis elegans β-Catenin proteins BAR-1, WRM-1 and HMP-2 make distinct protein interactions but retain functional redundancy in vivo, *Genetics* 159:159-172.



Figure 18. Kinetics of Gal-Screen® Assay. Gal-Screen® assays were performed with mammalian or yeast cells, with Reaction Buffer A and B. Assays were measured on a microplate luminometer. Data is presented as the percentage of the maximum signal.



Figure 19. Gal-Screen® Assay Detection of ψ 2BAG \propto Cells. Gal-Screen® assays were performed on ψ 2BAG \propto cells, which constitutively express β -Galactosidase. Assays were performed in several microplate formats, in media either with or without phenol red. Measurements were performed on a CCD microplate luminometer.

- Noueiry AO, Chen J, Ahlquist P (2000) A mutant allele of essential, general translation initiation factor DED1 selectively inhibits translation of a viral mRNA, *Proc Natl Acad Sci USA* 97(24):12985-1990.
- Palli SR, Kapitskaya MZ, Kumar MB, Cress DE (2003) Improved ecdysone receptor-based inducible gene regulation system, Eur J Biochem 270:1308-1315.
- 12. Parikh UM, Koontz DL, Chu CK, Schinazi RF, Mellors JW (2005) In vitro activity of structurally diverse nucleoside analogs against Human Immunodeficiency Virus type 1 with the K65R mutation in reverse transcriptase, Antimicrobial Agents and Chemotherapy 49(3):1139-1144.
- Saito K, Tomigahara Y, Ohe N, Isobe N, Nakatsuka I, Kaneko H (2000) Lack of significant estrogenic or antiestrogenic activity of pyrethroid insecticides in three in vitro assays based on classic Estrogen Receptor a-mediated mechanisms, *Toxicological Sciences* 57:54-60.
- Shields CM, Taylor R, Nazarenus T, Cheatle J, Hou A, Tapprich A, Haifley A, Atkin AL (2003) Saccharomyces cerevisiae Ats1p interacts with Nap1p, a cytoplasmic protein the controls bud morphogenesis, *Curr Genet* 44:184-194.
- 15. Varadarajan R, Sharma D, Chakraborty K, Patel M, Citron M, Sinha P, Yadav R, Rashid U, Kennedy S, Eckert D et al (2005) Characterization of gp120 and its single-chain derivatives, gp120-024012 and gp120-M9: Implications for targeting the CD4i epitope, Human Immuno Virus vaccine design. J Virology **79(3)**:1713-1723.
- Ward BM, Weisberg AS, Moss B (2003) Mapping and functional analysis of interaction sites within the cytoplasmic domains of the Vaccinia Virus A33R and A36R envelope proteins, *J Virol* 77(7):4113-4126.
- Wehrman TS, Casipit CL, Gewertz NM, Blau HM (2005) Enzymatic detection of protein translocation, Nature Methods 2(7):521-527.
- Yan YX, Boldt-Houle DM, Tillotson BP, Gee MA, D'Eon BJ, Chang XJ, Olesen CE, Palmer MA (2002) Cell-based high-throughput screening assay system for monitoring G protein-coupled receptor activation using beta-galactosidase enzyme complementation technology, *J Biomol Screen* 7(5):451-9.

Galacto-Star[™] Assay System One-step β–galactosidase reporter gene assay

Description

Galacto-*Star*^{**} Assay System is a chemiluminescent reporter assay system designed for the rapid, and sensitive detection of β -galactosidase in cell lysates. The chemiluminescent assay exhibits over three orders of magnitude greater sensitivity compared to colorimetric β -galactosidase assays.

Galacto-Star assay system includes Galacton-Star chemiluminescent β -galactosidase substrate. Automatic reagent injection is not required. Cell lysate is mixedwithReactionBuffercontainingGalacton-StarsubstrateandSapphire-II[™] enhancer.LightemissionfromGalacto-Starreachesmaximumin60-90minutes and remains constant for at least 1 hour (Figure 20). After incubation at room temperature, the signal is measured in a luminometer.

Advantages

Galacto-*Star* chemiluminescent assay for β -galactosidase is among the most sensitive reporter gene assays available (Figure 21). As few as 10 femtograms of β -galactosidase (20,000 molecules) are detectable [9]. High sensitivity makes this system ideal for detection of weak expression and for transfection normalization with other sensitive reporter gene assays. Colorimetric, fluorometric, and isotopic assays cannot rival the dynamic range of Galacto-*Star* Assay System, which spans femtogram to nanogram levels of protein. Unlike Galacto-Light[™] and Galacto-Light Plus[™] systems (see page 16) which require the sequential addition of substrate and enhancer, Galacto-*Star* Assay System is a one-step assay procedure, following lysate preparation, simplifying the detection procedure. The chemiluminescent assay is compatible with lysis buffers used with luciferase assays, making it ideal for transfection normalization.

Applications

The Galacto-Star system is formatted with a choice of lysis buffers for use in mammalian or yeast cells. The Galacto-Star assay system is used widely for traditional reporter gene assays in transfected mammalian cells [3, 4, 11], and in insect cells [16]. A wide variety of applications have been performed, including viral function assays with β -gal-encoding pseudovirions [17] and MAGI cells [6, 8], normalization of siRNA transfection [7], and as a reporter read-out for epitope recognition by an engineered CTL hybridoma cell line [15]. The Galacto-Star system has been used to assay tissue extracts of transgenic mice made with β -gal-tagged mouse embryonic stem cells [5, 12]. The system is also formatted for use with yeast cells, and is ideally suited for reporter gene assays in yeast [10], or the study of protein:protein interactions with the yeast two-hybrid system [2, 13]. Galacton-Star substrate has been used for reporter gene assays in bacterial cells with modified lysis reagents [18]. Two novel applications demonstrated have been a cell death assay, by measurement of β -gal reporter enzyme released into culture media [14], and a stop codon read-through assay using a constitutively-expressed β -gal-luciferase fusion construct [1].

The Galacto-Star assay system has wide application to assays that use β -gal reporter as a read-out for gene expression in many cell types and tissues from whole animals, or as a functional read-out for viral function, immune cell activation, cell death, and mRNA processing.



Figure 20. Detection of β -Galactosidase with Galacto-Star[™] Assay System.



Figure 21. Sensitivity of Galacto-Star™ Assay System.



luminometer plates and culture or treat as desired

Place plate in luminometer and measure β-galactosidase light emission

Figure 22. β-Galactosidase Reporter Gene Assay with Galacto-Star™ System.

measure β-galactosidase light emiss

SIZE

CAT#

Galacto-Star™ System Assay System, continued

One-Step β –galactosidase reporter gene assay

APPLIED BIOSYSTEMS ORDER INFO

For Research Use Only. Not for use in diagnostic procedures.

Product Configuration

Galacto-Star [™] System Standard Size	Galacto- <i>Star</i> ™ Mammalian Cell Reporter Gene Assay System	standard size	T1012
T1012 – Mammalian Cells T1019 – Yeast Cells Capacity: 200 single tube assays	Galacto- <i>Stai</i> ™ Mammalian Cell Reporter Gene Assay System	large size	T1014
Capacity: 600 assays with microplate format Contents:	Galacto- <i>Stai</i> ™ Mammalian Cell Reporter Gene Assay System	screening size	T1013
1.2 mL Galacton-Star ^a substrate 60 mL Reaction Buffer Diluent with Sapohire/II" enhancer	Galacto- <i>Star</i> ™ Yeast Cell Reporter Gene Assay System	standard size	T1019
70 mL Mammalian Lysis Solution or 80 mL 5X Z Yeast Lysis Buffer	Galacto- <i>Star</i> ™ Yeast Cell Reporter Gene Assay System	large size	T1021
Galacto- <i>Star</i> [™] System Large Size	Galacto- <i>Star</i> ™ Yeast Cell Reporter Gene Assay System	screening size	T1020
T1014 - Mammalian Cells	Galacto- <i>Star</i> ™ Reaction Buffer Diluent with Galacton- <i>Star</i> ® Chemiluminescent Substrate	180 mL 3.6 mL	T1056

T1014 – Marmanan Cells T1021 – Yeast Cells Capacity: 600 single tube assays Capacity: 1,800 assays with microplate format

Contents: • 3.6 mL Galacton-Star[®] substrate

180 mL Reaction Buffer Diluent

containing Sapphire-II[™] enhancer

- 210 ml Mammalian Lysis Solution
- or 240 mL 5X Z Yeast Lysis Buffer

Galacto-Star™ System Screening Size

T1013 – Mammalian Cells

T1020 – Yeast Cells Capacity: 15,000 assays with microplate format Contents:

- 30 mL Galacton-Star® substrate
- 1.5 L Reaction Buffer Diluent
- containing Sapphire-II[™] enhancer • 1.75 L Mammalian Lysis Solution
- or 2 L 5X Z Yeast Lysis Buffer
- β-Galactosidase (Galacto-Star™ Reporter Gene Assay System; for a complete reference list, please see the Applied Biosystems website.)
- Carnes J, Jacobson M, Leinwand L, Yarus M (2003) Stop codon suppression via inhibition of eRF1 expression, RNA 9:648-653.
- Daub M, Jockel J, Quack T, Weber CK, Schmitz F, Rapp UR, Wittinghofer A, Block C (1998) The RafC1 cysteine-rich domain contains multiple distinct regulatory epitopes which control Rasdependent Raf activation, *Mol Cell Bio* 18:6698-6710.
- De Bosscher K, Hill CS, Nicolas FJ (2004) Molecular and functional consequences of Smad4 C-terminal missense mutations in colorectal tumour cells, *Biochem J* 379:209-216.
- Fan S, Gao M, Meng Q, Laterra JJ, Symons MH, Coniglio S, Pestell RG, Goldberg ID, Rosen EM (2005) Role of NF-kappaB signaling in hepatocyte growth factor/scatter factor-mediated cell protection, Oncogene 24:1749-1766.
- Farhadi HF, Lepage P, Forghani R, Friedman HC, Orfali W, Jasmin L, Miller W, Hudson TJ, Peterson AC (2003) A combinatorial network of evolutionarily conserved myelin basic protein regulatory sequences confers distinct glial-specific phenotypes, J Neurosci 23(32):10214-10223.
- Fouts T, Godfrey K, Bobb K, Montefiori D, Hanson CV, Kalyanaraman VS, DeVico A, Pal R (2002) Crosslinked HIV-1 envelope-CD4 receptor complexes elicit broadly cross-reactive neutralizing antibodies in rhesus macaques, Proc Natl Acad Sci USA 99(18):11842-7.
- Keeton EK, Brown M (2005) Cell cycle progression stimulated by tamoxifen-bound estrogen receptor-alpha and promoter-specific effects in breast cancer cells deficient in N-CoR and SMRT, Mol Endocrinol 19(6):1543-1554.
- Kensinger RD, Catalone BJ, Krebs FC, Wigdahl B, Schengrund CL (2004) Novel polysulfated galactose-derivatized dendrimers as binding antagonists of human immunodeficiency virus type 1 infection, Antimicrob Agents Chemother 48(5):1614-1623.
- Martin CS, Olesen CEM, Liu B, Voyta JC, Shumway JL, Juo RR, Bronstein I (1997) Continuous sensitive detection of -galactosidase with a novel chemiluminescent 1,2-dioxetane, Bioluminescence and Chemiluminescence: Molecular Reporting with Photons England p 525-528.

- Olesnicky NS, Brown AJ, Dowell SJ, Casselton LA (1999) A constitutively active G-protein-coupled receptor causes mating self-compatibility in the mushroom Coprinus, *Embo J* 18:2756-2763.
- 11. Pittler SJ, Zhang Y, Chen S, Mears AJ, Zack DJ, Ren Z, Swain PK, Yao S, Swaroop A, White JB(2004) Functional analysis of the rod photoreceptor cGMP phosphodiesterase alpha-subunit gene promoter: Nrl and Crx are required for full transcriptional activity, J Biol Chem 279(19):19800-19807.
- Rokosh DG, Simpson PC (2002) Knockout of the alpha 1A/C-adrenergic receptor subtype: the alpha 1A/C is expressed in resistance arteries and is required to maintain arterial blood pressure, Proc Natl Acad Sci USA 99(14):9474-9479.
- Scharf KD, Heider H, Hohfeld I, Lyck R, Schmidt E, Nover L (1998) The tomato Hsf system: HsfA2 needs interaction with HsfA1 for efficient nuclear import and may be localized in cytoplasmic heat stress granules, *Mol Cell Biol* 18(4):2240-2251.
- 14. Schotte P, Denecker G, Van Den Broeke A, Vandenabeele P, Cornelis GR, Beyaert R (2004) Targeting Rac1 by the Yersinia effector protein YopE inhibits caspase-1-mediated maturation and release of interleukin-1beta, J Biol Chem 279(24):25134-25142.
- Serna A, Ramirez MC, Soukhanova A, Sigal LJ (2003) Cutting edge: efficient MHC class I crosspresentation during early vaccinia infection requires the transfer of proteasomal intermediates between antigen donor and presenting cells, J Immunol 171:5668-5672.
- Silver SJ, Davies EL, Doyon L, Rebay I (2003) Functional dissection of eyes absent reveals new modes of regulation within the retinal determination gene network, *Mol Cell Biol* 23(17):5989-5999.
- Simmons G, Rennekamp AJ, Chai N, Vandenberghe LH, Riley JL, Bates P (2003) Folate receptor alpha and caveolae are not required for Ebola virus glycoprotein-mediated viral infection, J Virol 77(24):13433-13438.
- Singh MP, Greenstein M (2005) A simple, rapid, sensitive method detecting homoserine lactone (HSL)-related compounds in microbial extracts, J Microbiol Methods 65(1):32-37.

Description

Galacto-Light[™] and Galacto-Light Plus[™] assay systems are designed for rapid, ultrasensitive β -galactosidase reporter gene assays. These chemiluminescent assays exhibit over three orders of magnitude greater sensitivity compared to colorimetric β -galactosidase assays and are performed in a fraction of the time required for assays for chloramphenicol acetyl transferase (CAT) and ELISA-based reporter assays.

Galacto-Light Assay Systems includes Galacton[®] chemiluminescent β -galactosidase substrate and requires the use of a luminometer equipped with an injector. Galacto-Light Plus Assay systems includes Galacton-Plus[®] substrate and may be used in tube luminometers without injectors. For use with a microplate luminometer, an automated injector is highly recommended for accurate results. Galacton-Plus substrate exhibits prolonged signal emission compared with Galacton substrate (Figure 23). The Reaction Buffer is designed to aid in discrimination of bacterial β -galactosidase reporter activity from endogenous mammalian activity. This enables sensitive detection even in cell lines with relatively high levels of endogenous β -galactosidase activity [3].

Advantages

Galacto-Light and Galacto-Light Plus Assay Systems for β -galactosidase are among the most sensitive reporter gene assays available. As few as 2 femtograms of β -galactosidase (2,600 molecules) are detectable with Galacto-Light Plus Assay System [3]. High sensitivity allows detection of weak expression. Colorimetric, fluorometric, and isotopic assays cannot rival the dynamic range of Galacto-Light and Galacto-Light Plus Assay Systems, which ranges from femtogram to nanogram levels of protein (Figure 24). Galacto-Light Plus Assay System shows up to a five-fold increase in signal-to-noise over Galacto-Light Assay System, making it more suitable for monitoring lower levels of reporter activity. The assays use a lysis solution compatible with luciferase assays, making them ideal for transfection normalization; both assays are performed rapidly in a luminometer with comparable levels of sensitivity. Both assays require addition of two reagents following lysate preparation; one reagent addition requires the use of a luminometer with an automatic injector for microplate assays. New users are recommended to use the Galacto-*Star*^{**} or Gal-Screen[®] Assay Systems, which provide simpler assay protocols without the need for injection, while providing identical high sensitivity detection.

Applications

The Galacto-Light and Galacto-Light Plus Assay Systems are widely used for traditional reporter gene assays in transfected mammalian cell lines in culture [14, 7, 19], primary culture cells [8], tissue extracts from transgenic mice [17], frog embryo extracts [10], and *Drosophila* embryo extracts [15]. A variety of applications have been performed, such as viral function assays with β -gal-encoding MAGI cells [16], and targeted gene expression for gene therapy [6]. These assay systems have also been utilized in yeast reporter gene assays, including *Schizosaccharomyces* [11] and *Candida* [23], and for the study of protein:protein interactions with the yeast two-hybrid system [1, 6] and DNA:protein interactions with the one-hybrid system [24]. Galacto-Light Plus assay system has been used for reporter gene assays in *Pseudomonas* bacterial cells [20].

Several novel applications have been performed with the Galacto-Light assay systems, including cytotoxicity, by measurement of β -gal reporter







Figure 24. Detection of β-Galactosidase with Galacto-Light[™] System.



Figure 25. β-Galactosidase Reporter Gene Assay with Galacto-Light Plus™ System.

Galacto-Light[™] and Galacto-Light Plus[™] Assay Systems, continued β-Galactosidase reporter gene assay

enzyme released into culture media [12], and an RNA trans-splicing assay performed in primary fetal fibroblasts by measuring reconstitution of β -gal activity from partial transcripts [13]. These assays provide highly sensitive detection for β -galactosidase complementation assays, and have been used for intracellular detection of protein-protein interactions [22] and cell fusion [21].

The Galacto-Light[™] and Galacto-Light Plus[™] assay systems have wide application to assays that use the β -gal reporter enzyme as a functional read-out, enabling highly sensitive detection in many different types of cells and organisms.

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Galacto-Light™ Reporter Gene Assay System	standard size	T1006
Galacto-Light™ Reporter Gene Assay System	large size	T1010
Galacto-Light™ Reporter Gene Assay System	screening size	T1008
Galacto-Light Plus™ Reporter Gene Assay System	standard size	T1007
Galacto-Light Plus™ Reporter Gene Assay System	large size	T1011
Galacto-Light Plus™ Reporter Gene Assay System	screening size	T1009
Galacto-Light™ Reaction Buffer Diluent with Galacton® Chemiluminescent Substrate	120 mL 1.2 mL	T1054
Galacto-Light™ Reaction Buffer Diluent with Galacton-Plus® Chemiluminescent Substrate	120 mL 1.2 mL	T1055
Light Emission Accelerator (Galacto-Light™)	210 mL	T2084
Light Emission Accelerator-II (Galacto-Light Plus™)	210 mL	T2222

For Research Use Only. Not for use in diagnostic procedures

β-Galactosidase (Galacto-Light™ Reporter Gene Assay System; for a complete reference list, please see the Applied Biosystems website.)

- Dunoyer P, Thomas C, Harrison S, Revers F, Maule A (2004) A cysteine-rich plant protein potentiates Potyvirus movement through an interaction with the virus genome-linked protein VPg, J Virol 78(5):2301-2309.
- Fulton R, Van Ness B (1993) Luminescent reporter gene assays for luciferase and b-galactosidase using a liquid scintillation counter, *BioTechniques* 14:762-763.
- Jain VK, Magrath IT (1991) A chemiluminescent assay for the quantitation of b-galactosidase in the femtogram range: application to quantitation of b-galactosidase in lacZ-transfected cells, Anal Biochem 199:119-124.
- Lee SW, Trapnell BC, Rade JJ, Virmani R, Dichek DA (1993) In vivo adenoviral vector-mediated gene transfer into balloon-injured rat carotid arteries, *Circulation Research* 73:797-807.
- Moessler H, Mericskay M, Li Z, Nagl S, Paulin D, Small JV (1996) The SM22 promoter directs tissue-specific expression in arterial but not venous or visceral smooth muscle cells in transgenic mice, *Development* 122:2415-2425.
- Ozturk-Winder F, Renner M, Klein D, Muller M, Salmons B, Gunzburg WH (2002) The murine whey acidic protein promoter directs expression to human mammary tumors after retroviral transduction, Cancer Gene Therapy 9:421-431.
- Petersen RK, Jorgensen C, Rustan AC, Froyland L, Muller-Decker K, Furstenberger G, Berge RK, Kristiansen K, Madsen L (2003) Arachidonic acid-dependent inhibition of adipocyte differentiation requires PKA activity and is associated with sustained expression of cyclooxygenases, *Journal* of Lipid Research 44:2320-2330.
- Saller RM, Ozturk F, Salmons B, Gunzburg WH (1998) Construction and characterization of a hybrid mouse mammary tumor virus/murine leukemia virus-based retroviral vector, J Virol 72(2):1699-1703.
- Shaper NL, Harduin-Lepers A, Shaper JH (1994) Male germ cell expression of murine b4-galactosyltransferase, J Biol Chem 40:25165-25171
- Gove C, Walmsley M, Nijjar S, Bertwistle D, Guille M, Partington G, Bomford A, Patient R (1997) Over-expression of GATA-6 in Xenopus embryos blocks differentiation of heart precursors, *EMBO J* 16(2):355-368.
- Remacle JE, Albrecht G, Brys R, Braus GH, Huylebroeck D (1997) Three classes of mammalian transcription activation domain stimulate transcription in Schizosaccharomyces pombe, EMBO J 16(18):5722-5729.
- Schafer H, Schafer A, Kiderlen AF, Masihi KN, Burger R (1997) A highly sensitive cytotoxicity assay based on the release of reporter enzymes, from stably transfected cell lines, *J Immunol Methods* **204**:89-98.

Product Configuration

Galacto-Light[™] or Galacto-Light Plus[™] System Standard Size

T1006 – Galacto-Light[™] System

T1007 – Galacto-Light Plus[™] System Capacity: 200 single tube assays

- Capacity: 600 assays with microplate format
- Contents:
- 0.4 mL Galacton[®] (Galacto-Light[™]) or Galacton-Plus[®] substrate (Galacto-Light Plus[™])
- 40 mL Reaction Buffer Diluent
- 70 mL Light Emission Accelerator (Galacto-Light[™]) or Light Emission Accelerator-II (Galacto-Light Plus[™])
 70 mL Lysis Solution

Galacto-Light[™] or Galacto-Light Plus[™] System Large Size

T1010 – Galacto-Light[™] System

T1011 – Galacto-Light Plus[™] System

Capacity: 600 single tube assays

Capacity: 1800 assays with microplate format Contents:

- 1.2 mL Galacton[®] (Galacto-Light[™]) or Galacton-Plus[®] substrate (Galacto-Light Plus[™])
- 1.2 mL Galactor (Galactoright 7 of Galactorin lus' substrate (Galactoright 1 lus'
 120 mL Reaction Buffer Diluent
- 210 mL Light Emission Accelerator (Galacto-Light™ System) or Light Emission Accelerator-II (Galacto-Light Plus™)
- 210 mL Lysis Solution

Galacto-Light[™] or Galacto-Light Plus[™] System Screening Size

T1008 – Galacto-Light[™] System

T1009 – Galacto-Light Plus[™] System Capacity: 15.000 assays with microplate format

- Contents:
- 10 mL Galacton[®] (Galacto-Light[™]) or Galacton-Plus[®] substrate (Galacto-Light Plus[™])
- 1 L Reaction Buffer Diluent
- 1.75 L Light Emission Accelerator (Galacto-Light[™]) or Light Emission Accelerator-II (Galacto-Light Plus[™])
- 1.75 L Lysis Solution

β-Galactosidase (Galacto-Light Plus™ Reporter Gene Assay System; for a complete reference list, please see the Applied Biosystems website.)

- Duan D, Yue Y, Engelhardt JF (2001) Expanding AAV packaging capacity with trans-splicing or overlapping vectors: A quantitative comparison, *Molecular Therapy* 4(4):383-391.
- Graslund T, Li X, Magnenat L, Popkov M, Barbas CF (2005) Exploring strategies for the design of artificial transcription factors, J Biol Chem 280(5):3707-3714.
- Liu QX, Jindra M, Ueda H, Hiromi Y, Hirose S (2003) Drosophila MBF1 is a co-activator for Tracheae Defective and contributes to the formation of tracheal and nervous systems, *Development* 130:719-728.
- Neurath AR, Strick N, Li YY (2002) Anti-HIV-1 activity of anionic polymers: a comparative study of candidate microbicides, *BMC Infectious Diseases* 2:27 (www.biomedcentral.com/1471-2334/2/27).
- Ryan AJ, Fisher K, Thomas CP, Mallampalli RK (2004) Transcriptional repression of the CTP: phosphocholine cytidyltransferase gene by sphingosine, Biochem J 382:741-750.
- Schumacher S, Laass K, Kant S, Shi Y, Visel A, Gruber AD, Kotlyarov A, Gaestel M (2004) Scaffolding by ERK3 regulates MK5 in development, *EMBO J* 23:4770-4779.
- Treeck O, Diedrich K, Ortmann O (2003) The activation of an extracellular signal-regulated kinase by oestradiol interferes with the effects of trastuzumab on HER2 signalling in endometrial adenocarcinoma cell lines, *European Journal of Cancer* **39**:1302-1309.
- Yarwood JM, Volper EM, Greenberg EP (2005) Delays in Pseudomonas aeruginosa quorum-controlled gene expression are conditional, Proc Natl Acad Sci USA 102(25):9008-9013.
- Charlton CA, Mohler WA, Radice GL, Hynes RO, Blau HM (1997) Fusion competence of myoblasts rendered genetically null for N-cadherin in culture, J Cell Biol 138:331-336.
- Rossi F, Charlton CA, Blau HM (1997) Monitoring protein-protein interactions in intact eukaryotic cells by beta-galactosidase complementation, Proc Natl Acad Sci USA 94:8405-8410.
- 23. Stoldt VR, Sonneborn A, Leuker CE, Ernst JF (1997) Efg1p, an essential regulator of morphogenesis of the human pathogen Candida albicans, is a member of a conserved class of bHLH proteins regulating morphogenetic processes in fungi, *EMBO J* 16(8):1982-1991.
- Wolf SS, Roder K, Schweizer M (1997) Construction of a reporter plasmid that allows expression libraries to be exploited for the one-hybrid system, *Biotechniques* 20(4):568-574.

Reporter Gene Assays and Reagents

Luc-Screen[®] Assay System Extended-glow firefly luciferase reporter gene assay

Description

The Luc-Screen[®] assay system with extended-glow light emission is designed for sensitive detection of firefly luciferase reporter enzyme, especially for high throughput screening assays. Luciferase is an ideal reporter due to the high sensitivity of detection and the absence of endogenous luciferase activity in mammalian cells. Luc-Screen conveniently couples in-well cell lysis in the presence of culture medium with a high sensitivity assay that exhibits extended-glow light emission kinetics. Light signal can be measured between 10 minutes and several hours after adding assay reagents. Luc-Screen system is designed for maximum assay flexibility in a high-throughput format and can be used in luminometers without automatic injectors.

Luc-Screen system is formulated to provide a convenient and easy-to-use firefly luciferase assay that is optimized for use in high-throughput screening. Cells are seeded into opaque white tissue culture (TC)-treated microplates (not supplied) or clear-bottom/opaque white side TC-treated microplates, if desired. The system has two reagents that are added to cells in culture medium (with or without phenol red) in microplate wells. Although the presence of phenol red causes some decrease in signal intensity, assay

sensitivity remains unaffected. Cell lysis occurs during an initial 10 minute incubation; within this time, light signal reaches plateau. Light emission persists with a signal half-life of 4-5 hours, providing flexibility in the time between reagent addition and measurement.

Advantages

Luciferase is one of the most sensitive reporter enzymes available. The Luc-Screen[®] system detects fewer than 50 femtograms of pure enzyme in culture medium samples (Figure 26). The high sensitivity of Luc-Screen[®] is complemented by a wide dynamic range. A linear signal is obtained with the Luc-Screen assay from 50 femtograms to 100 nanograms of pure enzyme in culture medium, a dynamic range of six orders of magnitude.

Applications

The Luc-Screen reporter gene assay system is ideal for high throughput firefly luciferase reporter gene expression assays in mammalian cells, and has been used for gene expression assays [1], compound screening [3] and large-scale promoter function assays [2].







Figure 27. Luc-Screen® Assay: Forskolin Induction of pCRE-Luc Transfected Cells. Cells were transfected with a luciferase construct under the control of a cAMP-responsive promoter, and were then seeded into a 96-well plate. Several wells were treated with forskolin to induce cAMP and stimulate luciferase induction in a screening-like assay. Measurements were performed on a CCD microplate luminometer.



Seed cells into tissue culture (TC) - treated

luminometer plates and culture or treat as desired

Add Buffer 1 and Buffer 2

Incubate 10-60 min



Place plate in luminometer and measure luciferase light emission

Figure 28. Luc-Screen[®] Reporter Gene Assay

Luc-Screen® Assay System, continued

Extended-glow firefly luciferase reporter gene assay

Product Configuration

Luc-Screen^{*} reagents come as two buffers which are added to culture medium in a 2:1:1 (culture medium : Buffer 1 : Buffer 2) ratio. Luc-Screen system is available in several sizes so that customers can match their needs with an appropriate kit.

	ASSAYS PER KIT	BUFFER 1	BUFFER 2
T1035	200	2 x 5 mL	2 x 5 mL
T1033	1,000	2 x 25 mL	2 x 25 mL
T1036	5,000	1 x 250 mL	1 x 250 mL
T1034	10,000	2 x 250 mL	2 x 250 mL

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Luc-Screen® Extended-Glow Firefly Luciferase Assay System	200 assays	T1035
Luc-Screen® Extended-Glow Firefly Luciferase Assay System	1,000 assays	T1033
Luc-Screen® Extended-Glow Firefly Luciferase Assay System	5,000 assays	T1036
Luc-Screen® Extended-Glow Firefly Luciferase Assay System	10,000 assays	T1034

For Research Use Only. Not for use in diagnostic procedures.

Luciferase (Luc-Screen® Reporter Gene Assay System; (For a complete reference list, please see the Applied Biosystems website.)

- Ahmed CM, Wills KN, Sugarman BJ, Johnson DE, Ramachandra M, Nagabhushan TL, Howe JA (2001) Selective expression of nonsecreted interferon by an adenoviral vector confers antiproliferative and antiviral properties and causes reduction of tumor growth in nude mice, *Journal of Interferon and Cytokine Research* 21:399-408.
- Coleman SL, Buckland PR, Hoogendoorn B, Guy C, Smith K, O'Donovan MC (2002) Experimental analysis of the annotation of promoters in the public database, *Human Molecular Genetics* 11(16):1817-1821.
- Wang X, Miyake H, Okamoto M, Saito M, Fujisawa J, Tanaka Y, Izumo S, Baba M (2002) Inhibition of the Tax-dependent Human T-Lymphotrpic Virus Type I replication in persistently infected cells by the fluoroquinolone derivative K-37, *Molecular Pharmacology* **61(6)**:1359-1365.

Reporter Gene Assays and Reagents

Luciferase Assay System

Luciferase reporter gene assay system

Description

Luciferase Assay System is a bioluminescent assay system designed for rapid, sensitive detection of firefly luciferase expressed by transfected cells. The kit incorporates the substrate luciferin and proprietary reagents to enhance light emission. This enhanced luciferase/luciferin reaction produces a light signal that decays with a half-life of approximately 5 minutes. The assay is compatible with the lysis buffer included in Tropix β -galactosidase assay kits, making this assay ideal for co-transfections.

Cell lysate is mixed with Substrate A, which contains reagents necessary for the luciferase reaction. Light signal from the luciferase enzyme present in the extract is measured immediately after the injection of Substrate B, containing luciferin.

Advantages

The high sensitivity and absence of endogenous luciferase activity in the majority of cell types makes luciferase an excellent reporter enzyme. The wide dynamic range of the assay enables accurate measurement of luciferase concentration from the femtogram to nanogram range (Figure 29). The wide dynamic range of seven orders of magnitude cannot be rivaled by fluorimetric or colorimetric assays for other reporter proteins.

Applications

The Luciferase Assay System provides a traditional flash kinetics luciferase assay system for measuring gene expression from a firefly luciferase reporter construct in mammalian cells [1, 2]. This assay system requires the use of a luminometer with injection capabilities, and the reaction kinetics provide a very rapid assay read-out.



Figure 29. Detection Limit and Dynamic Range of Luciferase Assay Kit.

Product Configuration

Luciferase Assay System Standard Size

T1000

Capacity: 200 assays Contents: • 20 mL of Luciferase Assay System Substrate A^{*} • 20 mL of Luciferase Assay System Substrate B^{*}

• 20 THE OF LUCIER ASE ASSAY SYSTEM SUBSTRATE B

70 mL of Luciferase Assay System Lysis Solution

Luciferase Assay System Large Size

T1002

Capacity: 600 assays

Contents: • 3 x 20 mL of Luciferase Assay System Substrate A*

3 x 20 mL of Luciferase Assay System Substrate B*

• 210 mL of Luciferase Assay System Lysis Solution

Luciferase Assay System Screening Size

T1001

Capacity: 5,000 assays

Contents:

5 x 100 mL of Luciferase Assay System Substrate A*

5 x 100 mL of Luciferase Assay System Substrate B*

1.75 L of Luciferase Assay System Lysis Solution

*Supplied lyophilized

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Luciferase Reporter Gene Assay System	standard size	T1000
Luciferase Reporter Gene Assay System	large size	T1002
Luciferase Reporter Gene Assay System	screening size	T1001

For Research Use Only. Not for use in diagnostic procedures.

Luciferase (Luciferase Assay System; for a complete reference list, please see the Applied Biosystems website.)

- Poser S, Impey S, Xia Z, Storm DR (2003) Brain-derived neurotrophic factor protection of cortical neurons from serum withdrawal-induced apoptosis is inhibited by cAMP J Neuroscience 23(11):4420-4427.
- Singh RP, Dhawan P, Golden C, Kapoor GS, Mehta KD (1999) One-way cross-talk between p38MAPK and p42/44MAPK, J Biol Chem 274(28):19593-19600



Figure 30. Luciferase Reporter Gene Assay.

Dual-Light[®] Assay System

Combined firefly luciferase and β -Galactosidase reporter gene assay

Description

Dual-Light* luminescent reporter gene assay is designed for the rapid and sensitive detection of firefly luciferase and β -galactosidase in the same sample. The use of two reporter genes, one as an experimental reporter and the other as a constitutively-expressed transfection control reporter, is very widely used and is often necessary to accurately quantitate activity from experimental reporter constructs. This assay enables the measurement of firefly luciferase and β -galactosidase in a single aliquot of cell lysate. Luciferase is typically used as the experimental reporter, and β -galactosidase is typically quantitated from a co-transfected constitutive expression vector to determine transfection efficiency. The luciferase measurement is then normalized to the β -galactosidase measurement. Light signal from each enzymatic reaction is measured sequentially in a luminometer with automatic injectors.

First, luciferase reporter enzyme activity is quantitated with an enhanced luciferase reaction. Following a 30-60 minute incubation and addition of a light emission accelerator, β -galactosidase reporter enzyme activity is determined with Galacton-Plus[®] substrate.

Advantages

Both reporter enzyme measurements are combined into a single sequential assay protocol using only one aliquot of extract for greater convenience and precision. The entire assay is completed in less than one hour. The wide dynamic range (Figure 31) of this dual assay enables accurate measurement of firefly luciferase and β -galactosidase concentrations over seven orders of magnitude, from the femtogram to nanogram range [9].

Applications

Dual-Light^{*} reporter gene assay system has been very widely used for reporter quantitation/transfection normalization from transiently transfected mammalian cell lines [1, 5, 6, 7, 8], as well as transfected primary cells [2, 3, 4, 10]. In addition, it has been used with a modified lysis buffer to quantitate luciferase and β -galactosidase activities from a novel reporter fusion construct in yeast cells [11].



Figure 31. Detection of Firefly Luciferase and $\beta\mbox{-}Galactosidase$ with Dual-Light* System.



Figure 32. Firefly Luciferase and β-Galactosidase Reporter Gene Assays with Dual-Light[®] System.

Dual-Light[®] Assay System, continued

Combined firefly luciferase and β -Galactosidase reporter gene assay

Product Configuration

Dual-Light® System Standard Size

T1003

Capacity: 200 combined assays Contents:

- 0.2 mL Galacton-Plus[®] substrate
- 5 ml. Buffer A*
- 22 mL Buffer B*
- 25 mL Light Emission Accelerator-II
- 70 mL Lysis Solution

Dual-Light® System Large Size

T1005

Capacity: 600 combined assays

- Contents:
- 0.6 mL Galacton-Plus[®] substrate
- 3 x 5 mL Buffer A*
- 3 x 22 mL Buffer B*
- 75 mL Light Emission Accelerator-II
- 210 mL Lysis Solution

Dual-Light® System Screening Size

T1004

Capacity: 4,000 combined assays Contents:

- 4 mL Galacton-Plus® substrate
- 20 x 5 mL Buffer A*
- 20 x 22 mL Buffer B*
- 500 mL Light Emission Accelerator-II
- 1.4 L Lysis Solution
- *Supplied lyophilized

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Dual-Light® Reporter Gene Assay System	standard size	T1003
Dual-Light® Reporter Gene Assay System	large size	T1005
Dual-Light® Reporter Gene Assay System	screening size	T1004
Dual-Light® Buffer A	5 mL	T2176
Dual-Light® Buffer B vith Galacton-Plus® Substrate	22 mL/0.2 mL	T1053
Light Emission Accelerator-II	210 mL	T2222

For Research Use Only. Not for use in diagnostic procedures.

 β -Galactosidase/Luciferase (Dual-Light® Reporter Gene Assay System; for a complete reference list, please see the Applied Biosystems website.)

- Babb R, Bowen BR (2003) SDP1 is a peroxisome-proliferator-activated receptor g2 co-activator that binds through its SCAN domain, *Biochem J* **370**:719-727.
- Bourcier T, Sukhova G, Libby P (1997) The nuclear factor k-B signaling pathway participates in dysregulation of vascular smooth muscle cells in vitro and in human atherosclerosis, J Biol Chem 272(25):15817-15824.
- Brown, AM, Lemke G (1997) Multiple regulatory elements control transcription of the peripheral myelin protein zero gene, J Biol Chem 272(46):28939-28947.
- Brown JD, DiChiara MR, Anderson KR, Gimbrone MA Jr, Topper JN (1999) MEKK-1, a component of the stress (stress-activated protein kinase/c-Jun N-terminal kinase) pathway, can selectively activate Smad2-mediated transcriptional activation in endothelial cells, J Biol Chem 274:8797-8805.
- Crowe DL, Chandraratna RAS (2004) A retinoid X receptor (RXR)-selective retinoid reveals that RXR-a is potentially a therpautic target in breast cancer cell lines, and that it potentiates antiproliferative and apoptotic responses to peroxisome proliferator-activated receptor ligands, *Breast Cancer Research* 6:R546R555 (DOI 10.1186/bcr913).
- Fedele M, Pentimalli F, Baldassarre G, Battista S, Klein-Szanto AJ, Kenyon L, Visone R, De Martino I, Ciarmiello A, Arra C et al (2005) Transgenic mice overexpressing the wild-type form of the HMGA1 gene develop mixed growth hormone/prolactin cell pituitary adenomas and natural killer lymphomas, Oncogene 24:3427:3435.
- Figueroa C, Vojtek AB (2003) Akt negatively regulates translation of the ternary complex factor Elk-1, Oncogene 22:5554-5561.
- Hollenberg AN, Susulic VS, Madura JP, Zhang B, Moller DE, Tontonoz P, Sarraf P, Spiegelman BM, Lowell BB (1997) Functional antagonism between CCAAT/Enhancer binding protein-a and peroxisome proliferator-activated receptor-g on the leptin promoter, J Biol Chem 272(8):5283-5290.
- Martin CS, Wight PA, Dobretsova A, Bronstein I.(1996) Dual luminescence-based reporter gene assay for luciferase and b-galactosidase, *BioTechniques* 21(3):520-524.
- Takemoto, Sun MJ, Hiroki J, Shimokawa H, Liao JK (2002) Rho-kinase mediates hypoxia-induced downregulation of endothelial nitric oxide synthase, *Circulation* **106**:57-62.
- Williams I, Richardson J, Starkey A, Stansfield I (2004) Genome-wide prediction of stop codon readthrough during translation in the yeast Saccharomyces cerevisiae, *Nucl Acids Res* 32(22):6605-6616.

Phospha-Light<u>[™] Assay System</u>

Secreted placental alkaline phosphatase (SEAP) reporter gene assay

Description

Phospha-Light[™] assay system is a chemiluminescent reporter gene assay for the sensitive detection of secreted placental alkaline phosphatase (SEAP). SEAP is a reporter protein that is secreted into the cell culture media and detected by testing aliquots of media [2]. The SEAP gene product is a truncated form of human placental alkaline phosphatase. Detection of nonsecreted placental alkaline phosphatase (PLAP) is also possible [16].

A combination of proprietary Tropix® reagents provides highly sensitive detection with the Phospha-Light Assay System. The system incorporates CSPD® high performance alkaline phosphatase substrate, Emerald™ luminescence enhancer, and a unique buffer system designed to specifically inhibit endogenous non-placental alkaline phosphatase activity.

A sample of cell culture media is transferred to a microfuge tube and diluted with Phospha-Light dilution buffer. After heating 30 minutes at 65°C, the sample is transferred to a luminometer tube or microplate. Phospha-Light assay buffer containing differential alkaline phosphatase inhibitors (to achieve high sensitivity in media containing non-placental alkaline phosphatase) and reaction buffer containing CSPD substrate and Emerald enhancer are then added. After a 20-minute incubation at room temperature, the chemiluminescent signal is measured in a luminometer. The unique properties of human placental alkaline phosphatase, including its heat stability and resistance to L-homoarginine, enable it to be distinguished from endogenous non-placental alkaline phosphatase activity.

Light emission kinetics observed with Phospha-Light Assay System reagents are shown (Figure 34). The persistent glow of the chemiluminescent signal permits the use of simple luminometers without injectors or scintillation counters for measurement.

Advantages

Secreted reporter proteins eliminate the need for cell lysis. A population of cells can be monitored over time and remains intact for further experimentation. The Phospha-Light assay is one of the easiest and fastest methods for optimizing transfection efficiency. This chemiluminescent method allows detection of fewer than 10 femtograms of SEAP (Figure 33), which is three orders of magnitude more sensitive than colorimetric detection of SEAP. The wide linear range of six orders of magnitude of the Phospha-Light assay enables accurate intra-assay comparisons without measuring several sample dilutions.

Applications

The Phospha-Light reporter gene assay system has been used widely for reporter gene assays to measure gene expression in established cell lines [4] and in transfected primary cells [15, 21], including as a gene knockdown/ RNA interference read-out [5]. The Phospha-Light reporter gene assay has been used for a wide variety of viral functional assays, including viral gene expression assays [8], viral replication [9], viral fusogenicity [1], virus neutralization and viral-mediated cell-cell fusion [10], and viral infectivity [14]. Use of the SEAP reporter protein is very enabling for in vivo reporter gene assays, by assaying serum samples from transgenic, transfected or viral vector-infected animals. The Phospha-Light reporter gene assay system has been used to measure SEAP levels in sera from transgenic or transfected whole animals, including mouse [18], rat [17], marmoset [7], monkey [13] and pig sera [12], and in chicken egg allantoic fluid [22]. The mouse SEAP protein (mSEAP) has recently been developed for improved SEAP protein stability in transgenic mice, and the Phospha-Light[™] system has been used for sensitive detection of mSEAP [18].











Figure 35. Secreted Alkaline Phosphatase Reporter Gene Assay with Phospha-Light[™] System.

Reporter Gene Assays and Reagents

Phospha-Light[™] Assay System, continued

Secreted placental alkaline phosphatase (SEAP) reporter gene assay

In addition to reporter gene (gene expression) applications, the Phospha-Light assay system is used to measure SEAP as a functional reporter for receptor-ligand binding assays with a SEAP-ligand chimera [20], proteasemediated secretion [19], and for secretion pathway activity [11]. Finally, the Phospha-Light[™] assay system has also been used for the cellular measurement of non-placental alkaline phosphatase as a biomarker [4].

Sensitive detection of SEAP or PLAP reporter protein with the Phospha-Light assay system enables a large number of applications in many areas of life science research, including gene expression, viral function assays, vaccine development, development of viral vectors and gene delivery methods for gene therapy, in vivo gene expression monitoring and novel cellular functional assays.

Product Configuration

Phospha-Light[™] System Standard Size

T1015

Capacity: 200 single tube assays (triplicate assays of 66 samples or duplicate assays of 83 samples) Capacity: 400 assays with microplate format (triplicate assays of 133 samples or duplicate assays of 166 samples)

Contents

- 1.0 mL CSPD[®] chemiluminescent substrate
- 19 mL Phospha-Light[™] Reaction Buffer Diluent with Emerald[™] enhancer
- 20 mL Phospha-Light[™] Assay Buffer
- 5 mL 5X Dilution Buffer
- 50 µL Positive Control Placental Alkaline Phosphatase

Phospha-Light[™] System Large Size

T1017

Capacity: 600 single tube assays (triplicate assays of 200 samples or duplicate assays of 250 samples) Capacity: 1,200 assays with microplate format (triplicate assays of 400 samples or duplicate assays of 500 samples)

Contents:

- 3.0 mL CSPD[®] chemiluminescent substrate
- 57 mL Phospha-Light[™] Reaction Buffer Diluent with Emerald[™] enhancer
- 60 mL Phospha-Light[™] Assay Buffer
- 15 mL 5X Dilution Buffer
- 50 µL Positive Control Placental Alkaline Phosphatase

Phospha-Light[™] System Screening Size

T1016

Capacity: 10,000 assays with microplate format (triplicate assays of 3,333 samples or duplicate assays of 4,166 samples)

Contents:

- 25 mL CSPD® chemiluminescent substrate
- 475 mL Phospha-Light[™] Reaction Buffer Diluent with Emerald[™] enhancer
- 500 mL Phospha-Light[™] Assay Buffer
- 125 mL 5X Dilution Buffer
- 425 µL Positive Control Placental Alkaline Phosphatase

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Phospha-Light™ Reporter Gene Assay System	standard size	T1015
Phospha-Light™ Reporter Gene Assay System	large size	T1017
Phospha-Light™ Reporter Gene Assay System	screening size	T1016

For Research Use Only. Not for use in diagnostic procedures.

Secreted Placental Alkaline Phosphatase (Phospha-Light™ Reporter Gene Assay System; for a complete reference list, please see the Applied Biosystems website.)

- Alexander L, Illyinskii PO, Lang SM, Means RE, Lifson J, Mansfield K, Desrosiers RC (2003) Determinanants of increased replicative capacity of serially passaged Simian Immunodeficiency Virus with nef deleted in Rhesus monkeys, J Virol 77(12):6823-6835.
- Berger J, Hauber J, Hauber R, Geiger R, Cullen BR (1988) Secreted placental alkaline phosphatase: a powerful new quantitative indicator of gene expression in eukaryotic cells. Gene 66: 1-10.
- Bronstein I, Fortin JJ, Voyta JC, Juo RR, Edwards B, Olesen CE, Lijam N, Kricka LJ (1994a) Chemiluminescent reporter gene assays: Sensitive detection of the GUS and SEAP gene products, *BioTechniques* 17:172-178.
- Brown MA, Zhao Q, Baker KA, Naik C, Chen C, Pukac L, Singh M, Tsareva T, Parice Y, Mahoney A et al (2005) Crystal structure of BMP-9 and functional interactions with pro-region and receptors, *J Biol Chem* 280(26):25111-25118.
- Cao HB, Wang A, Martin B, Koehler DR, Zeitlin PL, Tanawell AK, Hu J (2005) Down-regulation of IL-8 expression in human airway epithelial cells through helper-dependent adenoviral-mediated RNA interference, *Cell Research* 15(2):111-119.
- Cullen B, Malim M (1992) Secreted placental alkaline phosphatase as a eukaryotic reporter gene, Methods Enzymol 216:362-368.
- Duboise M, Guo J, Czajak S, Lee H, Veazey R, Desrosiers RC, Jung JU (1998) A role for Herpesvirus Saimiri orf14 in transformation and persistent infection, J Virol 72(8):6770-6776.
- Hobbs WE, Brough DE, Kovesdi I, DeLuca NA (2001) Efficient activation of viral genomes by levels of Herpes Simplex Virus ICPO insufficient to affect cellular gene expression or cell survival, *J Virol* **75(7)**:3391-3403.
- Hwang DR, Tsai YC, Lee JC, Huang KK, Lin RK, Ho CH, Chiou JM, Lin YT, Hsu JT, Yeh CT (2004) Inhibition of Hepatitis C virus replication by arsenic trioxide, *Antimicrobial Agents and Chemotherapy* 48(8):2876-2882.
- Johnson WE, Morgan J, Reitter J, Puffer BA, Czajak S, Doms RW, Desrosiers RC (2002) A replication-competent, neutralization-sensitive variant of Simian Immunodeficiency Virus lacking 100 amino acids of envelope, J Virol 76(5):2075-2086.
- 11.Kagan JC, Stein MP, Pypaert M, Roy CR (2004) Legionella subvert the functions of Rab1 and Sec22b to create a replicative organelle, J Exp Med 199(9):1201-1211.
- Khan AS, Smith LC, Abruzzese RV, Cummings KK, Pope MA, Brown PA, Draghia-Akli R (2003) Optimization of electroporation parameters for the intramuscular delivery of plasmids in pigs, DNA and Cell Biology 22(12):807-814.
- 13 Latta-Mahieu M, Rolland M, Caillet C, Wang M, Kennel P, Mahfouz I, Loquet I, Dedieu JF, Mahfoudi A, Trannoy E, Thuillier V (2002) Gene transfer of a chimeric trans-activator is immunogenic and results in short-lived transgene expression, *Human Gene Therapy* **13**:1611-1620.
- Pohlmann S, Krumbiegel M, Kirchhoff F (1999) Coreceptor usage of BOB/GPR15 and Bonzo/ STRL33 by primary isolates of human immunodeficiency virus type 1, J Gen Virol 80:1241-1251.
- Poser S, ImpeyS, Xia Z, Storm DR (2003) Brain-derived neurotrophic factor protection of cortical neurons from serum withdrawal-induced apoptosis is inhibited by cAMP, J Neuroscience 23(11):4420-4427.
- O'Connor KL, Culp LA (1994) Quantitation of two histochemical markers in the same extract using chemiluminescent substrates, *Biotechniques* 17(3):502-509.
- Riera M, Chillon M, Aran JM, Cruzado JM, Torras J, Grinyo JM, Fillat C (2004) Intramuscular SP1017-formulated DNA electrotransfer enhances transgene expression and distributes hHGF to different rat tissues, *Journal of Gene Medicine* 6:111-118.
- Rubenstrunk A, Orsini C, Mahfoudi A, Scherman D (2003) Transcriptional activation of the metallothionein I gene by electric pulses in vivo: Basis for the development of a new gene switch system, *Journal of Gene Medicine* 5:773-783.
- Sakai J, Rawson RB, Espenshade PJ, Cheng D, Seegmiller AC, Goldstein JL, Brown MS (1998) Molecular identification of the sterol-regulated luminal protease that cleaves SREBPs and controls lipid composition of animal cells, *Molecular Cell* 2:505-514.
- Zabeau L, Defeau D, Van der Heyden J, Iserentant H, Vandekerckhove J, Tavernier J (2004) Functional analysis of leptin receptor activation using a Janus kinase/signal transducer and activator of transcription complementation assay, *Molecular Endocrinology* 18(1):150-161.
- 21. Zhang J, Ou J, Bashmakov Y, Horton JD, Brown MS, Goldstein JL (2001) Insulin inhibits transcription of IRS-2 gene in rat liver through an insulin response element (IRE) that resembles IREs of other insulin-repressed genes, *Proc Natl Acad Sci USA* 98(7):3756-3761.
- Zhao H, Peeters BPH (2003) Recombinant Newcastle Disease Virus as a viral vector: Effect of genomic location of foreign gene on gene expression and virus replication, J Gen Virol 84:781-788.

Immunodetection Products

Introduction to Chemiluminescent ELISA.26cAMP-Screen® System28cAMP-Screen Direct® System30ELISA-Light™ System32
Introduction to
Chemiluminescent Membrane-Based Immunodetection
Western-Star [™] System
Western-Light [™] System

Introduction to Chemiluminescent ELISAs

ELISA

Enzyme-linked immunosorbent assays (ELISAs) are used for sensitive analyte detection. ELISAs can be formatted in several configurations on a variety of solid supports.

A direct sandwich ELISA is often used for the detection of larger molecules with multiple antigenic sites, usually a protein. In this format (Figure 36, Panel A), a solid support is coated with a capture antibody specific for the protein of interest; the capture antibody forms an immunocomplex with the antigen from the sample; a detector antibody covalently derivatized with alkaline phosphatase, specific for a second antigenic site on the captured protein, is then added. After washing, substrate with enhancer is added. The enzyme-generated signal is proportional to the concentration of antigen.

A variation of this assay is used to screen hybridomas for the production of monoclonal antibodies (Figure 36, Panel B). In this case, the solid phase is coated with antigen, hybridoma culture supernatant containing a monoclonal antibody is added, and then an enzyme-labeled secondary antibody specific for the monoclonal antibody is used to detect the antigen-captured monoclonal antibody.

Competitive immunoassays can be set up in two modes (Figure 36, Panel C and D). In these assays, a competition occurs between bound and free antigen, or labeled and unlabeled antigen, for available antibody binding sites. This type of assay results in an inverse standard curve. A decrease in signal is observed for increasing concentrations of antigen.

Chemiluminescent Immunoassays

The use of 1,2-dioxetane enzyme substrates permits the ultrasensitive detection of analytes by ELISA. The most common enzyme used in conjunction with 1,2-dioxetane substrates for ELISA applications is alkaline phosphatase. Tropix® 1,2-dioxetane substrates have been used for detecting thyroidstimulating hormone (TSH) with a commercially available ELISA employing an alkaline phosphatase label [2]. The results obtained with CSPD® substrate with Sapphire-II™ enhancer (Figure 37) show a significant improvement in signal-to-background performance at all concentrations of TSH compared to those obtained with the colorimetric substrate, *p*-nitrophenyl phosphate (pNPP). This benefit can be expected when a colorimetric direct sandwich ELISA is converted to 1,2-dioxetane/enhancer chemiluminescence. Sandwich immunoassay formats with 1,2-dioxetane substrates have been used for calculation of antigen-antibody binding constants [10], and quantitation of animal and human proteins from plasma and tissue extracts [9,15,18].

Competitive ELISAs can also be highly sensitive with the incorporation of 1,2-dioxetane substrates [1, 7, 13]. Because the standard curve in a competitive ELISA exhibits maximum signal intensity at the lowest concentrations of analyte, it may be necessary to adjust reagent concentrations to optimize detection of low analyte concentrations. The sensitivity of chemiluminescent detection permits the use of lower concentrations of capture antibody and competing antigen.

CDP-Star[®] substrate with Sapphire-II enhancer or Emerald-II[™] enhancer has become widely used for immunoassay protein detection applications such as detection of viral antigens [14] and plasma proteins [9].



Figure 36. Immunoassay Schemes Employing an Alkaline Phosphatase Label for Detection with CSPD® Substrate.

Introduction, continued

Additional Applications

Whole-cell ELISA

Immunoassay detection of surface antigens on whole cells has been demonstrated with 1,2-dioxetane chemiluminescent detection [19]. β -galactosidase enzyme conjugates can also be used with Galacton-*Star** substrate with Sapphire-II[®] enhancer for chemiluminescent immunoassay detection, particularly for whole cell ELISA applications [12] that may exhibit high levels of cellular alkaline phosphatase.

Protein Detection Applications

Anti-phosphopeptide immunoassays with CSPD® or CDP-Star® substrates and Sapphire-II or Emerald-II[™] enhancers have been developed for quantitation of several protein kinase activities. These include PKA, PKC, CAM-KII, receptor interacting protein and src kinases [16], WaaP protein tyrosine kinase and sugar kinase [22] and p38 kinase [8]. In addition a receptor binding assay of a neurotrophic factor to a tyrosine kinase receptor has been demonstrated [20].

Nucleic Acid and Nucleic Acid-Protein Interaction Detection Applications

DNA probe hybridization assays, DNA:protein interaction assays, and DNA aptamer binding assays are often formatted in microplate wells or on other solid phases. AP-labeled probes, hapten-labeled probes or antibodies to DNA:RNA duplexes [3] can be used to detect hybridization or binding with AP-conjugated detection reagents and 1,2-dioxetane substrates. Chemiluminescent enzyme-linked oligonucleotide assay (ELONA) has been used to quantitate DNA aptamer binding to protein [6]. CDP-*Star* substrate is used in microplate-based assay systems for detection of viral RNA or DNA by immunodetection [17, 21]. In addition, detection of chemical:DNA adducts in mammalian tissue extracts has been demonstrated using CDP-*Star* substrate with Emerald-II enhancer [5].



Figure 37. Comparison of Chemiluminescent and Colorimetric Detection of TSH.

REFERENCES

- Ashihara, Y, H Saruta, S Ando, Y Kikuchi and Y Kasahara (1994). Sensitive chemiluminescent EIA for FT3 using CSPD® or AMPPD®, p. 321-324. In Bioluminescence and Chemiluminescence: Fundamentals and Applied Aspects. AK Campbell, LJ Kricka and PE Stanley, editors. John Wiley & Sons, Chichester, England.
- Bronstein, I, JC Voyta, GHG Thorpe, LJ Kricka and G Armstrong. (1989). Chemiluminescent assay of alkaline phosphatase applied in an ultrasensitive enzyme immunoassay of thyrotropin. *Clin. Chem.* 35:1,441-1,446.
- Carpenter, WR, TE Schutzbank, VJ Tevere, KR Tocyloski, N Dattagupta and KK Yeung (1993). A transcriptionally amplified DNA probe assay with ligatable probes and immunochemical detection. *Clin. Chem.* **39**:1,934-1,938.
- Charizopoulou, N, S Jansen, M Dorsch, F Stanke, JR Dorin, H-J Hedrich and B Tummler (2004). Instability of the insertional mutation in CftrTgH(neoim)Hgu cystic fibrosis mouse model. BMC Genetics (http://www.biomedcentral.com/1471-2156/5/6)
- Divi, RL, FA Beland, PP Fu, LS Von Tungeln, B Schoket, JE Camara, M Ghei, N Rothman, R Sinha and MC Poirier (2002). Highly sensitive chemiluminescence assay for benzo[a]pyrene-DNA adducts: validation by comparison with other methods, and use in human biomonitoring. *Carcinogenesis* 23(12):2043-2049.
- Drolet, DW, L Moon-McDermott and TS Romig (1996). An enzyme-linked oligonucleotide assay. Nature Biotech 14:1021-1025.
- Fimbel, S, H Dechaud, C Grenot, L Tabard, F Claustrat, R Bador and M Pugeat (1995). Use of non-radioactive labels for half-life measurement of sex hormone-binding globulin in the rabbit. Steroids 60:686-692.
- Fiscella, M, JW Perry, B Teng, M Bloom, C Zhang, K Leung, L Pukac, K Florence, A Concepcion, B Liu, Y Meng, C Chen, EC Elgin, P Kanakaraj, TE Kaufmann, J Porter, R Cibotti, Y Mei, J Zhou, G Chen, V Roschke, G Komatsoulis, B Mansfield, S Ruben, I Sanyal and T-S Migone (2003). TIP, a T-cell factor identified using high-throughput screening increases survival in a graft-versus-host disease model. Nature Biotech 21:302-307.
- Forhead, AJ, L Thomas, J Crabtree, N Hoggard, DS Gardner, DA Giussani and AL Fowden (2002). Plasma leptin concentration in fetal sheep during late gestation: Ontogeny and effect of glucocorticoids. *Endocrinology* 143(4):1166-1173.
- Ge, L, A Lupas, S Peraldi-Roux, S Spada and A Pluckthun (1995). A mouse Ig kappa domain of very unusual framework structure loses function when converted to the consensus. *J Biol Chem* 270:12446-12451.
- Heider, H and C Schroeder (1997). Focus luminescence assay: macroscopically visualized foci of human cytomegalovirus and varicella zoster virus infection. J Virological Methods 66:311-316.
- Ho, HH, MT Gilbert, DR Nussenzveig and MC Gershengorn (1999). Glycosylation is important for binding to human calcitonin receptors. *Biochemistry* 38:1866-1872.
- Jordan, T, L Walus, A Velickovic, T Last, S Doctrow and H Liu (1996). A competitive chemiluminescent enzyme-linked immunosorbent assay for the determination of RMP-7 in human blood. J Pharm Biom Anal 14:1653-1662.
- Kimura, T, A Rokuhara, A Matsumoto, S Yagi, E Tanaka, K Kiyosawa and N Maki (2003). New enzyme immunoassay for detection of Hepatitis B virus core antigen (HBcAg) and realtiaon between levels of HBcAg and HBV DNA. J Clin Microbiol 41(5):1901-1906.
- Legris, F, J Martel-Pelletier, J-P Pelletier, R Colman and A Adam (1994). An ultrasensitive chemiluminoenzyme immunoasasy for the quantification of human tissue kininogens: application to synovial membrane and cartilage. J Immunol Methods 168:111-121.
- Lehel, C, S Daniel-Issakani, M Brasseur and B Strulovici (1997). A chemiluminescent microtiter plate assay for sensitive detection of protein kinase activity. Anal Biochem 244:340-346.
- 17. Neisters, HGM, M Krajden, L Cork, M de Medina, M Hill, E Fries and ADME Osterhaus (2000). A multicenter study evaluation of the Digene Hybrid Capture II signal amplification technique for detection of Hepatitis B virus DNA in serum samples and testing of EUROHEP standards. J Clin Microbiol **38(6)**:2150-2155.
- Nishizono, I, S Lida, N Suzuki, H Kawada, H Murakami, Y Ashihara and M Okada (1991). Rapid and sensitive chemiluminescent enzyme immunoassay for measuring tumor markers. *Clin Chem* 37:1639-1644.
- Roffman, E and N Frenkel (1991). Chemiluminescent quantitation of lymphocyte surface antigens. J Immunolog Methods 138:129-131.
- Sanicola, M, C Hession, D Worley, P Carmillo, C Ehrenfels, L Walus, S Robinson, G Jaworski, H Wei, R Tizard, A Whitty, RB Pepinsky and RL Cate (1997). Glial cell line-derived neurotrophic factor-dependent RET activation can be mediated by two different cell-surface accessory proteins. Proc Natl Acad Sci USA 94:6238-6243.
- Schmitz, G and A Dotzauer (1998). Proof of hepatitis A virus negative-sense RNA by RNA/DNAhybrid detection: A method for specific detection of both viral negative- and positive-strand RNA species. Nuc Acids Res 26(22):5230-5232.
- Zhao, X, CQ Wenzel and JS Lam (2002). Nonradiolabeling assay for Waap, an essential sugar kinase involved in biosynthesis of core lipopolysaccharide of Pseudomonas aeruginosa. J Biol Chem 277(7):4722-4730.

Immunodetection Products

cAMP-Screen[®] System Cyclic AMP immunoassay

Description

cAMP-Screen[®] Immunoassay System enables ultrasensitive determination of cyclic AMP (cAMP) levels in cell lysates. This competitive immunoassay is formatted with maximum flexibility to permit either manual assay or automated high-throughput screening. cAMP-Screen assay utilizes the highly sensitive chemiluminescent alkaline phosphatase (AP) substrate CSPD[®] with Sapphire-II[™] luminescence enhancer. The Ready-to-Use substrate/enhancer reagent generates sustained-glow light emission that is measured 30 minutes after addition.

cAMP-Screen immunoassay is formulated to be compatible with automated high-throughput screening instrumentation (Figure 38). Cells are seeded into plates, cultured, and treated with test compounds as desired. Cell lysates are prepared in either the presence or absence of culture media. Lysates are incubated with a cAMP-AP conjugate and an anti-cAMP antibody in a coated microplate; the resulting immune complexes are captured in the plate. In samples without cAMP, all of the cAMP-AP conjugate is captured on the coated surface, resulting in a high signal. In the presence of cAMP, the amount of cAMP-AP conjugate captured decreases as a result of competition for binding with unlabeled cAMP, causing a reduced signal (see Figure 40); signal reduction is proportional to the amount of cAMP present in the cell lysate. After washing to remove unbound cAMP-AP, the chemiluminescent substrate is added, and the resulting glow signal is measured in a luminometer.

Advantages

Chemiluminescent detection with the cAMP-Screen assay provides the highest sensitivity of any commercially available cAMP assay. As few as 60 femtomoles of cAMP can be detected with cAMP-Screen 96-well system. cAMP-Screen system is available in both 96-well and 384-well formats. Both formats have a wide assay dynamic range with the 96-well version able to detect cAMP concentrations from 0.06 to 6,000 picomoles and the 384-well version able to detect cAMP concentrations from 0.2 to 200 picomoles (without the need for sample dilution or manipulations such as acetylation). This is especially helpful in cell-based assays, when measuring G_{s} - or G_{i} -coupled agonist stimulation and/or inhibition.

Assay precision is very high for cAMP-Screen immunoassay. Intra-assay precision for duplicate samples is typically 3% or less. To ensure optimal plate %CVs, it is imperative that the substrate/enhancer be allowed to reach the maximum glow signal. Once the substrate/enhancer reaches the glow signal, the plate can be read for hours with little or no degradation of the signal. This is useful in screening where several plates are compared to each other. In addition, the assay exhibits exceptionally low cross-reactivity with other adenosine-containing or cyclic nucleotides.

Applications

The cAMP-Screen assay system is designed for quantitation of cellular cAMP for functional assays of receptor activation. cAMP-Screen has been used with established cell lines for functional measurements with endogenous receptors, cell lines with exogenously expressed ligand receptors on the cell surface, primary cell cultures, and tissues [1] in response to treatment with the appropriate ligands. The cAMP-Screen assay system has been used for receptor characterization [5, 7], orphan receptor ligand identification [4], and the characterization of novel chimeric receptors [3]. In addition, cAMP-Screen assay system can be used for high throughput screening assays [2, 8] for compounds which stimulate or interfere with these signal transduction pathways.



Figure 38. cAMP-Screen® Assay: Effect of Forskolin and an Antagonist on cAMP Levels. cAMP-Screen assays were performed using the indicated concentration of forskolin to induce cAMP in the absence or presence of a forskolin antagonist. Measurements were performed on a microplate luminometer.



Figure 39. cAMP Immunoassay with cAMP-Screen® System.

cAMP-Screen[®] System, continued Cyclic AMP immunoassay

Product Configurations

cAMP-Screen® Immunoassay System is available in either a 96-well or 384-well plate format. Each kit format is available in several sizes to provide the ideal-size kit for evaluation, R&D, or screening.

Component	T1500	T1501	T1502	T1504
Plate Format	96-well	384-well	96-well	384-well
Assays per Kit	192	768	960	19,200
Assay/Lysis Buffer	25 mL	65 mL	2 x 65 mL	2 L
cAMP Standard	2 mL	5 mL	2 x 5 mL	100 mL
Anti-cAMP Antibody	14 mL	20 mL	2 x 35 mL	500 mL
cAMP-AP Conjugate	100 mL	250 mL	2 x 250 mL	5 mL
Conjugate Dilution Buffer	10 mL	25 mL	2 x 25 mL	500 mL
Wash Buffer	500 mL (1X)	1 L (1X)	2 x 1 L (1X)	5 x 1 L (5X)
CSPD [®] /Sapphire-II [™] RTU	25 mL	25 mL	2 x 65 mL	650 mL Substrate/Enhancer Solution
Precoated Plates	2 plates	2 plates	10 plates	50 plates

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
cAMP-Screen® 96-Well Immunoassay System	192 assays	T1500
cAMP-Screen® 384-Well Immunoassay System	960 assays	T1501
cAMP-Screen® 96-Well Immunoassay System	768 assays	T1502
cAMP-Screen® 384-Well Immunoassay System	19,200 assays	T1504
Assay/Lysis Buffer	65 mL	T2327
1X Wash Buffer	1 L	T2337
5X Wash Buffer	1 L	T2356

For Research Use Only. Not for use in diagnostic procedures.

REFERENCES

(For a complete cAMP-Screen[®] Immunoassay System reference list, please see the Applied Biosystems website.)

- Andre P, SM Delaney, T LaRocca, D Vincent, F DeGuzman, M Jurek, B Koller, DR Phillips and PB Conley (2003). P2Y12 regulates platelet adhesion/activation, thrombus growth, and thrombus stability in injured arteries. J Clin Invest 112:398406.
- Chiulli, AC, K Trompeter and M Palmer (2000). A novel high throughput chemiluminescent assay for the measurement of cellular cyclic adenosine monophosphate levels. J Biomol Screening 5(4):239–247.
- Gupte, J, G Cutler, JL Chen and H Tian (2004). Elucidation of signaling properties of vasopressin receptor-related receptor 1 by using the chimeric receptor approach. Proc Nat Acad Sci USA 101(6):1508-1513.
- He W, FJP Miao, DCH Lin, RT Schwandner, Z Wang, J Gao, JL Chen, H Tian and L Ling (2004). Citric acid cycle intermediates as ligands for orphan G-protein-coupled receptors. *Nature* 429:188-193.
- Inbe, H, S Watanabe, M Miyawaki, E Tanabe and JA Encinas (2004). Identification and characterization of a cell-surface receptor, P2Y15, for AMP and adenosine. J Biol Chem 279(19):19790-19799.
- Kolachala, V, V Asamoah, L Wang, S Srinivasan, D Merlin and SV Sitaraman (2005). Interferon-γ down-regulates adenosine 2b receptor-mediated signaling and short circuit current in the intestinal epithelia by inhibiting the expression of adenylate cyclase. J Biol Chem 280(6):4048-4057.
- Nickolls, SA, MI Cismowski, X Wang, M Wolff, PJ Conlon and RA Maki (2003). Molecular determinants of melanocortin 4 receptor ligand binding and MC4/MC3 receptor selectivity. J Pharmacol Exp Therapeutics 304(3):1217-1227.
- Vater, A, F Jarosch, K Buchner and S Klussmann (2003). Short bioactive Spiegelmers to migraine-associated calcitonin gene-related peptide rapidly identified by a novel approach: Tailored-SELEX. Nuc Acids Res 31(21):e130 (DOI: 10.1093/nar/gng130.





Figure 40. Dynamic Range of cAMP-Screen® Assay.

cAMP-Screen Direct[®] System Cyclic AMP immunoassay

Description

cAMP-Screen Direct[®] Immunoassay System enables ultrasensitive quantitation of cyclic AMP (cAMP) levels in cell lysates. With the cAMP-Screen Direct assay system, cells are plated directly in the antibody pre-coated cAMP-Screen Direct assay plates. The pre-coated plates supplied are clearbottom, enabling examination of cells prior to start of cAMP quantitation assay. This assay system eliminates the need for a separate cell culture plate. This competitive immunoassay is formatted with maximum flexibility to permit either manual assay or automated high-throughput screening. cAMP-Screen Direct assay system utilizes the highly sensitive chemiluminescent alkaline phosphatase (AP) substrate CSPD[®] with Sapphire-II[™] luminescence enhancer. The Ready-to-Use substrate/enhancer reagent generates sustained-glow light emission that is measured 30 minutes after addition.

cAMP-Screen Direct immunoassay system is formatted to be compatible with automated high-throughput screening instrumentation (Figure 41). Cells are seeded into pre-coated assay plates (clear-bottom, opaque white wells), cultured, and treated with test compounds as desired. Cell lysates are then prepared in either the presence or absence of culture media, and additional cAMP-Screen Direct system reagents are added directly to the lysate in the culture/assay plate. Lysates are incubated with cAMP-AP conjugate and anti-cAMP antibody, and the resulting immune complexes are captured on the antibody-coated surface of the plate. In samples without cAMP, all of the cAMP-AP conjugate is captured on the coated surface, resulting in a high signal. In the presence of cAMP, the amount of cAMP-AP conjugate captured decreases, as a result of competition for binding with unlabeled cAMP causing a reduced signal (Figure 41); signal is inversely proportional to the amount of cAMP present in the cell lysate. After washing to remove unbound cAMP-AP, the chemiluminescent substrate/enhancer is added, and the resulting signal is measured in a luminometer without reagent injection.

Advantages

The cAMP-Screen Direct Immunoassay System offers several advantages for high throughput screening applications compared to the original cAMP-Screen^{*} assay system. With cAMP-Screen Direct system, cells are grown directly in the cAMP quantitation assay plate without subsequent lysate transfer, eliminating the need for a separate culture plate, thereby reducing total assay costs, simplifying automation requirements for assay execution, and offering increased intra-assay accuracy. Chemiluminescent detection employed by cAMP-Screen Direct assay provides identical sensitivity to the cAMP-Screen system, the highest sensitivity of any commercially available cAMP assay. As with the cAMP-Screen Direct 96-well system. The cAMP-Screen Direct system is available in both 96-well and 384-well formats. Both the cAMP-Screen and cAMP-Screen Direct assay formats have a wide assay dynamic range with the 96-well version able to detect cAMP concentrations from 0.06 to 6,000 picomoles and the 384-well version able to detect cAMP concentrations from 0.2 to 200 picomoles (without the need for sample dilution or manipulations such as acetylation).

Assay precision is very high for cAMP-Screen Direct Immunoassay System. Intra-assay precision for duplicate samples is typically 3% or less. In addition, the assay exhibits exceptionally low cross-reactivity with other adenosinecontaining or cyclic nucleotides.

For cells requiring specialized growth surfaces, we do not recommend use of cAMP-Screen Direct microplates, or making surface modifications of the cAMP-Screen Direct microplates. In this case, the appropriate surface-modified plates should be used for cell growth, followed by lysate preparation and use of cAMP-Screen assay system.

Applications

The cAMP-Screen Direct assay system is designed for quantitation of cellular cAMP for functional assays of receptor activation. cAMP-Screen Direct system is designed for use with established cell lines for functional measurements with endogenous receptors, cell lines with exogenously expressed ligand receptors, and primary cell cultures [1] in response to treatment of cells with the appropriate ligands. cAMP-Screen Direct system is particularly suited for high throughput screening applications for identification of compounds which stimulate or interfere with receptor function and signal transduction.



Figure 41. Sensitivity of the cAMP-Screen Direct[®] System. cAMP assay performed with the cAMP-Screen Direct system. cAMP standards were used in a cAMP-Screen Direct assay plate following growth of different densities of HEK293 cells in assay wells. Detection sensitivity of exogenously added cAMP is unchanged following growth of cells in assay plate. Signal intensity differences result from basal cellular levels of cAMP. Signal was measured with the TR717^m microplate luminometer.



culture as desired

Add assay reagents

n and incubate 30 min

Add substrate solution



Place plate in luminometer and measure AP light emission

Figure 42. cAMP Immunoassay with cAMP-Screen Direct® System.

cAMP-Screen Direct[®] System, continued Cyclic AMP immunoassay

Product Configurations

cAMP-Screen Direct^{*} Immunoassay System is available in either a 96-well or 384-well plate format. Each kit format is available in several sizes to provide the ideal-size kit for evaluation, R&D, or screening.

Component	T1505	T1506	T1507	T1508
Plate Format	96-well	384-well	96-well	384-well
Assays per Kit	192	768	960	19,200
Assay/Lysis Buffer	25 mL	65 mL	2 x 65 mL	2 L
cAMP Standard	2 mL	5 mL	2 x 5 mL	100 mL
Anti-cAMP Antibody	14 mL	20 mL	2 x 35 mL	500 mL
cAMP-AP Conjugate	100 mL	250 mL	2 x 250 mL	5 mL
Conjugate Dilution Buffer	10 mL	25 mL	2 x 25 mL	500 mL
Wash Buffer	500 mL (1X)	1 L (1X)	2 x 1 L (1X)	5 x 1 L (5X)
CSPD®/Sapphire-II™ RTU	25 mL	25 mL	2 x 65 mL	650 mLSubstrate/Enhancer Solution
Precoated Plates	2 plates	2 plates	10 plates	50 plates

SIZE	CAT#
192 assays	T1505
768 assays	T1506
960 assays	T1507
19,200 assays	T1508
65 mL	T2327
1 L	T2337
1 L	T2356
	192 assays 768 assays 960 assays 19,200 assays 65 mL 1 L

For Research Use Only. Not for use in diagnostic procedures.

REFERENCES

 Zhong, H, L Belardinelli, T Maa, I Feoktistov, I Biaggioni and D Zeng (2004). A2B adenosine receptors increase cytokine release by bronchial smooth muscle cells. Am J Respir Cell Mol Biol 30:118-125.







Figure 43. Quantitation of cAMP Levels. Comparison of the cAMP-Screen Direct[®] System with the cAMP-Screen[®] System (page 28) for quantitation of cellular cAMP levels. For the cAMP-Screen Direct system, SK-N-MC cells were cultured for four days in cAMP-Screen Direct assay plates and then treated with compounds. The entire assay was performed in the same microplate. For the cAMP-Screen[®] system, SK-N-MC cells were cultured for four days in a standard microplate and treated with compounds. Cell lysate was transferred to a cAMP-Screen System plate and assayed. (A) NPY-mediated inhibition of isoproterenol (10 µM)-stimulated cAMP production. (B) NPY-mediated inhibition

ELISA-Light[™] System

Immunoassay

Description

ELISA-Light[™] Immunoassay System is used for rapid ultrasensitive antigen detection in enzyme immunoassays employing an alkaline phosphatase label. ELISA-Light system offers an excellent combination of high sensitivity and throughput for assays measured in a luminometer.

ELISA-Light System significantly extends the limits of detection of most immunoassays compared to colorimetric and isotopic methods. The kit contains reagents optimized for chemiluminescent ELISAs employing alkaline phosphatase labels. The kit is available with a choice of Ready-to-Use chemiluminescent substrate formulations containing CSPD* or CDP-*Star** substrates with a luminescence enhancer (Sapphire-II[™] or Emerald-II[™] enhancers). I-Block[™] blocking reagent, a highly purified casein screened to effectively prevent the nonspecific binding of alkaline phosphatase labeled reagents, antibodies and analytes to solid surfaces, is also included.

Applied Biosystems offers several Ready-to-Use substrate/enhancer formulations optimized for specific assays or luminometers. Formulations containing Sapphire-II enhancer are suitable for most assays performed in tube or microplate luminometers. Formulations containing Emerald-II enhancer are recommended for less-sensitive instruments that require more intense luminescent signals to detect low levels of enzyme. The specific formulation that will provide the highest sensitivity for a specific assay can vary. Therefore an ELISA sampler kit containing all four substrate/enhancer formulations is available to facilitate assay development.

Applications

The ELISA-Light system is used for a variety of protein detection assays, including standard immunoassay formats (see background section), phosphopeptide immunoassays for protein kinase activity quantitation [8,16,22], receptor binding assays [20], and viral foci imaging with immunodetection [11].

ELISA-Light system is also useful for nucleic acid detection applications in microplate-based detection assays. The system reagents are used for highly sensitive DNA probe capture assays, including quantitative detection of labeled PCR products [4] and RNA:DNA hybrids [3,17,21] and immunodetection of chemical-DNA adducts [5].

Product Configuration

ELISA-Light[™] System Standard Size

Capacity: 500-1.300 assays

Contents:

• 100 mL of Ready-to-Use chemiluminescent substrate with enhancer

100 mL of 10X Assay Buffer
7.5 g of I-Block[™] blocking reagent

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#			
ELISA-Light [™] Immunoassay System					
Sampler Kit with 25 mL of each formulation listed below (100 mL total)	sampler size	T1022			
with 100 mL CSPD® 0.4 mM Ready-to-Use with Sapphire-II™	standard size	T1023			
with 100 mL CSPD® 0.4 mM Ready-to-Use with Emerald-II™	standard size	T1024			
with 100 mL CDP- <i>Star</i> ® 0.4 mM Ready-to-Use with Sapphire-II™	standard size	T1025			
with 100 mL CDP- <i>Star</i> ® 0.4 mM Ready-to-Use with Emerald-II™	standard size	T1026			
ELISA-Light™ System Accessories available separately, see page 44					
AP-Labeled Secondary Antibody Conjugates					
AP Conjugate, Goat anti-Rabbit IgG	100 µL	T2191			
AP Conjugate, Goat Anti-Mouse IgG and IgM	100 µL	T2192			
I-Block [™] Blocking Reagent	30 g	T2015			
Avidx-AP™ Streptavidin-Alkaline Phosphatase Conjugate	1 mL	T2016			

For Research Use Only. Not for use in diagnostic procedures.

REFERENCES

For a complete list of immunoassay applications with the ELISA-Light™ system and CSPD® and CDP-Star substrates please see page 27.



Plus CSPD[®] Substrate and Enhancer

Light is Produced Proportional to the Amount of Antigen

Figure 44. Chemiluminescent Sandwich Immunoassay.

Introduction to Chemiluminescent Membrane-based Immunodetection

Chemiluminescent Immunodetection

Western blotting is a widely used technique for the immunodetection proteins. Protein samples are separated by polyacrylamide gel electrophoresis and electrophoretically transferred to either nitrocellulose, polyvinylidene fluoride (PVDF), or nylon membranes. The resulting blot is incubated with an antigen-specific monoclonal or polyclonal antibody to form a specific antigen-antibody complex. This complex is usually detected with a speciesspecific secondary antibody, bearing an enzyme label.

Visualization of an enzyme-labeled secondary antibody can be accomplished with colorimetric or chemiluminescent methods. For colorimetric methods, enzyme-catalyzed reactions result in the localized deposition of a colored product on the membrane. The accumulating product often interferes with enzyme activity, limiting sensitivity, and is difficult to remove from the membrane, preventing reuse of the blot. Furthermore, the color usually fades and is difficult to photograph. In addition, colorimetric visualization methods are not quantitative. For chemiluminescent detection, a chemiluminescent substrate solution is added, and the resulting localized light emission is imaged either on X-ray film or with a CCD-based imaging platfrom.

Advantages of Chemiluminescence

Western-Light[™] and Western-*Star[™]* chemiluminescent Western blotting systems offer highly sensitive, rapid, nonisotopic protein detection [9,12]. Detection is performed on nitrocellulose, nylon, or PVDF membranes. With nitrocellulose membranes, an enhancing agent, Nitro-Block[™] or Nitro-Block-II[™] enhancer, is added to provide maximum sensitivity. Detection of alkaline phosphatase-labeled immunocomplexes is



Figure 45. Western-Light[™] and Western-Star[™] Systems Immunodetection Procedure.

Panel A. Nitrocellulose



Panel B. PVDF



(A) Human brain extract (10, 5, 2.5, 1.25 µg) was electrophoretically separated and transferred to Tropifluor[™] PVDF and nitrocellulose membrane.
(B) Blots were detected with the Western-Light[™] system and CSPD[®] substrate or the Western-Star[®] system and CDP-Star[®] substrate, with enhancer for nitrocellulose and without enhancer for PVDF. Blots were imaged on Kodak XAR-5 film for 30 sec, 20 min after substrate incubation.

performed with the CSPD[®] substrate (Western-Light System; see page 36) or CDP-*Star*[®] substrate (Western-*Star*[®] System; see page 34). Following incubation with substrate, blots are imaged immediately on X-ray film, with exposure times ranging from 30 seconds to 45 minutes with CSPD substrate, and from 1 second to 30 minutes with CDP-*Star* substrate, providing a hard copy of results. CDP-*Star* substrate provides a five- to ten-fold higher signal intensity than CSPD substrate, enabling much shorter exposure times [10]. Detection of alkaline phosphatase labels with 1,2-dioxetane substrates, unlike detection of horseradish peroxidase with enhanced luminol systems, produces a chemiluminescent signal with glow kinetics, enabling multiple exposures over many hours. Glow kinetics, combined with the high intensity signal generated with CDP-*Star* substrate, makes the Western-*Star* detection system ideal for chemiluminescent phosphor screen and CCD camera imaging [17,18].

Biotinylated Secondary Antibodies and Proteins

Immunoblotting can be performed with a biotinylated secondary antibody followed by streptavidin-alkaline phosphatase conjugate and CSPD or CDP-*Star* substrate, and often results in increased signal intensity with shorter exposure times compared to that obtained with alkaline phosphatase conjugated secondary antibody. Tropix* chemiluminescent substrates have been used for highly sensitive detection of biotin-labeled proteins [12]. Other applications include detection of immunoprecipitated biotin labeled proteins [11,13,15] and biotinylated peptide probes.

Detection of Phosphorylated Proteins

Detection of phosphotyrosine or phosphoserine and phosphothreonine-containing proteins on immunoblots has been widely used for the analysis of protein kinase substrates and activities [14,16]. Chemiluminescent immunoblot detection of phosphorylated proteins can be performed with the Western-*Star* or Western-Light chemiluminescent detection systems in conjunction with the appropriate anti-phospho-epit-ope primary antibody. A non-casein blocker should be substituted for the I-Block[™] blocking reagent, since it is casein and contains phosphorylated amino acids.

Detection System	Western- <i>Star</i> ™ see page 392	Western-Light™ see page 394
Key Attribute	More sensitive than Western-Light™	Established system with good sensitivity
Substrate	CDP-Star [®]	CSPD®
Detection Scheme	Alkaline Phosphatase-2° Ab Conjugate	Alkaline Phosphatase-2° Ab Conjugate
Film Exposure Time	1 sec - 30 min	2 min - 45 min
Special Features	Highest signal intensity— ideal for CCD camera and phosphor screen imaging	Original chemiluminescent Western system

Figure 47. Selection Guide for Western Blot Detection Systems. Applied Biosystems offers two Western blot detection systems.

REFERENCES

See Western-Light[™] and Western-*Star*[™] Immunodetection System sections (pages 35 and 37) for references.

Figure 46. Detection with CDP-Star® Substrate vs. CSPD® Substrate.

Immunodetection Products

Western-Star[™] System

Immunodetection

Description

Western-Star^{**} Immunodetection System is a highly sensitive chemiluminescent immunodetection system that provides speed, flexibility, control of film exposure time and the ability to use alternative imaging systems. Western-Star system incorporates CDP-Star^{*} substrate for detection of secondary antibody-alkaline phosphatase conjugates.

Advantages

Western-Star system generates a high intensity chemiluminescent signal that persists from hours to days, depending on the membrane type. CDP-Star substrate provides a five- to tenfold higher signal intensity than CSPD[®] substrate. Images can be generated immediately on X-ray or instant film to provide permanent, hard-copy results. In addition, the high signal intensity and long-lived signal generated with CDP-Star substrate provides an ideal quantitative detection system for chemiluminescent phosphor screen and CCD camera imaging systems.

Applications

Western-Star Immunodetection System is used for highly sensitive immunoblot detection of proteins, biotinylated proteins [1], and phosphoproteins in protein extracts from many sources, including cell cultures and tissues. These reagents are compatible with multiple types of membranes, including PVDF, nylon and nitrocellulose. Detection is performed directly with AP-labeled secondary antibodies, or with indirect detection of biotinylated antibodies or biotinylated proteins with a streptavidin-AP conjugate (Avidx-AP[™], see Reagents and Accessories section, page 44). Western-Star System has been used for confirmation of protein knockdown for siRNA gene expression regulation experiments [4].



Figure 48. Western-Star™ System vs. Enhanced Luminol Detection. Human brain extract was electrophoretically separated and transferred to PVDF. Blots were incubated with a monoclonal anti-actin, followed by Western-Star™ detection or with a horseradish peroxidase (HRP)-conjugated secondary antibody, followed by enhanced luminol substrate detection. Blots were imaged on Kodak XAR-5 Xray film. T = time after substrate incubation and E = exposure time.



Figure 49. Immunodetection with Western-Star[™] System.
Western-Star[™] System, continued

Immunodetect

Product Configuration

Western-Star[™] System Standard Size

T1048 with Goat Anti-Rabbit IgG AP Conjugate T1046

with Goat Anti-Mouse IgG+IgM AP Conjugate Capacity: 30 membrane blots (10 cm x 10 cm)

Contents:

• 100 mL CDP-Star® 0.25 mM Ready-to-Use substrate

- secondary antibody AP conjugate
- 30 g l-Block[™] blocking reagent
- 5 mL Nitro-Block-II $^{\scriptscriptstyle \rm M}$ chemiluminescence enhancer
- 150 mL 10X Assay Buffer concentrate
- 30 development folders (14 cm x 19 cm)

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Western- <i>Star</i> ™ System		
with Goat Anti-Mouse IgG+IgM AP Conjugate	standard size	T1046
with Goat Anti-Rabbit IgG AP Conjugate	standard size	T1048
Western- <i>Star</i> ™ System Accessories (availa	ible separately, see page	e 44)
Goat Anti-Rabbit IgG AP Conjugate	100 µL	T2191
Goat Anti-Mouse IgG+IgM AP Conjugate	100 µL	T2192
I-Block™ Blocking Reagent	30 g	T2015
Tropifluor™ PVDF Membrane (15 cm x 15 cm)– 5 Membranes/Pack	– 1 pack	T2234
Development Folders (14 cm x 19 cm)— 30 Folders/Pack	1 pack	T2258
or Research Lise Only. Not for use in diagnostic procedu	res	

For Research Use Only. Not for use in diagnostic procedures.

REFERENCES

(For a complete Western-Star[™] immunodetection system reference list, please see the Applied Biosystems website.)

- Araga, S, L Xu, K Nakashima, M Villain and JE Blalock (2000). A peptide vaccine that prevents experimental autoimmune myasthenia gravis by specifically blocking T cell help. FASEB J 14:185-196.
- Elmariah, SB, MA Crumling, TD Parsons and RJ Balice-Gordon (2004). Postsynaptic TrkB-mediated signaling modulates excitatory and inhibitory neurotransmitter receptor clustering at hippocampal synapses. J Neurosci 24(10):2380-2393.
- Fahrenkrog, B, U Sauder and U Aebi (2004). The S. cerevisiae HtrA-like protein Nma111p is a nuclear serine protease that mediates yeast apoptosis. J Cell Sci 117:115-126.
- Grunweller, A, C Gillen, VA Erdmann and J Kurreck (2003). Cellular uptake and localization of a Cy3-labeled siRNA specific for the serine/threonine kinase Pim-1. Oligonucleotides 13:345-352.
- Kobayashi, S, R Nantz, T Kitamura, R Higashikubo and N Horikoshi (2005). Combined inhibition of extracellular signal-regulated kinases and HSP90 sensitizes human colon carcinoma cells to ionizing radiation. Oncogene 24:3011-3019.
- Schaller, O, R Fatzer, M Stack, J Clark, W Cooley, K Biffiger, S Egli, M Doherr, M Vandevelde, D Heim, B Oesch and M Moser (1999). Validation of a western immunoblotting procedure for bovine PrPSc detection and its use as a rapid surveillance method for the diagnosis of bovine spongform encephalopathy (BSE). Acta Neuropathol **98**:437-443.
- Tian, J, A He, AG Lawrence, P Liu, N Watson, AJ Sinskey and J Stubbe (2005). Analysis of transient polyhydroxybutyrate production in Wautersia eutropha H16 by quantitative western analysis and transmission electron microscopy. J Bacteriol 187(11):3825-3832.
- Walz, A, S Park, JP Slovin, J Ludwig-Muller, YS Momonoki and JD Cohen (2002). A gene encoding a protein modified by the phytohormone indoleacetic acid. Proc Natl Acad Sci USA 99(3):1718-1723.

Immunodetection Products

Western-Light[™] System

Immunodetection

Description

Western-Light[™] System is a versatile non-isotopic, chemiluminescent immunodetection system that provides speed, sensitivity, flexibility, and control of film exposures. Western-Light System incorporates high-performance CSPD[®] substrate for detection of a secondary antibody-alkaline phosphatase conjugate.

Advantages

Western-Light System detection generates a chemiluminescent signal that persists from hours to days, depending on the type of membrane used. This enables an image to be generated on X-ray or instant film to provide a permanent hard copy of results. Multiple images may be rapidly and easily acquired with excellent control over film exposures. Detection of sub-pico-gram levels of protein is easily attainable (the ultimate sensitivity depends on primary antibody quality). New users of 1,2-dioxetane substrates for Western blot detection are recommended to use the Western-Star[™] immuno-detection system for maximum signal intensity and sensitivity.

Applications

Western-Light immunodetection system is used for highly sensitive immunoblot detection of proteins, biotinylated proteins, and phosphoproteins in protein extracts from many sources, including cell cultures and tissues. These reagents are compatible with multiple types of membranes, including PVDF, nylon and nitrocellulose. Detection is performed directly with AP-labeled secondary antibodies, or with indirect detection of biotinylated antibodies or biotinylated proteins with a streptavidin-AP conjugate (Avidx-AP[™], see Reagents and Accessories section, page 44).

Product Configuration

Western-Light[™] System Standard Size

T1047 with Goat Anti-Rabbit IgG AP Conjugate

T1045 with Goat Anti-Mouse IgG+IgM AP Conjugate Capacity: 30 membrane blots (10 cm x 10 cm) Contents:

- 100 mL CSPD® 0.25 mM Ready-to-Use substrate
- secondary antibody AP conjugate
- 30 g I-Block[™] blocking reagent
- 5 mL Nitro-Block[™] chemiluminescence enhancer
- 150 mL 10X Assay Buffer concentrate
- 30 development folders (14 cm x 19 cm)

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#
Western-Light [™] System		
with Goat Anti-Mouse IgG+IgM AP Conjugate	standard size	T1045
with Goat Anti-Rabbit IgG AP Conjugate	standard size	T1047
Western-Light™ System Accessories (available	e separately, see page	44)
Goat Anti-Rabbit IgG-AP Conjugate	100 µL	T2191
Goat Anti-Mouse IgG+IgM AP Conjugate	100 µL	T2192
I-Block™ Blocking Reagent	30 g	T2015
Tropifluor™ PVDF Membrane (15 cm x 15 cm)– 5 Membranes/Pack	– 1 pack	T2234
Development Folders (14 cm x 19 cm)— 30 Folders/Pack	1 pack	T2258



Figure 50. Immunodetection with Western-Light[™] System.

Western-Light[™] System, continued

Panel A. Tropifluor™ PVDF



15 7.5 3.8 1.9 0.9 0.5 0.2 0.1 .06 .03 .015 .007 .004 (ng)

Panel B. Nitrocellulose (with Nitro-Block™ enhancer)



(ng)

Figure 51. Western-Light[™] Immunodetection with Tropifluor[™] PVDF vs. Nitrocellulose. Purified human transferrin was electrophoretically separated and transferred to membrane, incubated with rabbit anti-transferrin polyclonal antibody, detected using Western-Light[™] system, and exposed immediately on X-ray film for 15 minutes. As shown, less than 3.6 pg of protein is detectable on both Tropifluor PVDF membrane (without Nitro-Block enhancer treatment) and nitrocellulose membrane (with Nitro-Block[™] enhancer treatment).

REFERENCES

(For a complete Western-Light[™] immunodetection system reference list, please see the Applied Biosystems website.)

- Bronstein, I, JC Voyta, OJ Murphy, L Bresnick and LJ Kricka (1992). Improved chemiluminescent western blotting procedure. *BioTechniques* 12:748-753.
- 10. Bronstein, I, CEM Olesen, CS Martin, G Schneider, B Edwards, A Sparks and JC Voyta (1994). Chemiluminescent detection of DNA and protein with CDP™ and CDP-Star® 1,2-dioxetane enzyme substrates, p. 269-272. In Bioluminescence and Chemiluminescence: Fundamentals and Applied Aspects. Campbell, AK, Kricka, LJ, and Stanley, PE, eds. John Wiley, Chichester, England.
- 11. Chang, AC, DR Salomon, S Wadsworth, M-JP Hong, CF Mojcik, S Otto, EM Shevach and JE Coligan (1995). a3b1 and a6b1 integrins mediate laminin/merosin binding and function as costimulatory molecules for human thymocyte proliferation. J Immunol 154:500-510.
- Gillespie, PG and AJ Hudspeth (1991). Chemiluminescence detection of proteins from single cells. Proc Natl Acad Sci USA 88:2563-2567.
- Gillespie, SKH and S Wasserman (1994). Dorsal, a Drosophila Rel-like protein, is phosphorylated upon activation of the transmembrane protein Toll. *Mol Cell Biol* 14:3559-3568.
- Kamps, MP and BM Sefton (1988). Identification of novel polypeptide substrates of the v-src, v-yes, v-fps, v-ros, and v-erb-B oncogenic tyrosine protein kinases utilizing antisera against phospho-tyrosine. Oncogene 2:305-315.
- Lantz, LM and KL Holmes (1995). An improved nonradioactive cell surface labeling technique for immunoprecipitation. *BioTechniques* 18:56-62.
- Morla, AO and JYJ Wang (1986). Protein tyrosine phosphorylation in the cell cycle of BALB/c3T3 fibroblasts. Proc Natl Acad Sci USA 83:8191-8195.
- Nguyen, Q, W Stubblebine, C Ragsdale, F Witney, MC Lorence and D Heffelfinger (1993). Detection of chemiluminescent western blot by reusable phosphor imaging screen. J NIH Res 5:79-80.
- Nguyen, Q and DM Heffelfinger (1995). Imaging and quantitation of chemiluminescence using photoexcitable storage phosphor screen. Anal. Biochem. 226:59-67.



4 Nucleic Acid Membrane-Based Detection Products

Introduction	39
Southern-Star™ System	41
Southern-Light [™] System	42

Introduction

Chemiluminescent Nucleic Acid Detection

Traditional nucleic acid detection methods employ radioisotopes such as ³²P, which have several disadvantages. Exposures often require several days and labeled probes must be prepared frequently because of radiolytic decay. Furthermore, the radiation safety precautions required are inconvenient and disposal costs are continually increasing.

Applied Biosystems has developed non-isotopic chemiluminescent technologies for detecting nucleic acids utilizing enzyme conjugates, such as alkaline phosphatase conjugates or alkaline phosphatase-labeled DNA probes. As shown, alkaline phosphatase activates the chemiluminescent CDP-*Star*^{*} substrate to produce light. A distinct band is formed since dephosphorylated CDP-*Star* substrate has a strong affinity for hydrophobic sites of the membrane. The localized persistent "glow" of light can be imaged on X-ray or instant film or with CCD camera-based imaging systems. Film exposures of 5 to 60 minutes are easily performed and sub-picogram quantities of labeled DNA can be detected.

Indirect Detection

Nucleic acid detection with CDP-Star or CSPD® chemiluminescent substrate is performed with either a direct or indirect alkaline phosphatase label with Southern-Star[™] or Southern-Light[™] systems. Indirect detection of biotin-labeled probes is accomplished with a streptavidin-alkaline phosphatase conjugate (Avidx-AP[™], see page 44). This conjugate is optimized to exhibit minimum nonspecific binding and superior sensitivity. Detection of other hapten labels, such as digoxigenin, fluorescein or DNP, with the appropriate antibody-alkaline phosphatase conjugate, can also be performed successfully.

Alkaline Phosphatase-oligonucleotide Conjugates

The use of oligonucleotides covalently derivatized with alkaline phosphatase is currently the most effective method for exploiting the inherent sensitivity of dioxetane chemiluminescence. This direct method offers high sensitivity, low background, short hybridizations, and relatively rapid film exposures [20]. A single copy gene can be detected in 0.25 µg of human genomic DNA using an alkaline phosphatase-derivatized oligonucleotide [6,17]. Use of CDP-*Star* substrate reduces exposure times from greater than 12 hours to less than 2 hours [17] in many applications.



Figure 52. Detection of Biotin-Labeled DNA with CDP-Star® Substrate.



Figure 53. Nucleic Acid Labeling Methods for Chemiluminescent Detection.

1,2-Dioxetane chemiluminescent substrates, including CDP-*Star* and CSPD substrates, have also been successfully utilized in sensitive, non-radioactive detection in a variety of membrane-based nucleic acid detection methods, including:

- Northern blotting [3,13,15,18]
- RNA knockdown detection [10]
- capped mRNA detection [14]
- DNA macroarray hybridization for DNA analysis [12] and gene expression analysis [19]
- direct detection of PCR-amplified DNA [2,9]
- SSCP analysis [21], differential display [1]
- plaque hybridization
- gel shift assays [8]
- RNase protection assays [11]
- RNA probe detection of RNA binding proteins [22]

1,2-Dioxetane chemiluminescent substrates have also been employed in a chemiluminescent in situ hybridization assay with DIG-labeled probes [16].

Nucleic Acid Membrane-Based Detection Products

Introduction, continued

Detection System	Key Attribute	Substrate	Detection Scheme	Film Exposure Time	Special Features
Southern-Star™ for detection of biotin- labeled probes see page 41	Most sensitive system for biotin detection	CDP-Star®	Streptavidin- AP Conjugate (AvidxAP™)	1 min–2 hr	Compared to ³² P, reduces overnight film exposures to less than 2 hr. Applications include: Southern and Northern blotting, colony/plaque screening, and direct detection of AP-oligo probes
Southern-Light™ for detection of biotin- labeled probes see page 42	Established system with good sensitivity	CSPD®	Streptavidin- AP Conjugate (Avidx-AP [™])	5 min–10 hr	Original system for chemiluminescent Southern blotting. Applications include Southern and Northern blotting, colony/ plaque screening, and direct detection of AP-oligo probes.

Figure 54. Selection Guide for Southern/Northern Blot Detection Systems. Applied Biosystems offers several Southern blotting systems. High sensitivity in non-radioactive Southern blots is achieved by efficient probe labeling and following proper blocking and wash steps as described in protocols. The chart below summarizes the differences and relative merits of each system. Southern-Star[®] system, which provides the highest signal intensity and sensitivity, is recommended to new users for all applications.

REFERENCES

(For a complete Southern-Light[™] and Southern-Star[™] nucleic acid detection reference lists, please see the Applied Biosystem website.)

- An, G, G Luo, RW Veltri and SM O'Hara (1996). Sensitive, nonradioactive differential display method using chemiluminescent detection. *BioTechniques* 20:342-346.
- Bettens, F., W.J. Pichler, and A.L. de Weck (1991). Incorporation of biotinylated nucleotides for the quantification of PCR-amplified HIV-1 DNA by chemiluminescence. *Eur J Clin Chem Clin Biochem* 29:685-688.
- Borst, A, MT Raimer, DW Warnock, CJ Morrison and BA Arthington-Skaggs (2005). Rapid acquisition of stable azole resistance by Candida glabrata isolates obtained before the clinical introduction of Fluconazole. Antimicrobial Agents and Chemotherapy 49(2):783-787.
- Bronstein, I, JC Voyta, KG Lazzari, OJ Murphy, B Edwards and LJ Kricka (1990). Rapid and sensitive detection of DNA in Southern blots with chemiluminescence. *BioTechniques* 8:310-313.
- Bronstein, I, CEM Olesen, CS Martin, G Schneider, B Edwards, A Sparks and JC Voyta (1994). Chemiluminescent detection of DNA and protein with CDP and CDP-Star 1,2-dioxetane enzyme substrates, p. 269-272. In AK Campbell et al. (eds.), *Bioluminescence and Chemiluminescence: Fundamentals and Applied Aspects*. John Wiley, Chichester, England.
- Cate, R, C Ehrenfels, M Wysk, R Tizard, JC Voyta, O Murphy III and I Bronstein (1991). Genomic Southern analysis with alkaline phosphatase conjugated oligonucleotide probes and the chemiluminescent substrate AMPPD. GATA 8:102-106.
- Davis, MA, DD Hancock, TE Besser and DR Call (2003). Evaluation of pulsed-field gel electrophoresis as a tool for determining the degree of genetic relatedness between strains of Escherichia coli 0157:H7. J Clin Microbiol 41(5):1843-1849.
- Düring, K (1993). Non-radioactive detection methods for nucleic acids separated by electrophoresis. J. Chromatography 618:105-131.
- Fong, T-T, DW Griffin and EK Lipp (2005). Molecular assays for targeting human and bovine enteric viruses in coastal waters and their application for library-independent source tracking. Applied and Environmental Microbiology 71(4):2070-2078.
- Grunweller, A, C Gillen, VA Erdmann and J Kurreck (2003). Cellular uptake and localization of a Cy3-labeled siRNA specific for the serine/threonine kinase Pim-1. *Oligonucleotides* 13:345-352.
- 11. Jaskula, JC (1996). Ultrasensitive nonisotopic mRNA detection. Applications in Molecular Biology (Biomedical Products suppl.) **November:**10-

- Jenkins, BD, GF Steward, SM Short, BB Ward and JP Zehr (2004). Fingerprinting diazotroph communities in the Chesapeake Bay by using a DNA microarray. Applied and Environmental Microbiology 70(3):1767-1776.
- Löw, R and T Rausch (1994). Sensitive, nonradioactive northern blots using alkaline transfer of total RNA and PCR-amplified biotinylated probes. *BioTechniques* 17:1026-1030.
- 14. Mizumoto, H, M Tatsuta, M Kaido, K Mise and T Okuno (2003). Cap-independent translational enhancement by the 3' untranslated region of Red Clover Necrotic Mosaic Virus RNA1. J Virol 77(22):12113-12121.
- Murata, M, Y Okimura, K lida, M Matsumoto, H Sowa, H Kaji, M Kojima, K Kangawa and K Chihara (2002). Ghrelin modulates the downstream molecules of insulin signaling in hepatoma cells. *J Biol Chem* 277(7):5667-5674.
- 16. Musiani, M, A Roda, M Zerbini, P Pasini, G Gentilomi, G Gallinella and S Venturoli (1996). Chemiluminescent in situ hybridization for the detection of cytomegalovirus DNA. Am J Pathology 148(4):1105-11
- Price, DC (1996). Chemiluminescent substrates for detection of restriction fragment length polymorphism. Science and Justice 36:275-282.
- Suar, M, JR van der Meer, K Lawlor, C Holliger and R Lal (2004). Dynamics of multiple lin gene expression in Sphingomonas paucimobilis B90A in response to different hexachlorocyclohexane isomers. Applied and Environmental Microbiology **70(11)**:6650-6656.
- Tadlock, L, Y Yamagiwa, C Marienfeld and T Patel (2003). Double-stranded RNA activates a p38 MAPK-dependent cell survival program in biliary epithelia. Am J Physiol Gastrointest Liver Physiol 284:G924-G932.
- Vary, CPH, M Carmody, R LeBlanc, T Hayes, C Rundell and L Keilson (1996). Allele-specific hybridization of lipoprotein lipase and factor-V Leiden missense mutations with direct label alkaline phosphatase-conjugated oligonucleotide probes. *Genetic Analysis: Biomolecular Engineering* 13:59-65.
- Weiss, N, I Eggersdorfer and C Keller (1996). Multiplex-PCR-based single-strand conformation polymorphism protocol for simultaneous analysis of up to five fragments of the low-density-lipoprotein receptor gene. *BioTechniques* 20:421-429.
- Werner, R, HP Mühlbach and MC Guitton (1995). Isolation of viroid-RNA-binding proteins from an expression library with nonradioactive-labeled RNA probes. *BioTechniques* 19:218-222.

Southern-Star[™] System

Nucleic acid detection system for biotin-labeled probes

Description

Southern-Star[™] nucleic acid detection system is a non-isotopic, chemiluminescence-based products for detection of biotin-labeled nucleic acids in Southern and Northern blotting with high sensitivity and ease of use. Southern-Star detection of biotin-labeled nucleic acid probes incorporates CDP-Star* substrate and Avidx-AP[™] streptavidin-alkaline phosphatase conjugate. Southern-Star system reagents are compatible with other labeling systems, including digoxigenin, fluorescein and DNP with the appropriate antibody-alkaline phosphatase conjugates, as well as alkaline phosphataselabeled oligonucleotides.

Advantages

CDP-Star substrates the chemiluminescent substrate included in Southern-Star system, generates a high intensity chemiluminescent signal that may persist for 2-3 days on nylon membrane, allowing multiple images to be acquired. CDP-Star substrate provides up to a ten-fold higher signal intensity than CSPD[®] substrate (included with the Southern-Light[™] system (Figure 55), and reaches maximum light level more rapidly. Images can be obtained immediately on X-ray or instant film to provide permanent, hard-copy records. In addition, the high signal intensity and long duration produced by CDP-Star substrate makes Southern-Star system detection ideal for chemiluminescent phosphor screen and CCD camera imaging systems.

Applications

Southern-*Star* nucleic acid detection system is used for highly sensitive detection of DNA for Southern blotting applications [7], including DNA detection and gene expression analysis on DNA microarrays [12,19] and direct detection of labeled DNA on membranes. It can also be used for additional nucleic acid detection applications, including colony and plaque screening, Northern blotting [3,15], RNA knockdown analysis [10], capped mRNA detection [14], RNase protection, and gel-shift assays.



Figure 55. Southern-Star[™] Nucleic Acid Detection with CDP-Star[®] vs. CSPD[®] Substrates. Yeast S. cerevisiae genomic DNA (EcoRI and BgI II digested) was electrophoretically separated and transferred to nylon membrane. Blots were hybridized with a single copy RPB 1 probe and detected with Avidx-AP[™] conjugate and 0.25 mM CSPD[®] (Panel A) or CDP-Star[®] (Panel B) substrate, according to the Southern-Star[™] system protocol. Blots were incubated with substrate for 75 min and then exposed to Xray film for 4 min.

Product Configuration

Southern-Star[™] System Standard Size

T1040

- Capacity: 30 membrane blots (10 cm x 10 cm)
- Contents: • 100 mL CDP-Star® 0.25 mM Ready-to-Use substrate
- 120 µL Avidx-AP[™] streptavidin-alkaline phosphatase conjugate
- 30 g l-Block[™] blocking reagent
- 150 mL 10X Assay Buffer concentrate
- 30 development folders (14 cm x 19 cm)

T1039

T1040 System with Tropilon-Plus™ Membrane

Southern-Star[™] System Small Size

T1042

- Capacity: 10 membrane blots (10 cm x 10 cm) Contents:
- 30 mL CDP-Star[®] 0.25 mM Ready-to-Use substrate
- 40 µL Avidx-AP[™] streptavidin-alkaline phosphatase conjugate
- 7.5 g l-Block[™] blocking reagent
- 50 mL 10X Assay Buffer concentrate
- 10 development folders (14 cm x 19 cm)

APPLIED BIOSYSTEMS ORDER INF	O SIZE	CAT#		
Southern-Star™ Nucleic Acid Detection System for Biotin-Labeled Probes				
with Tropilon-Plus™ Membrane standard size, 30 blots T1039				
	standard size, 30 blots	T1040		
	small size, 10 blots	T1042		
Southern- <i>Star</i> ™ System Accessories (a	Southern- <i>Star</i> ™ System Accessories (available separately, see page 44)			
I-Block [™] Blocking Reagent	30 g	T2015		
Avid <i>x</i> -AP [™] Streptavidin-AP Conjugate	1 mL	T2016		
Tropilon-Plus™ Nylon Membrane	1 roll	T2232		

(30 cm x 200 cm)



Nucleic Acid Membrane-Based Detection Products

Southern-Light[™] System

Nucleic acid detection system for biotin-labeled probes

Description

Southern-Light[™] nucleic acid detection system is a non-isotopic, chemiluminescence-based product for detection of biotin-labeled nucleic acids in Southern and Northern blotting with high sensitivity and ease of use. Southern-Light detection of biotin-labeled nucleic acid probes incorporates CSPD[®] substrate and Avidx-AP[™] streptavidin-alkaline phosphatase conjugate. Southern-Light system reagents are compatible with other labeling systems, including digoxigenin, fluorescein and DNP with the appropriate antibody-alkaline phosphatase conjugates, as well as alkaline phosphataselabeled oligonucleotides.

Advantages

The convenience, speed, sensitivity, and versatility of chemiluminescent Southern blotting makes radioactive methods obsolete. Each component of the Southern-Light system is optimized for maximum signal intensity and low background, including CSPD high-performance chemiluminescent alkaline phosphatase substrate, Avidx-AP streptavidin-alkaline phosphatase conjugate, and I-Block[™] blocking reagent. The Southern-Light system procedure is able to detect less than 1 pg of target DNA. Most exposures require less than 60 minutes and are performed with standard X-ray or instant film. Light emission persists for days enabling multiple re-exposures. Additionally, membranes are easily stripped and reprobed.

Applications

Southern-Light nucleic acid detection system is used for highly sensitive detection of DNA for Southern blotting applications and direct detection of labeled DNA on membranes. This system can also be used for additional nucleic acid detection applications, including colony and plaque screening, Northern blotting, RNase protection, and gel-shift assays. The Southern-Light system is also useful for detection of biotinylated proteins in non-isotopic immunoprecipitation procedures.



Figure 57. Nucleic Acid Detection with Fluorescein-Labeled Probe. Southern-Light™ system detection of *S. cerevisiae* single copy RNA polymerase gene RPB1 using a fluorescein-labeled probe. Lanes 1 through 5 contain 2 μg, 400 ng, 80 ng, 16 ng, and 3.2 ng total DNA, respectively.

Add Avidx-AP[™] Conjugate Add CSPD[®] Substrate CSPD[®] Light Light Light Light Light Light Light Film image F

Figure 58. Nucleic Acid Detection with Southern-Light[™] System.

Product Configuration

Southern-Light[™] System Standard Size

T1037

Capacity: 30 membrane blots (10 cm x 10 cm) Contents:

- 100 mL CSPD® 0.25 mM Ready-to-Use substrate
- 120 µL Avidx-AP[™] streptavidin-alkaline phosphatase conjugate
- 30 g I-Block[™] blocking reagent
- 150 mL 10X Assay Buffer concentrate
 30 development folders (14 cm x 19 cm)

T1038 T1037 System with Tropilon-Plus™ Membrane

Southern-Light[™] System Small Size

T1041

Capacity: 10 membrane blots (10 cm x 10 cm) Contents:

- 30 mL CSPD[®] 0.25 mM Ready-to-Use substrate
- 40 μL Avidx-AP $^{\rm \tiny m}$ streptavidin-alkaline phosphatase conjugate
- 7.5 g l-Block[™] blocking reagent
- 50 mL 10X Assay Buffer concentrate
- 10 development folders (14 cm x 19 cm)

APPLIED BIOSYSTEMS ORDER INFO) SIZE	CAT#	
Southern-Light [™] Nucleic Acid Detection System for Biotin-Labeled Probes			
	standard size, 30 blots	T1037	
with Tropilon-Plus™ Membrane	standard size, 30 blots	T1038	
	small size, 10 blots	T1041	
Southern-Light™ Accessories (availabl	e separately, see page 44)	
I-Block™ Blocking Reagent	30 g	T2015	
Avidx-AP [™] Streptavidin-AP Conjugate	1 mL	T2016	
Tropilon-Plus™ Nylon Membrane			
(30 cm x 200 cm)	1 roll	T2232	

5 Reagents and Accessories for Chemiluminescence

Reagents and Accessories for Chemiluminescence

AP-Labeled Secondary Antibody Conjugates

Secondary antibody-alkaline phosphatase conjugates are optimized for Western blotting and immunoassay procedures incorporating 1,2-dioxetane enzyme substrates. These labeled secondary antibodies provide minimum nonspecific binding and are recommended for applications employing Tropix chemiluminescent substrates. The conjugates are supplied at approx. 0.3 mg/mL in a buffered solution containing 50% glycerol. Goat anti-rabbit IgG (H+L) exhibits minimum cross-reactivity with human serum proteins. Goat anti-mouse IgG+IgM (H+L) exhibits minimum cross-reactivity with human, bovine, and horse serum proteins. The recommended initial working dilution for these labeled antibodies is 1:5,000; the best working dilution may need to be optimized for different uses.

Avidx-AP[™] Streptavidin-AP Conjugate

Avidx-AP[™] conjugate is a streptavidin-alkaline phosphatase conjugate developed for detection of biotinylated nucleic acids or proteins. The conjugate is tested and optimized for low background with chemiluminescent detection procedures. Conjugation methods developed at Applied Biosystems result in minimal nonspecific binding of conjugate to membranes. Avidx-AP conjugate is available as part of Southern-Light[™] and Southern-Star[™] detection systems and is also sold separately, for use with other applications, such as immunoassay detection with biotinylated antibodies, microplate detection assays for nucleic acids labeled with biotin, or other assays that incorporate a biotin label. The suggested working dilution of Avidx-AP conjugate is from 1:5,000 to 1:20,000.

DEA (Diethanolamine) Buffer Concentrate

DEA is the buffer of choice in applications employing concentrated Tropix^{*} substrates for the detection of alkaline phosphatase. DEA buffers provide favorable reaction kinetics and greater sensitivity compared to most buffer systems. DEA is supplied 99% pure (MW = 105.1 g/mol; density = 1.097 g/mL). The suggested working concentration is 0.1 M DEA (pH 10) with 1 mM MgCl₂.

I-Block[™] Protein-based Blocking Reagent

I-Block[™] reagent is a highly purified casein-based blocking reagent. The reagent provides superior blocking compared to both dried milk and BSA. Unlike other casein-based blocking reagents, I-Block reagent is essentially biotin-free. I-Block reagent is tested in assays using Tropix substrates and alkaline phosphatase conjugates, and is useful as a blocking reagent in membrane-based and immunoassay applications. The suggested working concentration is 0.2% (w/v) for detection of nucleic acids on neutral or positively-charged nylon membranes, and for immunoassays and protein detection on membranes (nitrocellulose, PVDF, or neutral nylon). For protein detection on positively-charged nylon membrane (such as Tropilon-Plus[™] membrane), a concentration of 3% is recommended. I-Block blocking solution is prepared in either Tris- or phosphate-buffered saline buffer with heating (40-50°C).

Tropifluor[™] Polyvinylidene Fluoride (PVDF) Membrane

Tropifluor[™] PVDF transfer membrane is a 0.45 µm pore-size membrane for protein blotting that provides the high binding (125 µg/cm²) and retention capacity inherent in PVDF (polyvinylidene fluoride) membranes. The low background and resulting superior signal-to-noise performance obtained with Tropifluor membrane enables picogram detection levels of proteins in immunoblotting applications. Tropifluor membrane is selected to yield low background levels with the use of Tropix substrates and conjugates, enabling ultrasensitive detection. PVDF membranes exhibit high mechani-

cal strength that enable multiple re-probing. Film exposures using Tropix substrates with Tropifluor membrane range from 1 second to 15 minutes. The combination of chemiluminescent detection with Western-Light[™] or Western-Star[™] immunoblotting detection systems and the high proteinbinding capacity of Tropifluor membrane provides the greatest sensitivity available in Western blotting.

Tropilon-Plus[™] Positively Charged Nylon Membrane

Tropilon-Plus[™] nylon membrane is a 0.45 µm pore-size, positively-charged membrane for nucleic acid blotting. The inherent positive charge offers increased binding capacity (>100 µg/cm²) and signal retention in nucleic acid detection experiments. In addition, Tropilon-Plus membrane has excellent lot-to-lot consistency with low within-lot binding capacity variances.

Tropilon-Plus membrane is optimized for chemiluminescent detection using Tropix substrates and alkaline phosphatase conjugates. Background levels associated with chemiluminescent methods can vary dramatically with other nylon membranes. Tropilon-Plus membrane provides a high-intensity signal with low levels of nonspecific background.

Tropilon-Plus membrane offers several advantages compared to neutral nylon membrane and other types of positively-charged nylon membranes. DNA binding capacity is higher and background signal is lower using the Tropilon-Plus membrane.

APPLIED BIOSYSTEMS ORDER INFO	SIZE	CAT#	
AP-Conjugate - Goat Anti-Rabbit IgG	100 µL	T2191	
AP-Conjugate - Goat Anti-Mouse IgG+IgM	100 µL	T2192	
Avidx-AP™ Streptavidin-AP Conjugate	1 mL	T2016	
DEA (Diethanolamine) Buffer Concentrate	120 µL	T2027	
I-Block™ Blocking Reagent	30 g	T2015	
Tropifluor™ PVDF Membrane (15 cm x 15 cm) - 5 Membranes/Pack	1 pack	T2234	
Tropilon-Plus™ Nylon Membrane (30 cm x 200 cm)	1 roll	T2232	
Development Folders (14 cm x 19 cm) - 30 Folders/Pack	1 pack	T2258	
Assay Buffer (10X) for Immunoassay, Western, and Southern Blotting	150 mL	T2187	

Learn more or order the products on this page at www.appliedbiosystems.com

North America

United States and Canada (English) Tel: 800.327.3002 or 650.638.5800 Fax: 650.638.5998

United States and Canada (French) Tel: 800.668.6913 Fax: 650.638.5875

Press "1" for Reagents and Consumables, including Ambion and Tropix brand products

Press "2" for Instruments and Software

Press "3" for Product Information

Press "4" for Service

Press "5" for Technical Support

Latin America

Argentina

Tel: +54.11.4854.7775 Fax: +54.11.4857.0884

Brazil

Tel: +55.11.5070.9600 Fax: +55.11.5070.9694

Mexico

Tel: +52.55.5535.3610 Fax: +52.55.5566.2308

Europe Austria

Tel: +43.1.867.35.75 Fax: +43.1.867.35.72

Belgium Tel: +0800.77074 Fax: +32.2.582.18.86

Denmark Tel: +45.45.58.60.00 Fax: +45.45.58.60.01

Finland

Tel: +358.9.69.37.94.27 Fax: +358.9.69.37.94.26

France

Tel: +33.1.69.59.85.85 Fax: +33.1.69.59.85.00

Germany

Tel: +49.6151.9670.0 Fax: +49.6151.9670.5599

Italy Tel: +39.039.83891 Fax: +39.039.838.9492

Netherlands

Tel: +0800.224.7253 Fax: +31.180.392.409

Norway

Tel: +47.23.16.25.75 Fax: +47.23.16.25.74

Portugal

Tel: +800206.639 (local only) Fax: +800.206.640

Spain

Tel: +34.91.484.6900 Fax: +34.91.806.1206

Sweden

Tel: +46.8.619.44.00 Fax: +46.8.619.44.01

Switzerland

Tel: +41.41.799.77.77 Fax: +41.41.790.06.76

United Kingdom

Applied Biosystems Tel: +44.1925.825650 Fax: +44.1925.282502

Ambion Tel: +44.1480.373.020 Fax: +44.1480.373.010

European Managed Territories

Africa Tel: +27.11.478.0411 Fax: +27.11.478.0349

Czech Republic, Slovakia Tel: +420.2.3536.5189 Fax: +420.2.3536.4314

Hungary Tel: +36.1.471.8989 Fax: +36.1.471.8980

Poland

Tel: +48.22.866.4010 Fax: +48.22.866.4020

Russian Federation and CIS Tel: +7.495.781.8191 Fax: +7.495.781.8192

South East Europe Tel: +385.1.3460.839 Fax: +385.1.3460.840

Baltics, Middle East, and West Asia Tel: +44.1925.282481 Fax: +44.1925.282509

Japan, Eastern Asia, China and Oceania

Australia Tel: +61.3.9730.8600 Fax: +61.3.9730.8799

China Beijing

Tel: +86.10.6410.6608 Fax: +86.10.6410.6617

Guangzhou Tel: +86.20.8760.9229

Shanghai Tel: +86.21.6473.6366

Hong Kong Tel: +852.2756.6928 Fax: +852.2756.6968

Japan

Tel: +81.3.5566.6100 Fax: +81.3.5566.6501

Korea Tel: +822.592.7238/9 Fax: +822.599.8713

New Zealand Tel: +0800.446.416

Singapore

Tel: +65.6896.2168 Fax: +65.6896.2147

Taiwan

Tel: +886.2.2358.2838 Fax: +886.2.2358.2839

Cell Biology Sales

US (West), Sales Specialist Jimmy Slack Tel: (858) 964-8152

US (East), Sales Specialist

Albert S. Faucon Tel: (617) 990-6960 Albert.Faucon@appliedbiosystems.com

James.Slack@appliedbiosystems.com

US, Field Applications Specialist

Dr. Carol Khodier Tel: (781) 280-5662 khodiecs@appliedbiosystems.com

Europe, Sr. Sales Manager

Damian Leyland Tel: +44-777-8556553 Damian.E.Leyland@eur.appliedbiosystems.com

Europe, Sales Specialist

Dr. Cäcilia Simon Tel: +49-173-3474608 Caecilia.Simon@eur.appliedbiosystems.com

Europe, Sr. Field Applications Specialist Scott Cribbes, PhD

Tel: +44-7836-618378 Scott.Cribbes@eur.appliedbiosystems.com

Asia Pacific (outside Japan) Sales Andrew Kyriazis Tel: +61 411 469 122 Andrew.Kyriazis@appliedbiosystems.com

US Technical Support

Patrick Moore Tel: (508)383-7958 moorepl@appliedbiosystems.com Information is subject to change without notice.

Applied Biosystems Corporate Headquarters

850 Lincoln Centre Drive Foster City, CA 94404 USA
 Toll-free Tel:
 800.327.3002

 Tel:
 650.638.5800

 Fax:
 650.638.5998