

QMS™ TEICOPLANIN APPLICATION
Beckman Coulter AU5800®



Catalog No. 0374645

Intended for the quantitative determination of teicoplanin in human serum or plasma on automated clinical chemistry analyzers as an aid in the management of patients receiving teicoplanin therapy.

For In Vitro Diagnostic Use Only
For Use Outside the US Only


Intended Use The information provided in this application sheet is intended as a supplement to the package insert. Refer to the package insert for information on intended use, reagent storage, reagent preparation, specimen collection, specimen preparation, specimen storage, quality control, and additional performance data.



Ordering Information

Item	Size	Catalog No.
QMS Teicoplanin Assay	R1 21 mL, R2 9 mL	0374645
QMS Teicoplanin Calibrator Set	6 levels, 1 x 1.0 mL ea.	0374652
QMS Teicoplanin Control Set	3 levels 1 x 2.0 mL ea.	0374660

To place an order or for technical service contact:

USA	In Europe
(800) 232-3342 Fax 510-979-5002	Tel: +49 (0)851-88 6890 Fax: +49 (0)851-88 68910

 **Microgenics Corporation, part of Thermo Fisher Scientific**
46500 Kato Road, Fremont, CA 94538 USA
U.S. Toll free: (800) 232-3342 / Tel: (510) 979-5000
U.S. Toll free fax: (888) 527-8001 / Fax: (510) 979-5420

  B·R·A·H·M·S GmbH, Neuendorfstrasse 25, 16761, Hennigsdorf, Germany
Tel: +49 (0)851-88 6890/ Fax: +49 (0)851-88 68910

Continued on next page

**Reagent
Storage**

Refer to the package insert for information on reagent storage.

**Procedure for
Analyzer**

Refer to the operator's manuals for information on analyzer operation.

**Results and
Data
Interpretation**

Results for samples will be printed in $\mu\text{g/mL}$.

Beckman Coulter AU5800 Parameters

QMS Teicoplanin Assay

Parameters		Specific Test Parameters												
General	LIH	ISE	HbA1c		Calculated Test	Range								
Test Name:		TCO	<	>	Type:	Serum	Operation	Yes						
Sample Volume	2.2	μL	Dilution	0	μL	OD Limit								
Pre-Dilution Rate	1	▽	Diluent Bottle	#	▽	Min.OD	-2.00	Max.OD	3.00					
Rgt. Volume	R1(R1-1)	162	μL	Dilution	0	μL	Reagent OD Limit							
	R1-2		μL	Dilution		μL	1 st .	Low	-2.00	High	3.00			
							Last	Low	-2.00	High	3.00			
	R2(R2-1)	54	μL	Dilution	0	μL	Dynamic Range		Low	3	High	100		
Common Rgt. Type	None		Name	None		Correlation Factor A	1		B	0				
Wavelength	Pri	700	▽nm	Sec.	None		▽nm	Factor for Maker	A		1	B	0	
Method	FIXED1		▽		Onboard Stability Period		#	Day	#	Hour				
Reaction Slope	+		▽		LIH Influence Check		#	▽						
Measuring Point1	1 st	12	Last		17	Lipemia		▽						
Measuring Point2	1 st		Last			Icterus		▽						
Linearity Limit						Hemolysis		▽						
Lag Time Check	No		▽											

Parameters		Specific Test Parameters						
General	LIH	ISE	HbA1c		Calculated Test	Range		
Test Name:		TCO	<	>	Type:	Serum	▽	
Value/Flag:		#	▽					
Level		Low	#	High	#			
Specific Ranges:		From		To				
	Sex	Year	Month	Year	Month	Low	High	
<input type="checkbox"/>	1.	#	▽	#	#	#	#	
<input type="checkbox"/>	2.	#	▽	#	#	#	#	
<input type="checkbox"/>	3.	#	▽	#	#	#	#	
<input type="checkbox"/>	4.	#	▽	#	#	#	#	
<input type="checkbox"/>	5.	#	▽	#	#	#	#	
<input type="checkbox"/>	6.	#	▽	#	#	#	#	
	7.	Standard demographics						
	8.	Not within expected values						
Panic Value	Low	#	High	#	Unit	μg/mL	Decimal Places	#

Beckman Coulter AU5800 Parameters

QMS Teicoplanin, *continued*

Parameters		Calibration Parameters			
Calibrators	Calibration Specific	STAT Table Calibration			
Test Name: # ▾	< ▾ > ▾	Type Serum ▾	Cuvette ▾		
<input type="checkbox"/> Use Serum Cal					
Calibration Type: 6AB ▾		Formula: EIA Type 1 ▾		Counts: 2	
< Calibrator Parameters >					
	Calibrator	OD	Conc	OD Range Low High	Slope - ▾
Point-1	# ▾		0	-2.0 3.0	Allowable Range Check <input type="checkbox"/> Reagent Blank <input type="checkbox"/> Calibration Advanced Calibration Operation No ▾ Interval (RB/ACAL) ▾ <input type="checkbox"/> Lot Calibration
Point-2	# ▾		5	-2.0 3.0	
Point-3	# ▾		10	-2.0 3.0	
Point-4	# ▾		25	-2.0 3.0	
Point-5	# ▾		50	-2.0 3.0	
Point-6	# ▾		100	-2.0 3.0	
Point-7	▾				
Point-8	▾				
Point-9	▾				
Point-10	▾				
<Point Cal. For Master Curve>		No. of Correction Points	Use Master Curve ▾		
	Calibrator	OD	Conc	OD Range Low High	Stability
Point-1	▾				Reagent Blank # Day # Hour
Point-2	▾				Calibration # Day # Hour
MB Type Factor		1-Point Calibration Point		- ▾	<input type="checkbox"/> with Conc-0

User defined

© 2017 Thermo Fisher Scientific Inc. All rights reserved.
 All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.
 AU Series Systems are trademarks of Beckman Coulter.

End