


Fermentas  
now sold as  
Thermo Scientific

Introducing Thermo Scientific  
Pierce Prestained and Unstained  
Protein Ladders for easy protein  
molecular weight determination.

Reach greater **highs** (and **lows**) with new protein ladder choices.



Fermentas  
now sold as  
Thermo Scientific

The Thermo Scientific Pierce Protein Research Product Portfolio now includes protein ladders made with Fermentas technology. The Pre-stained Protein Ladders are recommended for approximate determination of molecular weight, monitoring the progress of the electrophoresis run and the efficiency of the protein transfer to the membrane during Western blotting procedures. For precise determination of molecular weights of target proteins in any buffer system, we recommend using our unstained protein ladders.

**Highlights:**

- **Cost effective** – high quality without the high price
- **Performance** – sharp protein band resolution and consistent migration patterns provide easy molecular weight determination
- **Convenient** – protein ladders are ready to load with no heating or boiling required
- **Reliable** – exceptional lot-to-lot consistency and reproducibility

**new protein ladder choices**



# Reach new heights in gel electrophoresis

## Ordering Information

Product #	Description	Pkg. Size
26610	Pierce® Unstained Protein MW Marker	2 x 1mL
26612	Pierce Prestained Protein MW Marker	500µL
26614	PageRuler™ Unstained Protein Ladder	2 x 250µL
26616	PageRuler Prestained Protein Ladder	2 x 250µL
26617	PageRuler Prestained Protein Ladder	10 x 250µL
26619	PageRuler Plus Prestained Protein Ladder	2 x 250µL
26620	PageRuler Plus Prestained Protein Ladder	10 x 250µL
26634	Spectra® Multicolor Broad Range Protein Ladder	2 x 250µL
26623	Spectra Multicolor Broad Range Protein Ladder	10 x 250µL
26625	Spectra Multicolor High Range Protein Ladder	2 x 250µL
26628	Spectra Multicolor Low Range Protein Ladder	250µL
26630	PageRuler Unstained Broad Range Protein Ladder	2 x 250µL
26632	PageRuler Unstained Low Range Protein Ladder	2 x 250µL
26637	PageRuler Unstained High Range Protein Ladder	2 x 250µL


## Specialty Ladders

Product #	Description	Pkg. Size
26635	PageRuler Prestained NIR Protein Ladder	2 x 250µL
84785	SuperSignal® Molecular Weight Protein Ladder	250µL
84786	SuperSignal Enhanced Molecular Weight Protein Ladder	250µL




Visit [www.thermoscientific.com/protein-custom](http://www.thermoscientific.com/protein-custom) to request a quote, or contact our Bulk and Custom Sales department at [LVC.rockford@thermofisher.com](mailto:LVC.rockford@thermofisher.com).

[www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce)



Thermo Scientific Pierce Protein Ladders						Application							
						SDS PAGE				Western blotting			
Product Name	Product #	Former Fermentas Part #	# Proteins	MW Range	Reference Bands	Protein MW determination	Protein band visualization	Monitoring electrophoresis run	Coomassie/silver/fluorescent staining	NIR detection	Monitoring protein transfer	Chemiluminescent band visualization	
Colorimetric	Unstained	PageRuler Unstained Broad Range Protein Ladder	26630	SM1881	11	5-250kDa 20, 50 and 100kDa	Best	—	—	Best	—	—	Good†
		PageRuler Unstained Protein Ladder	26614	SM0661	14	10-200kDa 50kDa	Good	—	—	Good	—	—	Good†
		PageRuler Unstained Low Range Protein Ladder	26632	SM1891	8	3.4-100kDa 25kDa	Best	—	—	Best	—	—	Good†
		PageRuler Unstained High Range Protein Ladder	26637	—	8	60-250kDa 150kDa	Best	—	—	Best	—	—	Good†
		Pierce Unstained Protein Molecular Weight Marker	26610	SM0431	7	14.4-116kDa —	Good	—	—	Good	—	—	—
	Prestained	PageRuler Prestained Protein Ladder	26616, 26617	SM0671, SM0672	10	10-170kDa Green 10kDa; Orange 70kDa	Good	Good	Good	—	Good‡	Good	—
		PageRuler Plus Prestained Protein Ladder	26619, 26620	SM1811, SM1812	9	10-250kDa Green 10kDa; Orange 25 and 70kDa	Good	Good	Good	—	Good‡	Good	—
		Spectra Multicolor Broad Range Protein Ladder	26634, 26623	SM1841, SM1842	10	10-260kDa Green 10 and 50kDa; Orange 40, 70 and 260kDa; Pink 140kDa	Good	Best	Best	—	Good‡	Best	—
		Spectra Multicolor Low Range Protein Ladder	26628	SM1861	6	1.7-40kDa Green 10kDa; Orange 40kDa	Good	Best	Best	—	Good‡	Best	—
		Spectra Multicolor High Range Protein Ladder	26625	SM1851	8	40-300kDa Green 50kDa; Orange 70 and 300kDa	Good	Best	Best	—	Good‡	Best	—
Pierce Prestained Protein Molecular Weight Marker	26612	SM0441	6	20-120kDa —	Good	Good	Good	—	Good	Good	—		
Chemi. Lum.	SuperSignal Molecular Weight Protein Ladder	84785, 84786	—	8	20-150kDa —	Good	—	—	Good	—	Good	Best	
NIR	PageRuler Prestained NIR Protein Ladder	26635	—	10	11-250kDa 55kDa	Good	Good	Good	—	Best	Good	—	

† with Strep-Tactin® HRP or AP Conjugates  
‡ Blue and green prestained bands fluoresce in NIR region



Migration Patterns of PageRuler Prestained Protein Ladder

Thermo Scientific PageRuler Prestained Protein Ladder														
Gel type	Tris-Glycine						Tris-Acetate <sup>††</sup>		Bis-Tris <sup>††</sup>					
Gel concentration	4-20%	8-16%	10-20%	8%	10%	12%	15%	3-8%	7%	4-12%	10%		12%	
Running buffer	Tris-Glycine						Tris-Acetate		MOPS	MES	MOPS	MES	MOPS	MES
Apparent Molecular Sizes (kDa)														
% length of gel	10	170	170	170	170	170	170	170	130	130	150	140	140	140
20	130	130	100	130	130	100	70	70	115	115	140	115	115	115
30	100	100	70	100	100	55	40	55	80	80	80	80	80	80
40	70	70	55	70	55	40	35	40	65	65	65	65	65	65
50	55	55	40	55	40	35	25	35	50	50	50	50	50	50
60	40	40	35	40	35	25	25	25	40	40	40	40	40	40
70	35	35	25	35	25	25	15	15	30	30	30	30	30	30
80	25	25	15	25	15	15	10	10	25	25	25	25	25	25
90	15	15	10	15	10	10	10	10	15	15	15	15	15	15
100	10	10	10	10	10	10	10	10	10	10	10	10	10	10

Migration Patterns of PageRuler Plus Prestained Protein Ladder

Thermo Scientific PageRuler Plus Prestained Protein Ladder															
Gel type	Tris-Glycine						Tris-Acetate <sup>††</sup>		Bis-Tris <sup>††</sup>						
Gel concentration	4-20%	8-16%	10-20%	8%	10%	12%	15%	3-8%	7%	4-12%	10%		12%		
Running buffer	Tris-Glycine						Tris-Acetate		MOPS	MES	MOPS	MES	MOPS	MES	
Apparent Molecular Sizes (kDa)															
% length of gel	10	250	250	250	250	250	250	250	130	130	205	185	190	185	190
20	130	130	100	130	130	100	70	70	115	115	190	115	115	115	115
30	100	100	70	100	100	55	40	55	80	80	80	80	80	80	80
40	70	70	55	70	55	35	25	35	65	65	65	65	65	65	65
50	55	55	40	55	40	35	25	25	50	50	50	50	50	50	50
60	40	40	35	40	35	25	25	25	40	40	40	40	40	40	40
70	35	35	25	35	25	25	15	15	30	30	30	30	30	30	30
80	25	25	15	25	15	15	10	10	25	25	25	25	25	25	25
90	15	15	10	15	10	10	10	10	15	15	15	15	15	15	15
100	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10

Migration Patterns of Spectra Multicolor Broad Range Protein Ladder

Thermo Scientific Spectra Multicolor Broad Range Protein Ladder															
Gel type	Tris-Glycine						Tris-Acetate <sup>††</sup>		Bis-Tris <sup>††</sup>						
Gel concentration	4-20%	8-16%	10-20%	8%	10%	12%	15%	3-8%	7%	4-12%	10%		12%		
Running buffer	Tris-Glycine						Tris-Acetate		MOPS	MES	MOPS	MES	MOPS	MES	
Apparent Molecular Sizes (kDa)															
% length of gel	10	260	260	260	260	260	260	260	140	140	225	225	235	225	235
20	140	140	100	140	140	100	70	70	115	120	120	115	120	115	120
30	100	100	70	100	100	50	40	50	80	80	80	80	80	80	80
40	70	70	50	70	50	40	35	40	65	65	65	65	65	65	65
50	50	50	40	50	40	35	25	35	50	50	50	50	50	50	50
60	40	40	35	40	35	25	25	25	40	40	40	40	40	40	40
70	35	35	25	35	25	25	15	15	30	30	30	30	30	30	30
80	25	25	15	25	15	15	10	10	25	25	25	25	25	25	25
90	15	15	10	15	10	10	10	10	15	15	15	15	15	15	15
100	10	10	10	10	10	10	10	10	10	10	10	10	10	10	10

Migration Patterns of Spectra Multicolor High Range and Low Range Protein Ladder

Thermo Scientific Spectra Multicolor High Range Protein Ladder									
Gel type	Tris-Glycine					Tris-Acetate <sup>††</sup>		Bis-Tris <sup>††</sup>	
Gel concentration	4-12%	4%	6%	8%	10%	4-20%	3-8%	7%	4-12%
Running buffer	Tris-Glycine					Tris-Acetate		MOPS	
Apparent Molecular Sizes (kDa)									
% length of gel	10	300	300	300	300	300	270	270	270
20	250	250	250	250	250	250	205	205	185
30	180	180	180	180	180	180	150	150	140
40	130	130	130	130	130	130	120	120	115
50	100	100	100	100	100	100	85	85	80
60	70	70	70	70	70	70	65	65	65
70	50	50	50	50	50	50	50	50	50
80	40	40	40	40	40	40	40	40	40
90	100	100	100	100	100	100	100	100	100
100	100	100	100	100	100	100	100	100	100


Thermo Scientific Spectra Multicolor Low Range Protein Ladder												
Gel type	Tris-Tricine											
Gel concentration	4-20%	10-20%	14%	16%	18%	14-SP %	16-SP %	18-SP %	10-20% <sup>††</sup>	16% <sup>††</sup>	10-20% <sup>†††</sup>	
Running buffer	Tris-Tricine SDS											
Apparent Molecular Weights (kDa)												
% length of gel	10	40	40	40	40	40	40	40	40	40	40	40
20	25	25	25	25	25	25	25	25	25	25	25	25
30	15	15	15	15	15	15	15	15	15	15	15	15
40	10	10	10	10	10	10	10	10	10	10	10	10
50	5	5	5	5	5	5	5	5	5	5	5	5
60	2	2	2	2	2	2	2	2	2	2	2	2
70	10	10	10	10	10	10	10	10	10	10	10	10
80	5	5	5	5	5	5	5	5	5	5	5	5
90	2	2	2	2	2	2	2	2	2	2	2	2
100	100	100	100	100	100	100	100	100	100	100	100	100

Note: The apparent molecular size of each protein (kDa) has been determined by calibration against an unstained protein ladder in each electrophoresis condition.  
 †† – Migration patterns were determined using respective NuPAGE® precast gels.

SP – small peptide SDS-polyacrylamide gel with higher degree of crosslinking (c~5% instead of usual ~3%).  
 †† – migration patterns were determined using respective Novex® precast gels.  
 ††† – migration patterns were determined using respective Ready Gel® precast gel.


### Thermo Scientific Pierce Protein Molecular Weight Marker

The Thermo Scientific Pierce Prestained Protein Molecular Weight Marker is designed to monitor protein migration during SDS-polyacrylamide gel electrophoresis and monitor protein transfer onto membranes. Pierce Unstained Protein Molecular Weight Marker produces sharp bands on SDS-polyacrylamide gel following staining with Coomassie-based or silver stains and is recommended for accurate sizing of proteins on SDS PAGE and Western blots.




### Thermo Scientific PageRuler Prestained Protein Ladders

The Thermo Scientific PageRuler and PageRuler Plus Prestained Protein Ladders are comprised of highly purified recombinant proteins bound to three different chromophores for easy protein molecular weight determination. The green and orange reference bands in the PageRuler Protein Ladders provide easy orientation of the protein bands.




### Thermo Scientific PageRuler Unstained Protein Ladders

Thermo Scientific PageRuler Unstained Protein Ladders are mixtures of recombinant, highly purified proteins designed for accurate sizing of proteins by SDS-PAGE. The ladders resolve into clearly identifiable sharp bands when analyzed by SDS-PAGE and stained with Coomassie Blue or silver stain. Each ladder contains protein bands of greater intensity that serve as reference bands. Proteins can be detected on Western blots by staining with Ponceau S, or Coomassie Blue (except for 3.4kDa and 5kDa polypeptides) or by using Strep-Tactin conjugates.




### Thermo Scientific Spectra Multicolor Protein Ladders

Four different chromophores bound to proteins in Thermo Scientific Spectra Multicolor Broad Range Protein Ladders produce a brightly colored ladder with an easy-to-remember pattern. The Spectra Multicolor High Range and Low Range Protein Ladders are designed for large and small protein analysis, respectively.




### Thermo Scientific SuperSignal Molecular Weight Protein Ladders

Thermo Scientific SuperSignal Molecular Weight Protein Ladders contain a ready-to-use stabilized mixture of eight recombinant proteins ranging in size from 20 to 150kDa. These recombinant proteins bind antibodies used in the Western blot through an IgG binding site. The protein markers can then be visualized either using appropriate substrates for enzyme-labeled antibodies or via fluorescent dye-labeled antibodies. Use SuperSignal Enhanced Molecular Weight Protein Ladders with mouse monoclonals.



### Thermo Scientific PageRuler Prestained NIR Protein Ladder

The PageRuler Prestained NIR Protein Ladder is a mixture of 10 proteins (11 to 250kDa) that are blue-stained and fluor-labeled for near-IR fluorescent visualization and protein sizing. The protein MW markers in this ladder resolve into sharp bands when analyzed by SDS-PAGE and are labeled with a fluorescent dye for visualization with instruments equipped with near-infrared (NIR) fluorescence detection. The bands are also directly visible because the proteins are prestained blue.



© 2012 Thermo Fisher Scientific Inc. All rights reserved. These products are supplied for laboratory or manufacturing applications only. Prices listed herein were accurate at the time of printing. Visit our website for up-to-date prices. Facebook is a registered trademark of Facebook, Inc. Novex and NuPAGE are registered trademarks of Invitrogen. Ready Gel is a registered trademark of Bio Rad. Strep-Tactin is a registered trademark of IBA GmbH. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

#### Life Science Research

**Belgium and Europe, the Middle East and Africa Distributors**  
Tel: +32 53 85 71 84

**France**  
Tel: 0 800 50 82 15

**The Netherlands**  
Tel: 076 50 31 880

**Germany**  
Tel: 0228 9125650

**United Kingdom**  
Tel: 0800 252 185

**Switzerland**  
Tel: 0800 56 31 40

Email: perbio.euromarketing@thermofisher.com  
www.thermoscientific.com/perbio

For other regions, visit  
[www.thermoscientific.com/piercedistributors](http://www.thermoscientific.com/piercedistributors)

#### United States

Tel: 815-968-0747 or 800-874-3723  
Customer Assistance E-mail:  
Pierce.CS@thermofisher.com  
[www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce)

