



# BL21-AI™ Competent Cells—maximum expression with tightest regulation

BL21-AI™ One Shot® Chemically Competent *E. coli* give you:

- Tightest T7 regulation
- High yields of protein
- Convenient One Shot® format

## Tightly regulated, highly inducible

BL21-AI™ One Shot® Competent *E. coli* is an arabinose-inducible strain designed to give you the maximum protein expression with the tightest regulation available from a T7 expression system. It's the only strain that gives you both tight regulation and high yields, making it great for high-level expression of


toxic proteins. Because it has the tightly regulated arabinose-inducible *araBAD* promoter upstream of the T7 RNA polymerase gene, you can use it with any T7 promoter-based vector. All that in the convenient, cost-effective One Shot® format.

## Turn on expression

The results dramatically demonstrate the effectiveness of BL21-AI™ *E. coli*. To show tight regulation and high inducibility, the *lacZ* gene was cloned into the T7-based pET-201 and expressed in BL21-AI™ and BL21(DE3)pLysS strains. In the uninduced state, BL21-AI™ *E. coli* yielded the lowest levels of β-galac-


tosidase expression, whereas in the induced state BL21-AI™ gave the highest levels of expression (Figure 1). In addition to the lowest basal expression, BL21-AI™ *E. coli* yield higher protein levels than the BL21(DE3)pLysS strain (Figure 2).

**Figure 1 - Regulation of β-galactosidase expression in BL21-AI™ and BL21(DE3)pLysS**



BL21-AI™ *E. coli* offers lower basal (uninduced) and higher induced levels of T7 RNA polymerase mediated gene expression than BL21(DE3)pLysS. Cells were induced by adding arabinose (0.2% w/v) and IPTG (1 mM) for BL21-AI™, and IPTG (1 mM) for BL21(DE3)pLysS.

**Figure 2 - Expression of a toxic protein in BL21-AI™ and BL21(DE3)pLysS**



Expression of a toxic protein is enhanced in BL21-AI™ when compared to BL21(DE3)pLysS. A model toxic protein, elastin-like peptide-TEV protease fusion (ELP-TEV), was expressed from pET11b. BL21-AI™ *E. coli* was either uninduced (U) or induced (I) with 0.1% arabinose and 1 mM IPTG. BL21(DE3)pLysS was either un-induced (U) or induced (I) with 1 mM IPTG.

## Give it a shot

With BL21-AI™ One Shot® Chemically Competent *E. coli*, you get a complete kit that you can depend on for top performance. Each One Shot® kit contains 21 individual tubes of competent cells, 6 ml of S.O.C. medium, and pUC18 positive control. You only thaw the cells you need. No efficiency-zapping freeze-thaw cycles. No money wasted on unused cells. Perform your transformation directly in the tube and get great results every time.



## Take control

Use BL21-AI™ with any T7 promoter-containing vector for tightly controlled, high-level protein expression. You'll get

the regulation and yield your experiment requires. Call and order today.

### Product

BL21-AI™ One Shot® Chemically Competent *E. coli*

### Efficiency\*

$> 1 \times 10^8$

### Quantity

20 x 50  $\mu$ l

### Cat. no.

C6070-03

\* Transformation efficiency is calculated as cfu/ $\mu$ g of supercoiled pUC18 plasmid.

Genotype: F<sup>-</sup> *ompT hsdS<sub>B</sub>(r<sub>B</sub><sup>-</sup> m<sub>B</sub><sup>-</sup>) gal dcm araB::T7RNAP-tetA*.



710-011097 031402