







Loading dye	Cat. #	Size, ml	Composition	Features	Applications	Migration of dyes (1% agarose, TAE or TBE buffers)	Picture of tracking dyes*
<b>6X DNA Loading Dye</b>	R0611	5x1	<b>6X Solution</b> <ul style="list-style-type: none"> <li>10 mM Tris</li> <li>0.03% bromophenol blue</li> <li>0.03% xylene cyanol FF</li> <li>60% glycerol</li> <li>60 mM EDTA (pH 7.6 adjusted with NaOH)</li> </ul>	<ul style="list-style-type: none"> <li>Two-color tracking of DNA migration during electrophoresis.</li> <li>Dyes do not interfere with UV visualization of DNA fragments.</li> <li>EDTA inhibits metal dependent nucleases.</li> </ul>	<ul style="list-style-type: none"> <li>Preparation of DNA for loading on agarose or polyacrylamide gels.</li> </ul>	<b>Xylene cyanol FF:</b> TAE: 4160 bp TBE: 3030 bp  <b>Bromophenol blue:</b> TAE: 370 bp TBE: 220 bp	
<b>6X MassRuler DNA Loading Dye</b>	R0621	5x1	<b>6X Solution</b> <ul style="list-style-type: none"> <li>10 mM Tris</li> <li>0.03% bromophenol blue</li> <li>60% glycerol</li> <li>60 mM EDTA (pH 7.6, adjusted with NaOH)</li> </ul>	<ul style="list-style-type: none"> <li>One-color tracking of DNA migration during electrophoresis.</li> <li>Dyes do not interfere with UV visualization of DNA fragments.</li> <li>EDTA inhibits metal dependent nucleases.</li> </ul>	<ul style="list-style-type: none"> <li>Analysis of large DNA molecules.</li> <li>Preparation of DNA for loading on agarose or polyacrylamide gels.</li> </ul>	<b>Bromophenol blue:</b> TAE: 370 bp TBE: 220 bp	
<b>6X Orange DNA Loading Dye</b>	R0631	5x1	<b>6X Solution</b> <ul style="list-style-type: none"> <li>10 mM Tris</li> <li>0.15% orange G</li> <li>0.03% xylene cyanol FF</li> <li>60% glycerol</li> <li>60 mM EDTA (pH 7.6, adjusted with NaOH)</li> </ul>	<ul style="list-style-type: none"> <li>Two-color tracking of DNA migration during electrophoresis.</li> <li>Dyes do not interfere with UV visualization of DNA fragments.</li> <li>EDTA inhibits metal dependent nucleases.</li> </ul>	<ul style="list-style-type: none"> <li>Analysis of small DNA molecules.</li> <li>Preparation of DNA for loading on agarose or polyacrylamide gels.</li> </ul>	<b>Xylene cyanol FF:</b> TAE: 4160 bp TBE: 3030 bp  <b>Orange G:</b> TAE/TBE: <50 bp	
<b>6X TriTrack DNA Loading Dye</b>	R1161	5x1	<b>6X Solution</b> <ul style="list-style-type: none"> <li>10 mM Tris</li> <li>0.03% bromophenol blue</li> <li>0.03% xylene cyanol FF</li> <li>0.15% orange G</li> <li>60% glycerol</li> <li>60 mM EDTA (pH 7.6, adjusted with NaOH)</li> </ul>	<ul style="list-style-type: none"> <li>Three-color tracking of DNA migration during electrophoresis.</li> <li>Dyes do not interfere with UV visualization of DNA fragments.</li> <li>EDTA inhibits metal dependent nucleases.</li> </ul>	<ul style="list-style-type: none"> <li>Preparation of DNA for loading on agarose or polyacrylamide gels.</li> </ul>	<b>Xylene cyanol FF:</b> TAE: 4160 bp TBE: 3030 bp <b>Bromophenol blue:</b> TAE: 370 bp TBE: 220 bp <b>Orange G:</b> TAE/TBE: <50 bp	
<b>6X DNA Loading Dye &amp; SDS Solution</b>	R1151	5x1	<b>6X Solution</b> <ul style="list-style-type: none"> <li>0.03% bromophenol blue</li> <li>0.03% xylene cyanol FF</li> <li>60% glycerol</li> <li>1% SDS</li> <li>100 mM EDTA (pH 7.6, adjusted with Tris)</li> </ul>	<ul style="list-style-type: none"> <li>1% SDS eliminates DNA-protein interactions, prevents appearance of additional bands due to annealing of DNA molecules with cohesive ends.</li> <li>EDTA inhibits metal-dependent nucleases.</li> </ul>	<ul style="list-style-type: none"> <li>Analysis of DNA samples containing high amounts of DNA binding proteins.</li> <li>Kinetic experiments.</li> <li>DNA agarose gel analysis following restriction digestion, ligation or dephosphorylation reactions.</li> </ul>	<b>Xylene cyanol FF:</b> TAE: 4160 bp TBE: 3030 bp <b>Bromophenol blue:</b> TAE: 370 bp TBE: 220 bp	
<b>2X RNA Loading Dye</b>	R0641	1	<b>2X Solution</b> <ul style="list-style-type: none"> <li>95% formamide</li> <li>0.025% SDS</li> <li>0.025% bromophenol blue</li> <li>0.025% xylene cyanol FF</li> <li>0.025% ethidium bromide</li> <li>0.5 mM EDTA</li> </ul>	<ul style="list-style-type: none"> <li>Two-color tracking of DNA and RNA fragment migration during electrophoresis.</li> <li>Dyes do not interfere with UV visualization of DNA or RNA fragments.</li> <li>Formamide based denaturation DNA and RNA fragments.</li> </ul>	<ul style="list-style-type: none"> <li>Preparation of DNA for loading on denaturing gels</li> <li>Preparation of RNA for loading on agarose or polyacrylamide gels.</li> </ul>	<b>Xylene cyanol FF:</b> TAE: 4160 bp TBE: 3030 bp <b>Bromophenol blue:</b> TAE: 370 bp TBE: 220 bp	

\* for more detailed information regarding the migration rates of dyes in agarose and polyacrylamide gels see Table 7.1 and Table 7.2 on p.374.