invitrogen



# Take a closer look at exosomes

#### Vesicle transport and cell signaling

The 2013 Nobel Prize<sup>™</sup> in Physiology or Medicine was jointly awarded to three scientists for their discovery of vesicle transport in cells, a clear indication of the importance of vesicles in physiology and their applied potential in diagnostics and therapy.

Exosomes and other extracellular vesicles are powerful mediators and strong biomarker candidates. There's an urgent need for the development of new methods for vesicle isolation and characterization.

To answer this need, we're continually expanding our product offerings. The products presented here can help simplify the overall study of exosomes, and are also applicable to translational research.

An up-to-date list of published scientific papers citing the use of these products can be found at **thermofisher.com/exosomes** 

#### "Exosomes: the next small thing"

We asked ten prominent scientists to share their thoughts on the field of exosome research. Based on these interviews, a fascinating six-part video documentary was produced. This miniseries tells the story of exosome research and its impact on important research areas such as cancer, immunology, stem cell research, and potential future therapeutic and diagnostic applications.

Watch the six-part video series at **thermofisher.com/exosomesdocumentary** 

#### **ExosomesTalk**

Join the conversation and post questions at **exosomestalk.com**. This forum is curated by scientists with substantial experience working with exosomes.



## **Exosome** isolation

The total exosome isolation reagents allow for efficient precipitation and high recovery of intact exosomes (Figure 1).

- Typically 15–20 minutes of hands-on time
- A lot easier and less tedious than ultracentrifugation

## **Exosome characterization**

Preenriched exosomes can easily be visualized by flow cytometry while bound to the surface of Dynabeads<sup>™</sup> magnetic beads, enabling the detection of specific exosomal markers.

- Clear and defined FSC/SSC for easier gating
- Typically less than 1 hour of hands-on time

## **Ordering information**

| Starting sample for total exosome isolation | Cat. No. |
|---|----------|
| Cell culture media (100 mL)                 | 4478359  |
| Blood serum (30 mL)                         | 4478360  |
| Blood plasma (20 mL)                        | 4484450  |
| Urine (50 mL)                               | 4484452  |
| Other body fluids (up to 30 mL)             | 4484453  |

By targeting classical exosome surface proteins,
Dynabeads™ magnetic separation technology allows you

• Scalable protocol with minimal hands-on time

(preenriched from cell culture) (Figure 2).

to easily pull out specific exosome subpopulations

• Enables you to obtain a highly pure exosome subset

## **Ordering information**

| Subpopulation analysis  | Quantity | Cat. No. |
|---|----------|----------|
| Human CD9 detection   | 100 rxns | 10620D   |
| Human CD81 detection  | 100 rxns | 10622D   |
| Human EpCam detection   | 100 rxns | 10624D   |
| Human CD63 (isolation/detection)  | 150 rxns | 10606D   |
| Flexible streptavidin-based system (use your own biotinylated antibody) | 150 rxns | 10608D   |

Light microscopy can be used to visualize exosomes. Specialized spin columns can be used to remove unincorporated dyes.

- Remove low molecular-weight contaminants
- Allows for buffer exchange and desalting

## **Ordering information**

| Subpopulation isolation   | Quantity | Cat. No. |
|---|----------|----------|
| Human CD9 isolation   | 50 preps | 10614D   |
| Human CD81 isolation  | 50 preps | 10616D   |
| Human EpCam isolation   | 50 preps | 10618D   |
| Human CD63 (isolation/detection)  | 30 preps | 10606D   |
| Flexible streptavidin-based system (use your own biotinylated antibody) | 30 preps | 10608D   |

## **Ordering information**

| Supporting products                              | Quantity   | Cat. No. |
|--|------------|----------|
| BODIPY™ TR Ceramide                              | 250 µg     | D7540    |
| SYTO™ RNASelect™ Green<br>Fluorescent Cell Stain | 100 μL     | S32703   |
| Exosome Spin Columns (for dye removal, MW 3000)  | 30 columns | 4484449  |

## Isolation of exosomal cargo

The Invitrogen™ Total Exosome RNA & Protein Isolation Kit allows for the isolation of total RNA—or both RNA and proteins—from your sample of preenriched exosomes.

- Highly pure total RNA (including small-RNA fraction)
- Recover protein and RNA from the same sample

## **Ordering information**

| Isolation of total RNA & proteins         | Quantity | Cat. No. |
|---|----------|----------|
| Total Exosome RNA & Protein Isolation Kit | 40 preps | 4478545  |

By using fast and gentle Dynabeads magnetic separation technology, you can isolate proteins and protein complexes from preenriched exosomes—typically in only 30 minutes.

- 10-50X concentration of exosomal proteins
- Helps to significantly minimize background

## **Ordering information**

| Immunoprecipitation products                    | Quantity | Cat. No. |
|---|----------|----------|
| Exosome Immunoprecipitation Reagent (Protein A) | 1 mL     | 10610D   |
| Exosome Immunoprecipitation Reagent (Protein G) | 1 mL     | 10612D   |

## Analysis of exosomal cargo

Specific monoclonal antibodies allow for the detection of cellular and exosomal antigens (Figure 2). Exosomal RNA and membrane components can also be labeled using fluorescent dyes, while unincorporated dye can be removed by using specific spin columns (Figure 3).


- Antibodies verified for western analysis
- Remove low molecular weight (<3,000 MW) contaminants</li>

## **Ordering information**

| Antibodies for western analysis         | Quantity        | Cat. No.              |
|---|-----------------|-----------------------|
| Anti-Human CD9                          | 0.2 mL          | 10626D                |
| Anti-Human CD63                         | 0.2 mL          | 10628D                |
| Anti-Human CD81                         | 0.2 mL          | 10630D                |
|   |                 |                       |
| Supporting products                     | Quantity        | Cat. No.              |
| Supporting products  BODIPY TR Ceramide | Quantity 250 μg | <b>Cat. No.</b> D7540 |
|   |                 |                       |

# invitrogen

## A Total Exosome Isolation Reagent



### **B** Ultracentrifugation




Figure 1. Analysis of exosomes recovered from HeLa cell media. (A) Exosomes recovered with the Invitrogen™ Total Exosome Isolation Reagent (from cell culture media) have a yield and size distribution comparable to (B) exosomes isolated following a traditional ultracentrifugation protocol with a sucrose gradient. Profiles as analyzed on a NanoSight™ LM10 instrument show all particles to be smaller than 300 nm; most are about 50–150 nm in size.





Figure 2. Analysis of exosomal markers CD81 and CD9. Exosomes were preenriched from SW480 cell cultures by ultracentrifugation. Immunoaffinity-based Dynabeads magnetic separation was used to further purify (A) CD81-positive or (B) CD9-positive exosomes from 15 μL preenriched samples (lanes 1–3). The isolated subpopulations were subjected to western blot analysis with antibodies against CD81 and CD9. M: molecular weight marker. C: Control with 7.5 μL preenriched exosomes without further isolation.








Figure 3. Uptake by HeLa cells of labeled exosomes. An Invitrogen™ FLoid™ Cell Imaging Station was used. Red: Invitrogen™ Alexa Fluor™ 594 phalloidin; blue: DAPI; green: Invitrogen™ SYTO™ RNASelect™ stain.

## Find out more at thermofisher.com/exosomes

