




Unleash the
value!

ThermoFisher
SCIENTIFIC

**GlycanAssure™: A High Throughput
& High Resolution Glycan Analysis Platform**

- Introduction
 - Glycosylation and Glycan Diversity
 - Importance of Protein Glycosylation
 - Glycan Labeling Chemistry
 - Current Glycan Analysis Challenges
- Thermo Fisher GlycanAssure Platform
 - Sample Prep
 - Instrumentation
 - Software
 - Multicap CE for Glycan Analysis – Literature
 - GlycanAssure Workflow
- Sample Data Sets
 - Magnetic Bead Based Workflow vs. Carbon Column
 - Consistency Across Varying Glycoprotein Inputs
 - Reproducibility of CE Separation
 - Variability Across Capillaries
 - Variability Across Instruments & Capillary Arrays
 - Glycan Spike Studies
 - Improved Glycan Separation Using Thermo Fisher Dyes
 - Thermo Fisher 3500 Vs. Competition
- Summary

Glycans - Introduction



- Glycosylation – Attachment of Glycans (Carbohydrates or Oligosaccharides) to proteins
- Glycans are made up of Monosaccharides
- More than 70% of biotherapeutics are glycosylated
- Glycans are important for protein stability, folding, assembly, signaling, etc.
- Glycosylation is a critical quality attribute (CQA) of biotherapeutics

Glycan Diversity Poses an Analytical Challenge


Oligomer	Composition	Peptide	Oligosaccharide
Dimer	AA / AB	1 / 2	11 / 20
Trimer	AAA / ABC	1 / 6	120 / 720
Tetramer	AAAA / ABCD	1 / 24	1424 / 34560
Pentamer	AAAAA / ABCDE	1 / 120	17872 / 2144640

Glycoconjugate Biosynthesis

- Glycans are not templated
- Glycan structure is determined by sequential glycosyltransferase action

More than 50% of Human Proteins are Glycosylated

- Glycosylation is Heterogeneous
 - Structures of Sugars Attached
 - Sites To Which They are attached
 - Anomericity



Many Conditions Affecting Humans Involve Glycans

Pregnancy


N-glycosylation of human serum transcortin and thyroxine-binding globulin

Aging

Outer arm galactosylation of human serum IgG

Alcoholism

Increase in relative number of serum asialo-transferrin glycoforms



Cancer

Unique oligosaccharides
Such as CA 125 (Ovarian)
CA 19-9 (Colon)

Liver Disease


Increase in number of sialic acid residues on fibrinogens

Tuberculosis

agalactosyl IgG

Ref: Dwek et al., Ann. Rev. Biochem., 57 (1988) 785-838

Biotherapeutics: Glycosylation is a CQA



Ensures glycans/glycosylation is well characterized in biotherapeutics

- **FDA (USA) and EMA (Europe) define the regulatory requirements**
- **Biotherapeutic manufactures are legally obligated to comply with these**

US FDA, November 3, 2010

Many complex biologics have sugar chains that impact their efficacy and safety, and industry experts have noted that the FDA has made several statements over the years suggesting a primary concern that biosimilars have the same sugar structure for interchangeability and those sugar structures do not change with time.

EMA Guidelines 2009

On the development, production, characterization and specifications for monoclonal antibodies and related products 'glycan structures should be characterized, and particular attention should be paid to their degree of of mannosylation, galatosylation, fucosylation and sialylation'.

European Medicines Agency
EMA/CHMP/BWP/157653/2007 (2009)

Draft Guidance for Analysis of Biosimilars

As monoclonal antibody molecules come off patent, biosimilars will start entering the market.

FDA (USA) and EMA (Europe) consider the degree of glycosylation to be a critical factor in determining the degree of “similarity” to the original approved drug

Similar: Additional analytical data or other studies are necessary to determine if observed differences are within an acceptable range to consider the proposed biosimilar. For example, the agency says that “glycosylation plays an important role in the PK of certain protein products. Manufacturing process conditions may impact glycosylation. Comparative PK and PD studies of the proposed biosimilar product and the reference product help resolve that some difference in **glycosylation identified in analytical studies would be within an acceptable range to consider the proposed biosimilar product to be highly similar to the reference product.**”

Reference:

Draft Guidance for Industry - Clinical Pharmacology Data to Support a Demonstration of Biosimilarity to a Reference Product

U.S. Department of Health and Human Services, Food and Drug Administration


Center for Drug Evaluation and Research (CDER) & Center for Biologics Evaluation and Research (CBER)

May 2014 Biosimilars

European Medicines Agency EMEA/CHMP/BWP/157653/2007 (2009)



Which Glycans are Analyzed?

- **N-linked glycans**
 - Attached to the amide side chain of Asparagine in peptide sequence: **Asn- AA_x-Ser/Thr**
- **O-linked glycans**
 - Attached to the hydroxyl group on Serine or Threonine residues
- Other forms (less studied)
 - Glycophosphatidylinositol anchors
 - Attached to protein C-terminus
 - C-glycosylation on tryptophan residues
 - S-linked glycosylation
 - Linked through S atom in Cys and Met residues



Asparagine (Asn)


N-Linked



Serine (Ser)


Threonine (Thr)

O-Linked




Our method (and almost all current CE/HPLC methods) analyze ONLY N-glycans

Glycan Labeling Chemistries



How N-Glycans are Analyzed?

IgG (mAb)



Current Glycan Analysis Challenges



- Single channel separation with low throughput (HPLC/UHPLC & single cap CE)
- Poor data quality of high throughput methods (Chip based instruments)
- Labor Intensive workflow (multiple pipetting steps, use of spin columns, etc.)
- Use of toxic Sodium cyanoborohydride chemistry
- Use of lengthy vacuum centrifugation steps
- Generic software taking long analysis time
- Non-integrated solution from multiple vendors
- Commercial sample prep kits with high cost per sample
- No validation support for result/solution

Thermo Fisher GlycanAssure



- Very high throughput (96 samples/7-9hrs)
- No compromise in CE separation time
- Easy magnetic bead based workflow
- Fewer pipetting steps & hands-on time
- Multiple fluorescent dyes for glycan labeling
- No use of Sodium cyanoborohydride
- No vacuum centrifugation steps
- Parallel analysis of 24 samples on AB DNA Analyzer
- Fit for use “app” style software
- Complete integrated solution with sample prep kits, CE instrument & software
- Reduced cost per sample

First Fully Integrated System Combining Throughput & Data Quality

Product Overview




APTS Kit
PN: A28676




Turquoise Kit
PN: A28678




Teal Kit
PN: A28677



3500 & 3500XL
PN: A30467 & A30556



Data Acquisition Software
PN: A30750



Data Analysis Software
PN: A30751

General

First in class glycan analysis platform that combines both high throughput and data quality

Specific*

- Faster
- Simpler
- High Quality
- Cost-effective
- Integrated solution

*Compared to other Commercial & Homebrew Products

GlycanAssure Sample Prep Kits



- Three fluorescent dyes with distinct properties
- Traditional APTS and two Thermo proprietary dyes
- Faster glycan labeling with Teal & Turquoise dyes (30min)
- Magnetic beads for glycan purification and excess dye removal
- No need for excess dye removal with Teal (96 samples in 7hrs)
- Each kit with PNGase-F enzyme and beads for 96 samples

9-Aminopyrene-1,4,6-trisulfonic acid (**APTS**)

GlycanAssure CE Instrumentation



Key Features:

- 8 & 24-capillary systems for medium & high throughput analysis
- Single-line, 505 nm, solid-state, long-life laser that utilizes a standard power supply and requires no heat-removal
- Powerful, integrated data acquisition and analysis software that provides real-time assessment of data
- Radio frequency identification (RFID) technology that tracks key consumables data and records administrative information
- Advanced multiplexing capabilities for glycan analysis with up to six unique dyes
- Unrivaled application flexibility—one array and one polymer are used for most applications
- Simple setup, operation, and maintenance



GlycanAssure CE Instrumentation




- Industry “Gold Standard” CE instrument
- Parallel analysis of 8 or 24 samples (**< 2min/sample**)
- Robust capillary array
- Calibration across capillaries using internal standard

GlycanAssure Software

- **Acquisition Dashboard**
 - Quick access to instrument status and data
- **Run Setup**
 - Simplified and Intuitive chevron based, step by step workflow
- **Reports**
 - Easily create reports with choice of templates
- **Library**
 - Stores historical information for plates, experiments, methods and reference
- **Analysis Dashboard**
 - Quick view of analyzed and un-analyzed samples. Allows to create data analysis projects
 - Quick access to Favorite Projects
- **Data Processing**
 - Quick comparison of analysis methods, alignment, normalization, and smoothing
 - Library of analysis methods and recommendations for glycan analysis
- **Analysis**
 - Simplified graphics interface for manual integrations



GlycanAssure: Data Acquisition Software



applied biosystems

Applied Biosystems GlycanAssure

Glycan Acquisition Software

Version 1.0

For Research Use Only. Not for use in diagnostic procedures. © 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. CO125724 0915

ThermoFisher
SCIENTIFIC

GlycanAssure: Data Acquisition Software

Experiment Status 🔄

Run Progress...
 Experiment Name: ---
 Status: Idle
 Estimated Completion Time: HH:MM:SS Monitor Run

Experiments
 Not Started **0/0**

Last Experiment
 09 Feb 2016 02:31:37

Experiments
 Last Week **2**

Experiments
 Last Month **20**

20160209 QCLabeld... 09 Feb 2016 | 02:31:37

20160205 Training AE 05 Feb 2016 | 03:57:55

SV P2 DAY 16 28 Jan 2016 | 10:35:05

SV P2 DAY 15 27 Jan 2016 | 11:32:37

SV P2 DAY 14 26 Jan 2016 | 10:2:41

Recent User Activity

New Plate Setup New Experiment Group

Instrument 🔄 ?

3500xL

Genetic Analyzer: Idle

Instrument Door: Close

Oven Status:

Laser: OFF Oven: OFF

EP: OFF Oven Door: Close

Detection Cell Temperature: **25.11 °C**

Oven Temperature: **30.23 °C**

Instrument Sensor Instrument Sensor

Pre-heat: Pre-heat

Spatial **Spectral**

Consumables 🔄 ?

Polymer Pouch

Polymer POP7

Lot Number 1509077

Part Number 4393714

Total Injections Allowed 50

Injections Done

Anode Buffer

Buffer Type ABC

Lot Number 1510077

Part Number 4393925

Total Runs Allowed 100

Injections Done

Cathode Buffer

Buffer Type CBC

Lot Number 1505284

Part Number 4408256

Total Runs Allowed 100

Injections Done

Capillary Array

50cm - 24 cap

Lot Number 49

Part Number 4404689

Total Runs Allowed 160

Injections Done

GlycanAssure: Data Analysis Software

applied
biosystems

Applied Biosystems™ GlycanAssure™ Data Analysis Software

Version 1.0

For Research Use Only. Not for use in diagnostic procedures. © 2016 Thermo Fisher Scientific Inc.
All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries
unless otherwise specified. CO020995 0116

ThermoFisher
SCIENTIFIC

GlycanAssure: Data Analysis Software

GlycanAssure™ Data Analysis v1.0 - Untitled [AFLP] - gm Is Logged In Database OPTIXE2

Data Analysis Software | Dashboard | Processing | Analysis | Report | Manage

1-20 of 40 Projects | Create New Project


Analyzed Projects: 39
Unanalyzed Projects: 1

Projects

Date Modified | Project Name | Advanced Search | Reset

 ABRF POSTE... 17-02-2016...	 ABRF POSTE... 17-02-2016...	 ABRF POSTE... 17-02-2016...	 abrf poste... 17-02-2016...	 ABRF POSTE... 17-02-2016...	 ABRF POSTE... 17-02-2016...	 PA3500 APT... 16-02-2016...	 Cinderella... 16-02-2016...
 Cinderella... 16-02-2016...	 CINDERELLA... 16-02-2016...	 CINDERELLA... 16-02-2016...	 Sugarcube ... 12-02-2016...	 Sugarcube ... 12-02-2016...	 Sugarcube ... 12-02-2016...	 Sugarcube ... 12-02-2016...	 Sugarcube ... 12-02-2016...
 Sugarcube ... 12-02-2016...	 Sugarcube ... 12-02-2016...	 sys val te... 11-02-2016...	 SYS VAL P2... 10-02-2016...				

GlycanAssure: Data Analysis Software



Reports of AB Multi Capillary CE for Glycan Analysis

Glycobiology vol. 11 no. 4 pp. 275–281, 2001

Ultrasensitive profiling and sequencing of N-linked oligosaccharides using standard DNA-sequencing equipment

Nico Callewaert, Steven Geysens, Francis Molemans, and Roland Contreras¹
Unit of Fundamental and Applied Molecular Biology, Department of Molecular Biology, Ghent University and Flanders Interuniversity Institute for Biotechnology, K.L.-Ledeganckstraat 35, 9000 Ghent, Belgium
Received on July 21, 2000; revised on November 16, 2000; accepted on November 16, 2000

(Suzu and D ioniza For phore glyca which

Display Settings: Abstract Send to:

J Proteome Res. 2010 Dec 3;9(12):6655-64. doi: 10.1021/pr100802f. Epub 2010 Nov 2.

Optimized workflow for preparation of APTS-labeled N-glycans allowing high-throughput analysis of human plasma glycomes using 48-channel multiplexed CGE-LIF.

Ruhaak LR¹, Hennig R, Huhn C, Borowiak M, Dolhain RJ, Deelder AM, Rapp E, Wührer M.

Author information

Abstract
High-throughput methods for oligosaccharide analysis are required when searching for glycan-based biomarkers. Next to mass spectrometry-based methods, which allow fast and reproducible analysis of such compounds, further separation-based techniques are needed, which allow for quantitative analysis. Here, an optimized sample preparation method for N-glycan-profiling by multiplexed capillary gel electrophoresis with laser-induced fluorescence detection (CGE-LIF) was developed, enabling high-throughput glycosylation analysis. First, glycans are released enzymatically from denatured plasma glycoproteins. Second, glycans are labeled with APTS using 2-picoline borane as a nontoxic and efficient reducing agent. Reaction conditions are optimized for a high labeling efficiency, short handling times, and only limited loss of sialic acids. Third, samples are subjected to hydrophilic interaction chromatography (HILIC) purification at the 96-well plate format. Subsequently, purified APTS-labeled N-glycans are analyzed

Electrophoresis 2006, 27, 1363–1367

Kay Vogel
Joachim Kuhn
Knut Kleesiek
Christian Götting

Institut für Laboratoriums- und Transfusionsmedizin, Herz- und Diabeteszentrum Nordrhein-Westfalen, Universitätsklinik der Ruhr-Universität Bochum, Bad Oeynhausen, Germany

Received August 8, 2005
Revised October 8, 2005
Accepted October 10, 2005

Short Communication

A novel ultra-sensitive method for the quantification of glycosaminoglycan disaccharides using an automated DNA sequencer

Analysis of glycosaminoglycans (GAGs) is a key to understanding alterations in extracellular matrix composition and cell membrane. In this report we describe the analysis of GAG disaccharide as an example of GAG disaccharide analysis using an automated DNA sequencer-assisted GAG disaccharide analysis. The presented methodology allows nanomolar quantification of sulfonic acid (APTS)-derived GAG disaccharides.

mAbs 6:1, 185–196; January/February 2014; © 2014 Landes Bioscience REPORT


High-throughput glycosylation analysis of therapeutic immunoglobulin G by capillary gel electrophoresis using a DNA analyzer

Dietmar Reusch,^{1,*} Markus Habeger,¹ Tobias Kailich,¹ Anna-Katharina Heidenreich,¹ Michael Kampe,¹ Patrick Bulau,¹ and Manfred Wührer^{2,3}

¹Pharma Biotech Development Penzberg; Roche Diagnostics GmbH; Penzberg, Germany; ²Center for Proteomics and Metabolomics; Leiden University Medical Center; Leiden, The Netherlands; ³Division of BioAnalytical Chemistry; Department of Chemistry and Pharmaceutical Sciences; VU University Amsterdam; Amsterdam, The Netherlands


Keywords: monoclonal antibody (mAb), IgG glycosylation, automation, multiplexed capillary electrophoresis, DNA analyzer, HILIC-UPLC, APTS labeling, LC-MS

GlycanAssure APTS Workflow




*Total 9hrs taking into account initial glycoprotein prep time and time between steps. Total time of the workflow is only 7hrs. CE takes 45min to analyze 24 samples

GlycanAssure: Teal Workflow



*Total 7hrs taking into account initial glycoprotein prep time and time between steps. Total time of the workflow is only 5hrs. CE takes 45min to analyze 24 samples


GlycanAssure: Turquoise Workflow



*Total 8hrs taking into account initial glycoprotein prep time and time between steps. Total time of the workflow is only 5.5hrs. CE takes 45min to analyze 24 samples

Magnetic Bead Workflow Vs Carbon Spin Columns

- Glycans purified by carbon column or magnetic beads in triplicate
- Glycans were labeled with APTS for separation on 3500xL




	Peak	1	2	3	4	5	6	7	8	9	10
Carbon column	Rel % Area	0.69%	3.11%	3.54%	0.20%	0.16%	0.16%	1.55%	0.29%	1.96%	0.68%
	CV%	5.12%	4.88%	5.75%	4.31%	5.38%	3.29%	3.80%	0.39%	2.52%	1.76%
Magnetic Beads	Rel % Area	0.68%	2.82%	3.09%	0.23%	0.16%	0.18%	1.66%	0.26%	1.94%	0.68%
	CV%	1.91%	3.63%	4.13%	1.69%	1.17%	1.52%	0.50%	2.28%	1.37%	1.12%

	Peak	11	12	13	14	15	16	17	18	19	20	21
Carbon column	Rel % Area	12.28%	6.49%	15.64%	0.51%	2.80%	18.15%	7.58%	5.25%	0.56%	15.95%	2.43%
	CV%	2.64%	1.76%	0.95%	18.63%	1.26%	1.19%	2.15%	1.41%	1.57%	1.68%	2.62%
Magnetic Beads	Rel % Area	12.77%	6.44%	15.76%	0.53%	2.66%	18.38%	7.67%	5.16%	0.54%	16.00%	2.39%
	CV%	0.88%	1.46%	1.07%	1.53%	0.85%	0.59%	1.00%	0.46%	2.69%	0.60%	0.64%

Consistent Glycan Quantitation


Peak #	10µg		50µg		100µg	
	%Area	%CV	%Area	%CV	%Area	%CV
1	0.90%	4.05%	0.72%	3.06%	0.64%	5.07%
2	0.61%	6.31%	0.45%	3.31%	0.39%	4.77%
3	1.64%	2.22%	1.69%	2.46%	1.66%	2.68%
4	0.42%	7.09%	0.42%	4.67%	0.38%	4.26%
5	8.54%	2.47%	8.31%	2.19%	8.18%	1.48%
6	0.76%	5.89%	0.55%	8.97%	0.42%	6.47%
7	22.51%	2.04%	23.07%	4.25%	23.71%	1.60%
8	0.63%	4.50%	0.71%	9.54%	0.71%	3.36%
9	3.73%	0.83%	3.82%	2.60%	3.85%	1.28%
10	23.52%	0.43%	23.33%	2.52%	23.84%	0.65%
11	9.98%	0.52%	10.16%	3.37%	9.91%	0.88%
12	5.37%	1.00%	5.46%	3.50%	5.42%	1.00%
13	0.43%	4.27%	0.44%	4.53%	0.43%	2.12%
14	19.58%	1.16%	19.53%	2.91%	19.20%	0.85%
15	1.38%	2.25%	1.33%	4.83%	1.27%	1.84%

Reproducibility of Glycan Separation



Variability Across Capillaries of an Array

- Mixture of glycans from purified human serum IgG and bovine Fetuin
- Labeled with APTS & run on 3500xL CE (24 capillaries)




	Average	STDEV	CV%
1	1.8%	0.03%	1.5%
2	3.9%	0.08%	2.1%
3	1.4%	0.02%	1.5%
4	0.6%	0.02%	3.1%
5	0.8%	0.02%	2.3%
6	5.8%	0.43%	7.4%
7	9.1%	0.19%	2.1%
8	2.3%	0.22%	9.8%
9	2.7%	0.02%	0.8%
10	2.7%	0.02%	0.7%

	Average	STDEV	CV%
11	1.7%	0.03%	2.0%
12	5.4%	0.08%	1.4%
13	0.7%	0.08%	11.3%
14	2.4%	0.04%	1.5%
15	0.3%	0.02%	4.8%
16	0.7%	0.01%	1.4%
17	8.2%	0.07%	0.9%
18	4.2%	0.04%	1.0%
19	11.0%	0.09%	0.8%
20	1.8%	0.03%	1.7%

	Average	STDEV	CV%
20	1.8%	0.03%	1.7%
21	11.7%	0.08%	0.7%
22	5.0%	0.05%	1.0%
23	3.6%	0.03%	0.8%
24	10.6%	0.09%	0.9%
25	1.7%	0.02%	1.1%

*Data analysis was done using GeneMapper software that won't allow manual integration


Variability Across Multiple CE Instruments & Arrays



	Peak	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
1 st 3500 32 Inj	Average	3.4%	0.8%	1.0%	10.5%	3.3%	3.6%	8.4%	4.4%	3.7%	8.0%	15.3%	1.0%	3.0%	23.5%	10.2%
	CV%	2.8%	2.2%	2.3%	1.1%	2.0%	2.8%	9.2%	1.1%	0.9%	0.5%	0.9%	3.6%	1.5%	1.0%	1.3%
2 nd 3500 24 Inj	Average	3.3%	0.8%	1.1%	9.7%	2.8%	3.7%	8.1%	4.5%	3.7%	8.2%	15.7%	1.0%	3.1%	24.0%	10.4%
	CV%	2.0%	1.9%	1.8%	1.3%	2.3%	1.7%	0.9%	0.6%	0.7%	0.7%	0.4%	3.3%	0.8%	0.4%	1.2%
3 rd 3500 24 Inj	Average	3.5%	0.8%	1.0%	11.1%	3.6%	3.6%	7.9%	4.4%	3.6%	7.9%	15.1%	1.0%	3.0%	23.3%	10.1%
	CV%	1.3%	3.0%	3.2%	1.2%	2.3%	1.7%	0.9%	1.1%	1.2%	0.7%	0.6%	6.1%	1.2%	0.7%	1.7%
Cap 1	Average	3.5%	0.8%	1.0%	11.1%	3.7%	3.7%	7.9%	4.4%	3.6%	7.9%	15.1%	1.0%	3.0%	23.3%	10.0%
	CV%	1.5%	3.8%	2.7%	0.9%	1.9%	1.3%	0.6%	1.0%	1.1%	0.7%	0.3%	3.0%	1.1%	0.5%	0.9%
Cap 2	Average	3.4%	0.8%	1.1%	10.9%	3.6%	3.8%	7.9%	4.4%	3.6%	7.9%	15.2%	1.1%	3.0%	23.3%	10.1%
	CV%	1.2%	2.4%	2.2%	0.9%	1.9%	1.8%	1.0%	0.9%	0.7%	0.5%	0.5%	1.8%	0.9%	0.4%	1.3%
Cap 3	Average	3.4%	0.8%	1.1%	11.0%	3.6%	3.8%	8.0%	4.4%	3.6%	7.9%	15.2%	1.1%	3.0%	23.3%	9.9%
	CV%	1.5%	3.8%	2.1%	1.2%	2.3%	1.2%	0.6%	0.6%	0.8%	0.6%	0.3%	2.8%	1.2%	0.5%	0.7%




Glycan Spike Recovery

- 15 diverse types of purified glycans labeled with APTS
- Mixed to create large difference in concentrations (50% to <1%)
- A2F (sialylated), Man7 (high mannose), and G2 (complex) glycans were increased serially by 2x




Increased signals of three specific glycans while the rest remain constant

Better Glycan Separation Using Thermo Fisher Dyes




Thermo Fisher 3500 Vs Leading Single Capillary CE


Single cap
CE



3500XL



3500XL



GlycanAssure: Summary

- **Easy sample prep**
 - Magnetic bead based sample prep
 - Hands-on-time <3 hrs for 96 samples
 - Less no. of pipetting steps
 - No use of Sodium cyano borohydride
 - No vacuum centrifugation steps
- **Throughput**
 - Sample prep & data of 96 samples in 7-9hrs
- **Resolution**
 - Sialyated glycans
 - Structural isomers
 - Fucose species
 - High Mannose species
- **Dye Labeling**
 - Multiple dyes with superior sensitivity
- **Sensitivity**
 - Low glycoprotein input
- **Low Cost of Analysis**
 - Robust instrument and capillaries with low running cost
- **Software**
 - Data Acquisition & Analysis software with novel features
- **Integration**
 - Full sample prep, hardware & software solution

