



Rapid In-Process Testing using MycoSEQ™ and ViralSEQ™ Assays


March 11, 2015

In-Process Testing of Cell Cultures in Biopharma Manufacturing



The Challenge

- Cell cultures and bioreactors can become contaminated with difficult to detect Mycoplasma or viruses
- The spread of an undetected contamination event can have a devastating impact
- Potentially high cost per test and long testing cycles as a result of contamination for traditional culture based tests are not practical for routine contamination monitoring
- Development and implementation of rapid methods enabling routine in-process testing require specialized expertise and experience



The Alternative

- Qualify and implement fully integrated, rapid qPCR testing as a risk mitigation tool.
 - Actionable results in hours.




The Solution

- MycoSEQ™ kits: rapid, sensitive qPCR Mycoplasma detection
- ViralSEQ™ kits: rapid, sensitive qPCR MMV and Vesivirus detection

These assays have been implemented globally for use in routine in-process cell culture monitoring.

The Threat Agents

Bacteria




Non-Mycoplasma

- Generally fast growing and easily detectable, often visually, as the Bioreactor becomes turbid and the media changes color
- One exception is *P. acnes*, which grows very slowly in both cell culture and in the traditional Bioburden test

Mycoplasma

- Grow very well in cell culture but are very difficult to detect using traditional culture methods
- Contamination events generally are often undetected until results of the 28-day culture test are analyzed
- Late detection can often result in contamination of subsequent lots either due to shared raw materials or from contaminated equipment


Virus



Viral

- Mouse Minute Virus (MMV) and Vesivirus were the causative agents in essentially all known CHO and NS0 viral contamination events
- They infect rodent cell culture very efficiently. Very difficult, perhaps impossible with Vesivirus, to detect using traditional culture test methods


What is in-process testing?




Potential Risks of Undetected Contamination

Impact to Subsequent Lots


Cell Culture Manufacturing Process



Cell Culture Manufacturing Process



Cell Culture Manufacturing Process



The Impact of Contamination Can Be Minimized

Contamination events can lead to



In-Process Testing allows you to minimize the impact

- Detect contamination at earliest possible point
 - Even more critical for hard-to-detect viruses using adventitious agent test
- Protect facility from the spread of contamination
 - Limit contamination to current lots
 - Avoid introduction of contaminant into downstream processes
 - Quarantine suspect raw materials
- Correct contamination quickly
 - Rapidly investigate and assess scope of contamination
 - Decontamination is much easier if isolated to single bioreactor

Critical Considerations

qPCR for In-Process Testing

Key Attribute	MycoSEQ and ViralSEQ Features
Intended Use	<ul style="list-style-type: none">• Assays and protocols designed for in-process testing
Sensitivity	<ul style="list-style-type: none">• Consistent across species: 1-10 Genome Copies/PCR reaction• Sensitivity not dependant on metabolic state of organism, 1 cell or viral particle contains 1 genome copy.
Specificity	<ul style="list-style-type: none">• Specific for organisms targeted• Proven specificity in multiple studies
Results Interpretation	<ul style="list-style-type: none">• Objective by using acceptance criteria• Cycle at Threshold (Ct): Input DNA quantity• Melting Temperature (Tm): Amplicon Size• Derivative Value (DV): Amplicon Quantity
Controls	<ul style="list-style-type: none">• PCR Inhibition Control• Negative Control• Positive Control, discriminatory• Internal Positive control for ViralSEQ™ assays
Support	<ul style="list-style-type: none">• Specialized Field Application Scientists: On-site training• Experience with validation design and regulatory support• Equipment Validation services: IQ/OQ• Results interpretation advice from team with deep domain knowledge.

3 Components of MycoSEQ™ Assay Sensitivity

Rapid sample extraction with optional automation



High Percentage nucleic acid recovery



Detection Using Highest Sensitivity qPCR Assay



**High Sensitivity
contaminant
Detection typically
in under 4 hours**

Rapid Protocols for the following:

- Bioreactors
- Cell Therapy
- Tissue Therapy
- Raw Materials
- Media
- Serum

PrepSEQ™ Magnetic Bead Nucleic Acid Preparation

- High percentage recovery
- Manual or Automated options

MycoSEQ™ or Viralseq™ Detection Assay

- 1-10 copy per PCR reaction sensitivity
- Consistent across species or strains

Assay Sensitivity as a Function of Test Sample Volume

For rapid in-process testing, we recommend a 100 uL test sample volume


Starting Sample Volume	Sample Prep Elution Vol.	Assay Sensitivity* (GC/Rxn)	Sample Volume (PCR reaction)	Sensitivity (CFU or GC/mL)
10 mL	100 uL	10	1 mL	10
1 mL	100 uL	10	100 uL	100
100 uL	100uL	10	10 uL	1000

*The assay sensitivity stated above is not the lower limit of detection of the assay. We recommend that validation experiments are done at 10 GC or CFU/mL in order to obtain consistent results.


Mouse Minute Virus Assay and Control Design

MMV Detection Assay


MMV viral target – FAM™ dye



Discriminatory PCR positive & extraction control – VIC™ dye





IPC control – NED™ dye



- Unique design includes detection of target, discriminatory control, and internal positive control (IPC)
- Assay targets conserved region of MMV genome, enabling detection of all known MMV strains
- The discriminatory positive control enables monitoring sample prep efficiency and RT-PCR reaction while minimizing the risk of a false positive result due to cross-contamination of samples
- The internal positive control present in every PCR reaction enables assessment of PCR efficiency and detection of potential inhibition in every reaction to help reduce risk of false negative results.

High Sensitivity, Broad Dynamic Range


MMV Detection Assay




- At least 7 orders of linear dynamic range
- Detection limit ~10 genome copies

Vesivirus Assay and Control Design

Vesivirus targets



Discriminatory RT-PCR positive control





Internal positive Control



- Unique design, 2 independent viral targets, helps reduce risk of false negative result due to genetic drift
- Assay targets conserved regions of Vesivirus genome, enabling detection of all known Vesivirus strains
- The discriminatory positive control enables monitoring sample prep efficiency and RT-PCR reaction while minimizing the risk of false positive result due to cross-contamination of samples by positive control template.
- The Internal Positive Control (IPC) present in every PCR reaction enables assessment of PCR efficiency and detection of potential inhibition in every reaction to help reduce risk of false negative results.

High Sensitivity, Broad Dynamic Range

Vesivirus Detection Assay




- At least 6 orders of linear dynamic range
- Detection limit ~10 genome copies



High Sensitivity, Broad Dynamic Range Positive Control Template

10 to 10⁶ copies of RNA per RT-PCR reaction

- Blue: Vesivirus probe
- Green: Discriminatory positive control probe
- Red: Internal positive control





- At least 6 orders of linear dynamic range
- Detection limit ~10 copies of control



Mycoplasma Assay Sensitivity

Mycoplasma arginini



- At least 6 orders of linear dynamic range
- LOD: 1 genome copy/ PCR reaction
- Amplicon melt analysis: Additional confirmation of positive result.

Broad Detection of *Mycoplasma* Species


Inclusion Panel (Partial)

<i>Acholeplasma granularum</i>	<i>Mycoplasma genitalium</i>	<i>Mycoplasma synoviae</i>
<i>Acholeplasma laidlawii</i>	<i>Mycoplasma gypis</i>	<i>Mycoplasma testudinis</i>
<i>Acholeplasma pleciae</i>	<i>Mycoplasma hominis</i>	<i>Mycoplasma timone</i>
<i>Mycoplasma alkalescens</i>	<i>Mycoplasma hyorhinis</i>	<i>Spiroplasma citri</i>
<i>Mycoplasma alvi</i>	<i>Mycoplasma imitans</i>	<i>Spiroplasma endosymbiont</i>
<i>Mycoplasma anseris</i>	<i>Mycoplasma indienne</i>	<i>Spiroplasma insolitum</i>
<i>Mycoplasma arginini</i>	<i>Mycoplasma lagogenitalium</i>	<i>Spiroplasma kunkelii</i>
<i>Mycoplasma auris</i>	<i>Mycoplasma lipofaciens</i>	<i>Spiroplasma melliferum</i>
<i>Mycoplasma buccale</i>	<i>Mycoplasma mobile</i>	<i>Spiroplasma mirum</i>
<i>Mycoplasma californicum</i>	<i>Mycoplasma molare</i>	<i>Spiroplasma phoeniceum</i>
<i>Mycoplasma canadense</i>	<i>Mycoplasma mycoides</i>	<i>Spiroplasma poulsonii</i>
<i>Mycoplasma capricolum</i>	<i>Mycoplasma neurolyticum</i>	<i>Spiroplasma sp.</i>
<i>Mycoplasma caviae</i>	<i>Mycoplasma orale</i>	<i>Mycoplasma bovirhinis</i>
<i>Mycoplasma collis</i>	<i>Mycoplasma phocidae</i>	<i>Mycoplasma bovis</i>
<i>Mycoplasma cricetuli</i>	<i>Mycoplasma pirum</i>	<i>Mycoplasma bovigenitalium</i>
<i>Mycoplasma equirhinis</i>	<i>Mycoplasma pneumoniae</i>	<i>Mycoplasma canis</i>
<i>Mycoplasma fermentans</i>	<i>Mycoplasma salivarium</i>	<i>Mycoplasma felis</i>
<i>Mycoplasma gallinaceum</i>	<i>Mycoplasma simbae</i>	<i>Mycoplasma fastidiosum</i>
<i>Mycoplasma gallisepticum</i>	<i>Mycoplasma sp.</i>	<i>Mycoplasma muris</i>
<i>Mycoplasma gateae</i>	<i>Mycoplasma spumans</i>	<i>Mycoplasma pulmonis</i>

- >90 species
- Related *Acholeplasma laidlawii* and *Spiroplasma citri*
- Common isolated species recommended for testing and validation (listed in red)

Discriminatory Positive / Extraction Control

- Positive control DNA made with *Mycoplasma* amplicon modified to have melting temperature (T_m) outside range of *Mycoplasma* amplicons
- Allows additional level of confirmation of positive test results.
- Higher T_m allows discrimination between true *Mycoplasma* and accidental contamination of test sample with positive control.
- Enables simple extraction control spiking of test samples.



Sample Preparation Automation

- Same chemistry as the well established PrepSEQ™ kits
- Fully-automated DNA extraction from broad sample types
- Enables high percentage nucleic acid recovery
- Low to medium throughput – 13 samples per 30-45 min
- Minimal hands-on time required
- Easy to use, set up and walk away
- Closed system: all reagents packaged in sealed cartridges




Simultaneous Extraction of *Mycoplasma*, MMV, and Vesivirus Nucleic Acid

Mycoplasma	MMV	Vesivirus
Bacteria (no cell wall)	Non-enveloped DNA virus	Non-enveloped RNA virus
Double-stranded DNA genome	Single-stranded DNA genome	Single stranded RNA genome
~ 1 Mbp	~ 5 kb	7 - 8 kb
SYBR™ Green assay	TaqMan™ assay	TaqMan™ assay
Culture >21 days (10 cfu/mL from Bioreactor harvest)	1 TCID ₅₀ = ~2x10 ⁴ genomes	Culture >21 days 10 TCID ₅₀ /mL


- *M. arginini*, MMV, and Vesivirus (FCV) were spiked into 100 uL test samples, either PBS or 10⁵ CHO cells.
- Nucleic Acid was extracted using PrepSEQ 1-2-3 extraction procedure.
- Nucleic acid eluted in 100 uL, 10 uL used per PCR reaction.
- Mycoplasma SYBR™ green assay, Vesivirus RT-PCR assay and MMV TaqMan™ assay were run on the same plate using a single thermocycling protocol.

Detection of MMV, Mycoplasma and Vesivirus extracted from the same test sample


MMV




FCV (Vesivirus)



Mycoplasma



- All three organisms were detected with high efficiency
- No difference in recovery and detection from either PBS or 10⁵ CHO cells



We Provide Comprehensive Support

- Training by Field Applications Scientist
- Instrument IQ/OQ services
- Sample prep optimization
- Validation study design
- Support with regulatory submissions including attendance at meetings and support answering agency questions
- Data Interpretation advice from scientists with extensive backgrounds in mammalian cell culture manufacturing and contamination events




Appendix

Impact of Viral Contamination

NEW YORK, June 16 (Reuters)

A major bio pharmaceutical company has **HALTED PRODUCTION** of two of its top-selling drugs after detecting a virus at one of their plants that U.S. regulators had already cited for deficiencies. The company, whose shares fell 4.2 percent after making the announcement on Tuesday, said current inventories of several key orphan disease treatments were not sufficient to meet global demand. The virus strain, Vesivirus 2117, has not been shown to cause human infection but is known to interfere with the growth of cells used to produce biologic drugs, the biotechnology company said.

PDA/FDA Adventitious Agents Conference, December 2010



“
*A negative result from the
28 day Points to Consider
Adventitious Agent test is
absolutely meaningless.*

—Ivar Klajvin, Genentech

- ”
- Presentations on Contamination Events:
 - Genentech, 2, MMV
 - Amgen, 3, MMV
 - Genzyme, 1, Vesivirus
 - Merrimack, 1, MMV
 - None of these contaminations were first detected by the 28 day Points to Consider Adventitious Agent test.
 - PCR was the key analytical tool used in all phases of these events / investigations / decontaminations.


Applications of Ct using the MycoSEQ™ and ViralSEQ™ Assays

The Cycle at Threshold (Ct) is a quantitative parameter representing the quantity of DNA present in the test sample at the start of the PCR process.

- Viability assessment
- Acceptance criteria for test and control samples
- Estimation of contamination level in positive samples
- Comparison of results between experiments
- GC/CFU or infectious titer ratios
- Sample preparation protocol optimization
- Validation Studies
 - Compare your Ct to expected values for each organism or species
 - Estimation of LLOD by extrapolation from values obtained at 10 GC or CFU/sample



Lower Limit of Detection Assessment Using Principles of qPCR



- For efficient, linear qPCR assays, 1 C_t difference represents a 2-fold difference in the starting quantity of the target DNA.
- The MycoSEQ™ assay is highly efficient and linear for all species tested.
- Example calculation of LLOD:
The mean C_t of *M.arginini* at 10 GC/PCR reaction from the previous examples = 31.
 - 5 GC/PCR reaction = 32
 - 2.5 GC/PCR reaction = 33
 - 1.25 GC/PCR reaction = 34
 - 0.6 GC/PCR reaction = 35


This approach to LLOD estimation is supported by experimental results

Results from External Validation Study

Level of Detection, Part 1, 10 GC/mL from 10 mL sample: Summary of Results

<i>Mycoplasma</i> Species (Type Strain)	Total Number Tests/Positive Reactions	% Positive	Mean C _T (n=24)	SD	CV (%)
<i>A. laidlawii</i> PG8 ^T	24/24	100	33.87	0.625	1.8
<i>M. arginini</i> G230 ^T	24/24	100	30.90	0.99	3.2
<i>M. fermentans</i> PG18 ^T	24/24	100	32.21	1.68	5.2
<i>M. hominis</i> PG21 ^T	24/24	100	29.53	0.86	2.9
<i>M. hyorhinae</i> BTS7 ^T	24/24	100	29.22	0.85	2.9
<i>M. orale</i> CH19299 ^T	24/24	100	31.85	1.81	5.7
<i>M. pneumoniae</i> FH ^T	24/24	100	33.03	0.73	2.2
<i>M. salivarium</i> PG20 ^T	24/24	100	31.14	0.87	2.8
<i>M. synoviae</i> WVU 1853 ^T	24/24	100	33.25	0.89	2.7
<i>S. citri</i> R8A2 ^T	24/24	100	32.79	1.65	5.0

Results Interpretation/Acceptance Criteria



Cycle at Threshold (Ct)

- A measure of target DNA level at the beginning of the PCR reaction.
- Acceptance criteria:
 - For Routine Testing: Ct less than or equal to 36.



Derivative Value (DV)

- A measure of specific amplicon quantity generated during PCR reaction.
- Acceptance Criteria: DV greater than 0.08.



Melting Temperature (Tm)

- A measure of amplicon size and base composition (known for *Mycoplasma* using this assay).
- Acceptance Criteria: Tm between 75 and 81 degrees C.

All three criteria must be met for a test sample to be positive for the presence of *Mycoplasma* DNA.


Is there a risk of a false positive result due to detection of DNA from non-viable organisms?

Contaminant Detection Assay Features

Highly Sensitive	<ul style="list-style-type: none">• <i>Mycoplasma</i>: 1-3 genome copies/PCR reaction• Vesivirus: 10 genome copies/PCR reaction• Mouse Minute Virus: 10 genome copies/PCR reaction
Specific	<ul style="list-style-type: none">• No detection of organisms, verified by testing panel of genetically related organisms
Comprehensive	<ul style="list-style-type: none">• Detection of all known species or strains
Quantitative PCR Technology	<ul style="list-style-type: none">• Enables viability assessment• Closed tube analysis, no handling amplified target DNA
Proprietary Discriminatory Positive/Extraction Controls	<ul style="list-style-type: none">• Helps minimize risk of false positive test result from accidental cross contamination with positive control.• Enables extraction/detection control analysis
Internal Positive Control (IPC)	<ul style="list-style-type: none">• For Viral assays, present in every PCR reaction. Enables assessment of PCR efficiency and detection of potential PCR inhibition in every reaction to help reduce risk of false negative results.
Sample Preparation	<ul style="list-style-type: none">• Single, rapid sample prep for isolation of Nucleic Acid from all organisms (DS-DNA, SS-DNA and RNA)• Versatile, protocols available for small or large test sample volumes• Fully Automated option: increase throughput, improve efficiencies, help reduce risk of analyst error

Mycoplasma

The relationship between a Colony Forming Unit (CFU) and a Genome Copy (GC)



Mycoplasma Cell

Ideally, 1 GC equals 1 CFU, but that must be verified. If the ratio changes because the sample or stock has cells that contain DNA, but do not grow and represent as a CFU....

The sensitivity of an NAT test may be misleading...

Legal Disclaimer

For Research Use Only. Not for use in diagnostic procedures.

- © 2015 Thermo Fisher Scientific Inc. All rights reserved. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.