

Lost Calibration Values on the GeneAmp® 9700


Calibration values can occasionally be lost from the GeneAmp® 9700 due to a variety of factors. If the calibration values have been lost, an error message will appear when attempting a run.

The following document will explain the steps required to obtain and re-enter the lost calibration values for your GeneAmp® 9700.




The following are directions for entering calibration values into a GeneAmp® 9700.

When calibration values have been lost an error message will be displayed.



To obtain and re-enter calibration values, perform the following procedure. Go to the “Details” display by pressing the F3 key to determine which calibration values have been lost.



The three types of calibration values that can be lost are the:

- block sensor values
- power amplifier values
- electronics values

In the example above the Block Sensor values have been lost. Once you have determined which values have been lost, you will need the serial number of both the top module and the bottom module. The top module serial can be obtained by pulling the yellow plastic lever in the back of the machine out to 90 degrees. The lever will pull out to 45 degrees easily, and then you need to pull it the rest of the way. This will allow the top to slide forward about ½ an inch. You then will flip the top module upside down and the serial number will be on the underside of the top module.

GeneAmp® 9700s with a serial number starting at xxxS905xxxxx for the top module or 805S905xxxxx for the base have the calibration values on the instrument. They can be found with the serial numbers. You should write down these values, which contain a slope and offset, before putting the instrument back together. These values will be used with the following instructions.

If the values on the calibration sticker are expressed in exponential form (ex. 2.075410e-4), convert them to decimal numbers to enter into the instrument.

(ex. 2.075410e-4 = .0002075410)

When entering calibration values, the instrument is limited to 7 characters. If your calibration value has more than 7 digits, enter only the first seven.

For customers in the United States and Canada, please send and email to technical support at <http://faqs.appliedbiosystems.com>. Please include in the email:

- a. the serial number of the machine
- b. your name
- c. phone number
- d. fax number
- e. email address
- f. your address

Once this information is received, your calibration values will be resent to you via e-mail or Fax.

For Non- US and Canadian customers, please contact your local support group. You may find the contact info at:

<http://www.appliedbiosystems.com/support/contact/internationalsupport.cfm>.

Once the calibration values have been obtained, continue with the following pages for directions on entering the calibration values.

Entering Block Sensor Calibration Values on the GeneAmp® PCR System 9700

Press F4(Util)

```

08:00 AM      09/25/00      25.0°C
GeneAmp® PCR System 9700
Version
Name: TC001      User: <<AB>>
Run  Create  Edit  Util  User
F1    F2    F3    F4    F5
    
```

Press F1(Diag)

```

Utilities
Diag - Instrument diagnostics
TmCalc - Calculates melting temp
Config - Instrument configuration
Diag  TmCalc  Config  More  Exit
F1    F2    F3    F4    F5
    
```

Hold down the **Stop** key and press F3(TempVer)

```

Diagnostics
Hard - Hardware Diagnostics
System - System Performance Tests
TempVer - Temperature Verification
Upgrad - Firmware Upgrade
Hard  System  TempVer  Upgrad  Exit
F1    F2    F3    F4    F5
    
```

Press the down arrow key twice to highlight **Calibrations**, then press F1(Run).

```

Service Functions
Service Info  Block Params  Final Test
Base Serial #  Manual Ctrl  Manual TEC
Calibrations  Port Mode  Other Diags
Init BRAM  Stability  Log
Run  Exit
F1    F2    F3    F4    F5
    
```

Select F2(Block)

```

Calibrations
Elect - Electronics calibration
Block - Block Sensor calibration
PwrAmp - Power Amplifier calibration
Print - Print the Base calibration
Elect  Block  PwrAmp  Print  Exit
F1    F2    F3    F4    F5
    
```

Enter the new calibration values.
When complete, Press F1(Accept).
Press F5(Exit) **four** times to return to the main menu.

```

Block Sensor Calibration
Resist  Temp °C  Slope  Offset
1098.319  30.720  +x.xxxxxx  -xxx.xxxx
Accept  Calib  +/-  Default  Cancel
F1      F2      F3      F4      F5
  
```

Entering Electronics Calibration Values for the GeneAmp® PCR System 9700

Press F4(Util)

```

08:00 AM  09/25/00  25.0°C
GeneAmp® PCR System 9700
Version
Name: TC001  User: <<AB>>
Run  Create  Edit  Util  User
F1   F2   F3   F4   F5
  
```

Press F1(Diag)

```

Utilities
Diag - Instrument diagnostics
TmCalc - Calculates melting temp
Config - Instrument configuration
Diag  TmCalc  Config  More  Exit
F1    F2    F3    F4    F5
  
```

Hold down the **Stop** key and press F3(TempVer)

```

Diagnostics
Hard - Hardware Diagnostics
System - System Performance Tests
TempVer - Temperature Verification
Upgrad - Firmware Upgrade
Hard  System  TempVer  Upgrad  Exit
F1    F2    F3    F4    F5
  
```

Press the down arrow key twice to highlight **Calibrations**, then press F1(Run).

```

Service Functions
Service Info  Block Params  Final Test
Base Serial #  Manual Ctrl  Manual TEC
Calibrations  Port Mode  Other Diags
Init BRAM  Stability  Log
Run  Exit
F1   F2   F3   F4   F5
  
```

Select F1(Elect)

```

Calibrations
Elect - Electronics calibration
Block - Block Sensor calibration
PwrAmp - Power Amplifier calibration
Print - Print the Base calibration
Elect  Block  PwrAmp  Prnt  Exit
F1    F2    F3    F4    F5
  
```

Enter the new calibration values. When complete, Press F1(Accept). Press F5(Exit) **four** times to return to the main menu.

Electronics Calibration				
	Resistance	Slope	Offset	
Block	1103.124	+x.xxxxxx	+xxxx.xxx	
Hsink	9531.847	+x.xxxx-xx	+x.xxxx-xx	
Cover	2070.292	+x.xxxx-xx	+x.xxxx-xx	
Accept	Calb	+/-	Default	Cancel
F1	F2	F3	F4	F5

Entering Power Amplifier Calibration Values on the GeneAmp® PCR System 9700

Press F4(Util)

08:00 AM	09/25/00	25.0°C		
GeneAmp® PCR System 9700				
Version				
Name: TC001	User: <<AB>>			
Run	Create	Edit	Util	User
F1	F2	F3	F4	F5

Press F1(Diag)

Utilities				
Diag	- Instrument diagnostics			
TmCalc	- Calculates melting temp			
Config	- Instrument configuration			
Diag	TmCalc	Config	More	Exit
F1	F2	F3	F4	F5

Hold down the **Stop** key and press F3(TempVer)

Diagnostics				
Hard	- Hardware Diagnostics			
System	- System Performance Tests			
TempVer	- Temperature Verification			
Upgrad	- Firmware Upgrade			
Hard	System	TempVer	Upgrad	Exit
F1	F2	F3	F4	F5

Press the down arrow key twice to highlight **Calibrations**, then press F1(Run).

Service Functions				
Service Info	Block Params	Final Test		
Base Serial #	Manual Ctrl	Manual TEC		
Calibrations	Port Mode	Other Diags		
Init BRAM	Stability	Log		
Run				Exit
F1	F2	F3	F4	F5

Select F3(PwrAmp)

Calibrations				
Elect	- Electronics calibration			
Block	- Block Sensor calibration			
PwrAmp	- Power Amplifier calibration			
Print	- Print the Base calibration			
Elect	Block	PwrAmp	Print	Exit
F1	F2	F3	F4	F5

Enter the new calibration values.
 When complete, Press
 F1(Accept). Press
 F5(Exit) **four** times to
 return to the main menu.

Iset	Imeas	+Set	+1.143e+02	+2.078s+03
+10	+10.3	-Set	+1.029e+02	+2.071e+03
+5	5.2	0Set	+1.097e+02	+2.093e+03
-5	-5.0			
-10	-9.7			Ver B
<input type="button" value="Accept"/>	<input type="button" value="More"/>	<input type="button" value="+/-"/>	<input type="button" value="Default"/>	<input type="button" value="Cancel"/>
F1	F2	F3	F4	F5

104GU16-01

4371000 Rev B