


Parainfluenza Virus Research Using a Multiplex Real-Time RT-PCR Method and the ViiA™ 7 Real-Time PCR System



LABORATORY FOR MEDICAL MICROBIOLOGY AND IMMUNOLOGY, ST. ELISABETH HOSPITAL, TILBURG, THE NETHERLANDS

Left to right: Hèlen van Raak, John Rossen, PhD, and Astrid van Steen

The Laboratory for Medical Microbiology and Immunology performs both microbiological and serological analyses for the mid of Brabant (The Netherlands) and for specific analyses outside of that region. The lab performs tests for the detection of microorganisms (including bacteria, fungi, yeast, viruses, and parasites) that cause infectious disease in humans. Approaches used to analyze these organisms include: direct identification and isolation of the organisms (culture, antigen tests); identification of antibodies against the organisms (serology); and molecular techniques to identify the organisms by targeting their DNA/RNA.

APPLICATION

Virus Research and Differentiation

TECHNOLOGIES

ViiA™ 7 Real-Time PCR System
7500 Fast Real-Time PCR System
TaqMan® Universal PCR Master Mix

Abstract

Real-time PCR-based molecular virology testing requires assays that are sensitive, specific, and that can distinguish between virus types. Here, researchers from the Laboratory for Medical Microbiology and Immunology, St. Elisabeth Hospital, Tilburg,

The Netherlands, used a multiplex reverse transcription PCR (m-RT-PCR) assay in a research-only study to detect and differentiate human parainfluenza viruses (HPIVs 1–4) in 30 samples derived from human nasopharyngeal washes. The assay was performed using either the new Applied Biosystems® ViiA™ 7 Real-Time PCR System or the Applied Biosystems® 7500 Fast Real-Time PCR System. The data indicate that the ViiA™ 7 Real-Time PCR System produces results equivalent to those of the 7500 Fast Real-Time PCR System and offers several new features that allow increased productivity.

Introduction

Human parainfluenza viruses (HPIVs) are nonsegmented RNA viruses that are members of the Paramyxoviridae family [1]. HPIV type 1 (HPIV-1) and HPIV-3 belong to the Respirivirus (former Paramyxovirus) genus, while HPIV-2 and HPIV-4 (consisting

of subtypes A and B [2]) belong to the Rubulavirus genus. These pathogens cause upper and lower respiratory tract illnesses including pneumonia, croup, and bronchiolitis in infants and children.

The availability of HPIV-specific detection assays is important because other respiratory pathogens cause similar illnesses. In recent years, molecular methods have resulted in sensitive, specific, and rapid detection of respiratory viruses [3]. In this research-only study, a multiplex real-time RT-PCR assay to detect and differentiate parainfluenza viruses (HPIVs 1–4) was tested on 30 human samples using the Applied Biosystems® 7500 Fast Real-Time PCR System. The assay was also performed using the new Applied Biosystems® ViiA™ 7 Real-Time PCR System to evaluate data reproducibility on the two systems.



Materials and Methods

Figure 1 shows the workflow used to analyze 30 samples for HPIVs.

Primer Design

Primer and probe sequences were selected using Primer Express® software (Applied Biosystems) from an alignment of nucleotide sequences of HPIV-1, HPIV-2, HPIV-3, and HPIV-4. Highly conserved regions were selected for each virus (Table 1). Parainfluenza virus 1, 2, and 3 primers and probes were designed on the hemagglutinin-neuraminidase gene [4]. Parainfluenza virus 4 primer and probe sequences were designed on the nucleocapsid gene [5].

Nucleic Acid Extraction

Nasopharyngeal washes were collected, and 200 µL of each sample was used to isolate DNA and RNA using the MagNA Pure LC Total Nucleic Acid Isolation Kit (Roche Diagnostics), as described by the manufacturer. Purified nucleic acid was eluted in a final volume of 100 µL.

Reverse Transcription

Following sample preparation, cDNA was synthesized using MultiScribe™ Reverse Transcriptase (Applied Biosystems) and random hexamers according to the manufacturer's instructions.

Multiplex Real-Time PCR

HPIV detection was performed in parallel using a multiplex real-time PCR assay specific for HPIVs 1–4. Assays were performed in triplicate. To help ensure reliable results, each plate included a negative isolation control and a positive amplification control. Real-time PCR was performed

on the 7500 Fast Real-Time PCR System using the standard instrument protocol. The experiments were repeated on the ViiA™ 7 Real-Time PCR System.

Results and Discussion

Researchers at the Laboratory for Medical Microbiology and Immunology currently use an Applied Biosystems® 7500 Fast Real-Time PCR System. This platform is widely used in virology and microbiology laboratories in the Netherlands and enables the use of several dyes, including Cy®5, for multiplex real-time PCR analysis. This multiplexing capability allows researchers to expand the number of targets amplified in a single reaction and increase their throughput. The Applied Biosystems® ViiA™ 7 Real-Time PCR System incorporates the new OptiFlex™ System, which offers a higher level of multiplexing capability with six decoupled excitation and emission filters for detection of up to 21 colors (Table 2).

Figure 2 shows the amplification of HPIVs 1–4 using the multiplex real-time RT-PCR assay and the ViiA™ 7 Real-Time PCR System. The multiplex strategy was efficient, enabling high throughput and fast results. In addition, this approach limits the complexity of the assay, which can ultimately reduce costs and complexity of the validation process.

Results obtained on the 7500 Fast Real-Time PCR System using a 96-well Fast block were reproducible on the ViiA™ 7 Real-Time PCR System using a 384-well block (Table 3). The difference in the threshold cycle (C_t) value observed is due to the differences in plate format and sample volume. The new ViiA™ 7 Real-Time PCR System has a user-friendly software interface that made plate setup and results analysis easy to learn. In addition, the software data analysis assigns quality control flags, which gave more confidence in the results interpretation.

| Features | ViiA™ 7 Real-Time PCR System |
|----------------------|--|
| Block Configurations | 96-well, Fast 96-well, 384-well, TaqMan® Array Micro Fluidic Cards |
| Run Time | 30 minutes expected (Fast 96-well) 35 minutes (384-well) |
| Resolution | Down to 1.5-fold change for singleplex reaction |
| Excitation Source | OptiFlex™ System with halogen lamp |
| Detection Channels | Decoupled—6 emission, 6 excitation |
| 21 CFR Part 11 | Optional software module |
| Remote Monitoring | Available to monitor up to 4 instruments in real time and the status of up to 15 instruments |
| Data Export Format | User configurable: *.xls, *.txt, and 7900 formats |



A photograph of the ViiA 7 Real-Time PCR System, a laboratory instrument with a blue and white color scheme and a touchscreen interface on top.

Figure 1. Multiplex Real-Time RT-PCR Assay Workflow.




Figure 2. Linear Amplification Plot Showing Detection of HPIVs 1–4 Using the Multiplex RT-PCR Assay on the ViiA™ 7 Real-Time PCR System. The differences in primer and probe concentrations used for each target explain the diverse amplification exponential trends.

Table 1. Primer and Probe Sequences for Target Genes. DFO = Dragon Fly Orange; YY = Yakima Yellow. BHQ1 and BHQ2 are two FRET-based quenchers (Black Hole Quencher).

| Parainfluenza Virus 1 | |
|-----------------------|--|
| PIV1-F2 | AAA AAC TTA GGG TTA AAG ACA ATC CA |
| PIV1-R2 | GCC AGA TGT RTG TCY TTC CTG CTG GT |
| PIV1-2-DFO-BHQ2 | DFO- CAAACGATGGCTGAAAAGGGA -BHQ2 |
| Parainfluenza Virus 2 | |
| PIV2-F2 | CCA TTT ACC TAA GTG ATG GAA |
| PIV2-R2a | CGT GGC ATA ATC TTC TTT TT |
| PIV2-R2b | TGT GGC ATA ATC TTC TTT CT |
| PIV2-2-YY-BHQ1 | YY- AATCGCAAAAGCTGTTCAGTCAC -BHQ1 |
| Parainfluenza Virus 3 | |
| PIV3-F2 | CAG GAA GCA TTG TRT CAT CTG T |
| PIV3-R2 | ATA GTG TGT AAT GCA GCT YGT |
| PIV3-2-FAM-BHQ1 | FAM- ACCCAGTCATAACTACTCAACAGCAAC -BHQ1 |
| Parainfluenza Virus 4 | |
| PIV4-F1 | CAA AYG ATC CAC AGC AAA GAT TC |
| PIV4-R1 | ATG TGG CCT GTA AGG AAA GCA |
| PIV4-1-Cy5-BHQ2 | Cy5- GTATCATCATCTGCCAAATCGGCAATTAACA -BHQ2 |

Table 2. Filter Combinations and Dye Compatibility on the ViiA™ 7 Real-Time PCR System.

| Channel | Dye Examples | Excitation Filter | Emission Filter |
|---------|--|-------------------|-----------------|
| 1 | FAM™ , SYBR® , SYTO® 9 (MeltDoctor™) , fluorescein, SYPRO® Orange | 470 ± 15 nm | 520 ± 15 nm |
| 2 | VIC® , JOE™ , TET™ , HEX™ | 520 ± 10 nm | 558 ± 12 nm |
| 3 | TAMRA™ , NED™ , BODIPY® TMR-X | 550 ± 10 nm | 586 ± 10 nm |
| 4 | ROX™ , Texas Red® | 580 ± 10 nm | 623 ± 14 nm |
| 5 | LIZ™ | 640 ± 10 nm | 682 ± 14 nm |
| 6 | Alexa Fluor®, Joda-4 | 662 ± 10 nm | 711 ± 12 nm |

Dyes indicated in **boldface** are Applied Biosystems calibration dyes.

Table 3. Comparison of Results Using the ViiA™ 7 Real-Time PCR System and the 7500 Fast Real-Time PCR System.

| Target Gene | 7500 Fast Real-Time PCR System— Positive Samples | ViiA™ 7 Real-Time PCR System— Positive Samples |
|-------------|---|---|
| HPIV-1 | 7 | 7 |
| HPIV-2 | 1 | 1 |
| HPIV-3 | 22 | 22 |
| HPIV-4 | 3 | 3 |

Conclusion

The multiplex real-time RT-PCR assay presented here provides a robust, sensitive, and easily automatable workflow to detect HPIVs 1–4. Results obtained using the 7500 Fast Real-Time PCR System were reproducible on the ViiA™ 7 Real-Time PCR System. In addition, the ViiA™ 7 System's user-friendly interface and data analysis features offered advantages to the workflow and additional confidence in the results. This multiplex assay may be a valuable approach for HPIV detection in clinical research settings.

Acknowledgments

We acknowledge Joris Verheyde, PhD, Applied Biosystems Field Application Scientist, and Harold Verbakel, Senior Technician of the Laboratory for Medical Microbiology and Immunology, St. Elisabeth Hospital, Tilburg, The Netherlands, for their generous contributions to this collaboration.

References

- Murphy FA, Fauquet CM, Bishop DHL et al. (ed.). (1995) Virus Taxonomy. Sixth Report of the International Committee on Taxonomy of Viruses. *Arch Virol Suppl.* 10:1–586.
- Canchola J, Vargosko AJ, Kim HW et al. (1964) Antigenic variation among new isolated strains of parainfluenza type 4 virus. *Am J Hyg* 79:357–364.
- Kuypers J, Wright N, Ferrenberg J et al. (2006) Comparison of real-time PCR assays with fluorescent-antibody assays for diagnosis of respiratory virus infections in children. *J Clin Microbiol* 44(7):2382–2388.
- Templeton KE, Scheltinga SA, Beersma MF et al. (2004) Rapid and sensitive method using multiplex real-time PCR for diagnosis of infections by influenza A and influenza B viruses, respiratory syncytial virus, and parainfluenza viruses 1, 2, 3, and 4. *J Clin Microbiol* 42(4):1564–1569.
- van de Pol AC, van Loon AM, Wolfs TF et al. (2007) Increased detection of respiratory syncytial virus, influenza viruses, parainfluenza viruses, and adenoviruses with real-time PCR in samples from patients with respiratory symptoms. *J Clin Microbiol* 45(7):2260–2262.

Life Technologies offers a breadth of products DNA | RNA | protein | cell culture | instruments

For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.

© 2010 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. TaqMan is a registered trademark of Roche Molecular Systems, Inc. Cy is a registered trademark of GE Healthcare. **CO16313 1110**

Headquarters

5791 Van Allen Way | Carlsbad, CA 92008 USA | Phone +1 760 603 7200 | Toll Free in North America 800 955 6288

www.lifetechnologies.com

