


invitrogen



GeneArt seamless cloning, assembly, and mutagenesis tools

Engineered for predictability

ThermoFisher
SCIENTIFIC

Cut, edit, create. Break through with GeneArt genetic assembly tools.

Invitrogen™ GeneArt™ products are developed to provide innovative, robust solutions to meet the vast needs of genetic engineers, synthetic biology researchers, and molecular biologists. We've made every effort to provide easy-to-implement workflow solutions to deliver fast, quality results.




Free online construct and primer design tool

Convenient and intuitive, the Invitrogen™ GeneArt™ design tool is accessible to all researchers—free of charge—and is built to ease your workflow and maximize downstream results.

Step by step, the online tool guides users through experimental design and ordering—designing oligos for assembly or mutagenesis of DNA molecules *in silico*.

Explore this one-of-a-kind online tool at thermofisher.com/GeneArt

Contents



Seamless cloning and assembly Flexible workflow solutions to streamline your research	4
High-order genetic assembly Works well with complicated cloning designs	5
Site-directed mutagenesis Conveniently mutate 1, 2, or 3 sites from a single plasmid	6
Ordering information	7

Seamless cloning and assembly

Flexible workflow solutions to streamline your research

Invitrogen™ GeneArt™ Seamless Cloning and Assembly Kits enable *in vitro* cloning of up to 4 DNA fragments simultaneously into virtually any linearized vector, typically in 30 minutes, without extra DNA sequences, restriction endonucleases, or ligation. With potential construct sizes of up to 40 kb, our kits offer researchers the flexibility and convenience to complete basic, standard, and advanced cloning and assembly protocols.

- **Flexible**—use any vector of your choice
- **Precise**—no scars; clone what you want, where you want it
- **Efficient**—>90% cloning efficiency
- **Convenient**—use our free online web tool to design oligos and assemble DNA molecules *in silico*
- **Fast**—fragment assembly typically in half a day

GeneArt Seamless Cloning and Assembly Kits use a proprietary enzyme mix to recognize and precisely assemble DNA fragments sharing end homology. End homology is created by PCR amplification using primers designed to generate the correct overlap between adjacent DNA fragments to be assembled.

Invitrogen™ GeneArt™ Seamless Cloning Enzyme Mix

- Chemically competent cell transformation
- Linear cloning vector included
- Good for high throughput (HTP)
- Maximum construct size of 13 kb

GeneArt Seamless Cloning and Assembly Kit

- Good cloning efficiency
- Chemically competent *E. coli* host and linear cloning vector included
- Maximum construct size of 13 kb

Invitrogen™ GeneArt™ Seamless PLUS Cloning and Assembly Kit

- Superior cloning efficiency
- Chemically competent *E. coli* host, media, linear cloning vector included
- Maximum construct size of 40 kb
- Good for HTP
- Conjugative gene transfer to most gram-negative bacteria



GeneArt seamless cloning and assembly workflow





Figure 1. Cloning efficiency of the GeneArt Seamless PLUS Cloning and Assembly Kit with 4 precloned DNA fragments of 5 kb each, compared to other manufacturer's kits.

High-order genetic assembly

In addition to seamless assembly and high efficiency, the vector-independent Invitrogen™ GeneArt™ High-Order Genetic Assembly System for *in vivo* assembly delivers on:

Speed—clone up to 10 DNA fragments simultaneously into a single vector (up to 110 kb total)

Convenience—assemble existing DNA fragments without restriction digestion or PCR amplification

The GeneArt High-Order Genetic Assembly System relies on yeast's ability to take up and recombine DNA fragments with high efficiency via transformation-associated recombination, greatly reducing *in vitro* handling of DNA and eliminating the need for restriction digestion and ligation.




Figure 2. Cloning efficiencies for 1, 3, and 10 DNA fragments, 10 kb each, amplified using Invitrogen™ Platinum™ PCR SuperMix High Fidelity and cloned into pYES1L using the GeneArt High-Order Genetic Assembly System.

Site-directed mutagenesis

Conveniently mutate 1, 2, or 3 sites from a single plasmid

Invitrogen™ GeneArt™ Site-Directed Mutagenesis Systems provide a simple and highly efficient method for *in vitro* site-directed mutagenesis. These systems can generate base substitutions, deletions, or insertions in plasmids purified from *E. coli*, without specialized vectors, host strains, or restriction sites—an ideal solution for routine or complex mutagenesis.

- **Flexible**—choose from single site or multisite directed mutagenesis products
- **Precise**—alter up to 25 nucleotides when only one site is mutated
- **Efficient**—mutagenesis efficiency >90% (using a 3 kb plasmid)
- **Convenient**—use our free online web tool to design oligos and mutate DNA molecules *in silico*
- **Fast**—time-to-results is typically less than 3 hours (using a 10 kb or smaller plasmid)
- **Flexible**—use any vector of your choice

GeneArt Site-Directed Mutagenesis System


- 1 site/plasmid
- Efficiency >90%

GeneArt Site-Directed Mutagenesis PLUS System


- Perform mutagenesis on 1, 2, or 3 sites—all in the same plasmid
- Efficiency >90%

The GeneArt Site-Directed Mutagenesis Systems utilize mutagenic oligonucleotide primers to generate mutations. The mutagenesis protocol is streamlined by combining DNA methylation and amplification steps into a single reaction, eliminating post-mutagenesis digestion and purification steps.

These systems deliver superior mutagenesis performance with a wide range of plasmid sizes.



Principle of method for multisite directed mutagenesis.



GeneArt Site-Directed Mutagenesis PLUS System vs. competitor.

Ordering information

Product	Quantity		Cat. No.																																
GeneArt Seamless Cloning and Assembly Systems																																			
GeneArt Seamless PLUS Cloning and Assembly Kit		20 reactions	A14603																																
<table><tr><th>Component</th><th>Size</th><th>Quantity</th></tr><tr><td>2X GeneArt Enzyme Mix</td><td>100 µL/tube</td><td>1 tube</td></tr><tr><td>Linear pYES7L Vector (50 ng/µL)</td><td>40 µL/tube</td><td>1 tube</td></tr><tr><td>One Shot DH10B T1 SA Cells</td><td>21 x 50 µL/tube</td><td>1 tube</td></tr><tr><td>Control lacZ Insert</td><td>5 µL/tube</td><td>1 tube</td></tr><tr><td>Linear pUC19L Vector (50 ng/µL)</td><td>40 µL/tube</td><td>1 tube</td></tr></table>	Component	Size	Quantity	2X GeneArt Enzyme Mix	100 µL/tube	1 tube	Linear pYES7L Vector (50 ng/µL)	40 µL/tube	1 tube	One Shot DH10B T1 SA Cells	21 x 50 µL/tube	1 tube	Control lacZ Insert	5 µL/tube	1 tube	Linear pUC19L Vector (50 ng/µL)	40 µL/tube	1 tube	<table><tr><th>Component</th><th>Size</th><th>Quantity</th></tr><tr><td>pUC19 Control (10 pg/µL)</td><td>10 µL/tube</td><td>1 tube</td></tr><tr><td>S.O.C. Medium</td><td>6 mL/bottle</td><td>1 bottle</td></tr><tr><td>Control Insert (50 ng/µL)</td><td>5 µL/tube</td><td>1 tube</td></tr><tr><td>Stbl3/pRK2013 Glycerol Stock*</td><td>1 mL/tube</td><td>1 tube</td></tr></table>	Component	Size	Quantity	pUC19 Control (10 pg/µL)	10 µL/tube	1 tube	S.O.C. Medium	6 mL/bottle	1 bottle	Control Insert (50 ng/µL)	5 µL/tube	1 tube	Stbl3/pRK2013 Glycerol Stock*	1 mL/tube	1 tube	
Component	Size	Quantity																																	
2X GeneArt Enzyme Mix	100 µL/tube	1 tube																																	
Linear pYES7L Vector (50 ng/µL)	40 µL/tube	1 tube																																	
One Shot DH10B T1 SA Cells	21 x 50 µL/tube	1 tube																																	
Control lacZ Insert	5 µL/tube	1 tube																																	
Linear pUC19L Vector (50 ng/µL)	40 µL/tube	1 tube																																	
Component	Size	Quantity																																	
pUC19 Control (10 pg/µL)	10 µL/tube	1 tube																																	
S.O.C. Medium	6 mL/bottle	1 bottle																																	
Control Insert (50 ng/µL)	5 µL/tube	1 tube																																	
Stbl3/pRK2013 Glycerol Stock*	1 mL/tube	1 tube																																	
GeneArt Seamless Cloning Enzyme Mix		20 reactions	A14606																																
<table><tr><th>Component</th><th>Size</th><th>Quantity</th></tr><tr><td>2X GeneArt Enzyme Mix</td><td>100 µL/tube</td><td>1 tube</td></tr><tr><td>Linear pUC19L Vector (50 ng/µL)</td><td>40 µL/tube</td><td>1 tube</td></tr><tr><td>Control Insert (50 ng/µL)</td><td>5 µL/tube</td><td>1 tube</td></tr></table>	Component	Size	Quantity	2X GeneArt Enzyme Mix	100 µL/tube	1 tube	Linear pUC19L Vector (50 ng/µL)	40 µL/tube	1 tube	Control Insert (50 ng/µL)	5 µL/tube	1 tube																							
Component	Size	Quantity																																	
2X GeneArt Enzyme Mix	100 µL/tube	1 tube																																	
Linear pUC19L Vector (50 ng/µL)	40 µL/tube	1 tube																																	
Control Insert (50 ng/µL)	5 µL/tube	1 tube																																	
GeneArt Seamless Cloning and Assembly Kit		20 reactions	A13288																																
<table><tr><th>Component</th><th>Size</th><th>Quantity</th></tr><tr><td>10X Enzyme Mix</td><td>45 µL/tube</td><td>1 tube</td></tr><tr><td>5X Enzyme Buffer</td><td>90 µL/tube</td><td>1 tube</td></tr><tr><td>Linear pUC19L Vector (50 ng/µL), 4 control reactions</td><td>8 µL/tube</td><td>1 tube</td></tr><tr><td>Control Insert (50 ng/µL)</td><td>5 µL/tube</td><td>1 tube</td></tr></table>	Component	Size	Quantity	10X Enzyme Mix	45 µL/tube	1 tube	5X Enzyme Buffer	90 µL/tube	1 tube	Linear pUC19L Vector (50 ng/µL), 4 control reactions	8 µL/tube	1 tube	Control Insert (50 ng/µL)	5 µL/tube	1 tube	<table><tr><th>Component</th><th>Size</th><th>Quantity</th></tr><tr><td>One Shot TOP10 Chemically Competent <i>E. coli</i></td><td>50 µL/tube</td><td>21 tubes</td></tr><tr><td>pUC19 Control (10 pg/µL)</td><td>10 µL/tube</td><td>1 tube</td></tr><tr><td>S.O.C. Medium</td><td>6 mL/bottle</td><td>1 bottle</td></tr></table>	Component	Size	Quantity	One Shot TOP10 Chemically Competent <i>E. coli</i>	50 µL/tube	21 tubes	pUC19 Control (10 pg/µL)	10 µL/tube	1 tube	S.O.C. Medium	6 mL/bottle	1 bottle							
Component	Size	Quantity																																	
10X Enzyme Mix	45 µL/tube	1 tube																																	
5X Enzyme Buffer	90 µL/tube	1 tube																																	
Linear pUC19L Vector (50 ng/µL), 4 control reactions	8 µL/tube	1 tube																																	
Control Insert (50 ng/µL)	5 µL/tube	1 tube																																	
Component	Size	Quantity																																	
One Shot TOP10 Chemically Competent <i>E. coli</i>	50 µL/tube	21 tubes																																	
pUC19 Control (10 pg/µL)	10 µL/tube	1 tube																																	
S.O.C. Medium	6 mL/bottle	1 bottle																																	
GeneArt Linear pUC19L Vector for Seamless Cloning		20 reactions	A13289																																
GeneArt High-Order Genetic Assembly Systems																																			
GeneArt High-Order Genetic Assembly System		10 reactions	A13285																																
GeneArt High-Order Genetic Assembly System (with Yeast Growth Medium)		10 reactions + 2 L medium	A13286																																
GeneArt High-Order Linear pYES1L Vector		10 reactions	A13287																																
CSM Medium for Mav203 Yeast Cells		2 L	A13292																																
GeneArt High-Order Vector Conversion Cassette		10 reactions	A13291																																
GeneArt Site-Directed Mutagenesis Systems																																			
GeneArt Site-Directed Mutagenesis PLUS System		10 reactions	A14604																																
<table><tr><th>Component</th><th>Size</th></tr><tr><td>One Shot MAX Efficiency DH5α T1</td><td>1 box</td></tr><tr><td>DNA Methylase (4 units/µL)</td><td>12 µL</td></tr><tr><td>200X SAM (S-adenosine methionine)</td><td>10 µL</td></tr><tr><td>2X GeneArt Enzyme Mix</td><td>100 µL</td></tr><tr><td>10X PCR Enhancer</td><td>100 µL</td></tr></table>	Component	Size	One Shot MAX Efficiency DH5α T1	1 box	DNA Methylase (4 units/µL)	12 µL	200X SAM (S-adenosine methionine)	10 µL	2X GeneArt Enzyme Mix	100 µL	10X PCR Enhancer	100 µL	<table><tr><th>Component</th><th>Size</th></tr><tr><td>PCR Water</td><td>1.8 mL</td></tr><tr><td>pMSDM-White Vector (20 ng/µL)</td><td>5 µL</td></tr><tr><td>0.5 M EDTA</td><td>500 µL</td></tr></table>	Component	Size	PCR Water	1.8 mL	pMSDM-White Vector (20 ng/µL)	5 µL	0.5 M EDTA	500 µL	<table><tr><th>Component</th><th>Size</th></tr><tr><td>Control Primer Mix 1 (10 µM)</td><td>25 µL</td></tr><tr><td>Control Primer Mix 2 (10 µM)</td><td>25 µL</td></tr><tr><td>Control Primer Mix 3 (10 µM)</td><td>25 µL</td></tr></table>	Component	Size	Control Primer Mix 1 (10 µM)	25 µL	Control Primer Mix 2 (10 µM)	25 µL	Control Primer Mix 3 (10 µM)	25 µL					
Component	Size																																		
One Shot MAX Efficiency DH5α T1	1 box																																		
DNA Methylase (4 units/µL)	12 µL																																		
200X SAM (S-adenosine methionine)	10 µL																																		
2X GeneArt Enzyme Mix	100 µL																																		
10X PCR Enhancer	100 µL																																		
Component	Size																																		
PCR Water	1.8 mL																																		
pMSDM-White Vector (20 ng/µL)	5 µL																																		
0.5 M EDTA	500 µL																																		
Component	Size																																		
Control Primer Mix 1 (10 µM)	25 µL																																		
Control Primer Mix 2 (10 µM)	25 µL																																		
Control Primer Mix 3 (10 µM)	25 µL																																		
GeneArt Site-Directed Mutagenesis System		16 reactions	A13282																																
<table><tr><th>Component</th><th>Size</th></tr><tr><td>DNA Methylase (4 units/µL)</td><td>20 µL</td></tr><tr><td>200X SAM (S-adenosine methionine)</td><td>10 µL</td></tr><tr><td>10X Enhancer</td><td>100 µL</td></tr><tr><td>0.5 MEDTA</td><td>500 µL</td></tr></table>	Component	Size	DNA Methylase (4 units/µL)	20 µL	200X SAM (S-adenosine methionine)	10 µL	10X Enhancer	100 µL	0.5 MEDTA	500 µL	<table><tr><th>Component</th><th>Size</th></tr><tr><td>pUC19WHITE Control Plasmid (20 ng/µL)</td><td>100 ng</td></tr><tr><td>Control Primer Mix (10 µM)</td><td>25 µL</td></tr><tr><td>PCR Water</td><td>1.8 mL</td></tr><tr><td>5X Reaction Buffer</td><td>90 µL</td></tr></table>	Component	Size	pUC19WHITE Control Plasmid (20 ng/µL)	100 ng	Control Primer Mix (10 µM)	25 µL	PCR Water	1.8 mL	5X Reaction Buffer	90 µL	<table><tr><th>Component</th><th>Size</th></tr><tr><td>10X Enzyme Mix</td><td>45 µL</td></tr><tr><td>One Shot MAX Efficiency DH5α T1</td><td>1 box</td></tr></table>	Component	Size	10X Enzyme Mix	45 µL	One Shot MAX Efficiency DH5α T1	1 box							
Component	Size																																		
DNA Methylase (4 units/µL)	20 µL																																		
200X SAM (S-adenosine methionine)	10 µL																																		
10X Enhancer	100 µL																																		
0.5 MEDTA	500 µL																																		
Component	Size																																		
pUC19WHITE Control Plasmid (20 ng/µL)	100 ng																																		
Control Primer Mix (10 µM)	25 µL																																		
PCR Water	1.8 mL																																		
5X Reaction Buffer	90 µL																																		
Component	Size																																		
10X Enzyme Mix	45 µL																																		
One Shot MAX Efficiency DH5α T1	1 box																																		

invitrogen

Find out more at thermofisher.com/geneart

For Research Use Only. Not for use in diagnostic procedures. © 2015 Thermo Fisher Scientific Inc. All rights reserved.
All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **CO020001 1215**

ThermoFisher
S C I E N T I F I C