



Molecular biology

Prepare for discovery

Molecular biology workflow solutions







Unlocking scientific discovery: making molecular biology accessible to all

From the beginning, our scientists have set out to commercialize tools to help advance molecular biology. In the early 1990s, Invitrogen™ SuperScript™ Reverse Transcriptase and TOPO™ cloning technology were introduced to the market. TOPO cloning technology was followed by many other products to speed up and simplify previously complex techniques. Over the decades, we have continued to develop new products and improve old favorites.

Helping our customers apply these tools to solve new research challenges is our purpose. Here we provide foundational information and introduce the latest tools to aid your research.

Find additional information at [thermofisher.com/molecularmasterpiece](https://www.thermofisher.com/molecularmasterpiece).

Contents



Sample preparation

Tools and methods	7
Automation platforms	10
Automation-compatible kits	11
DNA and RNA isolation	14



Reverse transcription

Considerations	17
Reagent selection	18
Genomic DNA removal	21
Primers	22



PCR

Thermal cycler considerations	25
Plastics essentials	28
PCR enzymes	32
Oligo design and selection	34



Electrophoresis

Workflow	38
E-Gel product selection	41
Electrophoresis reagents	42



Cloning

Technologies overview	45
Restriction enzyme cloning	46
PCR cloning	48
Cloning with synthetic DNA	50
Transformation	51



Isothermal amplification

Overview	55
LAMP	56
MDA-WGA	58
RCA	59
RPA	60




Resources

Educational resources	63
Mobile apps	64
Custom Commercial Supply	64
FAQs	65
Ordering information	67



Sample preparation






Nucleic acid isolation is a crucial first step in the molecular biology workflow, whether you are isolating genomic DNA (gDNA) or RNA. Selecting nucleic acid purification products that are optimized to help provide maximum yield, purity, and integrity from virtually any sample type and application is important for your research success.


Advance your research at [thermofisher.com/kingfisher](https://www.thermofisher.com/kingfisher)

Portfolio of solutions for your nucleic acid isolation



DNA type

For gDNA extraction, cfDNA isolation, plasmid isolation, and DNA cleanup



RNA type

For purification of total RNA, cell-free RNA (cfRNA), transcriptome RNA, messenger RNA (mRNA), microRNA (miRNA) and other small RNA, and sequence-specific RNA capture

Learn more at thermofisher.com/kingfisherkits



Understanding common nucleic acid isolation methods




Automated purification instruments: automated processing of magnetic particles in a microplate format (e.g., Thermo Scientific™ KingFisher™ purification systems)

Samples are processed by moving magnetic beads (not liquid). The system utilizes magnetic rods covered with a disposable, specially designed tip comb and plates. The instrument functions without any dispensing or aspiration parts or devices. Before the run, samples and reagents, including magnetic particles, are dispensed into plates according to default protocols that are installed on the instrument.

Benefits:

- Process 6–96 samples per run
- 24- or 96-well plates for different input volumes
- Easily edit, modify, or create new protocols
- All the benefits of magnetic beads (below)



Magnetic beads: 0.5–1.0 μm particles with a paramagnetic core and modified shell (e.g., Applied Biosystems™ MagMAX™ kits and Invitrogen™ Dynabeads™ magnetic beads)

Samples are lysed in solution and allowed to bind nucleic acid to magnetic particles based on specific surface modifications. Application of an external magnetic field rapidly collects the particles. Rounds of release, wash, and recapture enable purification of the desired nucleic acid.

Benefits:

- No risk of clogging
- Increased target capture efficiency
- Rapid collection and concentration of sample
- Specialized equipment not required
- Scalability



Helpful tip

Invitrogen™ Dynabeads™ magnetic beads are suitable for a variety of applications in addition to nucleic acid purification.



Visit thermofisher.com/dynabeads to learn about additional applications, including cell isolation, immunoprecipitation, and chromatin immunoprecipitation (ChIP).

Learn more at thermofisher.com/sampleprep





Spin columns: Glass fiber, derivatized silica, or ion exchange membrane in column (e.g., Thermo Scientific™ GeneJET™ and Invitrogen™ PureLink™ kits)

Samples are lysed and passed through the membrane using centrifugal or vacuum force. Wash and elution solutions are subsequently passed through the membrane, and the sample is collected into a tube by centrifugation.

Benefits:

- Convenience
- Ease of use
- Throughput flexibility
- Specialized equipment not required



Organic extraction: Phenol-chloroform solution (e.g., Invitrogen™ DNAzol™ and TRIzol™ Reagents)

After homogenizing the sample with TRIzol Reagent, chloroform is added and the mixture separates into a clear upper aqueous layer containing RNA, an interphase layer, and a pink lower organic layer containing the DNA and protein. RNA is precipitated from the upper aqueous layer with isopropanol. DNA is precipitated from the interphase and organic layers with ethanol. Protein is precipitated from the phenol-ethanol supernatant with isopropanol.

Benefits:

- Efficient lysis of cells and tissue
- Rapid denaturation of nucleases
- Stabilization of nucleic acids
- Great for fatty and cartilaginous samples

Automation platforms: Find a model that meets your needs

Optimize and automate your DNA and RNA, cell, and protein purification workflows with KingFisher systems. When used with compatible bead-based reagents such as MagMAX and Dynabeads products, these versatile instruments enable the automation of DNA, RNA, protein, and cell isolation procedures.

The Thermo Scientific™ KingFisher™ PlasmidPro™ Maxi Processor uses a cartridge-based system to automate the preparation of endotoxin-free plasmid suitable for transfection. Learn more and request a demo at thermofisher.com/kingfisher.



KingFisher instrument:	Duo Prime	Flex	Apex	Presto	PlasmidPro
Instrument size	Compact benchtop	Benchtop	Benchtop	Benchtop—integrates with robotic liquid handler	Benchtop
Throughput level	Low to medium	High	High	Ultrahigh	Low
Processing volume range	<ul style="list-style-type: none"> • 50–1,000 µL: 12-pin magnet head • 200–5,000 µL: 6-pin magnet head 	<ul style="list-style-type: none"> • 20–100 µL: skirted PCR plate • 20–200 µL: 96-well plate • 50–1,000 µL: 96 deep-well plate • 200–5,000 µL: 24 deep-well plate 	<ul style="list-style-type: none"> • 15–1,000 µL: 96 deep-well plate • 15–200 µL: 96-well KingFisher standard plate • 10–80 µL: 96-well PCR plate • 30–5,000 µL: 24 deep-well plate • 30–200 µL: 96 storage tubes • 200–1,000 µL: 24 storage tubes 	<ul style="list-style-type: none"> • 50–1,000 µL: 96 deep-well plate • 200–5,000 µL: 24 deep-well plate • 50–150 µL: KingFisher 96 plate 	<ul style="list-style-type: none"> • 150 mL overnight bacterial culture up to 6 OD or 10⁶ cells/mL in LB medium
Samples per run	6 or 12	24 or 96	24 or 96	24 or 96	One
Customizable protocols	Yes	Yes	Yes, with touchscreen or PC software	Yes	No
Heating and cooling	<ul style="list-style-type: none"> • 10°C to 75°C (plate row block A) • 4°C to 75°C (elution strip block) 	<ul style="list-style-type: none"> • From 5°C above ambient temperature to 115°C 	<ul style="list-style-type: none"> • From 4°C above ambient temperature to 100°C • Cooling down to 4°C 	<ul style="list-style-type: none"> • From 5°C above ambient temperature to 115°C 	NA
Ultraviolet lamp	8 W (up to 16 hr)	No	2 UV lamps, max 23 hr 59 min	No	No
Additional details	For Research Use Only. Not for use in diagnostic procedures.	For Laboratory Use	For Laboratory Use	For Laboratory Use	For Research Use Only. Not for use in diagnostic procedures. Designed for automated endotoxin-free maxi preps of plasmid DNA.

Typical run with a KingFisher instrument



Prep plates



Select program



Load plates



Press start



25–120 min**

Run time

** Can vary depending on application and instrument.

The graphic shows the expected run time for sample preparation when using a KingFisher instrument and prefilled plates. It takes less than 15 minutes to prep plates, 1 minute to select the program or protocol, and 1 minute to load plates. You can then press start and walk away for 25–120 minutes, depending on the sample and analyte type.



Resource

Use our selection tool to find the right magnetic bead-based kit for your automated sample preparation.

Find out more at

thermofisher.com/kingfisherkits

Learn more at thermofisher.com/kingfisher



Selecting the right DNA and RNA isolation kits for your downstream research

Cancer research: Applied Biosystems™ MagMAX™ cell-free nucleic acid kits on KingFisher instruments are excellent for liquid biopsy research. They are optimized specifically for enrichment of cfDNA and cell-free total nucleic acid (cfTNA), enabling highly reproducible recovery, flexible input volumes, and efficient protocols. Learn more about [sample preparation solutions for liquid biopsy](#).

Application	Automation-ready extraction kit and reagents	Cat. No.
Cancer research	MagMAX Cell-Free DNA Isolation Kit	A29319
	MagMAX Cell-Free Total Nucleic Acid Isolation Kit	A36716
	MagMAX <i>mirVana</i> Total RNA Isolation Kit	A27828
	MagMAX FFPE DNA/RNA Ultra Kit	A31881
	MagMAX DNA Multi-Sample Ultra 2.0 Kit	A36570
	MagMAX Sequential DNA/RNA Kit	A65309
	Dynabeads FlowComp Human CD3 Kit	11365D
	Dynabeads FlowComp Human CD4 Kit	11361D
	Dynabeads FlowComp Human CD8 Kit	11362D
	Dynabeads Untouched Human T Cells	11344D
	Dynabeads Untouched Human CD4 T Cells Kit	11346D
	Dynabeads CD15*	11137D
	Dynabeads CD4	11145D
	Dynabeads Epithelial Enrich	16102
Dynabeads MyOne CD45 Leukocyte Depletion	11170D	

* Available in 2 mL sizes with automation scripts. Learn more about [cell isolation products](#).

DNA research: Using the Applied Biosystems™ MagMAX™ DNA Multi-Sample Ultra 2.0 Kit and KingFisher instruments, you can isolate gDNA from 50 µL to 2 mL of whole blood, saliva, buffy coat, buccal swabs, or other biological samples. The resulting purified gDNA is ideal for many downstream molecular biology applications such as real-time PCR (qPCR), next-generation sequencing (NGS), microarray analysis, and other applications.

Application	Automation-ready extraction kit and reagents	Cat. No.
Genomics	MagMAX DNA Multi-Sample Ultra 2.0 Kit	A36570
	MagMAX Viral/Pathogen II Nucleic Acid Isolation Kit	A48383R
	MagMAX <i>mirVana</i> Total RNA Isolation Kit	A27828
NGS	Dynabeads Streptavidin for Target Enrichment	65606D
	MagMAX Pure Bind Beads	A58521



Go directly to our kit selection tool at thermofisher.com/kingfisherkits

Selecting the right DNA isolation kits for your downstream research

Infectious disease research: Applied Biosystems™ MagMAX™ viral/pathogen kits on KingFisher instruments provide a sensitive and simple method for nucleic acid extraction from samples containing viruses or other pathogens.

Learn more at thermofisher.com/mvpprime.

Application	Automation-ready extraction kit and reagents	Cat. No.
Infectious disease research	MagMAX Viral/Pathogen Nucleic Acid Isolation Kit for HIV-1 dried blood spots	A53770
	MagMAX Viral/Pathogen II Nucleic Acid Isolation Kit	A48383R
	MagMAX Microbiome Ultra Nucleic Acid Isolation Kit, with bead tubes	A42358
	MagMAX DNA Multi-Sample Ultra 2.0 Kit	A36570
	MagMAX Prime Viral/Pathogen NA Isolation Kit	A58145
	Dynabeads Intact Virus Enrichment	10700D



Plasmid purification: Choose from our wide range of high-performing, cost-effective Thermo Scientific™ and Invitrogen™ kits for plasmid DNA isolation designed to isolate plasmid DNA at the purity and scale you need.

Learn more at thermofisher.com/plasmid.

Application	Technologies for plasmid DNA isolation	Cat. No.
Plasmid purification	GeneJET Plasmid Miniprep Kit	K0502
	GeneJET Endo-Free Plasmid Maxiprep Kit	K0861
	PureLink Fast Low-Endotoxin Midi Plasmid Purification Kit	A35892
	PureLink HiPure Plasmid Filter Maxiprep Kit	K210016
	PureLink Expi Endotoxin-Free Maxi Plasmid Purification Kit	A33073
	MagMAX Pro HT NoSpin Plasmid Miniprep Kit	A58309
	KingFisher PlasmidPro Maxi Processor Endotoxin-Free Cartridge	A54072



Helpful tip

Select the best plasmid preparation kit for your application with the [Plasmid Purification Selection Tool](#).

Go directly to our kit selection tool at thermofisher.com/kingfisherkits



Selecting the right DNA isolation kits for your downstream research

Wastewater DNA and RNA for disease surveillance: Applied Biosystems™ MagMAX™ wastewater kits offer an efficient and simple method for extracting high-quality nucleic acids from wastewater, sewage, or sludge samples for disease surveillance workflows. The purified DNA and RNA are excellent for use in a variety of downstream applications such as quantitative PCR, digital PCR, and NGS.

Application	Automation-ready extraction kit and reagents	Cat. No.
Surveillance	MagMAX Wastewater Ultra Nucleic Acid Isolation Kit with Virus Enrichment	A52610
	MagMAX Viral/Pathogen II Nucleic Acid Isolation Kit	A48383R
	MagMAX <i>mir</i> Vana Total RNA Isolation Kit	A27828
	MagMAX Microbiome Ultra Nucleic Acid Isolation Kit, with bead tubes	A42358
	MagMAX Prime Viral/Pathogen NA Isolation Kit	A58145

Streptavidin: Utilization of Invitrogen™ Dynabeads™ streptavidin magnetic beads offers an exceptional balance of capacity and yield, reproducibility, purity, and cost for smaller-scale isolation of specific proteins (e.g., immunoprecipitation (IP)) and protein complexes (co-immunoprecipitation (co-IP)) and other applications.

Application	Automation-ready extraction kit and reagents	Cat. No.
IP and co-IP	Dynabeads Streptavidin for Target Enrichment	65606D
	Dynabeads M-270 Streptavidin	65305
	Dynabeads M-280 Streptavidin	11205D
	Dynabeads MyOne Streptavidin T1	65601
	Dynabeads MyOne Streptavidin C1	65001

Go directly to our kit selection tool at thermofisher.com/kingfisherkits



Selecting the right DNA and RNA isolation kits

For the quality and performance you need, a full suite of products for DNA and RNA isolation is available for a wide range of sample types, throughputs, and input quantities. To use our online kit selection guide, go to thermofisher.com/rnaselection.



Applied Biosystems™ and Invitrogen™ technologies for DNA and total RNA isolation

Capabilities	Process a large amount of tissue	Fast isolation of RNA from a variety of samples	High-throughput purification of RNA and DNA	Process cells for gene expression
Kits	TRIzol reagents	PureLink kits	MagMAX kits	Cells-to-C _T kits
Prep time	30–60 min	<20 min	45 min	≤10 min
Sample types	Most samples, particularly those more difficult to lyse	Bacteria, liquid, blood, cells, yeast, plants, tissue	Blood, plants, saliva, urine, stool, soil, plasma, serum*	Cultured cells
Starting material	100 mg of tissue or 10 ⁷ cells	10 ⁸ cells, 200 mg of animal tissue, 250 mg of plant tissue, 0.2 mL of blood, 5 x 10 ⁸ yeast, 10 ⁹ bacteria	Variable depending on sample	1–100,000 cells
Yield	10 ⁶ epithelial cells: 8–15 µg; 100 mg tobacco leaf: 73 µg (variable depending on sample)	Up to 350 µg	Variable depending on sample	NA
High throughput-compatible	No	Yes	Yes	Yes
Technology	Organic extraction	Silica membrane spin column/filter plate	Magnetic beads	Crude lysate

* Specialty kits with optimized chemistry are available for extraction of cell-free DNA/RNA, total RNA, gDNA, and total nucleic acid.




Helpful tip


If you are not ready to process your RNA sample, simply store it in Invitrogen™ RNA_{later}™ Stabilization Solution for use at a later time. Visit thermofisher.com/stabilizerna.

Learn more at thermofisher.com/rnapreps





Reverse transcription



Reverse transcription is the synthesis of complementary DNA (cDNA) from an RNA template using reverse transcriptase enzymes. The cDNA can be used as a template for PCR amplification, cDNA library construction, RNA sequencing, and more. Selecting the right reverse transcriptase is critical to detecting low-abundance RNAs in a sample and obtaining high yields of full-length cDNA.



Find technical resources on reverse transcription
at [thermofisher.com/rteeducation](https://www.thermofisher.com/rteeducation)



Considerations for selecting the right reverse transcriptase

Sensitivity, thermostability, processivity, and inhibitor tolerance of reverse transcriptases all affect the quantity and length of cDNA synthesized.

Sensitivity

The ability of a reverse transcriptase to generate cDNA from the least amount of input RNA is an important attribute when working with low-copy genes or difficult sample sources where RNA may have already degraded.

Thermostability


Thermostable reverse transcriptases allow reactions to occur at higher temperatures, which help denature RNA with complex secondary structures or high GC content, for generation of longer cDNA, higher cDNA yields, and better coverage of RNA populations in the cDNA.

Processivity

Processivity is the enzyme's ability to add consecutive nucleotides without releasing the template. Highly processive reverse transcriptases allow synthesis of longer cDNA strands in a shorter reaction time, and overall better efficiency in making full-length cDNA.

Inhibitor tolerance

Compounds that have inhibitory effects on reverse transcriptases are common in RNA samples even after purification. Their sources include reagents used for RNA isolation and contaminants carried over from biological samples. Reverse transcriptases resistant to common inhibitors help minimize inconsistent or suboptimal results in cDNA-based assays.



(A) Sensitivity, (B) thermostability, (C) processivity, and (D) inhibitor tolerance of reverse transcriptases can affect the quantity and length of cDNA. P, BR, Q, BL, T, and N represent enzymes from other suppliers.

Learn more at thermofisher.com/reverse-transcription

Reverse transcription reagent selection guide

Ability to optimize reaction components and conditions

Complete cDNA synthesis kit with all reaction components

Most convenient and fewest pipetting steps for RT-qPCR applications

Product format	Stand-alone enzyme	First-strand cDNA synthesis kit	First-strand cDNA synthesis master mix for RT-qPCR
Recommended product	Invitrogen™ SuperScript™ IV Reverse Transcriptase	Invitrogen™ SuperScript™ IV First-Strand Synthesis System	Invitrogen™ SuperScript™ IV VILO™ Master Mix
Applications	RT-PCR, RT-qPCR, cloning, cDNA library construction, RACE, RNA-Seq	RT-PCR, RT-qPCR, cloning, cDNA library construction, RACE, RNA-Seq	RT-qPCR
Total RNA input	1 pg–5 µg	1 pg–5 µg	0.01 pg–2.5 µg
Optimal reaction temperature	50–55°C	50–55°C	50–55°C
Reverse transcription time	10 min	10 min	10 min
High cDNA yield with challenging or degraded RNA	•	•	•
Includes PCR step			
Available formats	• Stand-alone enzyme	• cDNA synthesis kit • cDNA synthesis kit with ezDNase enzyme	• SuperScript IV VILO Master Mix • SuperScript IV VILO Master Mix with ezDNase enzyme



Did you know?

These products are free of OPE and NPE, making them safer for aquatic life.

Learn more at thermofisher.com/superscript



Most convenient
and fewest pipetting steps
for RT-PCR applications

Optimized template
switching for RACE and
RNA-Seq applications

cDNA synthesis and amplification
directly from intact single cells or
low amounts of total RNA

Go from mammalian cell lysate
to cDNA synthesis without
isolating RNA

	One-step RT-PCR kit	Template switching master mix	cDNA PreAmp kit	Direct RT kit
	Invitrogen™ SuperScript™ IV UniPrime™ One-Step RT-PCR System	Invitrogen™ SuperScript™ IV Template Switching RT Master Mix	Invitrogen™ SuperScript™ IV Single Cell/Low Input cDNA PreAmp Kit	Invitrogen™ SuperScript™ IV CellsDirect™ cDNA Synthesis Kit
	RT-PCR	RT-PCR, RACE, RNA-Seq	RT-qPCR, RNA-Seq	RT-PCR, RT-qPCR
	0.01 pg–1 µg	1–1,000 cells or 2 pg–10 ng	1–1,000 cells or 2 pg–10 ng	1–10,000 cells
	50–55°C	50°C	50°C	50–55°C
	10 min	30 min	30 min	10 min
	•	•		•
	•		•	
	<ul style="list-style-type: none"> • SuperScript IV UniPrime One-Step RT-PCR System (colored) • SuperScript IV UniPrime One-Step RT-PCR System (dye-free) 	<ul style="list-style-type: none"> • Invitrogen™ SuperScript™ IV Template Switching RT Master Mix 	<ul style="list-style-type: none"> • SuperScript IV Single Cell/Low Input cDNA PreAmp Kit 	<ul style="list-style-type: none"> • SuperScript IV CellsDirect cDNA Synthesis Kit



Did you know?

The standard enzyme format is incompatible with lyophilization because of the glycerol in the storage buffer. The lyo-ready (lyophilization-ready) format of SuperScript reverse transcriptases has a glycerol content below 0.1% and offers greater stability for lyophilized molecular assay kits.

Learn more at thermofisher.com/lyoreadyenzymes.



RNase inhibitors

To ensure your RNA template is protected and reverse transcription is successful, it's necessary to include RNase inhibitors. This step is crucial for preserving the integrity of the sample and achieving accurate and reproducible results.

Note: All SuperScript IV kit formats either contain an RNase inhibitor or are supplied with one.




Product	Invitrogen™ SUPERase•In™ RNase Inhibitor	Invitrogen™ RNaseOUT™ Recombinant Ribonuclease Inhibitor	Applied Biosystems™ RNase Inhibitor	Invitrogen™ Ambion™ RNase Inhibitor, cloned
Mechanism	Protein-based inhibitor of nonhuman origin that noncovalently binds to and inhibits the most common and troublesome RNases, including RNase A, B, C, 1, and T1	Noncompetitive inhibitor of pancreatic-type ribonucleases such as RNase A	50 kDa recombinant enzyme used to inhibit RNase activity	Recombinant human protein produced in <i>E. coli</i> and inhibitor of neutral pancreatic RNase A-type enzymes
Sizes available	<ul style="list-style-type: none"> • 2,500 units • 10,000 units 	<ul style="list-style-type: none"> • 5,000 units 	<ul style="list-style-type: none"> • 2,000 units 	<ul style="list-style-type: none"> • 1,000 units • 2,500 units • 10,000 units

Genomic DNA removal


RNA purification methods, including protocols with DNase digestion on column, often fail to completely remove gDNA. Amplification of contaminating gDNA can cause inaccurate results. Moreover, traditional gDNA decontamination protocols with DNase I include time-consuming DNase inactivation steps under conditions that can damage RNA and affect results.

Invitrogen™ SuperScript™ IV VILO™ Master Mix is available in a format with the novel dsDNA-specific Invitrogen™ ezDNase™ Enzyme, which enables efficient, fast, and gentle gDNA removal from RNA samples to help ensure high confidence in RT-PCR and RT-qPCR results.

SuperScript IV VILO Master Mix cDNA synthesis workflow with ezDNase Enzyme



Traditional cDNA synthesis workflow with DNase I



Learn more at thermofisher.com/4vilo



Reverse transcription primers

The priming strategy you choose for reverse transcription is important for cDNA synthesis efficiency, consistency, and yield. Each primer type has its benefits and drawbacks, depending on the individual target RNA.


For full-length first-strand cDNA synthesis, oligo(dT) primers are recommended because of their specificity for eukaryotic mRNA, and they allow many different targets to be studied from the same cDNA pool. Typically, oligo(dT) primers are strings of 12–20 deoxythymidines. We offer oligo(dT) in different lengths and formats for flexibility in your reverse transcription experiments.

For target mRNA containing strong transcriptional pauses, random primers are better suited because they anneal throughout the target molecules. They are also an ideal choice for nonpolyadenylated RNA, such as bacterial RNA.

For two-step RT-PCR, a mixture of oligo(dT) and random primers is often used to achieve the benefits of both primer types. SuperScript IV VILO master mix is provided with an optimized ratio of primers.



Two most common primers used in reverse transcription



Helpful tip

To avoid poly(A) slippage during priming, anchored oligo(dT) primers can be used to anneal to the 5' end of the poly(A) tail of mRNA and prevent priming within the poly(A) tail. Learn more about selection of primers for reverse transcription at [thermofisher.com/rteeducation](https://www.thermofisher.com/rteeducation).

Learn more at [thermofisher.com/rtpimers](https://www.thermofisher.com/rtpimers)



PCR



The polymerase chain reaction (PCR) amplifies DNA, allowing its detection from very low sample inputs. The process involves 3 steps that are repeated over multiple cycles:

- 1) **Denaturation:** double-stranded DNA is heated to separate the strands
- 2) **Annealing:** primers, short DNA molecules, bind to specific regions of the DNA
- 3) **Extension:** DNA polymerase extends the primers

Find technical and educational resources about PCR at [thermofisher.com/pcreducation](https://www.thermofisher.com/pcreducation).

Thermal cyclers

Thermal cyclers, which automate the heating and cooling cycles required to amplify DNA, play a critical role in the success of PCR. The following are things to consider when selecting a thermal cycler.

Precise temperature control

Thermal cyclers with precise temperature control enable you to quickly and accurately determine optimal annealing temperatures. Several block technologies, including gradient and Applied Biosystems™ VeriFlex™ Blocks temperature control, are available. A VeriFlex Block employs a separate heating and cooling element in each temperature zone, allowing better control and precision of temperatures. Learn more about the technology at thermofisher.com/veriflextechnology.

Reliability

Thermal cyclers should be able to withstand repeated use, environmental stress, and shipping conditions. Component reliability can be tested using robotic assemblies for repeated testing of frequently used instrument components such as the heated lid, touchscreens, and temperature cycling modules. Applied Biosystems™ thermal cyclers adhere to stringent reliability criteria, which are reported at thermofisher.com/thermalcyclerreliability.

Temperature accuracy

Thermal cycler temperature accuracy is a key factor in the success or failure of a PCR reaction. It is particularly important during annealing temperature optimization,

which requires both accuracy and consistency in the thermal cycler block. If the temperature set point of the instrument does not correspond to the actual temperature of the block, further temperature optimization could be required. Review a study of temperature accuracy in a number of models, available at thermofisher.com/thermalcycleraccuracy.

Features

A variety of Applied Biosystems thermal cyclers are available to fit your applications and budget. Certain features may be important to you, depending on your needs. If you perform PCR optimization frequently, you will likely benefit from an instrument with a VeriFlex Block. If you would like to run optimized assays on a new or different thermal cycler, you can save re-optimization time by using a simulation mode.

If you want remote access to your instrument, you will appreciate the convenience of cloud-enabled thermal cyclers. They allow you to design and share protocols, schedule an instrument, start or stop a run, and check run status from anywhere, on any mobile device or desktop computer.

Fleet control



If you manage multiple thermal cyclers and users, you may benefit from a single interface for viewing all instruments at a glance and setting custom permissions by instrument, user, and method. Learn more at thermofisher.com/fleetcontrol.



Helpful tip

Using the right PCR plastics for your application and instrument can improve the reliability of your PCR results. Go to thermofisher.com/findplastics to determine the right PCR plastics for you.

VeriFlex Block temperature zones



Learn more at thermofisher.com/thermalcyclers



Select the Applied Biosystems thermal cycler that's right for you



Questions for customers	ProFlex™ PCR System <ul style="list-style-type: none"> › Do you share the device with colleagues? › Do you expect your throughput needs to change? › Do you want to access your instrument remotely? 	VeritiPro™ Thermal Cycler <ul style="list-style-type: none"> › Do you perform a lot of optimizations? › Do you want to access your instrument remotely?
Key benefits	Ultimate flexibility and throughput	Ultimate performance
Max sample throughput	768 reactions (2 x 384-well block) 480,000 reactions (flat block)	384 reactions
Max block ramp rate	6.0°C/sec	6.0°C/sec
Temperature optimization	6-zone VeriFlex Block on 96-well system 2-zone VeriFlex Block on 3 x 32-well system	6-zone VeriFlex Block on 96-well system
Compatible with Fleet Control Software	Yes	Yes
3D tour	thermofisher.com/proflex3dtour	thermofisher.com/veritipro3dtour

Looking for a CE-IVD-labeled thermal cycler?
 Learn more about VeritiPro™ Dx Thermal Cycler at thermofisher.com/veritiprodx



Cloud-enabled



Automation-ready

MiniAmp™ Thermal Cycler

- > Do you want an instrument with just the features needed for routine PCR?
- > Do you want to access your instrument remotely?

Routine PCR, elevated

96 reactions

3.5°C/sec

3-zone VeriFlex Block on MiniAmp Plus model

Yes

Automated Thermal Cycler

- > Do you want to place your instrument on a robotic platform now or in the future?

Designed for easy robotic integration

384 reactions

3.5°C/sec

None

Yes

thermofisher.com/atc3dtour

Find out more at thermofisher.com/thermalcyclers



PCR and qPCR plastics, seals, and accessories

Since PCR is a sensitive detection method, PCR plastics must be of high quality and free of contaminants and inhibitors to enable optimal performance. Regardless of the plastics format, proper fit and uniform heat transfer during thermal cycling are essential.

Manufacturing quality control

Applied Biosystems™ PCR and qPCR plastic consumables are manufactured in world-class facilities dedicated to the production of high-quality molecular biology-grade plastics. After manufacturing, all plastics undergo strict quality control.

Integrity testing: Each well of every plate is visually inspected and tested for leaks. This thorough quality control screening verifies every well is intact to protect all reactions.

Evaporation testing: Samples are run through PCR to test sealing performance. Well liquid volumes are analyzed post-PCR to verify seal integrity. This helps ensure that every production lot conforms to strict tolerances.

Biological testing: All PCR plastics from Thermo Fisher Scientific are tested and guaranteed to be free of DNA, RNases, nuclease inhibitors, and PCR inhibitors. PCR certificates are available for convenience and documentation purposes, upon request.

Materials

Applied Biosystems™ MicroAmp™ optical microplates are made of polypropylene for optimal transfer of thermal energy for efficient PCR. A select medical-grade polypropylene is utilized for its exceptional biocompatibility and inert properties.


Applied Biosystems™ MicroAmp™ EnduraPlate™ microplates are constructed with a polycarbonate frame. This rigid design prevents damage from robotic grippers and provides better tolerance of rapid heating and cooling, while retaining thin-walled polypropylene wells for efficient heat transfer to the reaction mixture.

Thermo Scientific™ Sustain™ Series— reduce the environmental impact of your lab

Explore PCR plastics from the Sustain Series. Manufactured under ISO 9001 guidelines within a Class 100,000 cleanroom from a biobased medical-grade virgin polypropylene that is chemically and molecularly identical to existing Thermo Scientific™ PCR plastics, these products are certified by the International Sustainability and Carbon Certification (ISCC+) system. This production approach is more sustainable and is equivalent to 1.96x* lower CO₂ equivalent emissions per kg of biobased polypropylene resin compared to traditional fossil fuel-based polypropylene resin.

thermofisher.com/sustainplastics

Greener 
by design™



* Compared to equivalent fossil fuel-derived polypropylene resin emissions, cradle to gate (kg CO₂e).

Find out more at thermofisher.com/pcrplastics



Applied Biosystems PCR and qPCR plastics are validated and tested for reliability and optimal performance. They are “Engineer Approved” for use with all Applied Biosystems thermal cyclers and qPCR instruments, and are available in a variety of 32-, 48-, 96-, and 384-well plates; tube strips; single tubes; caps; and seals. The table below provides a detailed comparison of each product. Easily find the PCR and qPCR plastics compatible with your instrument using the online selection tool at thermofisher.com/findplastics.

	Small-scale experiments with a few samples	Daily experiments	Reliable and sustainable PCR plastics	Complete-workflow experiments—ideal for automation	Automation-compatible
PCR plastics product	Single tubes, strips, caps, adhesive film, and accessories	MicroAmp™ optical microplates	Sustain Series products	MicroAmp™ EnduraPlate™ optical microplates	MicroAmp™ EnduraPlate™ optical reaction plates
Formats	<ul style="list-style-type: none"> • Single tubes • Single tubes with caps • 8-strip tubes with caps • 12-strip caps 	<ul style="list-style-type: none"> • 32-well • 48-well Fast • 96-well* • 96-well Fast* • 384-well* 	<ul style="list-style-type: none"> • 96-well (multiple) • 384-well • Strip tubes • Strip tube caps • Strip tubes with attached caps 	<ul style="list-style-type: none"> • 96-well • 96-well Fast • 384-well • 96-well full-skirted 	<ul style="list-style-type: none"> • 96-well • 96-well Fast • 384-well
DNA-, RNase-, PCR inhibitor-free	Yes	Yes	Yes	Yes	Yes
Colors available	Clear, or mixed packs containing red, orange, blue, and green	Clear	Clear	Single-color packs (red, blue, green, yellow, or clear) and 5-plate sampler (one of each color)	Clear
Barcode available	No	Yes (1 or 2 sides)	No	Yes (3 sides)	Yes (3 sides)
Automation-compatible	No	Yes (for those with * above)	Yes	Yes	Yes



Did you know?

Low-profile plastics, also referred to as “Fast” tubes or plates, are generally required for fast (0.1 mL) thermal blocks. Fast plastics utilize lower volumes (0.1 mL) than the standard (0.2 mL) tubes or plates. The low profile minimizes the air space above the reaction, helping reduce the effects of evaporation and enhancing thermal conductivity. Learn more about PCR and qPCR plastics at thermofisher.com/pcrplastics-education.

PCR reagents

DNA polymerase is an essential component for PCR because of its key role in synthesizing new DNA strands. Because of the sensitive and specific nature of PCR, it is important to choose high-quality enzymes and reagents to produce optimal results. The following are things to consider when choosing PCR enzymes.

Specificity

Nonspecific amplification is one of the major hurdles in PCR, since it can drastically impact the yield and sensitivity of target amplification. One way to help reduce nonspecific amplification is through the use of a hot-start DNA polymerase, which utilizes an antibody or chemical modification so that the polymerase becomes active only at the high temperature of the denaturation step. In addition to improving specificity, a hot-start DNA polymerase increases yield and allows convenient room-temperature setup for high-throughput applications.

Thermostability

Since thermal cycling is a key feature of the conditions that enable the repetitive chain reaction of amplifying DNA, thermostability of the DNA polymerase to be used is also an important feature. Highly thermostable DNA polymerases are recommended for amplifying GC-rich or long templates that often require prolonged high-temperature reactions.

Fidelity

The fidelity, or proofreading capability, of a DNA polymerase is based on its 3' to 5' exonuclease activity, which corrects misincorporated nucleotides. This function is critical in applications such as cloning, sequencing, and site-directed mutagenesis, for accurate replication of DNA sequences.

Processivity

A DNA polymerase's processivity is defined as the number of nucleotides being incorporated in a single binding event. This property often reflects synthesis rate and speed, as well as affinity for its substrates. Therefore, highly processive DNA polymerases are beneficial to amplify challenging templates such as long, GC-rich, or inhibitor-containing DNA.

Primer annealing temperature

The primer annealing temperature of each DNA fragment to be amplified often needs optimization when designing a PCR protocol. To help simplify annealing and enable co-cycling of PCR assays, consider a DNA polymerase with a reaction buffer that allows a universal annealing temperature of 60°C for primers.



Did you know?

The residual bacterial DNA in recombinant PCR enzymes poses challenges in microbial genome analysis, such as accurately detecting bacterial strains by 16S rRNA gene sequences. To enable confidence and success in microbial PCR assays, choose PCR enzymes with controlled low levels of residual bacterial and human genomic DNA.

Find out more at

thermofisher.com/broad-range-pcr.



Helpful tip

Direct PCR is a way to help simplify PCR experiments, save time, and prevent sample loss in the workflow. Direct PCR allows you to amplify target sequences directly from the samples without the need to first isolate and purify the DNA.




Find out more at

thermofisher.com/direct-pcr.

Find out more at thermofisher.com/pcrenzymes



Antibody-based hot-start DNA polymerase and its activation in PCR for enhanced specificity



Helpful tip

One of the most common PCR troubleshooting issues is the presence of unwanted bands, or nonspecific amplification. To reduce nonspecific amplification:

- Use hot-start PCR
- Optimize annealing temperature
- Check primer design
- Prevent DNA cross-contamination
- Decrease template and/or primer concentration
- Optimize Mg^{2+} concentration



Did you know?

High-quality dNTPs are essential for successful PCR reactions. Thermo Fisher is one of the few primary manufacturers of nucleotides. Our high-quality nucleotides include dNTPs, NTPs, and modified nucleotides to support PCR and related applications. Building your own master mix? Learn more about our offerings and robust quality testing.

Find out more at thermofisher.com/dNTPs.



Choose the right PCR reagent for your research needs

A comprehensive portfolio of PCR enzymes and master mixes is available with the high performance and consistency you need. Start with the selection guide below to find the best enzyme for common PCR applications.

DNA polymerase	Invitrogen™ Platinum™ SuperFi™ II DNA Polymerase	Invitrogen™ Platinum™ II Taq Hot-Start DNA Polymerase	Applied Biosystems™ AmpliTaq Gold™ 360 DNA Polymerase	Invitrogen™ Platinum™ Direct PCR Universal Master Mix
PCR type	High-fidelity PCR	Hot-start PCR	Hot-start PCR	Direct PCR
Capabilities	Highly accurate amplicon sequences, universal primer annealing, robust amplification of difficult targets	Universal primer annealing, fast DNA synthesis, detection of low-abundance targets	Chemical hot start	Detection of target DNA without genomic DNA purification
Technical specifications				
Fidelity compared to Taq polymerase	>300x	1x	1x	1x
Target length	Up to 20 kb*	Up to 5 kb	Up to 5 kb	Up to 8 kb
Hot-start modification	Antibody-mediated	Antibody-mediated	Chemical modification	Antibody-mediated
Speed	15–30 sec/kb	15 sec/kb	60 sec/kb	20 sec/kb
Universal primer annealing	Yes	Yes	No	Yes
Inhibitor tolerance	Yes	Yes	No	Yes
Blunt or 3'-A end	Blunt	3'-A	3'-A	3'-A
Compatible with Applied Biosystems™ TaqMan™ probes	No	Yes	Yes	No
Certified low level of bacterial gDNA	Yes	Yes	Yes	No
Applications				
Cloning and subcloning	•			
Site-directed mutagenesis	•			
GC-rich amplification	•	•	•	•
Template generation for sequencing	•	•	•	•
High-throughput PCR	•	•		•
Long PCR (up to 20 kb)	•			
Genotyping	•	•	•	•
Amplification of samples with suboptimal purity	•	•		•
Colony PCR	•	•	•	•
Multiplex PCR	•	•	•	•
Fast PCR	•	•		•

* Amplification of up to 40 kb fragment sizes is possible, but may require additional optimization of reaction conditions and primer design.



Did you know?

These products are free of OPE and NPE, making them safer for aquatic life.

Learn more at thermofisher.com/pcrenzymes



Innovations for superior PCR

To help your research move faster, we have continually improved our PCR enzymes and reagents. For example, the latest Invitrogen™ Platinum™ DNA polymerases are designed with the following key innovative features.

More robust and versatile

Advanced enzymatic engineering and methodology provide DNA polymerases with fast cycling, high tolerance of PCR inhibitors, and efficient amplification of challenging DNA like GC-rich sequences. These features help you amplify DNA targets confidently with speed and simplicity.

Find out more at thermofisher.com/platinumenzymes.



Universal primer annealing

The innovative Invitrogen™ Platinum™ PCR buffers enable universal primer annealing at 60°C. This design allows you to co-cycle different PCR assays (instead of running them sequentially), drastically reducing tedious optimization steps and saving time.

Find out more at thermofisher.com/universalannealing.

Direct gel loading

The latest Invitrogen™ Platinum™ DNA polymerase master mixes are available in a green buffer format that allows direct gel loading and eliminates tedious steps of dye addition, helping to reduce pipetting errors. DNA migration is easily tracked with two dyes (blue and yellow) that are readily visible during electrophoresis (depicted in the lanes for 5 and 15 min in the figure to the right).



Custom DNA oligos

With a complete range of custom-synthesized oligonucleotide primers, probes, and genes, we recognize the need for quality, reliability, and convenience.


Invitrogen™ custom DNA oligos are synthesized on state-of-the-art automated systems to help increase performance, speed, and capacity. Available in a range of synthesis scales, purification options, and modifications, oligos are analyzed by mass spectrometry to help ensure the quality of all end products. This means you will receive high-quality custom DNA oligos quickly and efficiently.

Choose the right oligos and purification methods for your applications at [thermofisher.com/oligos](https://www.thermofisher.com/oligos).

Best practices for primer design

Good primer design is essential for a successful PCR assay. For design tips, review the infographic below or go to [thermofisher.com/primerdesign](https://www.thermofisher.com/primerdesign).

Length




Sequence

~~TCGCGATATATATATC~~
~~ATCGGGGGGTACTAG~~

Avoid low-complexity sequences in the template and primer, and runs of the same nucleotide

GC content (%)



Avoid self-dimers (complementarity of >3 bases)



Annealing temperature (T_m)



Avoid secondary structures in primer design



Find out more at [thermofisher.com/oligos](https://www.thermofisher.com/oligos)

Primer design made easy

Whether you're performing PCR, cloning, or capillary electrophoresis (CE) sequencing, take advantage of the benefits offered by our robust and easy-to-use Primer3-based Invitrogen™ OligoPerfect™ Designer.

- **Speed up**—design primers for up to 50 genes at the same time
- **Store your data**—ability to save your projects
- **Work smarter**—recognizes .txt and .fasta file types
- **Order with ease**—select and add primers directly to your online cart from the design tool*

Try the OligoPerfect Designer at thermofisher.com/oligoperfect or visit the oligo utility hub for our full suite of tools and calculators at thermofisher.com/oligotools.



Did you know?

Thermo Fisher Scientific offers scales beyond 10 μmol up to kilograms, and employs a team of manufacturing scientists dedicated to custom method development for unique modifications.

For information on large-scale and complex project capabilities, visit thermofisher.com/largescaleoligos.

Purification method	Description	Benefit	Applications	25 nmol	50 nmol	200 nmol	1 μmol	10 μmol
Desalt 25 nmol, 10–100 nt; 50 nmol, 5–100 nt	Oligos are processed through a normal-phase chromatography column, which removes salts but not truncated sequences	A salt-free DNA solution, ready to use; suitable for many PCR and sequencing applications without further purification	<ul style="list-style-type: none"> • Endpoint PCR • Isothermal sequencing • Fluorescent sequencing 	•	•	•	•	•
Cartridge 50 nmol–1 μmol, 7–55 nt	Based on reverse-phase chromatography; removes truncated sequences from the completed synthesis	Provides full-length sequences necessary for some applications	<ul style="list-style-type: none"> • Antisense oligos (ASO) • First-strand cDNA synthesis for generation of libraries 	NA	•	•	•	NA
HPLC ≥50 nmol, 10–55 nt; long oligo HPLC available	Reverse-phase high-performance liquid chromatography (HPLC) removes truncated sequences or unincorporated labels the same way as cartridge purification	Ensures highly purified primer required in some applications (≥85% full length)	<ul style="list-style-type: none"> • Fluorescent sequencing • Gel shift assays 	NA	•	•	•	•
PAGE ≥200 nmol, 7–100 nt	Polyacrylamide gel electrophoresis (PAGE) is a method used to differentiate full-length product from truncated sequences based on size and conformation	Provides the highest percentage of full-length oligos (≥85%) required for certain demanding applications such as mutagenesis or adapter production	<ul style="list-style-type: none"> • Microarrays • Amplified fragment length polymorphism (AFLP) analysis • PCR using oligos with critical 5' sequences (e.g., restriction endonuclease sites, RNA polymerase promoters) • Production of cloning adapters • Site-directed mutagenesis 	NA	NA	•	•	•


* Pricing and "Add to cart" feature may vary by geographic region. For questions, please contact your local office or distributor at thermofisher.com/contactus.



The background features a white space filled with various colorful shapes: a red bean-like shape at the top left, a blue ring-like structure on the left, a large blue cylindrical object with horizontal ridges at the top center, several blue and green spheres on the right, and a large orange textured shape at the bottom right. A green horizontal band spans the bottom of the page. Overlaid on this band and extending upwards are several vertical columns of horizontal bars in shades of blue, orange, and red, resembling a gel electrophoresis pattern.

Electrophoresis





Nucleic acid electrophoresis is a common technique in molecular biology used to separate, identify, quantify, and purify nucleic acids. Setting up electrophoresis involves a number of steps to achieve optimal separation and analysis of nucleic acid samples, such as gel preparation, ladder selection, sample visualization, and gel documentation.

Find technical resources on nucleic acid electrophoresis at [thermofisher.com/na-electrophoresis-education](https://www.thermofisher.com/na-electrophoresis-education)



Nucleic acid electrophoresis


Choosing the right tools for nucleic acid electrophoresis can significantly improve and accelerate results, enabling you to address downstream applications sooner.

Determining the appropriate gel type and gel concentration is an essential step that will help streamline the separation of nucleic acids. Learn more about convenient reagents for agarose gel electrophoresis, including hassle-free precast Invitrogen™ E-Gel™ Agarose Gels and pour-your-own Invitrogen™ UltraPure™ Agarose reagents, in this section.


If you need ...	Convenience, rapid results, and a safer workflow	High-quality reagents, a versatile workflow, and cost savings
Product	E-Gel Agarose Gels	UltraPure Agarose
Product format	Precast agarose cassettes	Powder
Buffer	Dry—none*	TBE or TAE
Protocol time (approx.)	18 min	120 min
Ready to use	Yes	No
Get more information at	thermofisher.com/egel	thermofisher.com/ultrapure

* Note: This is a dry precast electrophoresis system.

E-Gel precast agarose electrophoresis system



Traditional DNA electrophoresis workflow



Find out more at thermofisher.com/electrophoresis

Simplify electrophoresis with E-Gel precast agarose cassettes

E-Gel precast gels

Using precast agarose gels can simplify the nucleic acid electrophoresis workflow. Invitrogen™ E-Gel™ precast gels are self-contained and ready for use with the agarose, electrodes, and DNA stain packaged inside a disposable cassette. There are no gels to pour, buffers to make, staining or destaining steps to perform, or gel boxes to assemble. Just load your samples and run.

E-Gel precast gels offer excellent resolution and clarity in ≤ 18 minutes and are ideal for analyzing PCR products, restriction digests, plasmid preparations, and genotyping products. To help simplify cloning workflows, Invitrogen™ E-Gel™ CloneWell™ II gels use a double-comb design to enable recovery of purified DNA for downstream applications, without the need for additional purification kits or steps.


Find out more at thermofisher.com/egel.



E-Gel precast gels—faster and safer workflows


Restriction digest analysis:

- Fast and easy preparation of vectors for cloning experiments using FastDigest restriction enzymes and the Invitrogen™ E-Gel™ Power Snap Plus Electrophoresis System




RNA analysis:

- Optimization of loading conditions of RNA markers on an Invitrogen™ E-Gel™ system



Genotyping:

- Convenient and accurate visualization of DNA fragments, such as universal target fragments, for genotyping



Helpful tips

E-Gel precast gels are available in a variety of formats for routine and high-throughput applications, with different stains and agarose percentages (0.8%, 1%, 2%, and 4%). To find the right gel for your needs, see the selection guide at thermofisher.com/egelselection.

Choose Invitrogen™ E-Gel™ DNA ladders for precise electrophoresis band analysis with exceptional DNA fragment purity and quality, reduced dye masking, and improved ladder migration on E-Gel precast agarose gels. Find out more at thermofisher.com/egel-ladders.

Find out more at thermofisher.com/egel




E-Gel application examples

Fast and accurate analysis of RNA with E-Gel EX agarose gels


As RNA is the starting material for many workflows, its quality is important. Because RNA tends to form secondary structures and is prone to degradation, traditional methods of agarose gel electrophoresis are challenging. Traditional methods require additional time for buffer preparation, hand-pouring of gels, and prolonged run times. Because denaturing agents need to be added, gels are run in a fume hood.

With Invitrogen™ E-Gel™ EX Agarose Gels, several labor-intensive steps are removed from this process. The use of toxic denaturing reagents in the gel is no longer required. To denature the RNA, samples are treated with 50% formamide and heated at 70°C for 10 min, and then chilled on ice prior to loading on the E-Gel EX gel. The figure to the right shows clear, crisp bands from a diverse collection of RNA samples. Lanes 1–8 are samples, lanes M and 9 are ladders.



Genotyping detection of GMOs in soy products

Invitrogen™ E-Gel™ 48 and 96 Agarose Gels are useful for the rapid testing of larger numbers of samples. In this example, store-bought foods with unknown GMO status were tested. Multiplex PCR was used to amplify a known GMO sequence and an internal control (F1). The PCR products were separated and visualized on an E-Gel 48-well 2% agarose gel with Invitrogen™ SYBR™ Safe stain. Read the **application note** for more details. Two bands indicate the presence of genetic modification in many of the samples.



Helpful tips

Improve your nucleic acid electrophoresis. View our [webinar](#) Nucleic Acid Agarose Gel Electrophoresis: Techniques, Innovations, and Common Errors

Select your E-Gel electrophoresis device

Integrate your electrophoresis running and imaging into a single small device

Enjoy speed and convenience




Cloud-enabled
Ethernet connectivity



Which system is right for you?	Invitrogen™ E-Gel™ Power Snap Plus Electrophoresis System <ul style="list-style-type: none"> > Do you use electrophoresis often? > Would you like to store your gel images on internal servers or the cloud? > Do you want to perform quantitative analysis of your gels? > Do you value speed? > Would you like to store, share, and analyze gel images online? 	Invitrogen™ E-Gel™ Power Snap Electrophoresis System <ul style="list-style-type: none"> > Do you run less than 25 samples at a time? > Do you want a simple and fast solution for your electrophoresis? > Do you need to save bench space?
Key difference	Low- to high-throughput analysis	Low-throughput analysis
Memory	64 GB	32 GB
Connectivity	USB drive, ethernet, Wi-Fi, printer	USB drive
Sample throughput	Up to 96 samples	Up to 22 samples
Applications	Genotyping, fast PCR analysis, routine electrophoresis, and cloning	Routine electrophoresis and cloning
Software	Invitrogen™ iBright™ Analysis Software	NA
Find out more	thermofisher.com/powersnapplus	thermofisher.com/powersnap

Electrophoresis reagents

Choose high-quality agarose to pour your own gels. Select a ladder of the correct size range and an improved DNA stain for best detection and size estimation of DNA fragments.

DNA stains

Detection of nucleic acid samples in gels can be improved using fluorescent dyes that are safer or more sensitive than ethidium bromide. Invitrogen™ SYBR™ Safe and SYBR™ Gold stains provide greater safety or sensitivity with lower background fluorescence than the conventional ethidium bromide stain.

SYBR Safe stain is specifically formulated to be less hazardous than ethidium bromide and reduces your exposure to UV light.

Find out more at thermofisher.com/stains and thermofisher.com/sybrsafe.

UltraPure reagents for electrophoresis

Invitrogen™ UltraPure™ reagents are specifically formulated to help meet your nucleic acid analysis and purification needs. Invitrogen™ UltraPure™ Agarose and other reagents are made from highly pure biochemicals for maximum reliability and superior performance.

Find out more at thermofisher.com/ultrapure.

DNA ladders

Invitrogen™ DNA ladders are available in a wide variety of size ranges (10 to 48,502 bp) and formats for different applications. To create DNA ladders of exceptional quality, each fragment is purified individually using proprietary chromatography-based technology. Our DNA ladders are stable during prolonged storage at room temperature and after multiple freeze-thaw cycles.

Find out more at thermofisher.com/ladders.



Did you know?

Chromatographically purified nucleic acid fragments are considered the trusted standard for ladders since the technology provides high control over quality, banding pattern, intensity, and quantity for ladder composition.

Learn more at thermofisher.com/na-electrophoresis-education.


Fluorescent nucleic acid gel stains

	Safer detection	Ultimate detection
	SYBR Safe stain	SYBR Gold stain
Sensitivity (dsDNA)	Sensitive (>3 ng)	Ultrasensitive (>0.1 ng)
Less hazardous and more environmentally friendly	•	
Improved cloning efficiency	•	•



Cloning





Molecular cloning and synthetic DNA assembly are essential techniques in biotechnology for DNA sequence assembly. Restriction enzymes, TOPO PCR cloning, Gibson Assembly cloning, Golden Gate cloning, and Gateway cloning enable precise assembly of simple and complex constructs.

Find technical resources on molecular cloning at
thermofisher.com/cloningeducation



Cloning and gene synthesis

From restriction enzymes to gene synthesis, a large portfolio of tools and resources is available to help you obtain high-quality cloned DNA for your next discovery.

	Move fragments from one vector to another	Clone PCR fragments	Express genes in multiple systems	Create large inserts from multiple fragments		Assemble or clone multiple fragments	Sequenced, cloned genes
Method	Thermo Scientific™ FastDigest™ restriction enzymes	Invitrogen™ TOPO™ cloning	Invitrogen™ Gateway™ cloning	Invitrogen™ GeneArt™ seamless cloning and GeneArt™ Gibson Assembly™ cloning kits	Invitrogen™ GeneArt™ Type IIS assembly	Invitrogen™ GeneArt™ Strings™ DNA Fragments	Invitrogen™ GeneArt™ Gene Synthesis
Key benefits	<ul style="list-style-type: none"> Familiarity, flexibility, convenience, time savings Universal protocol and complete digestion in 5–15 minutes in one buffer 100% buffer compatibility with downstream applications Direct loading on gels 	<ul style="list-style-type: none"> >95% efficiency, 5-minute PCR cloning Available in a variety of formats and sizes 	<ul style="list-style-type: none"> High-throughput and high-efficiency shuttling among multiple expression vectors 	<ul style="list-style-type: none"> Seamless multifragment assembly by homologous recombination Directional cloning of up to 15 fragments Up to 95% efficiency and 15-minute cloning 	<ul style="list-style-type: none"> One-tube seamless multifragment assembly by simultaneous restriction digestion and ligation Directional cloning of up to 8 fragments, for up to 20 kb total Efficient for repetitive and very small sequences 	<ul style="list-style-type: none"> Synthesized DNA fragments ready to clone via the method of your choice No starting DNA required Pool sequence-verified 	<ul style="list-style-type: none"> Custom-cloned genes in your choice of vector Sequence-verified Can be optimized for a specific host for maximal protein expression No starting DNA required
Technology basics	<ul style="list-style-type: none"> Restriction digestion and ligation 	<ul style="list-style-type: none"> Topoisomerase-based, ligase-free cloning 	<ul style="list-style-type: none"> Single-step, directional, and site-specific DNA recombination Restriction enzyme- and ligase-free 	<ul style="list-style-type: none"> End-terminal homology recombination using overlapping sequences Transformation-associated recombination (TAR) in <i>Saccharomyces cerevisiae</i> 	<ul style="list-style-type: none"> Type IIS restriction and ligation in a single reaction 	<ul style="list-style-type: none"> Linear dsDNA assembled from pooled synthetic oligonucleotides 200–3,000 bp, also available in library format with randomized bases 	<ul style="list-style-type: none"> DNA of interest cloned in vector 100% sequence-verified with quality assurance documentation
Online tools available	FastDigest selection tool	TOPO PCR selection tool	Vector list	Invitrogen™ GeneArt™ Primer and Construct Design Tool		Invitrogen™ GeneArt™ Instant Designer for design and optimization	
Needs DNA source material (gene in plasmid, library, etc.)	•	•	•	•	•		
Use your own vector	•		*	•	•	•	•

* Vector needs to be converted with Invitrogen™ Gateway™ Vector Conversion System with One Shot™ *ccdB* Survival™ 2 T1^R Competent Cells.

Discover more at thermofisher.com/cloning



Restriction enzyme cloning

Found naturally in bacteria, restriction enzymes recognize and cleave specific DNA sequences. Cleavage results in “sticky” ends (5′ or 3′ protruding ends) or blunt ends that allow DNA inserts to be cloned into vectors with compatible ends. Star activity, buffer compatibility, and varying protocols for complete digestion are some common hurdles in restriction digestion.

FastDigest restriction enzymes

To simplify cloning, we offer FastDigest enzymes—an advanced line of restriction enzymes that share buffer compatibility with downstream modifying enzymes. Benefits include:

- Complete digestion in 5–15 min
- Universal buffer for multiple digestions for any combination of enzymes
- No sequential digestions and buffer changes
- 176 unique enzymes
- Direct loading of reaction mixture on gels



See our enzyme selection tool at thermofisher.com/fastdigest.

DNA/RNA-modifying enzymes

Complete your cloning reaction with modifying enzymes that have 100% activity in FastDigest and FastDigest Green buffers. These include:

- **Thermo Scientific™ T4 DNA Ligase**—used to join the ends of DNA fragments and vectors in a fast 10 min reaction
- **Thermo Scientific™ FastAP™ Thermosensitive Alkaline Phosphatase**—used for the dephosphorylation of cloning vectors to prevent recircularization. Fast 10 min reaction with complete inactivation in a 5 min reaction
- **Thermo Scientific™ T4 Polynucleotide Kinase**—phosphorylates PCR products, oligonucleotides, and other DNA prior to ligation reactions
- **Thermo Scientific™ DNA polymerases**—used for DNA blunting by filling in 5′-overhangs and second-strand synthesis of cDNA
- **Invitrogen™ CorrectASE™ enzyme**—utilized to enhance the accuracy of synthetic gene and DNA fragment assembly by removing mismatches

Visit thermofisher.com/modifyingenzymes.




Type IIS restriction enzymes

A specific group of restriction enzymes called Type IIS endonucleases cleave DNA outside of their recognition sequences. In combination with DNA ligase, Type IIS restriction enzymes are utilized to drive the insertion of one or several DNA fragments into a recipient vector without the inclusion of residual restriction enzyme sites and other unwanted DNA sequences at fragment junctions (scarless cloning). Commonly referred to as Golden Gate cloning, this method is useful for the construction of large DNA fragments.


Find FastDigest Type IIS enzymes at [thermofisher.com/fastdigesttypeiis](https://www.thermofisher.com/fastdigesttypeiis).

For GeneArt Type IIS Assembly Kits, go to [thermofisher.com/typeiis](https://www.thermofisher.com/typeiis).

Type IIP restriction enzyme



Type IIS restriction enzyme



Did you know?

Both Gateway and Type IIS restriction enzyme (Golden Gate) cloning technologies have been used for CRISPR-based genome editing in plants. View an article and references [here](#).


PCR cloning

PCR cloning is a method in which double-stranded DNA fragments amplified by PCR are ligated into a vector. With PCR amplification, this cloning technique requires much less starting material for the insert sequence and allows introduction of new restriction and/or recombination sites to the 5' end of the inserts.

TOPO cloning

TOPO PCR cloning technology was developed to help improve cloning efficiency, simplify protocol setup, and accommodate a wide range of PCR insert sizes. TOPO cloning vectors are linearized by the activity of topoisomerase I (which also has a ligase function) that is covalently bound to the 3' phosphate on each end (see figure below). This system enables the vectors to be joined to PCR inserts with compatible ends (with up to 95% efficiency), without the need for additional ligation steps, in 5 minutes.

Find out more at thermofisher.com/topo.



Quickly find your TOPO cloning kit with our interactive selection tool. Search by application, vector, or desired competent cells at thermofisher.com/topoguide.




Did you know?

The Invitrogen™ TOPO™ XL-2 Complete PCR Cloning Kit provides all the necessary elements for highly efficient cloning of extra-long PCR products from 1 to 13 kb. thermofisher.com/topoxl2

Gateway cloning

To shuttle a DNA fragment of interest (insert) among vectors, the Gateway cloning system offers site-specific, recombinase-based cloning. It maintains the insert's proper orientation and reading frame during shuttling using the Gateway vectors. Once a gene is cloned into an entry clone, you can then move the DNA fragment into one or more destination vectors simultaneously.

Find out more at thermofisher.com/gatewaycloning.



Gateway cloning system reactions. The scheme shows the four types of plasmids and enzyme mixes involved in Gateway cloning reactions. Red arrow represents the fragment of interest. Adapted from Katzen F (2007) *Expert Opin Drug Discov* 2(4):571–589.

GeneArt Gibson Assembly and seamless cloning kits


Several methods are available to assemble multiple PCR or synthesized DNA fragments. These include seamless cloning, Type IIS restriction enzyme (Golden Gate) cloning, and Gibson Assembly cloning. Using our kits, cloning can typically take place within half a day.

Invitrogen™ GeneArt™ Gibson Assembly™ kits allow for the simultaneous assembly of up to 15 very large DNA fragments (up to a total of 100 kb) to create precise and seamless constructs with no additional sequences, in highly efficient reactions. This cloning method circumvents the need for multiple rounds of restriction enzyme analysis and digestion, DNA end repair, dephosphorylation, ligation, enzyme inactivation, and cleanup, and is a powerful tool in synthetic biology.

GeneArt Gibson Assembly kits offer these benefits:

- Assembly of up to 15 fragments to build seamless clones
- Cloning efficiencies up to >95%
- Choice of complete kits with competent cells or master mixes

Invitrogen™ GeneArt™ Seamless Cloning and Assembly Kits allow the assembly of up to 4 DNA fragments. The Invitrogen™ GeneArt™ Seamless PLUS Cloning and Assembly Kit comes with two linearized vectors for expression in either *E. coli* or *Saccharomyces cerevisiae*.



Did you know?

The Gibson Assembly method has been referenced in thousands of peer-reviewed publications and is a powerful method that can be used to seamlessly construct synthetic and natural genes, genetic pathways, and entire genomes [1].

Comparison of seamless cloning kits

Invitrogen product	GeneArt™ Seamless Cloning and Assembly Kit	GeneArt™ Seamless PLUS Cloning and Assembly Kit	GeneArt™ Type IIS assembly kits	GeneArt™ Gibson Assembly™ HiFi and EX Kits
Number of fragments for simultaneous assembly	Up to 4	Up to 4	Up to 8	Up to 6 (HiFi) and 15 (EX)
Reaction time	30 min	15–60 min	5–60 min	15–60 min
Homologous overlaps	15 nt	15–80 nt	4–6 nt	>20 nt
Ligase in the master mix	No	No	Yes	Yes
Vector included	Yes (pUC19L)	Yes (pYES7L and pUC19L)	Yes (pType IIS recipient vector)	No
Mutagenesis protocol	Yes	Yes	No	Yes

1. Gibson DG et al. (2009) Enzymatic assembly of DNA molecules up to several hundred kilobases. *Nat Methods* 6(5):343–345.

Discover more at thermofisher.com/seamless



Cloning with synthetic DNA

Save time by ordering your gene cloned into a vector of interest. Your gene can also be delivered as a DNA fragment similar to a PCR product. The only starting material required is the sequence.

GeneArt gene synthesis and cloned genes


Using the Invitrogen™ GeneArt™ Custom Gene Synthesis service is a reliable and cost-effective way to obtain customized DNA constructs with 100% sequence accuracy. The GeneArt gene synthesis service offers a variety of options, including accelerated production, cloning into your vector of choice, and plasmid preparation.

- An industry leader for reliable delivery and performance
- Proprietary, peer-reviewed, and widely used Invitrogen™ GeneArt™ GeneOptimizer™ algorithm optimizes genes for maximum protein expression, at no additional cost
- High-quality production with ISO 9001:2015 certification
- Delivery in tube format with a minimum of 5 µg DNA

GeneArt Strings DNA fragments

Invitrogen™ GeneArt™ Strings™ DNA fragments are an economical and time-saving alternative to PCR. They are made by the same high-quality process used for gene synthesis, and they can be prepared quickly for cloning.

- Up to 3 kb available with fast production
- Can be optimized with the GeneArt GeneOptimizer algorithm
- Delivered dried, ready for resuspension and cloning with any method



Find out more at thermofisher.com/genesyntesis

Transformation (cont.)

Medium- and high-throughput transformation

Performing bacterial transformations one by one can be very time-consuming and create a bottleneck in your experimental workflow. There are times when medium- and high-throughput transformation options are desired. Invitrogen™ MultiShot™ chemically competent cells provide three flexible product formats to meet your throughput needs.

Find out more at thermofisher.com/multishot



StripWell format

- Medium-throughput option
- Twelve 8-tube strips
- Suitable for 1–96 transformations
- Five *E. coli* strains available

FlexPlate format

- High-throughput option
- 96-well plate separates into 12 x 8-well segments
- Manual and automated platform transformations
- Six *E. coli* strains available



96-well plate

- Highest-throughput option
- Five 96-well plates
- Available with the TOP10 strain
- Stable replication of high copy number plasmids



Bacterial growth media

Sterile, quality growth media and plates are essential for selection and growth of transformed cells. Premade solutions offer convenience and consistency for your cloning workflows. Featured products include:



One Shot LB Agar Plates

- **Invitrogen™ One Shot™ LB Agar Plates with or without Antibiotics**—ready-to-use pre-poured plates with a longer shelf life of 3.5 months
- **Gibco™ LB Broth (1X)** and **Gibco™ Terrific Broth**—popular growth media, available as convenient ready-to-use liquids
- **Invitrogen™ MagicMedia™ *E. coli* Expression Medium**—autoinduction medium for protein expression in *E. coli*

Other available media include Invitrogen™ S.O.C., 2-YT Broth, M9 Minimal Salts, and Gibco™ Bacto™ CD Supreme Fermentation Production Medium (FPM).

Visit thermofisher.com/growthmedia to learn more.



Did you know?

Invitrogen competent cells can be provided in custom configurations per your request. Large and custom volumes as well as multiple formats are at your fingertips. Simply email us at customorders@thermofisher.com.




Did you know?

One Shot LB Agar Plates are available with a 3.5-month shelf life. Learn more at thermofisher.com/growthmedia.



Isothermal amplification





Isothermal amplification is a technique that utilizes enzymes, typically strand-displacing polymerases, to amplify nucleic acid sequences at a constant temperature. This results in continuous and exponential amplification that is not limited by the constraint of thermal cycling.

Find technical and educational resources about isothermal amplification at [thermofisher.com/isothermal](https://www.thermofisher.com/isothermal)



Isothermal amplification overview

Isothermal nucleic acid amplification techniques (INAATs) are fast alternatives to PCR that enable exponential amplification of nucleic acids at constant temperatures. PCR requires repeated denaturation steps at 95°C to separate the DNA strands for the primers to bind. In contrast, INAAT uses DNA polymerases with strand-displacing properties. This eliminates the need for temperature cycling, resulting in reaction times that are as short as 15 min. There are many known INAATs, including loop-mediated amplification (LAMP), multiple displacement amplification (MDA), whole genome amplification (WGA), recombinase polymerase amplification (RPA), and more. Different methods use specific enzymes and reaction conditions.

Once assembled, isothermal reactions can be run on a heat block, thermal cycler, or real-time PCR instrument. Using a real-time PCR instrument also enables fluorescent detection of reaction products. Alternatively, positive results can be detected visually using DNA-binding dyes like Invitrogen™ SYBR™ Green Nucleic Acid Gel Stain. Many instruments, like microplate readers, can be also used to measure the dye-labeled reaction products. Reactions can be further analyzed by agarose gel electrophoresis.


Comparison of seamless cloning kits

Technology	Reaction temperature	Reaction time	Detection method
Loop-mediated isothermal amplification (LAMP)	60–65°C	15–60 min	Fluorescence, colorimetry, turbidity, lateral flow
Multiple displacement amplification (MDA)	30–40°C	60–180 min	Fluorescence, colorimetry
Whole genome amplification (WGA)	30–40°C	60–180 min	Fluorescence, colorimetry
Recombinase polymerase amplification (RPA)	37°C	30–60 min	Fluorescence, lateral flow
Rolling circle amplification (RCA)	30–65°C	60–90 min	Fluorescence, colorimetry, turbidity

Loop-mediated isothermal amplification (LAMP)

LAMP overview

One commonly used isothermal amplification technique is loop-mediated isothermal amplification (LAMP), which utilizes a set of four to six primers and a strand-displacing polymerase, such as Invitrogen™ Lyo-ready Bst DNA Polymerase, to amplify target DNA at a constant temperature. The LAMP reaction produces a large amount of DNA amplicons, with a characteristic ladder-like pattern that can be visualized by gel electrophoresis, or that can be detected by turbidity, fluorescence, or colorimetry.



The amplification process begins with the invasion of an inner primer into the target nucleic acid sequence, followed by extension via a strand-displacing DNA polymerase. As the extension proceeds, the first product is displaced, and an outer primer anneals to the newly synthesized strand, forming a self-hybridizing loop structure. This structure contains multiple sites for amplification initiation and serves as a seed for exponential LAMP reactions.

Find out more at thermofisher.com/lamp.

SuperScript IV RT-LAMP Master Mix

Invitrogen™ SuperScript™ IV RT-LAMP Master Mix is a reverse transcription loop-mediated isothermal amplification (RT-LAMP)-based solution for fast and simple detection of various pathogens, including influenza virus, measles virus, *S. enterica*, *S. aureus*, SARS-CoV-2, and others. Our master mix reagents provide maximum flexibility to help optimize and accelerate your pathogen research and surveillance.

Product highlights

- **Fast**—RNA and DNA target detection in as little as 5 minutes with evolved Bst DNA polymerase
- **Efficient**—one-step reaction for reverse transcription of RNA to cDNA with Invitrogen™ SuperScript™ IV Reverse Transcriptase
- **Sensitive**—greater sensitivity and specificity utilizing Invitrogen™ RNaseOUT™ Recombinant Ribonuclease Inhibitor and an optimized buffer
- **Simple**—streamlined workflow: single-tube format, only requires a 65°C heating block
- **Flexible**—several options for evaluating results, including real-time and endpoint detection methods



Find out more about SuperScript IV RT-LAMP Master Mix [here](https://thermofisher.com/lamp).

Lyo-ready Bst DNA Polymerase

Lyo-ready Bst DNA Polymerase is an engineered version of the Bst DNA polymerase large fragment, which shows a significantly faster reaction speed, increased sensitivity, and tolerance to inhibitors.

Lyo-ready Bst DNA Polymerase provides maximum flexibility to optimize your LAMP reaction and works with various types of pathogens, including human adenovirus, measles virus, SARS-CoV-2, and other pathogens.

Product highlights

- **Fast**—amplifies targets in as little as 10 minutes
- **Sensitive**—provides sensitivity down to 50 copies
- **Robust**—amplifies even from inhibitor-containing RNA/DNA samples
- **Flexible**—provides the ability to optimize your LAMP or RT-LAMP reaction



Find out more about Lyo-ready Bst DNA Polymerase [here](#)

Tips:


1. Did you know? The term “lyo-ready” refers to an enzyme that is provided in a liquid formulation without glycerol, making it compatible with microfluidics-based systems and various downstream applications such as lyophilization. Furthermore, it maintains the necessary stability and activity levels for direct enzymatic reactions.
2. To minimize nonspecific amplification in LAMP, follow these steps:
 - a. **Prevent cross-contamination:** Use uncontaminated reagents and maintain a clean work environment.
 - b. **Enhance primer design:** Optimize primer sequences for improved specificity.
 - c. **Optimize reaction conditions:** Adjust Lyo-Ready Bst DNA Polymerase amount and reaction time.

Rolling circle amplification (RCA) and whole genome amplification (WGA)

WGA is a popular technique used to amplify small amounts of DNA, as little as that of a single cell, to obtain large quantities of product for downstream applications like next-generation sequencing. This method is based on multiple displacement amplification (MDA). When the starting material is circular, RCA is utilized. Both MDA-WGA and RCA use a strand-displacing DNA polymerase. Thermo Fisher Scientific offers two types: the wild-type Thermo Scientific™ phi29 DNA Polymerase, and the engineered Thermo Scientific™ EquiPhi29™ DNA Polymerase, which has exceptional performance compared to phi29 DNA Polymerase.

Multiple displacement amplification–whole genome amplification (MDA-WGA)


The process of MDA-WGA starts with the addition of random hexamer primers to the DNA sample, which hybridize to the template and initiate DNA synthesis. A strand-displacing polymerase then synthesizes DNA strands along the template from the multiple priming positions.



Learn more about rolling circle amplification at thermofisher.com/rca

Rolling circle amplification (RCA)

The process of RCA uses a circular DNA template, such as a plasmid or circularized oligonucleotide. DNA polymerase continuously synthesizes new DNA strands as it moves around the template. This leads to exponential amplification and the generation of concatemers containing numerous tandem repeats that are complementary to the circular template.



EquiPhi29 DNA Polymerase

EquiPhi29 DNA Polymerase is a proprietary phi29 DNA Polymerase mutant developed through *in vitro* protein evolution. This enzyme is significantly improved over phi29 DNA Polymerase in protein thermostability, reaction speed, product yield, and amplification bias. At the same time, it retains all the benefits of the wild-type enzyme, including high processivity (>70 kb), strong strand-displacing activity, and 3' to 5' exonuclease (proofreading) activity. EquiPhi29 DNA Polymerase is also available as a convenient kit with all required components for RCA and MDA-WGA.

Product highlights

- **Fast**—amplifies target in less than 2 hours
- **Sensitive**—provides sensitivity down to 1 fg of DNA
- **High yield**—yields up to 17 µg amplified DNA
- **Variety of applications**—compatible with cell-free DNA enrichment, and cell-free protein expression
- **Availability**—available in lyo-ready format



Did you know?

These products are free of OPE and NPE, making them safer for aquatic life.




Find out more about EquiPhi29 DNA Polymerase at thermofisher.com/equiph29.



Recombinase polymerase amplification (RPA)

RPA is a sensitive, low-temperature (39–42°C) isothermal DNA amplification technique that is initiated with two primers and carried out using a mixture of recombinase enzymes, a single-stranded DNA-binding protein (SSB), and a strand-displacing DNA polymerase such as Bst DNA polymerase. To detect RNA targets, a reverse transcriptase can be incorporated into the RPA reaction. Amplified DNA can be detected using various detection methods such as endpoint detection, real-time analysis, lateral flow assays, and CRISPR-Cas9 systems.

As shown below, RPA uses a recombinase complex made of UvsX and UvsY proteins, which binds to the primers in the presence of ATP. This complex seeks a sequence on the DNA homologous to the primer, and promotes strand opening and invasion by the primer. The SSB stabilizes the resulting D-loop structure, allowing the DNA polymerase to extend the primer while displacing the DNA strand.



Lyo-ready RPA Kit

The Invitrogen™ Lyo-ready RPA Kit includes all required components in a glycerol-free format for isothermal RPA, enabling high reaction sensitivity, specificity, and tolerance to inhibitors. The Lyo-ready RPA Kit contains separate reagents that provide maximum flexibility to optimize the RPA reaction and work with various types of pathogens (RNA or DNA), including *S. aureus*, *P. aeruginosa*, influenza and measles viruses, and others. Components for RPA and RT-RPA reactions are also available in a stand-alone format.

Product highlights

- **Sensitive**—provides sensitivity down to as low as one copy of target
- **Specific**—detects and amplifies a single target in a mixture of RNA and DNA
- **Fast**—amplifies target in as little as 20 minutes
- **Robust**—amplifies even from inhibitor-containing RNA and DNA samples




Learn more about RPA at thermofisher.com/rpa.

Find out more about the Lyo-ready RPA Kit [here](#).



Resources




We have compiled educational resources, mobile apps, frequently asked questions, and other information to help you achieve success in your research.



Educational resources

Suitable for new and experienced molecular biologists alike, our free online educational resources are designed to help you review the basics, build your expertise, or discover our latest innovative technologies. Explore our educational resources in the following areas of molecular biology.

Reverse transcription




thermofisher.com/rteducation

PCR enzymes, plastics, and thermal cyclers




thermofisher.com/pcreducation

Nucleic acid electrophoresis



thermofisher.com/na-electrophoresis-education

Restriction enzymes, molecular cloning, and competent cells



thermofisher.com/cloningeducation



Resources

Webinars: Watch live and recorded webinars for in-depth understanding of molecular and synthetic biology techniques and tools to help elevate your research.

thermofisher.com/mbwebinars

Videos: Experience entertaining and visual learning with our educational videos on molecular biology techniques, how-to guides, tips and tricks, and more.

thermofisher.com/mbvideos

Application notes: Read white papers and application notes from our R&D scientists on our product innovations.

thermofisher.com/mbliterature

Online tools: Use our interactive online tools for PCR annealing temperature, restriction enzyme information, product selection, and more.

thermofisher.com/mbtools

Find out more at thermofisher.com/mbschool



Table of contents

Podcast: Speaking of Mol Bio



Listen now to Speaking of Mol Bio. This podcast highlights trending applications in science and the molecular biological aspects of those applications. Our hosts delve into deep discussions with CEOs, R&D scientists, researchers, and key opinion leaders across the globe.

This podcast helps scientifically curious people—from all scientific and nonscientific backgrounds—understand how modern molecular biology applications can help push the boundaries in medicine, science, drug discovery, and in the cure and treatment of diseases.

Listen now at thermofisher.com/molbiopodcast.

Mobile apps



DailyCalcs—science calculator

The DailyCalcs app turns your phone into a science calculator to help simplify everyday tasks in the lab. The app features eight calculators: molarity, dilution, formula weight, transfection, unit conversions, culture vessel data, media conversions, and specific productivity. Find it in your preferred mobile application store.



Instrument Connect—remote monitoring

Instrument Connect allows you to view instrument status, monitor or schedule a run, and more on any cloud-enabled instrument, including the Applied Biosystems™ ProFlex™, SimpliAmp™, and MiniAmp™ PCR instruments. Access at thermofisher.com/connect.



PCR Quest—match-3 lab game


Test your PCR knowledge with our lab game—PCR Quest—where you travel from lab to lab crushing the world's toughest diseases. Download at thermofisher.com/pcrquest.

Custom Commercial Supply

As a leading supplier of molecular biology products, we offer tailored solutions for companies developing new molecular assays. Whether you're just starting or need a specific solution, we have what you need. Work with an experienced supplier that knows both raw materials and new technologies. Our dedicated business team is here to provide value beyond just our products.

What do our OEM solutions mean to you?

- Customization of products and services
- Consultation, partnership, and expertise
- Negotiated business terms
- Warranties and indemnification
- Commercial-use rights and obligations
- Risk and liability management



Find out more at thermofisher.com/mdx

Choosing greener products?

Look for the leaf.

For customers seeking greener products, Thermo Fisher has created an easy identification system with our green leaf symbol. When you see the green leaf symbol on one of our products, it means that product meets one or more of our environmental sustainability criteria.



Greener by design™



Frequently asked questions

Below are some common questions and answers to help you start or troubleshoot molecular biology experiments.

Sample preparation

Which kit should I use to isolate nucleic acids from my sample?

Choosing the right product is fundamental to ensuring proper lysis of cells and tissue, as well as sufficient yield and quality of isolated nucleic acids. Look to our selection guides ([see pages 11–14](#)) to help you decide according to nucleic acid type, sample source, experimental throughput, and format as well as downstream applications.

What are the key steps to preventing RNA degradation?

The basic lab precautions listed below can help minimize RNA degradation and avoid experimental inconsistency and failure.

- Use nuclease-free pipette tips and tubes
- Use nuclease-free water and reagents
- Regularly decontaminate work surfaces
- Properly stabilize RNA sources before storage

For more tips and troubleshooting advice on sample prep, visit [thermofisher.com/rnabasics](https://www.thermofisher.com/rnabasics) and [thermofisher.com/napsupport](https://www.thermofisher.com/napsupport).

Reverse transcription

How do I improve the efficiency of cDNA synthesis when working with challenging samples (e.g., low-abundance, degraded, inhibitor-containing, or GC-rich RNA)?

When working with challenging RNA samples, select a reverse transcriptase that is highly sensitive, processive, thermostable, and resistant to common inhibitors, to help you obtain the highest cDNA yield ([see pages 17–19](#)).

What are the benefits of using random primers, oligo(dT) primers, gene-specific primers, or oligo(dT)/random mixed primers in reverse transcription?

- Random primers are good to use with degraded RNA, RNA with high secondary structure, nonpolyadenylated RNA, or prokaryotic RNA.
- Oligo(dT) primers are an optimal choice for synthesis of full-length cDNA from eukaryotic mRNA. Applications include cDNA cloning, cDNA library construction, and 3' rapid amplification of cDNA ends (3' RACE).
- Gene-specific primers are designed based on known sequences of the target RNA. These primers offer the most specific priming and are commonly used in one-step RT-PCR.
- A mixture of oligo(dT) and random primers is often used in two-step RT-PCR to achieve the benefits of each primer type ([see page 22](#)).

For more tips and troubleshooting advice on reverse transcription, visit [thermofisher.com/rteducation](https://www.thermofisher.com/rteducation).

PCR amplification

How can I optimize primer annealing for PCR?

Traditionally, gradient thermal cyclers have been used to simultaneously assess multiple temperatures around the theoretical annealing point. Compared to gradient thermal cyclers, instruments with VeriFlex technology allow more precise temperature control for faster optimization of primer annealing ([see pages 26–27](#)).

Tedious optimization steps may be circumvented using the novel Invitrogen™ Platinum™ DNA polymerases. Their innovative buffers enable specific annealing at 60°C for most primers when they are designed following general primer design rules ([see pages 30–33](#)).



Frequently asked questions (cont.)

What do I need to run fast PCR?

PCR amplicons shorter than 1 kb can be amplified in as little as 40 minutes using “fast” enzymes (high processivity; **see page 32**), “fast” plastics (low profile and ultra-thin walls; **see page 29**), and “fast” thermal cyclers (fast ramp rate; **see pages 26–27**).

How can I prevent sample evaporation during PCR?

Proper sealing of your reactions will help prevent evaporation during PCR.

- When using adhesive film to seal a plate, be sure to properly align the seal to cover all wells and press firmly along all edges of the plate using an applicator tool.
- When sealing a plate using cap strips, ensure that the cap strips are compatible with the plate and thermal cycler being used. Be sure to align cap strips with each well of the plate and place firmly across the plate for a secure fit.
- Use the applicator tool (Cat. No. 4333183 or 4330015) or other comparable sealing tools as needed.

For more tips and troubleshooting advice on PCR, visit [thermofisher.com/pcreducation](https://www.thermofisher.com/pcreducation) and [thermofisher.com/pcrsupport](https://www.thermofisher.com/pcrsupport).

Nucleic acid electrophoresis

Why is it important to choose the right ladder when using E-Gel precast agarose gels?

Accurate analysis of electrophoresis bands often depends on the DNA ladder chosen for your gel run. E-Gel DNA ladders are formulated with ready-to-use buffers unique for E-Gel precast agarose gels, and DNA standards designed for optimal separation (**see page 42**).

Are there safer alternatives to ethidium bromide for staining nucleic acids in gel electrophoresis?

SYBR Safe DNA Gel Stain is a safer alternative to ethidium bromide and is commonly used in gel electrophoresis. SYBR Safe DNA stain is not classified as hazardous waste or as a pollutant under US federal regulations (**see page 42**).

Do I need a buffer to run the E-Gel system?

No. The E-Gel electrophoresis system does not require electrophoresis buffers like TBE or TAE. E-Gel cassettes already contain everything you will need and are classified as dry electrophoresis.

For more tips and troubleshooting advice on nucleic acid electrophoresis, visit [thermofisher.com/na-electrophoresis-education](https://www.thermofisher.com/na-electrophoresis-education) and [thermofisher.com/na-electrophoresis-support](https://www.thermofisher.com/na-electrophoresis-support).

Cloning

Do you have a buffer compatibility chart for restriction enzymes?

All FastDigest restriction enzymes are 100% active in one universal FastDigest buffer (**see page 46**). Hence, there is no buffer compatibility chart for FastDigest restriction enzymes.

What is the main difference between GeneArt Strings DNA Fragments and GeneArt Gene Synthesis?

GeneArt Strings DNA Fragments are custom-made, linear DNA fragments that are uncloned and double-stranded. GeneArt Gene Synthesis is a service offered for chemical synthesis, cloning, and sequence verification of genetic sequences (**see page 50**).

What are some key considerations for choosing competent cells for my cloning applications?

Genotype, transformation efficiency, growth rate, and throughput format are important factors in choosing competent cells for cloning. The genotype of a cell strain may determine growth conditions and suitability for transformation with specific DNA types (**see page 51**).

For more tips and troubleshooting advice on cloning, visit [thermofisher.com/cloningeducation](https://www.thermofisher.com/cloningeducation) and [thermofisher.com/cloningsupport](https://www.thermofisher.com/cloningsupport).



Ordering information

	Quantity	Cat. No.
Nucleic acid isolation		
PureLink Quick Plasmid Miniprep Kit	50 preps	K210010
PureLink HiPure Plasmid Filter Midiprep Kit	25 preps	K210014
PureLink HiPure Plasmid Filter Maxiprep Kit	10 preps	K210006
PureLink <i>Pro</i> Quick96 Plasmid Purification Kit	4 x 96 preps	K211004A
PureLink Quick Gel Extraction Kit	50 preps	K210012
TRIzol Plus RNA Purification Kit	50 preps	12183555
PureLink RNA Mini Kit	10 preps	12183020
PureLink Genomic DNA Mini Kit	10 preps	K182000
PureLink <i>Pro</i> 96 Genomic DNA Purification Kit	4 x 96 preps	K182104A
PureLink <i>Pro</i> 96 Viral RNA/DNA Purification Kit	4 plates	12280096A
PureLink Viral RNA/DNA Mini Kit	50 preps	12280050
PureLink Genomic Plant DNA Purification Kit	50 preps	K183001
MagMAX Microbiome Ultra Nucleic Acid Isolation Kit, with bead tubes	100 preps	A42358
MagMAX Microbiome Ultra Nucleic Acid Isolation Kit, with bead plate	100 preps	A42357
MagMAX Cell-Free DNA Isolation Kit	1 kit	A29319
MagMAX Cell-Free Total Nucleic Acid Isolation	1 kit	A36716
MagMAX FFPE DNA/RNA Ultra Kit	1 kit	A31881
MagMAX <i>mirVana</i> Total RNA Isolation Kit	96 reactions	A27828
MagMAX Sequential DNA/RNA Kit	1 kit	A65309
MagMAX Prime Viral/Pathogen NA Isolation Kit	Up to 600 reactions	A58145
KingFisher PlasmidPro Maxi Processor Endotoxin-Free Cartridge	4 cartridges	A54072
PureLink PCR Purification Kit	50 preps	K310001
PureLink Quick Gel Extraction and PCR Purification Combo Kit	50 preps	K220001
PureLink Quick Gel Extraction Kit	50 preps	K210012
	5 mL	A58521
MagMAX Pure Bind Beads	50 mL	A58522
	250 mL	A58523
	40 reactions	A35374
TaqMan Fast Advanced Cells-to-C ₁ Kit	100 reactions	A35377
Dynabeads M-270 Streptavidin	2 mL	65305
Dynabeads MyOne Streptavidin C1	2 mL	65001

	Quantity	Cat. No.
KingFisher instruments – For Laboratory Use.		
KingFisher Apex Purification System with 96 PCR Head	1 system	5400910
KingFisher Apex Purification System with 96 Combi Head	1 system	5400920
KingFisher Apex Purification System with 96 Deep-Well Head	1 system	5400930
KingFisher Apex Purification System with 24 Combi Head	1 system	5400940
KingFisher Flex Purification System with 96 PCR Head	1 system	5400610
KingFisher Flex Purification System with 96 KingFisher Plate	1 system	5400620
KingFisher Flex Purification System with 24 Deep-Well Head	1 system	5400640
KingFisher Flex Purification System with 96 Deep-Well Head	1 system	5400630
KingFisher Presto Purification System with 24 Deep-Well Head	1 system	5400840
KingFisher Presto Purification System with 96 Deep-Well Head	1 system	5400830

KingFisher instruments – For Research Use Only. Not for use in diagnostic procedures.		
KingFisher Duo Prime Purification System	1 system	5400110
KingFisher PlasmidPro Maxi Processor	1 instrument	A66427







Ordering information (cont.)




	Quantity	Cat. No.
Reverse transcription		
SuperScript IV Reverse Transcriptase	2,000 units	18090010
	10,000 units	18090050
SuperScript IV First-Strand Synthesis System	50 reactions	18091050
	200 reactions	18091200
SuperScript IV VILO Master Mix	50 reactions	11756050
	500 reactions	11756500
SuperScript IV VILO Master Mix with ezDNase Enzyme	50 reactions	11766050
	500 reactions	11766500
SuperScript IV One-Step RT-PCR System	25 reactions	12594025
	100 reactions	12594100
SuperScript IV UniPrime One-Step RT-PCR System	100 reactions	12597100
SuperScript IV CellsDirect cDNA Synthesis Kit	50 reactions	11750150
	500 reactions	11750350
SuperScript IV CellsDirect Lysis Reagents	500 reactions	11750550
	48 reactions	11752048
SuperScript IV Single Cell/Low Input cDNA PreAmp Kit	96 reactions	11752096
	192 reactions	11752192
SuperScript IV Template Switching RT Master Mix	48 reactions	A65423
	96 reactions	A65424
	5 x 96 Reactions	A65425
SUPERase-In RNase Inhibitor (20 U/μL)	2,500 units	AM2694
	10,000 units	AM2696
RNaseOUT Recombinant Ribonuclease Inhibitor	5,000 units	10777019
RNase Inhibitor	2,000 units	N8080119
Ambion RNase Inhibitor, cloned, 40 U/μL	2,500 units	AM2682
	10,000 units	AM2684
Ribonuclease H	30 units	18021014
Random Hexamers (50 μM)	5 nmol	N8080127

	Quantity	Cat. No.
Random Primers	100 μL	48190011
Oligo(dT) ₁₂₋₁₈ Primer	50 μL	18418012
Oligo(dT) ₂₀ Primer	50 μL	18418020
DNase I, Amplification Grade	100 units	18068015
Isothermal amplification		
Lyo-ready Bst DNA Polymerase	1,200 units (6 U/μL)	A56655
	6,000 units (6 U/μL)	A56656
	1,200 units (40 U/μL)	A56657
SuperScript IV RT-LAMP Master Mix	100 reactions	A51801
	400 reactions	A51802
	1,000 reactions	A51803
EquiPhi29 DNA Polymerase	250 units	A39390
	1,000 units	A39391
	5,000 units	A39392
EquiPhi29 DNA Amplification Kit	100 reactions	A65393
	500 reactions	A65394
phi29 DNA Polymerase (10 U/μL)	250 units	EP0091
	1,000 units	EP0092
	5,000 units	EP0094
Lyo-ready RPA Kit	100 reactions	A72127
	500 reactions	A72128
Lyo-ready T4 UvsX Protein	120 μg	A72124
Lyo-ready T4 UvsY Protein	120 μg	A72125
Lyo-ready T4 Gene 32 Protein	1,600 μg	A72123
L.b. Cas12a Nuclease	200 pmol	A40001796
	2,000 pmol	A40001795



Ordering information (cont.)




	Quantity	Cat. No.
PCR		
Thermal cyclers		
ProFlex PCR System, 3 x 32-well	1 instrument	4484073
ProFlex PCR System, 96-well	1 instrument	4484075
VeritiPro Thermal Cycler, 96-well	1 instrument	A48141
MiniAmp Plus Thermal Cycler	1 instrument	A37835
MiniAmp Thermal Cycler	1 instrument	A37834
Automated Thermal Cycler (ATC), 96-well, laptop	1 instrument	A31486
Plastics		
MicroAmp EnduraPlate Optical 96-Well Fast Multicolor Reaction Plates with Barcode	5 plates	4483493
MicroAmp Optical Adhesive Film	100 covers	4311971
MicroAmp Optical 96-Well Reaction Plate	10 plates	N8010560
MicroAmp Optical 8-Cap Strips	300 strips	4323032
MicroAmp Fast Optical 96-Well Reaction Plate, 0.1 mL	10 plates	4346907
MicroAmp Fast Reaction Tube with Cap, 0.1 mL	1,000 tubes	4358297
MicroAmp EnduraPlate Optical 384-Well Multicolor Reaction Plates with Barcode	5 plates	4483316
MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plates with Barcode	20 plates	4483354
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	20 plates	A32811
MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL	125 strips	A30589
MicroAmp EnduraPlate Optical 96-Well Full-Skirted Plates with Barcode, clear	50 plates	A31728
 PCR Plate, 384-well, standard, Sustain Series	50 plates	AB1384SS
 PCR Plate, 96-well, non-skirted, Sustain Series	25 plates	AB0600SS
 PCR Plate, 96-well, segmented, semi-skirted, Sustain Series	25 plates	AB0900SS
 Tubes and Flat Caps, strips of 8, Sustain Series	250 strips	AB1182SS

	Quantity	Cat. No.
PCR enzymes		
Platinum II <i>Taq</i> Hot-Start DNA Polymerase	100 reactions	14966001
	500 reactions	14966005
Platinum II Hot-Start PCR Master Mix (2X)	50 reactions	14000012
	200 reactions	14000013
Platinum II Hot-Start Green PCR Master Mix (2X)	50 reactions	14001012
	200 reactions	14001013
 AmpliTaq Gold 360 DNA Polymerase	100 units	4398813
	250 units	4398823
 AmpliTaq Gold 360 Master Mix	1 mL	4398876
	5 mL	4398881
Platinum SuperFi II DNA Polymerase	100 units	12361010
	500 units	12361050
Platinum SuperFi II PCR Master Mix	100 reactions	12368010
	500 reactions	12368050
Platinum SuperFi II Green PCR Master Mix	100 reactions	12369010
	500 reactions	12369050
 Platinum Direct PCR Universal Master Mix	100 reactions	A44647100
	500 reactions	A44647500
dNTP Set (100 mM)	4 x 250 µL	10297018
	8 x 1.25 mL	10297117

Order custom DNA oligos at thermofisher.com/orderoligos



Ordering information (cont.)

	Quantity	Cat. No.
Nucleic acid separation and analysis		
 SYBR Safe DNA Gel Stain	400 µL	S33102
SYBR Gold Nucleic Acid Gel Stain	500 µL	S11494
UltraPure DNase/RNase-Free Distilled Water	500 mL	10977015
 UltraPure Agarose	100 g	16500100
 TrackIt 100 bp DNA Ladder	100 applications	10488058
UltraPure TAE Buffer, 10X	4 L	15558026
E-Gel Agarose Gels with SYBR Safe DNA Gel Stain	10 gels	A42135
E-Gel Double Comb Agarose Gels with SYBR Safe DNA Gel Stain	10 gels	A42348
E-Gel EX Double Comb Agarose	10 gels	A44889
E-Gel CloneWell II Agarose Gels with SYBR Safe DNA Gel Stain, 0.8%	10 gels	G661818
 E-Gel Agarose Gels with SYBR Safe DNA Gel Stain, 2%	2 x 10 gels	A45204
E-Gel EX Agarose Gels, 2%, with SYBR Gold DNA Stain	20 gels	G402022
 E-Gel 1 Kb Plus DNA Ladder	100 applications	10488090
 E-Gel Sample Loading Buffer, 1X	4 x 1.25 mL	10482055
E-Gel Power Snap Plus Electrophoresis System	1 system	G9301* G9311**
E-Gel Power Snap Plus Electrophoresis System Starter Kit, 48-well, 1%	1 kit	G9341* G9331**
E-Gel Power Snap Plus Electrophoresis System Starter Kit, 48-well, 2%	1 kit	G9342* G9332**
E-Gel Power Snap Plus Electrophoresis System Starter Kit, 96-well, 1%	1 kit	G9391* G9381**
E-Gel Power Snap Plus Electrophoresis System Starter Kit, 96-well, 2%	1 kit	G9392* G9382**
E-Gel Power Snap Electrophoresis System Starter Kit, EX 2%	1 kit	G8342ST
 E-Gel 48 Agarose Gels with SYBR Safe DNA Gel Stain, 2%	8 gels	G820802
 E-Gel 96 Agarose Gels with SYBR Safe DNA Gel Stain, 2%	8 gels	G720802
 E-Gel 48 Agarose Gels with SYBR Safe DNA Gel Stain, 4%	8 gels	G820804

	Quantity	Cat. No.
Cloning and gene synthesis		
FastDigest BamHI	800 reactions	FD0054
	2,500 reactions	FD0055
FastDigest BcuI	20 reactions	FD1253
	50 reactions	FD1254
FastDigest BshTI	20 reactions	FD1464
FastDigest DpnI	50 reactions	FD1703
	100 reactions	FD1704
FastDigest EcoRI	800 reactions	FD0274
	2,500 reactions	FD0275
FastDigest KpnI	300 reactions	FD0524
	20 reactions	FD0593
FastDigest NotI	50 reactions	FD0594
	150 reactions	FD0595
	250 reactions	FD0596
FastDigest Sall	200 reactions	FD0644
FastDigest XbaI	300 reactions	FD0684
	750 reactions	FD0685
FastDigest XhoI	400 reactions	FD0694
	1,200 reactions	FD0695
FastDigest Esp3I (BsmBI) (IIS class)	20 reactions	FD0454
FastDigest BpiI (BbsI) (IIS class)	20 reactions	FD1014
FastDigest Eco31I (BsaI) (IIS class)	50 reactions	FD0293
	100 reactions	FD0294
T4 DNA Ligase (5 U/µL)	200 units	EL0014
	1,000 units	EL0011
FastAP Thermosensitive Alkaline Phosphatase (1 U/µL)	300 units	EF0654
	1,000 units	EF0651

* North America, Europe, Middle East, and Africa.

** Asia Pacific, Japan, Latin America, and greater China.



Ordering information (cont.)

	Quantity	Cat. No.
T4 Polynucleotide Kinase (10 U/μL)	500 units	EK0031
	2,500 units	EK0032
DNA Polymerase I (10 U/μL)	500 units	EP0041
	2,500 units	EP0042
Klenow Fragment (10 U/μL)	300 units	EP0051
	1,500 units	EP0052
T4 DNA Polymerase (5 U/μL)	1 mL of 5X reaction buffer	EP0061
	2 x 1 mL of 5X reaction buffer	EP0062
T7 DNA Polymerase (10 U/μL)	300 units	EP0081
CorrectASE Enzyme	50 reactions	A14972
TOPO TA Cloning Kit for Subcloning, without competent cells	25 reactions	450641
Zero Blunt TOPO PCR Cloning Kit, without competent cells	25 reactions	450245
pENTR/D-TOPO Cloning Kit, with One Shot TOP10 Chemically Competent <i>E. coli</i>	20 reactions	K240020
One Shot TOP10 Chemically Competent <i>E. coli</i>	21 x 50 μL	C404003
One Shot Stbl3 Chemically Competent <i>E. coli</i>	21 x 50 μL/tube	C737303
MAX Efficiency DH5α Competent Cells	5 x 200 μL	18258012
ElectroMAX DH10B Cells	5 x 100 μL	18290015
MAX Efficiency Stbl2 Competent Cells	5 x 200 μL	10268019
MultiShot TOP10 Chemically Competent <i>E. coli</i>	5 plates	C40005
MultiShot StripWell TOP10 Chemically Competent <i>E. coli</i>	1 rack	C409601
MultiShot StripWell BL21 Star (DE3) Chemically Competent <i>E. coli</i>	1 rack	C609601
MultiShot FlexPlate TOP10 Chemically Competent <i>E. coli</i>	1 plate	C4081201
MultiShot FlexPlate DH5α T1R Chemically Competent <i>E. coli</i>	1 plate	C4481201
MultiShot FlexPlate Stbl3 Chemically Competent <i>E. coli</i>	1 plate	C7381201

	Quantity	Cat. No.
GeneArt Gibson Assembly HiFi Master Mix	50 reactions	A46628
GeneArt Gibson Assembly EX Master Mix	50 reactions	A46636
GeneArt Seamless Cloning and Assembly Enzyme Mix	20 reactions	A14606
GeneArt Type IIS Assembly Kit, AarI	10 reactions	A15916
GeneArt Type IIS Assembly Kit, BsaI	10 reactions	A15917
GeneArt Type IIS Assembly Kit, BbsI	10 reactions	A15918
GeneArt High-Order Genetic Assembly System	10 reactions	A13285
Gateway BP Clonase II Enzyme Mix	20 reactions	11789020
Gateway LR Clonase II Enzyme Mix	20 reactions	11791020
MultiSite Gateway Pro Plus	20 reactions	12537100
LR Clonase II Plus Enzyme	20 reactions	12538120
Gateway Vector Conversion System with One Shot <i>ccdB</i> Survival Cells	1 kit	11828029
PCR Cloning System with Gateway Technology with pDONR 221 and OmniMAX 2 Competent Cells	20 reactions	12535029
PCR Cloning System with Gateway Technology with pDONR/Zeo and OmniMAX 2 Competent Cells	20 reactions	12535037
Gateway pDONR 221 Vector	6 μg	12536017
pENTR/D-TOPO Cloning Kit, with One Shot TOP10 Chemically Competent <i>E. coli</i>	20 reactions	K240020
pCR 8/GW/TOPO TA Cloning Kit with One Shot TOP10 <i>E. coli</i>	20 reactions	K250020

Order custom cloned genes or DNA fragments at thermofisher.com/geneartdashboard.



 Find out more at thermofisher.com/molecularmasterpiece

For Research Use Only. Not for use in diagnostic procedures. © 2020–2023, 2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan and AmpliTaq Gold are trademarks of Roche Molecular Systems, Inc., used under permission and license. TRizol and DNAzol are trademarks of Molecular Research Center, Inc. Gibson Assembly is a trademark of Telesis Bio, Inc., used under permission and license.

BROC-9207299 0425