GRAM IODINE

REF	R40056, Gram Iodine.....	250 ml/Btl, EA
REF	R40057, Gram Iodine.....	250 ml/Btl, 5/Pk
REF	R40077, Gram Iodine.....	Gallon

1. INTENDED USE

A reagent/set of reagents used in a technique for the differentiation of Gram-negative and Gram-positive microorganisms.

The device is for professional use only and is not intended for self-testing

2. SUMMARY AND EXPLANATION

The Gram stain method was developed in 1884 by the Danish bacteriologist, Christian Gram to differentiate bacterial cells from infected tissue.¹ Later it was discovered that the bacterial cell wall composition was the key to the Gram stain and would differentiate bacteria into two groups based on cell color after staining. Since that time there have been many modifications of the original technique.^{2,3}

3. PRINCIPLE

The primary stain is crystal violet, a basic dye taken up by all bacteria due to its ability to rapidly permeate the cell wall. It stains the protoplast of bacteria purple. The potassium-iodine mixture is the mordant which complexes with the primary stain in the cell. In gram-positive cells, the crystal violet-iodine complex is trapped in the cell due to a decrease in cell wall permeability caused by alcohol dehydration. Gram-positive cells are stained purple. In gram-negative cells, the complex is removed by the decolorizer due to an increase in permeability caused by solubility of the lipids in alcohol. The counterstain used must be a contrasting color to the primary stain and, in this case, it is safranin which stains gram-negative cells red.

4. REAGENTS

Potassium Iodide (CAS 7681-11-0).....	6.6 g
Iodine (CAS 7553-56-2).....	3.3 g
Demineralized Water (CAS 7732-18-5).....	990.0 ml

5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet for additional information.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

6. STORAGE

Store product in its original container at 20-25°C until used.

7. PRODUCT DETERIORATION

This product should not be used if (1) the color has changed from yellow-red, (2) the expiration date has passed, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{4,5}

9. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Gram Crystal Violet (R40052), (7) Gram Decolorizer (R40054), (8) Gram Safranin (R40058), (9) Glass slides, (10) Bunsen burner or slide warmer, (11) Microscope, immersion oil, (12) Staining rack, (13) Methanol (optional).

10. PROCEDURE

Add Gram Iodine concentrate to diluent. To each 250 ml bottle of diluent, add one 2.5 ml tube of concentrate. To each gallon of diluent, add one 40 ml tube of concentrate. Mix well. Use Gram Iodine within 3 months after reconstitution.

- 10.1 Make a thin smear of the material for study on a glass slide. Allow the slide to air-dry completely or use a slide warmer. Fix the slide by passing 3 or 4 times through the flame of a Bunsen burner or flood the smear with methanol and allow it to evaporate.
- 10.2 Place the slide on a staining rack and overlay with Gram Crystal Violet for 1 minute.
- 10.3 Wash thoroughly with water and overlay with Gram Iodine mordant for 1 minute.

10.4 Flood with Gram Decolorizer until the solvent flows colorless from the slide (10-30 seconds).

10.5 Rinse with water and overlay with Gram Safranin for 30 seconds.

10.6 Rinse with water and allow to air-dry.

11. INTERPRETATION

Gram-positive organisms stain purple.

Gram-negative organisms stain red.

12. QUALITY CONTROL

All lot numbers of Gram Iodine have been tested and found to yield acceptable stain results as listed in the Interpretation section. Positive and negative control slides should be tested prior to use of new lot numbers of stains and decolorizing reagent and at least weekly thereafter. If aberrant quality control results are noted, patient results should not be reported.

13. LIMITATIONS

- 13.1 Test positive and negative controls with each new lot number of stains and weekly thereafter to verify reagent efficacy and ensure the decolorization procedure is performed correctly.⁵
- 13.2 To obtain the most reliable results, remove clinical specimen isolates from an 18-24 hour culture growing on nonselective medium.⁵
- 13.3 Do not overheat smears when fixing slides. Excessive heating will cause atypical staining.⁵
- 13.4 Gram-positive organisms contained in a specimen may appear gram-negative if the patient is on antimicrobial therapy.⁵
- 13.5 Precipitated crystal violet can appear as coccoid shapes or fungal elements, as well as other artifacts or background material.⁵

14. BIBLIOGRAPHY

1. Bartholomew, J.W. and T. Mittwer. 1952. Bacteriol. Rev. 16:1-29.
2. Mittwer, T., J.W. Bartholomew, and B.J. Kallman. 1950. Stain Technol. 25:169-179.
3. Bartholomew J.W. and T. Mittwer. 1951. Stain Technol. 26:231-240.
4. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO
5. Garcia, L.S. 2010. Clinical Microbiology Procedures Handbook. 3rd ed., American Society for Microbiology, Washington, D.C.
6. Mangels, J.I., M.E. Cox, and L.H. Lindberg. 1984. Diagn. Microbiol. Infect. Dis. 2:129-137.

15. SYMBOL LEGEND

REF	Catalogue Number
IVD	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
	Keep away from sunlight
UDI	Unique Device Indicator
EC REP	Authorized representative in the European Community
UK CA	UK Conformity Assessed
CE	European Conformity Assessment
	Manufacturer

12076 Santa Fe Trail Drive
Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
(800) 255-6730 International: (913) 888-0939

For technical information contact your local distributor.
CAS (Chemical Abstracts Service Registry No.)

Manufactured for Remel Inc.

©2022 Thermo Fisher Scientific Inc. All rights reserved.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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Tillverkad för Remel Inc.
©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.