

Sample	Mean O.D.	Standard Deviation	%CV
1	1.40	0.014	1.0
2	0.78	0.038	3.8
3	0.26	0.008	3.1
4	0.12	0.004	2.9

The intra-assay or within-run CV was evaluated by testing 24 wells with each of 4 positive specimens. The mean intra-assay CV was 4.23%.

Sample	Mean O.D.	Standard Deviation	%CV
1	1.62	0.047	2.9
2	0.91	0.053	5.8
3	0.53	0.019	3.5
4	0.38	0.018	4.7

CROSS-REACTIVITY

The ProSpecT Giardia EZ Microplate Assay has been tested with stool specimens found to be O&P positive for a number of faecal parasites. No cross reactivity was observed with any of the infectious agents listed below.

<i>Ascaris lumbricoides</i>	<i>Entamoeba histolytica</i>
<i>Blastocystis hominis</i>	<i>Enterobius vermicularis</i>
<i>Chilomastix mesnili</i>	Hookworm
<i>Cryptosporidium parvum</i>	<i>Hymenolepis nana</i>
<i>Dientamoeba fragilis</i>	<i>Iodamoeba butschlii</i>
<i>Endolimax nana</i>	<i>Strongyloides stercoralis</i>
<i>Entamoeba coli</i>	<i>Trichuris trichiura</i>
<i>Entamoeba hartmanni</i>	

14. BIBLIOGRAPHY

- Addiss, D.G., 1991.**
Evaluation of a commercially available Enzyme-Linked Immunoabsorbent Assay for *Giardia lamblia* Antigen in Stool.
J. Clin. Microbiol. 29(6):1137-1142.
- Black, R.E. et al., 1977.**
Giardiasis in day care centres: evidence of person to person transmission.
Paediatrics 60(4):486-491.
- Centre for Disease Control, Atlanta, GA, 1978.**
Intestinal parasite surveillance - United States, 1976.
Morbidity and Mortality Weekly Report. 27(20):167-168.
- Craun, G.F., 1979.**
Waterborne Giardiasis in the United States: A Review.
Am. J. Public Health. 69(8):817-819.
- Craun, G.F., 1986.**
Waterborne Giardiasis in the United States 1965-84.
Lancet I:513-514.
- Kean, B.H. et al., 1979.**
Epidemic of amoebiasis and giardiasis in a biased population.
Br. J. Ven. Dis. 55:375-378.
- Meyers, J.D. et al. 1977.**
Giardia lamblia infection in homosexual men.
Br. J. Ven. Dis. 53:54-55.
- Rendtorff, R.C., 1979.**
The experimental transmission of *Giardia lamblia* among volunteer subjects in waterborne transmission of Giardiasis.
(Jakubowski, W., and J.C. Hoff, eds.), pp. 64-8. EPA, U.S., Cincinnati, OH.
- Rosoff, J.D., and H.H. Stibbs. 1986.**
Isolation and identification of a *Giardia lamblia*-Specific Stool Antigen (GSA 65) Useful in Coprodiagnosis of Giardiasis.
J. Clin. Microbiol. 23(5):905-910.
- Rosoff, J.D., and H.H. Stibbs., 1986.**
Physical and chemical characterization of a *Giardia lamblia*-Specific Antigen useful in the coprodiagnosis of Giardiasis.
J. Clin. Microbiol. 24(6):1079-1083.
- Rosoff, J.D. et al., 1989.**
Stool diagnosis of Giardiasis using a commercially available enzyme immunoassay to detect *Giardia*-Specific Antigen 65 (GSA 65).
J. Clin. Microbiol. 27(9):1997-2002.
- Schieven, B.C. and Z. Hussain., 1990.**
Evaluation of an enzyme immunoassay test kit for diagnosing infections with *Giardia lamblia*.
Serodiag and Immunother. 4:109-113.
- Sonnad, S., L. Bahrami, P. O'Hanley, 1991.**
Compatibility Assessment of Three Common Transport Media Systems with the ProSpecT™ / *Giardia* Immunoassay.
Presented at the 1991 American Society for Microbiology Meeting, Session 20.
- Tijssen, P.,**
Practice and Theory of Enzyme Immunoassays.
Laboratory Techniques in Biochemistry and Molecular Biology, R.H. Burdon and P.H. van Knippenberg, eds., Elsevier, N.Y., 1985, pp.14-16.
- Visvesvara, G.S., 1982.**
Giardiasis in Children.
J. Pediatric Gastroenterol Nutr. 1(4):463-465.

- Wolfe, M.S., 1979.**
Managing the patient with giardiasis: clinical, diagnostic and therapeutic aspects.
Waterborne Transmission of Giardiasis. (W. Jakubowski, and J.C. Hoff, eds.), pp.39-52. EPA, U.S., Cincinnati, OH.

ProSpecT™ is a registered trademark



Oxoid Ltd Wade Road Basingstoke Hants, RG24 8PW
UK

For technical assistance please contact your local distributor.

IFU X7590A, Revised January 2014

Printed in the UK