

The results of these evaluations are presented below:

		Acid Fast	
		+	-
ProSpecT	+	81	0
Cryptosporidium	-	0	133
		81	133 214

Sensitivity 81/81 = 100% (95.5 - 100%)

Specificity 133/133 = 100% (97.3 - 100%)

Numbers in parenthesis are 95% confidence intervals.

Clinical studies were conducted to evaluate the performance of the ProSpecT Cryptosporidium Microplate Assay. Specimens were obtained from hospital labs and the CDC. Patient populations represented in the specimen pool were symptomatic patients in normal prevalence populations, symptomatic patients in a high prevalence (HIV positive) population and asymptomatic patients from a day care population. Specimens were submitted unpreserved or preserved in 10% formalin or SAF. Samples were tested for *Cryptosporidium* by either acid-fast (AF) or immunofluorescent staining methods (IFA). A total of 212 specimens were tested; 134 were positive for *Cryptosporidium* specific antigen (CSA) and 78 were negative. The results with the ProSpecT Cryptosporidium Microplate Assay are presented below:

		Acid Fast	
		+	-
ProSpecT	+	130	0
Cryptosporidium	-	4	78
		134	78 212

Sensitivity 130/134 = 97% (92.5 - 99.2%)

Specificity 78/78 = 100% (95.4 - 100%)

Numbers in parenthesis are 95% confidence intervals.

A prospective trial was conducted at a large metropolitan hospital. All samples submitted for acid-fast staining for *Cryptosporidium* over a period of 4 months were included in the study. Samples were unpreserved and frozen at -20°C prior to testing with the ProSpecT Cryptosporidium Microplate Assay. The results of initial testing and the resolved data are presented below. Data was resolved by repeat testing of the 14 AF negative/ CSA positive samples. Six of the fourteen were reproducibly positive for CSA. Specific inhibition studies with antibody to CSA showed greater than 50% inhibition in all 6 samples. These 6 samples are considered to be true positives in the resolved data.

		Acid Fast		Resolved	
		+	-	+	-
ProSpecT	+	28	14	34	8
Cryptosporidium	-	1	335	1	335
		29	349	35	343 378

Sensitivity 28/29 = 97% (82.2 - 99.9%) 34/35 = 97% (85.1 - 99.9%)

Specificity 335/349 = 96% (93.4 - 97.8%) 335/343 = 98% (95.5% - 99.0%)

Numbers in parentheses are 95% confidence intervals.

ANALYTICAL SENSITIVITY

The ProSpecT Cryptosporidium Microplate Assay detects approximately 20 nanograms/ml of CSA.

REPRODUCIBILITY

The inter-assay or run-to-run coefficient of variation (CV) of the ProSpecT Cryptosporidium Microplate Assay was evaluated by selecting 10 positive specimens with varying optical density readings. Each sample was tested in 10 wells per day for five days. The mean inter-assay CV was 10.6%.

The intra-assay or within-run CV was evaluated by testing 24 wells with each of 5 positive specimens. The mean intra-assay CV was 2.52%.

CROSS-REACTIVITY

The ProSpecT Cryptosporidium Microplate Assay has been tested with stool specimens found to be O&P positive for a number of faecal parasites. No cross-reactivity was observed with any of the infectious agents listed below.

<i>Ascaris lumbricoides</i> (2)	<i>Giardia lamblia</i> (5)
<i>Blastocystis hominis</i> (4)	<i>Hymenolepis nana</i> (2)
<i>Chilomastix mesnili</i> (1)	<i>Iodamoeba butschlii</i> (2)
<i>Dientamoeba fragilis</i> (4)	<i>Isospora belli</i> (2)
<i>Endolimax nana</i> (3)	<i>Strongyloides stercoralis</i> (2)
<i>Entamoeba coli</i> (6)	<i>Taenia solium</i> (1)
<i>Entamoeba hartmanni</i> (2)	<i>Trichuris trichiura</i> (1)
<i>Entamoeba histolytica</i> (5)	

The numbers in parentheses indicate the number of specimens tested.

Cryptosporidiosis in Patients with AIDS.

J. Infect. Dis. 155:150.

- 11. Stehr-Green, J.K. et al., 1987.**
Shedding of Oocysts in immunocompetent individuals infected with *Cryptosporidium*.
Am. J. Trop. Med. 1 Hyg 36(2):338-342.
- 12. Taylor, J.P. et al., 1985.**
Cryptosporidiosis outbreak in a day-care centre.
Am. J. Dis. Child. 139:1023-1025.
- 13. Ungar, BLP. 1990.**
Enzyme-linked immunoassay for detection of *Cryptosporidium* antigens in faecal specimens.
J. Clin. Microbiol., 28:2491.

ProSpecT™ is a registered trademark of Remel Inc.



Oxoid Ltd Wade Road Basingstoke Hants, RG24 8PW UK

For technical assistance please contact your local distributor.

IFU X7591A, Revised April 2012

Printed in the UK

14 BIBLIOGRAPHY

- 1. Alpert, G. et al., 1986.**
Outbreak of Cryptosporidiosis in a Day-Care Centre.
Pediatrics 77(2):152-157.
- 2. Anon. 1984.**
Cryptosporidiosis among children attending day-care centres - Georgia, Pennsylvania, Michigan, California, New Mexico.
Morbidity Mortality Weekly Rep. 33(42):599-601.
- 3. Arrowood, M.J. and C.R. Sterling, 1989.**
Comparison of Conventional Staining Methods and Monoclonal Antibody-based Methods for *Cryptosporidium* Oocyst Detection.
J. Clin. Microbiol. 27(7):1490-1495.
- 4. Chapman, P.A., B.A. Rush and J. McLauchlin, 1990.**
An enzyme immunoassay for detecting *Cryptosporidium* in faecal and environmental samples.
J. Med. Microbiol. 32:233-237.
- 5. D'Antonio, R.G. et al., 1985.**
A waterborne outbreak of cryptosporidiosis in normal hosts.
Ann. Intern. Med. 103:886.
- 6. Dubey, J.P., C.A. Speer and R. Fayer, eds., 1990.**
Cryptosporidiosis of Man and Animals.
CRC Press.
- 7. Gallaher, M.M. et al., 1989.**
Cryptosporidiosis and Surface Water.
Am. J. Public Health 79(1):39-42.
- 8. Hayes, E.B. et al., 1989.**
Large Community Outbreak of Cryptosporidiosis Due to Contamination of a Filtered Public Water Supply.
N. Eng. J. Med. 320(21):1372-1376.
- 9. Leech, J.H., M.A. Sande and R.K. Root, eds., 1988.**
Parasitic Infections.
Churchill Livingstone.
- 10. Navin, T.R. and A.M. Hardy, 1984.**