

respiratory syncytial virus in clinical specimens by viral culture, direct and indirect immunofluorescence, and enzyme immunoassay. *J. Clin. Microbiol.* **26**: 588-591.

19. **Bartholoma N.Y. and Forbes B.A. (1989)** Successful use of shell vial centrifugation and 16 to 18 hour immunofluorescent staining for the detection of Influenza A and B in clinical specimens. *Am J Clin Pathol* **92**: 487-490.

20. **Krech T., Gerhard-Fsadni D., Hofmann N., Miller S.M. (1985)** In-terference of Staphylococcus aureus in the detection of *Chlamydia trachomatis* by monoclonal antibodies. *The Lancet*, 1161-1162.

21. **Tyeryar F.J. (1983)** Report of a workshop on Respiratory Syncytial Virus and Parainfluenza Viruses *Journal of Infectious Diseases* **148**: 588-598.`

22. **Caul E.O. (1986)** Personal communication (on file at Oxoid (Ely) Ltd).

23. **Chin D.Y., Mosley W.H., Poland J.D., Rush D., Belden E.A., Johnson O. (1963)** Epidemiologic Studies of Type B influenza in 1961-1962. *American journal of Public Health* **53**: 1068-1074

24. **McLean, D.M., Bannatyre, R.M. and Givan, K.F. (1967)** Myxovirus dissemination by air. *Canadian medical Association Journal* **89**: 1257-59



X7852A Marts 2017

 OXOID Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

For teknisk assistance skal du kontakte din lokale forhandler.

Brugsanvisninger og versioner på andre sprog kan downloades fra www.oxoid.com/IFU