

# Easy, rapid delivery of blood culture samples



## The Thermo Scientific VersaTREK REDOX Transporter for pneumatic tube systems

The Thermo Scientific™ VersaTREK™ REDOX™ Transporter is easy to use and provides excellent shock absorption for blood culture bottles sent through pneumatic tube systems. Made from a sturdy polypropylene resin and able to withstand the strain of continuous use, the device can be easily cleaned with common disinfection agents.

The REDOX Transporter offers:



Preformed bottle sections to easily accommodate either the REDOX 40mL or 80mL bottle, eliminating the need for two separate carriers. Bottle sets travel together!



Simple locking clasp closures for ease of use



Compatibility with most 3-, 4- and 6-inch pneumatic tube carriers\*



\*6428-30 dimensions: 10.9" L x 3.3" W x 2.7" H  
6428-33 dimensions: 10.8" L x 2.7" W x 2.8" H

# REDOX Transporter instructions:

1

**Examine Transporter for signs of wear. If excessive wear is noted, please discontinue use and obtain new Transporter.**

2

**Place bottle set into the REDOX Transporter.**

3

**Close the REDOX Transporter lid over the bottles, and press the two latches downward until the transporter is closed.**

4

**Place the REDOX Transporter into a large biohazard bag.**

5

**Place the REDOX Transporter into a pneumatic tube carrier.**

When utilizing a six-inch pneumatic tube carrier, it is highly recommended that some type of foam insert be used with the REDOX Transporter to prevent the Transporter from moving within the carrier during transport.

6

**Close the carrier. The REDOX Transporter is now ready to be sent through the pneumatic tube system.**



The REDOX Transporter (Part No. 6428-30 or 6428-33) is sold 10/box. If you would like to order this product, please contact your Thermo Fisher Scientific Microbiology Account Manager or Customer Service at 800.871.8909.

Pictures are for instructional use of the product and are not an exact representation of the product received.

©2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

**Contact Information:**

microbiology@thermofisher.com

USA +800 255 6730

International +44 (0) 1256 841144

LT2011B

July 2019

**ThermoFisher**  
SCIENTIFIC