


Phenom ParticleX AM Desktop SEM

Multi-purpose desktop SEM delivering purity at microscale





Automatically detect particles, satellites, and conglomerates in additive manufacturing powder.

A growing number of manufacturing companies are establishing scanning electron microscopy (SEM) systems in-house. This trend, from outsourcing to in-house analysis, is growing and the benefits, such as the ability to perform a broad range of automated desktop analyses, chemical classification and verification according to specific norms are clear.

Timely and accurate quality control are prerequisites for today's manufacturing. The Thermo Scientific™ Phenom™ ParticleX Desktop SEM is a versatile solution for high-quality analysis in-house. It gives you the ability to carry out speedy analysis, verification and classification of materials, supporting your production with fast, accurate and trusted data. The system is automated and offers multiple sample analysis, making testing and classification up to 10 times faster. Outsourcing typically takes up to 10 working days, whereas the Phenom ParticleX Desktop SEM gives you certainty within one day. The system is simple to operate and fast to learn, opening up the use of particle and material analysis to a wider group of users in-house. In addition to eliminating the need to outsource, the ease-of use and automation of the Phenom ParticleX Desktop SEM allow you to offload sample analysis from other SEMs in your laboratory.

The Phenom ParticleX Desktop SEM not only provides high quality SEM analysis, it is also designed to perform particle analysis of metal powders at the microscale for the additive industry, confirming that the product fulfills the quality requirements. All now made possible in-house and on your desktop.

Perception software

Automated additive powder analysis

Intuitive user interface limits training time

Day and night SEM operation enhances efficiency

Automated reporting with dedicated report templates

Reclassification and requantification of results

Suitable for all types of additive powder

High speed analysis with accurate and repeatable results

Versatile output formats: PDF, Word, Excel, CSV

Additive powder analysis

The Phenom ParticleX Desktop SEM for additive manufacturing enables automated scanning electron microscopy EDS. It is the proven solution for monitoring the three most critical characteristics of metal powders for powder-bed and powder-fed additive manufacturing processes: particle size distributions, individual particle morphology and foreign particles.


The Phenom ParticleX Desktop SEM for additive manufacturing measures various size and shape parameters, such as minimum and maximum diameter, perimeter, aspect ratio, roughness, and feret diameter. All of which can be displayed with 10%, 50%, 90% values (e.g. d10, d50, d90).

The supplied analysis recipes allow you to differentiate between satellite particles, spherical particles and deformed particles. Particles size distributions can be plotted as number-based or volume-based.

General usage

The Phenom ParticleX Desktop SEM features a chamber which includes an accurate and fast motorized stage that allows analysis of samples of up to 100 mm x 100 mm. In spite of this larger sample size, a proprietary loading shuttle keeps the loading cycle to an industry-leading 40 seconds or less. In practice this improves the throughput factors higher than other SEM systems.

The user interface is based on the proven ease-of-use



Specific parameters can be individually set.

technology applied in the successful Phenom Desktop SEM products. The interface enables both existing and new users to quickly become familiar with the system with a minimum of training.


The standard detector in the Phenom ParticleX Desktop SEM is a four-segment backscattered electron detector (BSD) that yields sharp images and provides chemical contrast information together with a fully integrated energy dispersive X-ray (EDS) system for elemental analysis. A secondary electron detector (SED) for surface sensitive imaging is optional.

Elemental analysis is provided by EDS technology, which allows users to analyze the chemical composition of their samples. Detailed chemical composition can be obtained from a micro volume via a spot analysis. Elemental distribution can be visualized with the elemental mapping option.


Elemental mapping and line scan

For the user, it is simply click and go to work with the elemental mapping and line scan functionality of the Phenom ParticleX Desktop SEM. The elemental mapping functionality visualizes the distribution of elements throughout the sample, and selected elements can be mapped at a user-specified pixel resolution and acquisition time.

The real-time mapping algorithm shows live build-up of



Deformed particles can easily be identified.




Poor powder characteristics can lead to 3D-printed failures.

the selected elements. The line scan functionality shows the quantified element distribution in a line plot. This is especially useful for coatings, paints and other applications with multiple layers for analyzing edges, coatings, and cross sections. Results of both the elemental mapping and line scan functionality can be easily exported by using an automated report template.

Secondary electron detector

A secondary electron detector (SED) is optionally available on the Phenom ParticleX Desktop SEM. The SED collects low-energy electrons from the top surface layer of the sample, making it the perfect choice to reveal detailed sample surface information. The SED can be of great use for applications where topography and morphology are important. This is often the case when studying microstructures, fibers or particles.

The integrated EDS allows chemical typing of each individual particle, allowing you to readily identify any foreign particulate in your powder from previous printing cycles.



Plot powder size distributions as volume-based or number-based.

Imaging	
Detection modes	
Light	Magnification range: 3–16x
Electron	Magnification range: 160–200,000x
Illumination	
Light	Bright field / dark field modes
Electron	<ul style="list-style-type: none"> Long lifetime thermionic source (CeB₆) Multiple beam currents
Acceleration voltages	<ul style="list-style-type: none"> Default: 5 kV, 10 kV and 15 kV Advanced mode: adjustable range between 4.8 kV and 20.5 kV imaging and analysis mode Perception: 15 kV
Vacuum levels	Low - medium - high
Resolution	<10 nm
Detector	
Standard	Backscattered electron detector, energy dispersive spectroscopy detector
Optional	Secondary electron detector
Digital image detection	
Light optical	Proprietary high-resolution color navigation camera, single-shot
Electron optical	High-sensitivity backscattered electron detector (compositional and topographical modes)
Image formats	
JPEG, TIFF, PNG	
Image resolution options	
960 x 600, 1920 x 1200, 3840 x 2400 and 7680 x 4800 pixels	
Data storage	
USB flash drive, Network, workstation with SSD	
Sample stage	
Computer-controlled motorized X and Y	

EDS	
Detector type	<ul style="list-style-type: none"> Silicon Drift Detector (SDD) Thermoelectrically cooled (LN₂ free)
Detector active area	25 mm ² or 70 mm ²
X-ray window	Ultra thin silicon nitride (Si ₃ N ₄) window allowing detection of elements B to Am
Energy resolution	Mn Kα ≤132 eV
Processing capabilities	Multi-channel analyzer with 2048 channels at 10 eV/ch
Max. input count rate	300,000 cps
Hardware integration	Fully embedded
Software	
<ul style="list-style-type: none"> Integrated in Phenom user interface Integrated column and stage control Auto-peak ID Iterative strip peak deconvolution Confidence of analysis indicator Export functions: CSV, JPG, TIFF, ELID, EMSA 	
Report	
DOCX format	
Elemental Mapping & Line Scan	
Elemental mapping	
Element selection	Individual user-specified maps, plus backscatter image and miximage
Backscatter image and mix-range	
Selected area	Any size, rectangular
Mapping resolution range	32 x 20 to 960 x 600 pixels
Pixel dwell time range	1–500 ms
Line scan	
Line scan resolution range	16–512 pixels
Line scan dwell time range	10–500 ms
Report	
DOCX format	
SED	
Detector type	Everhart Thornley

System

Dimensions & weight

Imaging module	316(w) x 587(d) x 625(h) mm, 75 kg
Diaphragm vacuum pump	145(w) x 220(d) x 213(h) mm, 4.5 kg
Power supply	260(w) x 260(d) x 85(h) mm, 2.3 kg
Monitor (24")	531(w) x 180(d) x 511(h) mm, 5.6 kg
Workstation	<ul style="list-style-type: none"> Powerful workstation including SSD storage 93(w) x 293(d) x 290(h) mm, 5.6 kg

Sample size

- Max. 100 mm x 100 mm (up to 36 x 12 mm pin stubs)
- Max. 40 mm height (optional up to 65 mm)

Scan area

100 mm x 100 mm

Sample loading time

Light optical	<5 s
Electron optical	<60 s

Site requirements

Ambient conditions

Temperature	15°C ~ 30°C (59°F ~ 86°F)
Humidity	Between 20% and 80% RH
Power	Single phase AC 100–240 Volt, 50/60 Hz, 163 W average, 348 W max

Recommended table size

150 x 75 cm, load rating of 150 kg

Sample holders and inserts

- Manual-Z sample holder
- 49-stub insert (optional)

 Learn more at thermofisher.com/phenom-particle-x-am

thermoscientific

For research use only. Not for use in diagnostic procedures. For current certifications, visit thermofisher.com/certifications

© 2024 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. DS0294-EN-05-2024