

## Helios 5 Laser System

### Fastest high-quality subsurface and 3D characterization at millimeter scale with nanometer resolution

The Thermo Scientific Helios 5 Laser System delivers unmatched capabilities for extreme large-volume 3D analysis, high-throughput sample preparation and nanoprototyping. Featuring an innovative, fully integrated femtosecond laser, it offers the fastest material removal rate with the highest cut face quality.

The Thermo Scientific™ Helios™ 5 Laser System is part of the fifth generation of the industry-leading Helios family. It combines the best-in-class Thermo Scientific Elstar™ SEM column, for ultra-high-resolution imaging and best analytical capabilities, with a high-throughput Tomahawk HT focused ion beam (FIB) column, for superior performance at all operating conditions, and a femtosecond laser, enabling *in situ* ablation at unprecedented material removal rates for fast millimeter-scale characterization at nanometer resolution.

#### Statistically relevant subsurface and 3D characterization

The femtosecond laser provides the capability to cut materials orders of magnitude faster than a FIB. For many materials, a large cross-section of hundreds of microns can be created within 5 minutes. Featuring a different removal mechanism (ablation vs ion sputtering), the laser can easily process challenging materials, such as non-conductive or ion beam-sensitive samples; for example: glass, ceramics, hard and soft polymers, biomaterials, graphite, etc. Moreover, a protective coating to reduce charge or to mitigate the curtain artifacts is not needed, further saving time and increasing the total throughput.

The Helios 5 Laser System is designed with all three beams (SEM/FIB/laser) pointing at a single coincident point, enabling very short laser/SEM switching times and no vacuum transfers. This provides high-precision *in situ* feedback, making it a unique solution for fully automated, ultra-large-volume 3D characterization using laser tomography.

#### Key Benefits

**Fastest material removal** for millimeter-scale cross sections with a material removal rate up to 15,000x faster than a typical Ga FIB

**Statistically relevant subsurface and 3D data analysis** by acquiring much larger volumes within a shorter amount of time


**Accurate and repeatable cut placement** with triple beam coincidence on the sample

**Fast characterization of deep subsurface features** via extraction of subsurface TEM lamella or chunks for 3D analysis

**High-throughput processing of challenging materials** such as non-conductive or ion beam-sensitive

**Fast and easy characterization of air-sensitive samples** without the need to transfer samples between different instruments for imaging and cross-sectioning

**All capabilities of proven Helios 5 UC platform**, including high-quality TEM and APT sample preparation and extreme high-resolution imaging capabilities




Laser tomography of a carbon fiber-reinforced composite in epoxy.

## High-quality, large-volume analysis

The ultra-short duration of the femtosecond laser pulses introduces almost no artifacts while delivering much higher site-specific accuracy compared to traditional mechanical polishing techniques, where, for example, delamination is often an issue. The femtosecond laser induces no heat impact, no microcracks, and no melting effects. In most cases, femtosecond laser-milled surfaces are clean enough for direct SEM imaging, and the quality is often sufficient for surface-sensitive EBSD mapping. When it is necessary to even further improve the surface quality, a FIB polishing procedure can be used to reveal ultra-fine features.


Built on the proven Helios 5 UC platform, this instrument incorporates a suite of state-of-the-art technologies and provides the highest performance for common DualBeam™ use cases, such as high-resolution S/TEM and atom probe tomography (APT) sample preparation and extreme high-resolution SEM imaging with the most precise materials contrast.



Schematics of the Helios 5 Laser System. Three beams converging at a single coincident point to allow fast switching between SEM imaging and laser processing, provide accurate and repeatable cut placement, increase throughput, and enable 3D characterization.

## Helios 5 Laser System enables new applications


- Large-volume cross-sectioning for rapid failure analysis
- Fast access to buried subsurface layers often not accessible by (P)FIB
- Fast, precise micromachining of complex shapes, tensile rods, chunks, microCT specimens, etc.
- Millimeter-scale 3D laser serial sectioning, including analytical characterization with EDS or EBSD
- Easy and fast characterization of air-sensitive samples (e.g., batteries) with femtosecond laser processing and high-resolution SEM imaging in the same vacuum chamber
- Fast, high-quality processing of challenging materials (e.g., carbon-based, ion beam-sensitive, and non-conductive)
- Extended range of correlative microscopy with deep subsurface sample extraction




Comparison of representative 3D volumes acquired within the same amount of time with FIB, plasma FIB, and femtosecond laser.

## Overview of Helios 5 Laser System applications

### Metals (conductive)



Steel produced by additive manufacturing.




3D quantitative analysis of Li-ion battery graphite anode.

### Semiconductors (conductive/non-conductive)




GPU stack.




3D reconstruction of a solder bump.

### Polymers (non-conductive, electron beam-sensitive)




Nylon 6.6 (textile).




Carbon fiber-reinforced epoxy composite.

### Glass/ceramics (non-conductive, thermally sensitive, shock-sensitive)



Packaged solid-state battery.



Glass fibers.



Description	FIB	PFIB	Femtosecond Laser
Max volume*	40x40x40 $\mu\text{m}^3$	200x200x200 $\mu\text{m}^3$	2,000x2,000x1,000 $\mu\text{m}^3$
Max current	65 nA	2.5 $\mu\text{A}$	~1 mA (ion-beam-current equivalent)
Slicing current**	10 nA	180 nA	74 $\mu\text{A}$
Spot size	1–4 nm	15 nm–15 $\mu\text{m}$	15 $\mu\text{m}$

\* Typical max volumes acquired within the same amount of time.

\*\* 3D slicing currents used in practice for many materials.

### Femtosecond laser specifications

<b>Laser integration</b>	Fully integrated in the chamber with the same coincident point of all 3 beams (SEM/FIB/laser), enabling accurate and repeatable cut placement and 3D characterization.	
<b>First harmonic</b>	Wavelength	1,030 nm (IR)
	Pulse duration	<280 fs
<b>Second harmonic</b>	Wavelength	515 nm (green)
	Pulse duration	<300 fs
<b>Optics</b>	Coincident point	WD = 4 mm (same as SEM/FIB)
	Objective lens	Variable (motorized)
	Polarization	Horizontal/vertical
<b>Repetition rate</b>	1 kHz–1 MHz	
<b>Beam position accuracy</b>	<250 nm	
<b>Protective shutter</b>	Automated SEM/FIB protective shutter	
<b>Software</b>	Laser control software	
	Laser 3D serial sectioning workflow	
	Laser 3D serial sectioning workflow with EBSD	
	Laser Scripting*	
<b>Safety</b>	Interlocked laser enclosure (Class 1 laser safety)	

\* With optional AutoScript 4 Software.

Learn more at [thermofisher.com/TriBeam](https://thermofisher.com/TriBeam)

thermo scientific

For research use only. Not for use in diagnostic procedures. For current certifications, visit [thermofisher.com/certifications](https://thermofisher.com/certifications)

© 2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. DS0375-EN-04-2022