

XRF analysis

Thermo Scientific UniQuant Software

Program for quantitative XRF analysis using Thermo Scientific ARL WDXRF Spectrometers

Thermo Scientific™ UniQuant Software is a complete analysis package for standardless semi-quantitative to quantitative XRF analysis using intensities measured by a sequential X-ray spectrometer. It represents a unique method in XRF analysis.


Introduction to UniQuant Software

The program is highly effective for analyzing samples for which no standards are available. As its name suggests, UniQuant Software unifies all types of samples into one analytical method. It is unique in that respect. Sample preparation is usually minimal or not required at all. Samples may be of very different nature, sizes, and shapes. Elements from fluorine (F) up to americium (Am)

(or their oxide compounds) are analyzed in very diverse samples like a piece of glass, a screw, metal drillings, lubricating oil, loose fly ash powder, polymers, phosphoric acid, thin layers on a substrate, soil, paint, the year rings of trees and in general those samples for which no standards are available.

Elements Be to O can also be analyzed in some applications providing that the required crystals are installed on the instrument goniometer. The reporting is in weight % along with an estimated error for each element.

The UniQuant Software application is integrated within Thermo Scientific™ OXSAS™ XRF Analysis Software as an option and is available for Thermo Scientific XRF sequential instruments.



UniQuant Software uses peak to peak acquisition allowing for best limits of detection.

Major features

- UniQuant Software is calibrated as standard for the analysis of 76 elements (F to Am).
- A new calibration (daughter calibration) can be easily derived from the original calibration (parent calibration), thus enabling to make calibrations which are specific to given matrices in order to further improve the accuracy of analysis using only a few international or in-house standards
- All calibrations are maintained over time by the setting-up samples (SUS) program by means of intensity correction directly after measurements. This program is based on 'pure' metal samples, very stable, almost indestructible, and easy to clean. One can compare the intensities from new measured samples with those measured years ago. This maintenance feature is unique and a huge time saver.
- More flexibility in the definition of compounds and possibility to use two XRF results to deduce compound concentrations through calculation or to calculate a surplus as element or oxide. This leads to more accurate results, for example for sulfides and stearates.
- As specific filters can be used to modify the X-ray excitation, any such tube filter can be specified per group of XRF lines.
- Ultra-light elements Be to O can also be analyzed in some applications providing that the necessary crystals are installed on the instrument. These elements require a matrix specific calibration, i.e., the calibration is no longer universal for these ultra-light elements.

Table 1 shows which 79 elements may be analyzed by UniQuant Software. UniQuant Software calculates absorption by H and Li. Argon has been included because it may be found in materials made under an Argon atmosphere.

Nature of samples

As stated above, the unknown samples may take a great variety of physical forms such as a:

- Solid disk of metal or a synthetic material
- Multi-element monolayer on a substrate
- Small piece of solid sample placed on a supporting film
- Pressed powder that may include a binder
- Very small amount of powder on a supporting film (usually in helium atmosphere)
- Solid solution of a mineral (a glass bead)
- Liquid sample from a small drop to a full cup (only in spectrometers that allow analysis of liquids in helium atmosphere)
- Filter aerosol sample

Analysis time

A totally unknown sample may be measured by the prescribed measuring channels (122 spectral positions) for determination of 76 elements (F to Am). The spectrometer time then is about 14 minutes. Ultra-light elements can be included in special applications with specific calibration and crystals. When included the spectrometer time is about 18 minutes.

Samples belonging to a known family may be measured by a smaller subset of measuring channels. For example, the analysis of routine waste disposal samples may be limited to 55 measuring channels with extra-long measuring times for relevant traces. The spectrometer time can then be as low as 5 minutes.

H																	<i>He</i>
Li	Be											B	C	N	O	F	<i>Ne</i>
Na	Mg											Al	Si	p	S	Cl	Ar
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	<i>Kr</i>
Rb	Sr	Y	Zr	Nb	Mo	<i>Tc</i>	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	<i>Xe</i>
Cs	Ba	L	Hf	Ta	w	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	<i>Po</i>	<i>At</i>	<i>Rn</i>
<i>Fr</i>	<i>Ra</i>	A															
		L	La	Ce	Pr	Nd	<i>Pm</i>	Sm	Eu	Gd	Tb	Dy	Ho	Er	Tm	Yb	Lu
		A	<i>Ac</i>	Th	<i>Pa</i>	U	<i>Np</i>	Pu	Am	<i>Cm</i>	<i>Bk</i>	<i>Cf</i>	<i>Es</i>	<i>Fm</i>	<i>Md</i>	<i>No</i>	<i>Lr</i>

Bold= analyzed elements / *Italic* = not known by UniQuant

Table 1: Elements analyzed.

Accuracy for majors and minors

The accuracy for concentrations higher than 1 weight % primarily depends on the physical nature of the sample. The errors are smallest for thick, full-area homogeneous samples and are quite acceptable for less favorable physical conditions. As a rule of thumb, the standard error in weight % for a major or minor constituent is equal to

$$\text{StdErr in Wt \%} = K \cdot \sqrt{\frac{\%C}{100}}$$

where $K=0.05$ to 0.15 depending on the element and the physical nature of the sample. UniQuant Software is intended to cover the widest possible concentration ranges while using one single set of calibration data. Here we are not thinking about a wide range of alloys or of oxide samples. The range that we mean includes samples like oils, polymers, beads, thin layers and all types of alloys.

For specific applications, where very high accuracy is required, UniQuant Software may use specially calibrated data sets, for example one for alloys and one for fused beads or glass. Then international or the client's own standards are used to firm up the calibration. This way of working may lead to the same high accuracy as with conventional analysis using regression analysis of standards.

Although using specialized data sets has not been the primary philosophy behind UniQuant Software, its application allows replacing conventional methods by the UniQuant Software with far less specialized analytical programs. Several UniQuant Software users have indeed done so.

Trace analysis

Precision and limits of detection

Precision (reproducibility) of the analysis of a given sample specimen depends only on counting statistics. For each analyzed element, UniQuant Software reports the standard deviation (sigma) in ppm which takes into account counting precision and the corrections that have been applied for background and line overlaps.

For large (full area) samples that are not highly diluted, the sigma's are surprisingly small, for example 1 or 2 ppm for measuring times of 4 to 10 seconds per analytical line. The sigma is smallest with lighter matrix samples, for higher atomic numbers, and with longer measuring times.

Trace elements (with $Z>20$) in heavier matrices can be well determined from 20 ppm onward. For light matrix samples like polymers, this value is 5 ppm or even lower. Thanks to peak-by-peak acquisition the counting time per element is much longer than in the case of a scan (40 to 100 times longer). This provides the best possible limits of detection.

Accuracy

The accuracy for traces is depending on the quality of the corrections made for:

- Background—well done by UniQuant Software
- Spectrometer's spectral impurities—well done by UniQuant Software
- Spectral line overlaps—uniquely solved by UniQuant Software. Very important!
- Matrix effects—solved by FP (fundamental parameters)
- Physical effects—UniQuant Software has unique ways to compensate for certain physical effects

Thanks to the counting times used which are multiples of 2 seconds, spinning of the sample during analysis (one turn every 2 seconds) can be activated to even out inhomogeneities of the sample like grinding grooves on metals.

Interactivity

UniQuant Software has been designed for a maximum of interactivity. A pre-condition is speed of calculation. The user interface is designed for a minimum of keystrokes or mouse operations.

The need for fast interactivity is illustrated by the following example:

UniQuant Software evaluates a totally unknown powder sample in a first calculation (5 seconds) for which the analyst assumes that this mineral sample consists of oxides. The results, however, show a very high content of sulfur. The analyst concludes that the sample is a sulfide ore. This means that elements like Pb, Zn, Fe, Mo are present as sulfides while elements like Si and Al occur as oxides. The original assumption assumed most elements to be as oxides, even sulfur. The sum of concentrations ended up at a level higher than 100%. Now, the analyst just changes oxides to sulfides in the input of the general data table and starts a second calculation. All this in a matter of just a few seconds.

On the other hand, UniQuant Software can be used on routine unknown samples that do not require advanced treatment by an analyst. Once the analyst has refined the proper calculation for a given sample type, measurement of other samples of the same type can be easily taken over by any operator.

Thin layer samples

Mono-layer samples

UniQuant Software can calculate the sample mass along with its associated standard deviation. At the same time the composition of the layer is calculated. If the layer is on a substrate with elements that are also present in the layer, UniQuant Software will take their effect into account.

Calculation of multiple samples

In order to short-cut a lot of work in case of multiple samples, a batch mode exists. Samples are “tagged” in a list and the “process” is started. It may be used for evaluating a suite of similar unknown samples or calibration samples as well.

Reporting

The results of calculations are in terms of % and ppm and can be viewed on screen or printed as a comprehensive form that appeals to the analyst. This form can be saved on a disk in various formats like DOC, XLS, or PDF or re-directed to a specified file on disk.

In addition, a report can be printed in a standard form (in the order of atomic number) that is intended for the analyst's client where elements (and oxides) are presented in the order of atomic number Z, or listed in descending order of concentration, or in alphabetic order of element names.

Each reported concentration is accompanied by a 'StdErr'. The practical confidence interval is ± 2 StdErr., but can be modified to 1 StdErr or 3 StdErr to remove false positives.

UniQuant Software as an OXSAS Software option

The UniQuant Software option includes:

- A license for the use of the UniQuant Software option within the OXSAS Software
- An extensive online description with amongst others:
 - Instructions for setting up the spectrometer
 - Instructions for initial calibration and maintenance
 - Handling of unknown samples
 - The theory about XRF and UniQuant Software
- A set of 10 calibration samples for initial setting-up and for maintenance
- Some small hardware, such as centering rings, to facilitate sample presentation
- An optional set of 53 pure element samples is available. The total set of samples is used for factory calibration of UniQuant Software.

UniQuant Software is calibrated at the factory for elements fluorine (Z=9) to americium (Z=95) on the goniometer of the Thermo Scientific™ ARL™ PERFORM'X™ XRF Spectrometer Series, Thermo Scientific™ ARL™ X900 XRF Spectrometer Series, or Thermo Scientific™ ARL™ OPTIM'X™ XRF Spectrometers ensuring a constant high quality of analytical results. On special request, ultra-light elements can be calibrated for specific matrices if the appropriate crystals are installed.

Element	Certificate	UniQuant
Cu	87.9	86.1
Pb	1.23	1.6
Sn	5.7	6.0
Zn	4.29	4.7
Mn	<0.005	<
Al	<0.005	<
Fe	0.05	0.064
Ni	0.66	0.82
P	0.032	0.026
As	0.012	<
Si	<0.005	<
Sb	0.07	0.065
S	0.035	0.088

CDA 922 (copper alloy).

Element	Certificate	UniQuant
Ni	53.7	52.5
Cr	19.17	19.0
Mo	9.96	9.9
Co	10.59	11.0
Al	1.50	1.5
Ti	3.19	3.3
Cu	0.026	<
Fe	1.47	1.6
Si	0.07	0.041
Mn	0.02	<
P	0.004	<
S	<0.002	<
Nb	0.050	0.042
V	0.020	0.017

RENE 41 (nickel alloy).

Element	Certificate	UniQuant	
		Solid	Drillings
Si	0.20	0.18	0.13
P	0.030	0.033	0.015
S	0.024	0.025	0.079
V	0.94	0.90	0.86
Cr	5.17	5.5	5.5
Mn	0.21	0.23	0.23
Fe	(59.8)	59.9	60.4
Co	10.2	10.2	10.3
Mo	1.07	1.1	1.3
W	22.4	20.3	20.2

SS484 tool steel (solid and drillings).

Compound	Certificate	UniQuant	
		Powder	Fusion
SiO ₂	19.4	20.9	17.2
Al ₂ O ₃	3.1	3.3	2.2
FeO+Fe ₂ O ₃	16.6+2.0	18.2	19.6
MnO	18.6	17.8	18.9
CaO	32.6	30.6	30.7
MgO	8.0	6.9	8.8
P ₂ O ₅	0.47	0.44	0.44
V ₂ O ₅	0.11	0.11	0.12
TiO ₂	0.53	0.51	0.49
Cr ₂ O ₃	0.2	0.22	0.23
S	0.24	0.47	0.14
Na ₂ O	0.1	<	0.28
K ₂ O	0.10	0.08	0.09

N.B. Fusion is 15:1 Dilution

S1 (basic slag).

Element	Certificate	UniQuant
Ni	47.5	48.07
Cr	21.9	21.61
Fe	18.35	18.39
Mo	8.6	8.29
Co	1.53	1.53
Mn	0.63	0.6
Si	0.52	0.56
W	0.46	0.39
Cu	0.12	0.15
Al	0.11	0.12
Nb	0.1	0.09
V	0.08	0.08
P	0.017	0.02
Ti	0.011	0.01

Hast X (nickel alloy).

Both simple analysis results and advanced user screen are available to the analyst.

Compound	Certificate	UniQuant
Fe	1.20	1.25
Zn	0.12	0.055
Pb	0.14	0.10
Cu	0.018	0.008
Cr	0.05	0.075
Ni	0.017	0.013
SiO ₂	N.A.	19.4
Al ₂ O ₃	N.A.	2.9
CaO	N.A.	5.8
MgO	N.A.	1.1
Na ₂ O	N.A.	1.8
S	<0.002	<

Rest C,H,N Calculated by UniQuant = 64.3%

N.A. = not available

N.B. Benefits of Matrix Identification

Environmental sample.

Element	Certificate	UniQuant
Fe	56.4	55.89
Ni	24.96	25.14
Cr	14.06	14.18
Ti	2.15	2.25
Mo	1.08	1.03
V	0.26	0.264
Nb	0.23	0.231
Co	0.23	0.221
Al	0.19	0.186
Si	0.17	0.145
Mn	0.13	0.115
Cu	0.08	0.087
P	0.011	0.014
S	0.005	0.012

Steel A286.

Compound	Certificate	UniQuant
SiO ₂	75.7	74.5
Al ₂ O ₃	12.08	11.6
Fe ₂ O ₃	2.00	2.4
CaO	0.78	0.95
Na ₂ O	3.36	3.7
K ₂ O	4.99	5.7
F	0.42	(0.57)
MgO	(0.06)	0.06

Nim G granite (majors).

Compound	Certificate	UniQuant
Ba	(120)	110
Ce	198	220
Mn	160	210
Nb	53	110
Pb	40	<
Rb	325	470
Th	52	<
Ti	540	980
Y	147	230
Zn	50	<
Zr	300	540

Nim G granite (traces ppm).

Element	Certificate	UniQuant
Al	160	180
Si	480	390
Cl	55	87
Ca	140	150
Ti	49	56

Rest C,H,N Calculated by UniQuant = 99.83%

Polypropylene (ppm Levels).

Compound	Certificate	UniQuant
CaO	62.14	62.18
SiO ₂	21.24	20.83
Fe ₂ O ₃	4.40	4.33
MgO	4.02	4.08
Al ₂ O ₃	3.68	3.64
SO ₃	2.22	2.92
K ₂ O	0.83	0.91
Na O	0.38	0.58
TiO ₂	0.20	0.2
SrO	0.037	0.032
P ₂ O ₅	0.1	0.087
Mn ₂ O ₃	0.12	0.12
F	(0.05)	—
Cl	(0.02)	0.027
ZnO	(0.03)	0.025


Cement NIST 1885.

Element	Certificate	UniQuant
Si	19.81	20.06
Al	20.76	20.72
Ti	1.43	1.36
Fe	0.831	0.763
K	0.513	0.55
Mg	0.113	0.105
Ca	0.0249	0.049
Zr	(0.05)	0.048
Na	0.0492	0.046
P	(0.02)	0.025
V	N.A.	0.027
S	N.A.	0.045
Cr	0.0227	0.021
Sr	0.0084	0.0076
Zn	(0.0087)	0.0071
Mn	0.0047	0.0039

Manual input: Lol = 13.3%
N.A. = not available

Flint clay NIST 97b.

Irregular shaped small samples




Polypropylene insert to locate sample

UniQuant results on small and irregular samples

3 drillings: 25.8 mg BCS 380 aluminium alloy			3 drillings: 100.3 mg BCS 364 leaded bronze		
Element	Certificate %	UniQuant %	Element	Certificate %	UniQuant %
Al	94.61	93.3	Cu	78.6	80.6
Si	2	2.64	Sn	9.3	7.8
Fe	1.15	1.18	Pb	9.2	11.4
Ni	0.91	0.94	Ni	0.28	0.33
Cu	0.9	0.83	Sb	0.18	0.16
Ti	0.22	0.28	Zn	0.13	0.118
Mg	0.18	0.08			
ZnO	(0.03)	0.025			

Other sample types



Loose powders or drilling Liquids

UniQuant typical results

Drillings			Lube oil		
Element	Certificate %	UniQuant %	Element	Certificate %	UniQuant %
Si	0.2	0.13	Ca	0.12	0.15
V	0.94	0.86	Mg	0.10	0.13
Cr	5.17	5.5	Zn	0.142	0.14
Mn	0.21	0.23	P	0.125	0.13
Fe	(59.7)	60.4	S	0.63	0.65
Co	10.2	10.3			
Mo	1.07	1.3			
W	22.4	20.2			

MThe rest 98.74 (C + H + N) is calculated by UniQuant.
UniQuant is used when no specific standards are available.

Liquids, loose powders and small particles are usually measured in helium atmosphere with a spectrometer with X-ray tube located below the sample, i.e., ARL PERFORM'X and ARL OPTIM'X Spectrometers.

For more information about UniQuant Software please visit thermofisher.com/uniquant