remel™

EN

BactiCard™ E. coli

REF R21108

25

1. INTENDED USE

BactiCard™ E. coli is a rapid, qualitative, enzymatic test for the presumptive identification of *Escherichia coli* grown on plated media by the detection of tryptophanase and β -glucuronidase. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

2. SUMMARY AND EXPLANATION

E. coli is the most commonly recovered bacterial species in the clinical laboratory with disease processes ranging from urinary tract infections to gram-negative sepsis.¹ In 1981, Godsey et al. found 4-methylumbelliferyl- β -D-glucuronide (MUG) a sensitive and selective substrate for detection of β -glucuronidase activity which is present in most *E. coli* (97%).² Trepeta and Edberg used the MUG test in conjunction with oxidase, indole, and lactose fermentation as a method for rapid, cost-effective identification of *E. coli*.³

3. PRINCIPLE

BactiCard E. coli consists of two separate test circles for MUG and indole (IND). In the presence of the enzyme β -glucuronidase, 4-methyl-umbelliferone is released from the substrate, 4-methylumbelliferyl- β -D-glucuronide. This compound (4-methylumbelliferone) is highly fluorescent and can be detected by longwave ultraviolet light. In the indole test, tryptophan and *p*-dimethylaminobenzaldehyde serve as substrates for detection of the enzyme, tryptophanase. This enzyme attacks the tryptophan molecule producing indole which reacts with *p*-dimethylaminobenzaldehyde to produce a red color.

4. REAGENTS AND MATERIALS SUPPLIED

BactiCard E. coli Cards

BactiCard E. coli Test Cards (25 – 5 sets of 5)
Reactive Ingredients:
4-Methylumbelliferyl- β -D-glucuronide
Tryptophan and *p*-Dimethylaminobenzaldehyde

BactiCard E. coli Rehydrating Fluid

Rehydrating Fluid (1 bottle, 6.0 ml) with preservative

BactiCard E. coli MUG Reagent

MUG Reagent (1 bottle, 6.0 ml) 4% Sodium Hydroxide

Applicator Sticks

Applicator Sticks (1 vial)

Instructions for Use (IFU)

5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully.

Caution!

- MUG Reagent may cause irritation to skin, eyes, and respiratory system.
- Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Contents

Sodium hydroxide 1310-73-2

Oxalic acid dihydrate 6153-56-6

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-352-323-3500 (Call Collect)

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

DANGER

US & EU

US & EU

H290	May be corrosive to metals.
H302	Harmful if swallowed.
H312	Harmful in contact with skin.
H315	Causes skin irritation.
H318	Causes serious eye damage.
P264	Wash face, hands and any exposed skin thoroughly after handling
P280	Wear protective gloves/protective clothing/eye protection/face protection
P260	Do not breathe dust/fume/gas/mist/vapors/spray
P234	Keep only in original container
P314	Get medical attention/advice if you feel unwell
P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P332+P313	If skin irritation occurs: Get medical advice/attention
P305+P351 +P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P360	Take off contaminated clothing and wash before reuse
P390	Absorb spillage to prevent material damage
P301+P312	IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell
P403+P233	Store in a well-ventilated place. Keep container tightly closed
P406	Store in corrosive resistant polypropylene container with a resistant inliner
P402	Store in a dry place
P310	Immediately call a POISON CENTER or doctor/physician.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P501	Dispose of contents/container to an approved waste disposal plant

6. STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 2-8°C until used. Do not freeze or overheat. Allow product to equilibrate to room temperature before use. Do not incubate prior to use. Once opened, maintain unused test cards at 2-8°C in the foil pouch provided. Seal to protect from moisture.

7. PRODUCT DETERIORATION

This product should not be used if (1) the color has changed, (2) the expiration date has passed, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{1,4}

9. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Oxidase reagent, (7) Gram stain reagents, (8) Microscope slides, (9) Microscope, and (10) Longwave ultraviolet light (Woods Lamp).

10. PROCEDURE

Test isolates acceptable for testing with BactiCard E. coli are oxidase-negative, gram-negative bacilli. Remove test isolates from pure, 18-72 hour cultures growing on blood agar.

- Label a Test Card with the specimen identification.
- Inoculate only one isolate per Test Card.

- Add one (1) drop of Rehydrating Fluid to the MUG test circle, only. Do not oversaturate the test area.
- Inoculate each test circle with a visible inoculum of the test isolate using an applicator stick (provided).
- Incubate at room temperature (RT) for 15 minutes.
- Observe the IND test circle for development of a red color, indicating a positive test.
- Add one (1) drop of MUG Reagent to the MUG test circle.
- Observe the MUG test circle in a darkened room with a longwave ultraviolet light for bright-blue fluorescence, indicating a positive test.

11. INTERPRETATION

Indole Test:

Positive Test - Color change to red

Negative Test - No color change

MUG Test:

Positive Test - Bright blue fluorescence

Negative Test - No fluorescence

12. EXPECTED VALUES

Organism	Indole	MUG	Oxidase	Lactose
<i>Escherichia coli</i>	+	+ ^a	–	+
<i>Shigella</i> Gp A, B, C	V	–	–	–
<i>Shigella sonnei</i>	–	–	–	–
<i>Klebsiella pneumoniae</i>	–	–	–	+
<i>Klebsiella oxytoca</i>	+	–	–	+
<i>Enterobacter</i>	–	–	–	+
<i>Serratia marcescens</i>	V	–	–	–
<i>Yersinia</i>	V	–	–	–
<i>Salmonella</i>	–	–	–	– ^b
<i>Citrobacter freundii</i>	–	–	–	V
<i>Citrobacter koseri</i>	+	–	–	V

+, most strains positive; –, most strains negative; V, variable

^a Most strains of *E. coli* O157:H7 are MUG-negative.⁵

^b Rare lactose-positive strains occur.⁴

13. QUALITY CONTROL

All lot numbers of BactiCard E. coli have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory procedures. If aberrant quality control results are noted, patient results should not be reported.

Quality control testing should be run with each shipment and new lot number received.

		RESULTS	
CONTROL	INCUBATION	IND	MUG
<i>Escherichia coli</i> ATCC® 25922	Ambient, 15 min. @ RT	+	+
<i>Klebsiella pneumoniae</i> ATCC® 27736	Ambient, 15 min. @ RT	–	–

14. PERFORMANCE CHARACTERISTICS⁷

Remel BactiCard E. coli was evaluated with 243 strains of gram-negative bacilli representing 15 genera and showed greater than 98% agreement with expected values.

15. LIMITATIONS

- This test is only part of the overall scheme for identification of *E. coli*. Additional testing may be required for definitive identification. Consult appropriate references for further instructions.⁴
- Serological testing may be required to differentiate certain closely related strains of *E. coli* and *Shigella*.¹
- Most strains of *E. coli* O157:H7 are MUG-negative. Enterohemorrhagic *E. coli* is associated with bloody diarrhea and abdominal cramps. Complications may occur in the very young and elderly, including hemolytic uremic syndrome (HUS) and thrombotic thrombocytopenic purpura. Latex agglutination tests

- may be used in conjunction with BactiCard E. coli for rapid and accurate detection of *E. coli* O157:H7.⁵
- False-negative MUG reactions may occur with inoculum removed from MacConkey Agar or Triple Sugar Iron Agar.⁵
 - Some strains of *Salmonella*, *Shigella*, *Staphylococcus*, and *Streptococcus* possess the enzyme β -glucuronidase and may produce a positive MUG test. A positive test for lactose fermentation can be used to rule out *Salmonella* and *Shigella* and a Gram stain will differentiate gram-positive cocci from *E. coli*.⁶

16. BIBLIOGRAPHY

- Winn, W., S. Allen, J. William, E. Koneman, G. Procop, P. Schreckenberger, and G. Woods. 2006. Koneman’s Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott Williams and Wilkins, Baltimore, MD.
- Godsey, J.H., M.R. Matteo, D. Shen, G. Tolman, and J.R. Gohlke. 1981. J. Clin. Microbiol. 13:483-490.
- Trepeta, R.W. and S.C. Edberg. 1984. J. Clin. Microbiol. 19:172-174.
- Jorgensen J.H., M.A. Pfaller et al. 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.
- Raj, P. 1993. Clin. Microbiol. Newsl. 15:89-93.
- Kilian, M. and P. Bulow. 1976. Acta. Pathol. Microbiol. Scand. 84:245-251.
- Data on file.

17. PACKAGING

REF	R21108	 25
------------	--------------	--

18. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Use By (Expiration Date)
	Importer
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Mark
	Manufacturer

©2022 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

 Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

www.thermofisher.com/microbiology
 Tel: (800) 255-6730 • International: (913) 888-0939

www.oxid.com/IFU
 Europe 00800 135 79 135 • US 1 855 2360 190
 CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU21108	March 2023 Addition of Importer symbol to Section 18

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©2022 Thermo Fisher Scientific Inc. Tutti i diritti riservati.
Tutti gli altri marchi sono di proprietà di Thermo Fisher Scientific Inc. e delle sue consociate.