remel

BactiDrop™ Dobell & O'Connor Iodine EN

REF R21510

INTENDED USE

BactiDrop Dobell & O'Connor Iodine is a stain used for the microscopic examination of fecal samples for the detection of intestinal protozoa. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having protozoal infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

SUMMARY AND EXPLANATION

This iodine solution was formulated by Dobell and O'Connor in 1921 for the staining of protozoan cysts in fecal specimens.¹ They found that a weak iodine solution did not coagulate fecal particles that could obscure or entrap organisms. It has also been reported that a weak iodine solution does not adversely affect the refractile nature of stained organisms, such as concentrated iodine does.²

PRINCIPLE

The iodine mount is useful for the examination of cysts and eggs, but trophozoites are killed. The chromatin material of amebic cysts stain dark brown against the yellow-brown cytoplasm, the nuclear structures are differentiated and the glycogen masses stain a mahogany brown.²

REAGENTS (CLASSICAL FORMULA)*

Potassium Iodide (CAS 7681-11-0)	20.0 g
Iodine (CAS 7553-56-2)	10.0 g
Demineralized Water (CAS 7732-18-5)	1000.0 ml

*Adjusted as required to meet performance standards.

PRECAUTIONS

CAUTION! May cause allergic skin reaction. May cause eye, skin, and respiratory tract irritation.

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet on company website for detailed information on chemical reagents.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device

STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat. Protect from light.

PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed.

SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{3,4}

MATERIALS PROVIDED

50 BactiDrop Dobell & O'Connor Iodine – glass ampule enclosed in a protective plastic sheath with a dropper dispenser tip (0.75/ampule)
1 Ampule crusher (reusable)

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Applicator sticks, plain and cotton-tipped, (4) Incubator, slide warmer, (5) Quality control organisms, (6) Glass slides, coverslips, (7) Microscope with a calibrated ocular micrometer, (8) Disposable glass or plastic pipettes.

PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion. Discard remaining portion of partially used ampule at end of workday.

Test Procedure:

Dispense stain as needed for procedures involving examination of intestinal protozoa.

QUALITY CONTROL

All lot numbers of BactiDrop Dobell & O'Connor Iodine have been tested using the following quality control organisms and have been found to be

acceptable. Testing with a known, fresh, reference specimen should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL

Giardia intestinalis
ATCC® 30888

RESULTS

Intact protozoa, mahogany brown in color



LIMITATIONS

1. Iodine preparations fade quickly if exposed to light.³
2. The karyosome of *Iodamoeba butschlii* does not stain with iodine.
3. The glycogen of older cysts may not stain with iodine.²
4. Specimens collected from patients receiving barium sulfate, mineral oil, bismuth, non-absorbable antidiarrheal preparations, antimarials, and certain antibiotics are unsuitable for examination.²
5. A minimum of three fecal specimens is recommended for routine examination for parasites prior to treatment.²
6. Protozoan organisms should not be identified on the basis of wet mount (with few exceptions). Permanent stained smears should be examined to confirm the identification.²

BIBLIOGRAPHY

1. Dobell, C. and F.W. O'Connor. 1921. Intestinal Protozoa of Man. William Wood, New York.
2. Clinical and Laboratory Standards Institute CLSI. 2005. Procedures for the Recovery and Identification of Parasites from the Intestinal Tract; Approved Guideline. 2nd ed. M28-A2. CLSI, Wayne, PA.
3. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C.
4. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.

Symbol Legend

	Catalog Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
	Use By (Expiration Date)
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Authorized European Representative
	Manufacturer
	UK Conformity Assessment
	European Conformity Assessment

Remel Inc.
12076 Santa Fe Trail Drive
Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
(800) 255-6730
International: (913) 888-0939

For technical information contact your local distributor.
©2022 Thermo Fisher Scientific Inc. All rights reserved
BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
CAS (Chemical Abstracts Service Registry No.)

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, EE. UU.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus filiales.