---

# remel

## BactiDrop™ Indole (Ehrlich's)

(English)

### INTENDED USE

Remel BactiDrop Indole (Ehrlich's) is a reagent recommended for use in qualitative procedures for the detection of indole production by microorganisms grown in indole broth, especially nonfermenters and anaerobes.

### SUMMARY AND EXPLANATION

Vracko and Sherris reported that the indole test could be performed using *p*-dimethylaminobenzaldehyde.<sup>1</sup> When performing the indole test, Ewing and Davis<sup>2</sup> recommend Ehrlich's reagent for use with nonfermenters such as *Flavobacterium*<sup>3</sup> (*Chryseobacterium*) and with anaerobes.<sup>4</sup>

### PRINCIPLE

Intracellular enzymes, collectively called "tryptophanase", mediate the production of indole by hydrolytic activity against the amino acid tryptophan. Indole combines with dimethylaminobenzaldehyde to form a red compound. The reaction occurs by a condensation process formed by an acid splitting of the protein. When using Ehrlich's reagent, small amounts of indole produced from the breakdown of tryptophan by tryptophanase must first be extracted with xylene before being detectable.<sup>5</sup> Optimal results are dependent on the medium in which the indole test is performed. Use of an inadequate peptone broth or glucose-containing medium may cause aberrant results.<sup>6</sup>

### REAGENTS (CLASSICAL FORMULA)\*

*p*-Dimethylaminobenzaldehyde (CAS 100-10-7) .8.68 g  
Hydrochloric Acid, Conc. (CAS 7647-01-0) ..... 174.0 ml  
Ethanol, 95% (CAS 64-17-5)..... 826.0 ml

\*Adjusted as required to meet performance standards.

### PRECAUTIONS

**WARNING!** Causes irritation and possible burns by all routes of exposure. Flammable liquid and vapor. This substance has caused adverse reproductive and fetal effects in humans. May cause central nervous system depression. May cause liver, kidney and heart damage. Repeated or prolonged exposure may cause erosion of exposed teeth. Corrosive to metal.

This product is for *In Vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully. Refer to Material Safety Data Sheet for additional information.

### STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat.

### PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed. Discard remaining portion of partially used ampule at end of workday.

### SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.<sup>7,8</sup>

### MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Xylene, (7) Indole broth.

### PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion.

### Test Procedure:

1. Inoculate a tube of indole broth with the test organism and incubate at 35-37°C for 18-24 hours.
2. After incubation, transfer 2 ml of broth to a separate tube.
3. Add 1 ml of xylene to the aliquot, shake the tube, and allow to stand so that the solvent rises to the surface.
4. Dispense 5 drops of BactiDrop Indole (Ehrlich's) down the side of the tube.
5. Observe for the development of a red ring just below the solvent layer.

### INTERPRETATION

Positive Test - Development of a red color just below the solvent layer

Negative Test - No color development

**NOTE:** If the 24 hour culture is negative, reincubate the remaining portion for an additional 24 hours and then retest.

### QUALITY CONTROL

All lot numbers of BactiDrop Indole (Ehrlich's) have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

**CONTROL***Bacteroides thetaiotaomicron*  
ATCC® 29741*Elizabethkingia meningoseptica*  
ATCC® 13253*Clostridium sordelii*  
ATCC® 9714*Prevotella melaninogenica*  
ATCC® 25845*Pseudomonas aeruginosa*  
ATCC® 27853*Veillonella parvula*  
ATCC® 10790**RESULTS**

Positive

Positive

Positive

Negative

Negative

Negative

**LIMITATIONS**

1. Peptone broths (other than tryptophan broth) should be qualified for use in the indole test by testing with a known positive indole-producing organism.<sup>6</sup>
2. Medium containing glucose should not be used for indole detection.<sup>6</sup>
3. Some organisms form indole but break it down as rapidly as it is produced. A false negative reaction may occur, especially among the *Clostridium*.<sup>6</sup>
4. Weakly false positive reactions may occur if a mixed culture of indole-positive and indole-negative organisms is used. Test only pure cultures.<sup>6</sup>
5. Ehrlich's reagent is not specific for indole. By using a solvent, such as xylene to first extract indole, the red color obtained in the solvent layer only occurs with indole and skatole.<sup>6</sup>
6. The optimal pH for tryptophanase activity is one that is slightly alkaline (pH 7.4 to 7.8). A decrease in pH (toward acidity) may result in a false negative or weakly false positive reaction.<sup>6</sup>
7. Cultures to be tested for indole production must be incubated aerobically. A decrease in oxygen tension decreases indole production.<sup>6</sup>

**BIBLIOGRAPHY**

1. Vracko, R. and J.C. Sherris. 1963. Am. J. Clin. Pathol. 39:429-432.
2. Ewing, W.H. and B.R. Davis. 1970. Media and Tests for the Differentiation of *Enterobacteriaceae*. CDC, Atlanta, GA.
3. King, E.O. 1967. The Identification of Unusual Pathogenic Gram-Negative Bacteria. CDC, Atlanta, GA.
4. Holdeman, L.V., Cato, E.P., and W.E.C. Moore. 1977. Anaerobe Laboratory Manual. Virginia Polytechnic Institute, Blacksburg, VA.
5. Baron, E.J. and S.M. Tenover. 1990. Bailey and Scott's Diagnostic Microbiology. 8<sup>th</sup> ed. The C.V. Mosby Co., St. Louis, MO.

6. McFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3<sup>rd</sup> ed. Lippincott Williams & Wilkins, Philadelphia, PA.
7. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10<sup>th</sup> ed. ASM Press, Washington, D.C.
8. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12<sup>th</sup> ed. Mosby Elsevier, St. Louis, MO.

**PACKAGING**BactiDrop Indole (Ehrlich's) (0.75 ml/Ampule):  
REF 21520 ..... 50 Ampules/Pk**Symbol Legend**

<b>REF</b>	Catalog Number
<b>IVD</b>	In Vitro Diagnostic Medical Device
<b>LAB</b>	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
<b>LOT</b>	Batch Code (Lot Number)
	Use By (Expiration Date)
<b>EC REP</b>	Authorized European Representative
	Manufacturer

	Remel Inc. 12076 Santa Fe Trail Drive Lenexa, KS 66215, USA www.thermofisher.com (800) 255-6730 International: (913) 888-0939
---	--

For technical information contact your local distributor.

BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries.

ATCC is a registered trademark of American Type Culture Collection.  
CAS (Chemical Abstracts Service Registry No.)

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