remel Rapid™ NH System

REF R8311001 20

1. INTENDED USE

The Rapid™ NH System is a qualitative micromethod using enzyme reactions to identify clinical isolates grown on agar of *Neisseria* species, *Haemophilus* species, *Moraxella* species and related microorganisms. The device is used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

A complete listing of the organisms addressed by Rapid NH System is provided in the Rapid NH Differential Chart.

2. SUMMARY AND EXPLANATION

Organisms belonging to the family Neisseriaceae are characterized as Gram-negative cocci, occurring in pairs or masses, or Gram-negative plump rods (often coccobacillary) in pairs or short chains. There are four genera within the family: *Neisseria*, *Moraxella*, *Acinetobacter*, and *Kingella*.¹ The natural habitat of these organisms is mucous membranes and only two species, *N. gonorrhoeae* and *N. meningitidis*, are considered to be primary pathogens.² Most other Neisseriaceae isolated from human infections have been classified as opportunistic pathogens. Because of this distinction among Neisseriaceae relative to human infection, the primary interest of the clinical laboratory has been the identification and confirmation of gonococcal and meningococcal isolates and the differentiation of these species from other Neisseriaceae.

Rapid NH System has been designed to definitively identify *N. gonorrhoeae*, *N. meningitidis*, and *Moraxella catarrhalis* and to differentiate these organisms from other species of *Neisseria*, *Moraxella*, and *Kingella*.^{2,7}

Species of the genus *Haemophilus* are obligate parasites that are associated with the respiratory tract of man and animals. *Haemophilus influenzae* is the etiologic agent of a variety of human infections, including chronic respiratory infection and meningitis. Other species are implicated in venereal disease and conjunctivitis. The differentiation of pathogenic *Haemophilus* from *Haemophilus* species which constitute normal flora is important laboratory information. The Rapid NH System will identify and differentiate *Haemophilus* spp., as well as biochemically type *H. influenzae* and *Haemophilus parainfluenzae*.^{2,8}

Rapid NH Panels are disposable plastic trays with 10 reaction cavities, which contain dehydrated reactants. The Panel allows simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in Rapid Inoculation Fluid is used as the inoculum which rehydrates and initiates test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison with the probability values in the Differential Chart (Table 4), or by use of the Rapid ERIC™ software.

3. PRINCIPLE

The tests used in the Rapid NH System are based upon the microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional tests and single-substrate chromogenic tests, described in Table 1.

4. REAGENTS

Rapid Inoculation Fluid	(R8325102, supplied separately)
(1 ml/Tube)	
KCl	6.0 g
CaCl ₂	0.5 g
Deminerlized Water	1000.0 ml
Rapid Nitrate A Reagent	(R8309003, supplied separately)
(15 ml/Bottle)	
Sulfanilic Acid	8.0 g
Glacial Acetic Acid	280.0 ml
Deminerlized Water	720.0 ml
Rapid Nitrate B Reagent	(R8309004, supplied separately)
(15 ml/Bottle)	
N,N-dimethyl-1-naphthylamine	6.0 g
Glacial Acetic Acid	280.0 ml
Deminerlized Water	720.0 ml

Table 1. Principles and Components of the Rapid NH System

Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography#
Before Reagent Addition:					
1	PRO	Proline <i>p</i> -nitroanilide	0.1%	Hydrolysis of the colorless amide substrate by specific enzymes releases yellow <i>p</i> -nitrophenol.	1-3, 7-10
2	GGT	γ -Glutamyl <i>p</i> -nitroanilide	0.12%		
3	ONPG	<i>o</i> -Nitrophenyl, β , D-galactoside	0.25%	Hydrolysis of the colorless glycoside substrate releases yellow <i>o</i> -nitrophenol.	1, 11
4	GLU	Glucose	2.0%	Utilization of the sugar substrate produces acid products which lower the pH and change the indicator.	1, 11
5	SUC	Sucrose	2.0%		
6	EST	Fatty Acid Ester	0.5%	Hydrolysis of the fatty acid ester produces acid products which lower the pH and change the indicator.	1
7	RES	Resazurin	0.1%	Hydrolysis of resazurin to resorufin results in a color change.	8
8	PO ₄	<i>p</i> -Nitrophenyl phosphate	0.1%	Hydrolysis of the colorless phosphoester releases yellow <i>p</i> -nitrophenol.	12
9	ORN	Ornithine	0.8%	Hydrolysis of ornithine produces basic products which raise the pH and change the indicator.	4, 6, 13
10	URE	Urea	0.36%	Hydrolysis of urea produces basic products which raise the pH and change the indicator.	6, 13
After Reagent Addition:					
8	NO ₂	Nitrite	1.2%	Reduction of nitrite to nitrogenous products is detected by the absence of the ability to diazotize nitrate reagents.	1, 2, 6
9	NO ₃	Nitrate	0.3%	Reduction of nitrate to nitrite is detected by the ability to diazotize nitrate reagents.	6, 13, 14
10	IND	Tryptophane	0.16%	Utilization of tryptophane results in the formation of indole which is detected with Rapid Spot Indole Reagent.	6, 13, 14

Rapid Spot Indole Reagent (R8309002, supplied separately)
(15 ml/Bottle)
p-Dimethylaminocinnamaldehyde 10.0 g
Hydrochloric Acid 100.0 ml
Deminerlized Water 900.0 ml

5. PRECAUTIONS AND WARNINGS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after use. Directions should be read and followed carefully.

Non-disposable apparatus should be sterilised by any appropriate procedure after use, although the preferred method is to autoclave for 15 minutes at 121°C; disposables should be autoclaved or incinerated. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated area swabbed with a standard bacterial disinfectant or 70% alcohol. Do NOT use sodium hypochlorite. Materials used to clean spills, including gloves, should be disposed of as biohazardous waste.

Do not use reagents beyond the printed expiration dates.

Do not use if there is any evidence of contamination or other signs of deterioration.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. In the event of malfunction do not use device.

Caution!

- Rapid Nitrate A Reagent, Rapid Nitrate B Reagent, and Rapid Spot Indole Reagent may cause irritation to skin, eyes, and respiratory system.
- Refer to Safety Data Sheet, available on company website, and product labeling for information on potentially hazardous components, for detailed information on reagent chemicals.

6. STORAGE

Rapid NH System, Spot Indole, and Nitrate A and B Reagents should be stored in their original containers at 2-8°C until use. Allow products to equilibrate to room temperature before use. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. Rapid Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

7. PRODUCT DETERIORATION

This product should not be used if (1) the expiration date has passed, (2) the plastic tray is broken or the lid is compromised, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{2,16,17}

9. MATERIALS SUPPLIED

20 Rapid NH Panels
20 report forms
2 chipboard incubation trays
Instructions for Use
1 color guide

10. MATERIALS REQUIRED BUT NOT SUPPLIED

Loop sterilization device
Inoculating loop, swabs, collection containers
Incubators, alternative environmental systems
Supplemental media
Quality control organisms
Gram stain reagents
Microscope slides
Oxidase reagent
Cotton swabs
Rapid Inoculation Fluid, 1 ml (R8325102)
McFarland #3 turbidity standard (R20413) or equivalent
Pipettes
Rapid Spot Indole Reagent (R8309002)
Rapid Nitrate A Reagent (R8309003)
Rapid Nitrate B Reagent (R8309004)
ERIC (Electronic Rapid Compendium, R8323600) (optional).

11. PROCEDURE

There are two alternative procedures for the Rapid™ NH System: the 1 hour procedure and the general procedure.

The 1-hour procedure is only applicable to suspected gonococci


obtained from urogenital specimens isolated on selective agars. The general procedure should be used for Neisseriaceae from all other body sites and isolated on all other media. *Haemophilus* and other bacteria should be tested using the general procedure.

Inoculum Preparation:

- Test organisms must be grown in pure culture and examined by Gram stain and oxidase test prior to use in the system.
Note: The cellular morphology and Gram stain characteristics should be carefully observed since coccobacillary rods may resemble diplococci in smears.
- Test organisms may be removed from a variety of nonselective and selective agar growth media. The following types of media are recommended:
Nonselective Media: Chocolate Agar; Tryptic Soy Agar with 5% Sheep Blood.
Selective Media: Thayer-Martin Agar; New York City Agar.
Notes:
 - When using the 1-hour procedure, only selective agars can be used.
 - Cultures used for inoculum preparation should preferably be 18-24 hours old. Slow-growing isolates may be tested using 48-hour cultures.
 - The use of media other than those recommended may compromise test performance.
- Using a cotton swab or inoculating loop, suspend sufficient growth from the agar plate culture in Rapid Inoculation Fluid (1 ml) to achieve a visual turbidity approximately equal to a McFarland turbidity standard number 3 or equivalent.
Notes:
 - Suspensions significantly less turbid than a McFarland standard number 3 will result in aberrant reactions.
 - Bacterial suspensions that are slightly more turbid than a McFarland standard number 3 will not affect test performance and are recommended for stock cultures, quality control strains, and the 1-hour procedure.
 - Suspensions should be mixed thoroughly and vortexed if required.
 - Suspensions should be used within 15 minutes of preparation.
- An agar plate may be inoculated for purity and any additional testing that may be required using a loopful of the test suspension from the inoculation fluid tube. Incubate the plate for at least 18-24 hours at 35-37°C.

Inoculation of Rapid NH Panels:

- Peel back the lid of the panel over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
- Using a pipette, gently transfer the entire contents of the Inoculation Fluid tube into the upper right-hand corner of the panel. Reseal the inoculation port of the panel by pressing the peel-back tab back in place.
- After adding the test suspension, and while keeping the panel on a level surface, tilt the panel back away from the test cavities at approximately a 45° (see below).

- While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.
- While maintaining a level, horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction

Table 2. Interpretation of Rapid NH Panel Tests*


Cavity #	Test Code	Reagent	Reaction		Comments
			Positive	Negative	
Before Reagent Addition:					
1	PRO	None	Yellow	Clear or tan	Any development of a distinct yellow color should be scored as positive.
2	GGT				
3	ONPG				
4	GLU				
5	SUC	None	Yellow, gold, or yellow-orange	Red, red-orange, or orange	Only a distinct yellow, gold, or yellow-orange should be scored as positive. All other colors should be scored as negative.
6	EST	None	Yellow, gold, or yellow-orange	Red, red-orange, or orange	Note: A red layer may form on the top of the cavity. Gently shake the panel or stir with an applicator stick before scoring.
7	RES	None	Pink	Purple, blue, or violet	Only the development of a significant pink color should be scored as positive. All other colors should be scored as negative.
8	PO ₄	None	Yellow	Clear, tan, straw, or very pale yellow	Only the development of a significant yellow color throughout the cavity should be scored as positive.
9	ORN	None	Red, violet, or purple	Yellow or orange	
10	URE				
After Reagent Addition:					
9	NO ₃	Rapid Nitrate A Rapid Nitrate B	Red or orange	Yellow	Any development of a red or orange color should be scored as positive
10	IND	Rapid Spot Indole	Brown or black	Orange or red	Any development of a brown or black color should be scored as positive.
8**	NO ₂	Rapid Nitrate A Rapid Nitrate B	Clear, tan, or straw	Pink or red	Any development of red or pink color should be scored as negative.

*NOTE: Panels should be read by looking down through the reaction wells against a white background.

**NO₂ Test: Perform the NO₂ test when the PRO test (cavity 1) is the only positive test and the test isolate is a Gram negative coccus (suspect *Neisseria* sp.). Negative test color development may be slow. Allow a full 5 minutes for color development before reading and recording NO₂ reactions.

cavities until the inoculum flows along the baffles into the reaction cavities (see below). This should evacuate all of the inoculum from the rear portion of the panel.

Note: If the panel is tilted too quickly, air may be trapped at the test cavity junction, restricting fluid movement.

- Return the panel to a level position. If necessary, gently tap the panel on the bench top to remove any air trapped in the cavities.

Notes:

- Examine the test cavities which should appear bubble-free and uniformly filled. Slight irregularities in test cavity fills are acceptable and will not affect test performance. If the panel is grossly misfilled, a new panel should be inoculated and the misfilled panel discarded.
- Complete the inoculation of each panel receiving inoculation fluid before inoculating additional panels.
- Do not allow the inoculum to rest in the back portion of the panel for prolonged periods without completing the procedure.

Incubation of Rapid NH Panels:

When using the 1 hour procedure, incubate inoculated panels at 35-37°C in a non-CO₂ incubator for 1 hour. When using the general procedure, incubate inoculated panels at 35-37°C in a non-CO₂ incubator for 4 hours. For ease of handling, panels may be incubated in the chipboard incubation trays provided with the kit.

Scoring of Rapid NH Panels:

Rapid NH panels contain 10 reaction cavities that provide 12 test scores and, if required, a 13th test score (NO₂). Test cavities 8 through 10 are bifunctional, containing two separate tests in the same cavity. Bifunctional tests are first scored before the addition of reagent providing the first test result, and then the same cavity is scored again after the addition of reagent to provide the second test result. Bifunctional test cavities are indicated with the first test above the bar and the second test below the bar. The Nitrite test (cavity 8), which is only required as indicated below in step 5, is designated with a box drawn around the reagent-requiring test.

Rapid NH Panel Test Location

Cavity #	1	2	3	4	5	6	7	8	9	10
Test code	PRO	GGT	ONPG	GLU	SUC	EST	RES	PO ₄	ORN	URE
								NO ₂	NO ₂	IND

- While firmly holding the Rapid NH panel on the benchtop, peel off the label lid over the reaction cavities by pulling the lower right hand tab up and to the left.
- Without the addition of any reagents, read and score cavities 1 (PRO) through 10 (URE) from left to right using the interpretation guide presented in Table 2. Panels should be read by looking down through the reaction wells against a white background. Record test scores in the appropriate boxes on the report form using the test code above the bar for bifunctional tests.
- Add the following reagents to the cavities indicated:
 - Add 2 drops of Rapid Nitrate A Reagent to cavity 9 (NO₃).
 - Add 2 drops of Rapid Nitrate B Reagent to cavity 9 (NO₂).
 - Add 2 drops of Rapid Spot Indole Reagent to cavity 10 (IND).
- While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.
- While maintaining a level, horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction

Note: Negative test color development may be slow. Allow a full five minutes before scoring as positive.

Table 3. Quality Control Chart for RapID NH Panels

Organism	PRO	GGT	ONPG	GLU	SUC	EST	RES	PO4	ORN	URE	NO3	IND
<i>Haemophilus influenzae</i> Biotype I ^a ATCC® 9006	–	–	–	+	–	–	–	+	+	+	+	+
<i>Aggregatibacter aphrophilus</i> ATCC® 7901	–	–	+	+	+	–	V	+	–	–	V	–
<i>Aggregatibacter aphrophilus</i> ATCC® 49146	V	+	+	+	+	–	V	+	–	–	+	–
<i>Oligella urethralis</i> ATCC® 17960	V	+	–	–	–	–	+	–	–	–	–	–
<i>Moraxella catarrhalis</i> ATCC® 8176	+	–	–	–	–	+	V	–	V	–	V	–

+, positive; –, negative; V, variable; (–), usually negative; (+), usually positive

^a Key indicator strains demonstrate acceptable performance of the most labile substrate in the system and reactivity in a significant number of wells, according to Clinical and Laboratory Standards Institute recommendations for streamlined quality control.²⁶

6. Reference the microcode obtained on the report form in ERIC for the identification.

12. RESULTS AND RANGE OF EXPECTED VALUES

The RapID NH Differential Chart (Table 4) and the *Haemophilus* Biotype Chart (Table 5) illustrate the expected results for RapID NH System. Differential Chart results are expressed as a series of positive percentages for each system test. This information statistically supports the use of each test and provides the basis, through numerical coding of digital test results, for a probabilistic approach to the identification of the test isolate.

Identifications are made using individual test scores from RapID NH panels in conjunction with other laboratory information (i.e., Gram stain, oxidase, growth on differential or selective media) to produce a pattern that statistically resembles known reactivity for taxa recorded in the RapID NH System database. These patterns are compared through the use of RapID NH Differential Chart (Table 4), or by derivation of a microcode and the use of ERIC.

13. QUALITY CONTROL

All lot numbers of RapID NH System have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported. Table 3 lists expected results for the selected battery of test organisms.

Notes:

- The quality control of RapID Reagents is accomplished by obtaining the expected reactions for tests requiring the addition of the reagents (cavities 8-10).
- Organisms which have been repeatedly transferred on agar media for prolonged periods may produce aberrant results.
- Quality control strains should be stored frozen or lyophilized. Prior to use, quality control strains should be transferred 2-3 times from storage on an agar medium that is recommended for use with RapID NH System.
- Formulations, additives, and ingredients of culture media vary from manufacturer to manufacturer and may vary from batch to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve quality control discrepancies.

14. LIMITATIONS

- The use of RapID NH System and the interpretation of results requires the knowledge of a competent laboratorian who is trained in general microbiological methods and who judiciously makes use of training, experience, specimen information, and other pertinent procedures before reporting the identification obtained using RapID NH System.
- Specimen source, oxidase reaction, Gram stain characteristics, and growth on selective agars should be considered when using RapID NH System.
- RapID NH System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.

Table 4 - RapID NH Differential Chart (see Section 12)

Organism	PRO	GGT	ONPG	GLU	SUC	EST	RES	PO ₄	ORN	URE	NO ₃	IND
<i>Aggregatibacter actinomycetemcomitans</i> ^a	46	98	1	95	16	59	39	58	1	0	79	0
<i>Aggregatibacter aphrophilus</i> ^b	2	92	96	98	98	12	90	99	0	0	90	0
<i>Aggregatibacter segnis</i> ^c	0	0	90	52	33	0	24	98	0	0	96	0
<i>Branhamella catarrhalis</i> ¹	39	0	0	0	0	99	92	0	8	0	63	0
<i>Cardiobacterium hominis</i>	12	96	0	79	63	0	93	0	0	0	0	99
<i>Eikenella corrodens</i>	86	0	9	0	0	0	23	0	97	0	93	0
<i>Gardnerella vaginalis</i>	99	0	77	91	44	1	12	5	0	0	0	0
<i>Haemophilus ducreyi</i>	0	0	0	8	0	0	0	79	0	0	96	0
<i>Haemophilus haemolyticus</i>	0	0	0	97	0	0	13	99	0	99	96	90
<i>Haemophilus influenzae</i> ^c	4	0	0	98	0	2	9	99	50	50	96	50
<i>Haemophilus parahaemolyticus</i>	0	0	88	98	95	0	12	99	34	99	99	0
<i>Haemophilus parainfluenzae</i>	2	0	93	98	96	0	12	99	50	50	96	50
<i>Kingella denitrificans</i>	93	0	0	91	0	5	91	0	0	0	94	0
<i>Kingella kingae</i>	19	0	0	76	0	0	57	98	0	0	0	0
<i>Moraxella atlantae</i>	0	22	0	0	0	8	90	96	0	0	0	0
<i>Moraxella lacunata</i>	0	0	0	0	0	0	96	5	0	0	88	0
<i>Moraxella nonliquefaciens</i>	8	0	0	0	0	0	94	0	0	0	96	0
<i>Moraxella osloensis</i>	4	39	0	0	0	0	93	27	0	0	9	0
<i>Neisseria cinerea</i>	99	0	0	0	0	0	68	0	0	0	0	0
<i>Neisseria elongata subsp. nitroreducens</i> ^d	79	0	0	0	0	0	91	0	0	0	89	0
<i>Neisseria flavescens</i>	92	9	0	0	0	0	99	96	0	0	0	0
<i>Neisseria gonorrhoeae</i>	98h	0	0	87	0	1	4	0	0	0	0	0
<i>Neisseria lactamica</i>	96	0	99	97	2	5	0	18	0	0	0	0
<i>Neisseria meningitidis</i>	93	98i	0	96	0	0	62	8	2	0	0	0
<i>Neisseria mucosa</i>	96	2	0	95	92	0	92	1	0	0	95	0
<i>Neisseria sicca/subflava</i>	97	14	0	95	96	2	96	9	0	0	0	0
<i>Neisseria weaveri/elongata</i> ^a	99	0	0	0	0	0	96	0	0	0	0	0
<i>Oligella ureolytica</i>	12	99	0	0	0	0	98	9	0	99	0	0
<i>Oligella urethralis</i>	39	99	0	0	0	0	99	0	0	0	0	0
<i>Pasteurella multocida</i>	0	0	0	97	96	0	46	99	94	0	98	99
<i>Psychrobacter phenylpyruvicus</i> ^f	0	7	0	0	0	18	95	5	31	95	62	0
<i>Suttonella indologenes</i> ^g	7	92	0	96	84	31	98	90	0	0	0	99

^aPreviously designated *Haemophilus actinomycetemcomitans*.

^bPreviously designated *Haemophilus aphrophilus*.

¹Includes biogroup *aegyptius*.

^dPreviously designated CDC group M-6.

^ePreviously designated CDC group M-5 (*Neisseria weaveri*).

^fPreviously designated *Moraxella phenylpyruvica*.

^gPreviously designated *Kingella indologenes*.

^hPRO-negative strains of *Neisseria gonorrhoeae* have been reported.¹⁸

ⁱGGT-negative strains of *Neisseria meningitidis* have been reported.²⁸

^jPreviously designated *Haemophilus segnis*

^kPreviously designated *Moraxella catarrhalis*

- Westley, J.R., P.J. Anderson, V.A. Close, B. Halpern, and E.M. Lederberg. 1967. Appl. Microbiol. 15:822-825.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.
- Saginur, R., B. Clecner, J. Portnoy, and J. Mendelson. 1982. J. Clin. Microbiol. 15:475-477.
- Takahashi, H., H. Tanaka, H. Inouyr, T. Kuroki, Y. Watanabe, S. Yamai, and H. Watanabe. 2002. J. Clin. Microbiol. 40:3035-3037.
- Koneman, E.W., S.D. Allen, W.M. Janda, P.C. Schreckenberger, and W.C. Winn. 1997. Color Atlas and Textbook of Diagnostic Microbiology. 5th ed. Lippincott Williams & Wilkins, Philadelphia, PA.

17. PACKAGING

REF R8311001 RapID NH System..... 20 Tests/Kit

18. SYMBOL LEGEND

REF	Catalogue Number
IVD	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
UDI	Unique Device Identifier
EC REP	Authorized representative in the European Community
UK CA	UK Conformity Assessed
CE	European Conformity Assessment
	Manufacturer

RapID™ and ERIC™ are trademarks of Thermo Fisher Scientific and its subsidiaries.

ATCC® is a registered trademark of American Type Culture Collection.

Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

www.thermofisher.com/microbiology
Tel: (800) 255-6730 • International: (913) 888-0939

www.oxid.com/IFU
Europe +800 135 79 135 • US 1 855 2360 190
CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU8311001	August 2023 Updated to meet IVDR requirements

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