remel

RapID™ STR System

EN

REF R8311003..... 20

1. INTENDED USE

The RapID™ STR is a qualitative micromethod using enzyme reactions to identify clinical isolates of *Streptococcus* species and other related Gram-positive bacteria grown on agar. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

The RapID STR System is intended to aid in the identification of Lancefield groups A, B, C, D, and G streptococci, viridans streptococci, and *Streptococcus pneumoniae*, *Enterococcus* spp., *Aerococcus* spp., *Gemella* spp., *Leuconostoc* spp., *Pediococcus* spp., *Weisella confusa*, and *Listeria monocytogenes*.¹⁻¹⁰ A complete listing of the organisms addressed by the RapID STR System is provided in the RapID STR Differential Chart.

2. SUMMARY AND EXPLANATION

RapID STR System is comprised of (1) RapID STR Panels and (2) RapID STR Reagent. RapID STR Panels are disposable plastic trays with 10 reaction cavities, which contain dehydrated reactants. The panel allows simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in RapID Inoculation Fluid is used as the inoculum which rehydrates and initiates test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison with the probability values in the Differential Chart (Table 4), or by use of the RapID ERIC™ software.

3. PRINCIPLE

The tests used in RapID STR System are based on microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional tests and single-substrate chromogenic tests, described in Table 1.

4. REAGENTS

RapID STR Reagent (provided with kit) (15 ml/Bottle)

Reactive ingredient per liter:

o-Dianisidine..... 0.9 ml

RapID Inoculation Fluid

(R8325102, supplied separately) (1 ml/Tube)

KCl 6.0 g

CaCl₂ 0.5 g

Demineralized Water..... 1000.0 ml

5. PRECAUTIONS AND WARNINGS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after use. Directions should be read and followed carefully.

Non-disposable apparatus should be sterilised by any appropriate procedure after use, although the preferred method is to autoclave for 15 minutes at 121°C; disposables should be autoclaved or incinerated. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated area swabbed with a standard bacterial disinfectant or 70% alcohol. Do NOT use sodium hypochlorite. Materials used to clean spills, including gloves, should be disposed of as biohazardous waste.

Do not use reagents beyond the printed expiration dates.

Do not use if there is any evidence of contamination or other signs of deterioration.

DANGER

US & EU

US & EU

US ONLY

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. In the event of malfunction do not use device.

Caution!

1. RapID STR Reagent is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May cause cancer, impair fertility, or cause harm to unborn child. Danger of serious irreversible effects.
2. Refer to Safety Data Sheet, available on company website, and product labeling for information on potentially hazardous components, for detailed information on reagent chemicals.

Composition / information on ingredients

Methyl alcohol 67-56-1

Acetic acid 64-19-7

2-Methoxyethanol 109-86-4

[1,1'-Biphenyl]-4,4'-bis(diazonium), 3,3'-dimethoxy-, (T-4)-tetrachlorozincate(2-) (1:1) Fast Blue B salt 14263-94-6

6. HAZARDS NOT OTHERWISE CLASSIFIED (HNOC)

Poison, may be fatal or cause blindness if swallowed. Vapor harmful. Cannot be made non-poisonous. WARNING! This product contains a chemical known in the State of California to cause birth defects or other reproductive harm.

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-352-323-3500 (Call Collect)

7. STORAGE

The RapID STR System should be stored in its original container at 2-8°C until use. Allow product to equilibrate to room temperature before use. DO NOT interchange reagents among different RapID systems. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. RapID Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

8. PRODUCT DETERIORATION

This product should not be used if (1) the expiration date has passed, (2) the plastic tray is broken or the lid is compromised, or (3) there are other signs of deterioration.

9. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{20,21}

10. MATERIALS SUPPLIED

- 20 RapID STR Panels
- 20 Report forms
- RapID STR Reagent (one plastic dropper-bottle containing reagent sufficient for 20 panels)
- 2 chipboard incubation trays
- Instructions for Use (IFU)
- 1 color guide

11. MATERIALS REQUIRED BUT NOT SUPPLIED

- Loop sterilization device
- Inoculating loop, swabs, collection containers
- Incubators, alternative environmental systems
- Supplemental media
- Quality control organisms
- Gram stain reagents
- Microscope slides
- Cotton swabs
- RapID Inoculation Fluid-1 ml (R8325102)
- McFarland #1 turbidity standard or equivalent (R20411)
- Pipettes
- ERIC (Electronic RapID Compendium, R8323600) (Optional).

12. CONTENTS SYMBOLS

STR Panels	STR Panels
Report Forms	RapID Report Forms
STR Reagent	STR Reagent
Incubation Trays	Incubation Trays

13. PROCEDURE

Inoculum Preparation:

1. Test organisms must be grown in pure culture, examined by Gram stain, and tested for hemolysis prior to use in this system.

Note: Hemolysis is enhanced by incubation anaerobically or in 5-7% CO₂.

2. Test organisms may be removed from nonselective agar growth media. The following types of media are recommended:

Tryptic Soy Agar (TSA) with or without 5% Sheep Blood; Chocolate Agar.

Notes:

- Some media containing or supplemented with mono- or disaccharides are not recommended since they may suppress glycolytic activity and reduce test selectivity.
- Plates used for inoculum preparation should preferably be 18-24 hours old. Slow growing isolates may be tested using 48-hour plates.
- The use of media other than those recommended may compromise test performance.
- 3. Using a cotton swab or inoculating loop, suspend sufficient growth from the agar plate in RapID Inoculation Fluid (1 ml) to achieve a visual turbidity equal to a #1 McFarland turbidity standard or equivalent.

Notes:

- Suspensions significantly less turbid than a #1 McFarland standard may result in aberrant reactions.
- Bacterial suspensions that are slightly more turbid than a #1 McFarland standard will not affect test performance and are recommended for stock cultures and quality control strains. However, bacterial suspensions significantly more turbid than a #1 McFarland standard may compromise test performance.

Table 1. Principles and Components of the RapID STR System

Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
Before Reagent Addition:					
1	ARG	L-arginine	2.0%	Hydrolysis of arginine releases basic products which raise the pH and change the indicator.	11-13
2	ESC	Esculetin	0.5%	Hydrolysis of glucoside releases esculetin which reacts with ferric ion forming a black compound.	12
3	MNL	Mannitol	1.5%	Utilization of the carbohydrate substrate produces acidic products which lower the pH and change the indicator.	1, 2, 11
4	SBL	Sorbitol	1.5%		
5	RAF	Raffinose	1.2%		
6	INU	Inulin	1.5%		
7	GAL	<i>p</i> -Nitrophenyl- β , β -galactoside	0.1%	Hydrolysis of colorless <i>p</i> -nitrophenyl-substituted glycoside or phosphoester releases yellow <i>p</i> -nitrophenol.	14-16
8	GLU	<i>p</i> -Nitrophenyl- β , β -glucoside	0.1%		
9	NAG	<i>p</i> -Nitrophenyl-n-acetyl- β , β -glucosaminide	0.1%		
10	PO ₄	<i>p</i> -Nitrophenyl phosphate	0.2%		

After Reagent Addition:


7	TYR	Tyrosine β -naphthylamide	0.05%	Hydrolysis of arylamide releases free β -naphthylamine which is detected with RapID STR Reagent.	15, 17-19
8	HPR	Hydroxyproline β -naphthylamide	0.08%		
9	LYS	Lysine β -naphthylamide	0.08%		
10	PYR	Pyrrolidine β -naphthylamide	0.1%		

- Suspensions should be mixed thoroughly and vortexed if required.
- Suspensions should be used within 15 minutes of preparation.


4. An agar plate may be inoculated for purity and any additional testing that may be required using a loopful of the test suspension from the inoculation fluid tube. Incubate the plate for 18-24 hours at 35-37°C.

Inoculation of RapID STR Panels:

1. Peel back the lid of the panel over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
2. Using a pipette, gently transfer the entire contents of the Inoculation Fluid tube into the upper right-hand corner of the panel. Reseal the inoculation port of the panel by pressing the peel-back tab back in place.
3. After adding the test suspension and while keeping the panel on a level surface, tilt the panel back away from the reaction cavities at approximately 45° (see below).

4. While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.

5. While maintaining a level, horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction cavities until the inoculum flows along the baffles into the reaction cavities (see below). This should evacuate all of the inoculum from the rear portion of the panel.

Note: If the panel is tilted too quickly, air may be trapped at the test cavity junction, restricting fluid movement.

6. Return the panel to a level position. If necessary, gently

Table 2. Interpretation of RapID STR System Tests* (see Color Guide for examples)

Cavity #	Test Code	Reagent	Reaction		Comments

Table 3. Quality Control Chart for RapID STR Panels

Organism	ARG	ESC	MNL	SBL	RAF	INU	GAL	GLU	NAG	PO ₄	TYR	HPR	LYS	PYR
<i>Enterococcus faecalis</i> ATCC® 29212	+	+	+	+	-	-	-	+	+	(-)	V	-	+	+
<i>Enterococcus hirae</i> ATCC 49479 (formerly <i>E. durans</i>)	+	+	-	-	-	-	+	-	+	-	V	+*	(-)	+
<i>Streptococcus gallolyticus</i> ^{a,b} ATCC® 9809	-	+	+	-	+	+	+	+	V	-	-	-	+	-
<i>Streptococcus pyogenes</i> ^a ATCC® 19615	+	-	-	-	-	-	-	(+)	-	+	+	-	+	+

+, positive; -, negative; V, variable; (-), usually negative; (+), usually positive

^a Key indicator strains demonstrate acceptable performance of the most labile substrate in the system and reactivity in a significant number of wells, according to Clinical and Laboratory Standards Institute recommendations for streamlined quality control.²⁶ ^b Previously *Streptococcus bovis*

*Note: *Enterococcus durans* may yield a very weak positive reaction in the HPR cavity. *Gemella morbillorum* ATCC® 27824 can also be used as a quality control strain for the HPR reaction. However, *E. durans* should still be used for quality control with the GLU and LYS cavities.

form in ERIC (Electronic RapID Compendium) for the identification.

14. RESULTS AND RANGE OF EXPECTED VALUES

The RapID STR Differential Chart (Table 4) illustrates the expected results for RapID STR System. Differential Chart results are expressed as a series of positive percentages for each system test. This information statistically supports the use of each test and provides the basis, through numerical coding of digital test results, for a probabilistic approach to the identification of the test isolate.

Identifications are made using individual test scores from RapID STR panels in conjunction with other laboratory information (i.e. Gram stain, hemolysis, colonial morphology, growth on differential or selective media) to produce a pattern that statistically resembles known reactivity for taxa recorded in the RapID STR System database. These patterns are compared through the use of the RapID STR Differential Chart, or by derivation of a microcode and the use of ERIC.

15. QUALITY CONTROL

All lot numbers of RapID STR System have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported. Expected results for selected quality control organisms are listed in Table 3.

Notes

- RapID STR Reagent quality control is accomplished by obtaining the expected reactions for tests requiring the addition of the reagent (cav. 7-10).
- Organisms which have been repeatedly transferred on agar media for prolonged periods may provide aberrant results.
- Quality control strains should be stored frozen or lyophilized. Prior to use, quality control strains should be transferred 2-3 times from storage on an agar medium that is recommended for use with RapID STR System.
- Formulations, additives, and ingredients of culture media vary from manufacturer to manufacturer and may vary from batch to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve quality control discrepancies.

16. LIMITATIONS

- The use of RapID STR System and interpretation of results requires the knowledge of a competent microbiologist, familiar with laboratory procedures, who is trained in general microbiological methods and who judiciously makes use of training, experience, specimen information, and other pertinent procedures before reporting the identification obtained using this system.
- Characteristics such as Gram stain reaction, hemolysis, and cellular morphology must be considered when using the RapID STR System.
- RapID STR System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.
- RapID STR System is designed for use with the taxa listed in the RapID STR Differential Chart. The use of organisms not specifically listed may lead to misidentifications.
- Expected values listed for RapID STR System tests may differ from conventional test results or previously reported information.
- The accuracy of RapID STR System is based on the statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single test found in the RapID STR System to establish the identification of a test isolate is subject to the error inherent in that test alone.

17. PERFORMANCE CHARACTERISTICS

The RapID STR System performance characteristics have been established by laboratory testing of reference and stock cultures at Remel and by clinical evaluations using fresh clinical and stock isolates.²²⁻²⁵

18. BIBLIOGRAPHY

- Lennette, E.H., A. Balows, W.J. Hausler, Jr., and J.P. Truant. 1980. Manual of Clinical Microbiology. 3rd ed. ASM, Washington, D.C.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- Poole, P.M. and G. Wilson. 1976. J. Clin. Microbiol. 29:740-745.
- Facklam, R.R. 1977. J. Clin. Microbiol. 5:184-201.
- Ruoff, K.L. and L.J. Kunz. 1982. J. Clin. Microbiol. 15:920-925.
- Berlutti, F., M.C. Thaller, S. Schippa, F. Pantanella, and R. Pompei. 1993. Int. J. Syst. Bacteriol. 43:63-68.
- Collins, M.D., D. Jones, J.A. Farlow, R. Klipper-Balz, and K.H. Schleifer. 1984. Int. J. Syst. Bacteriol. 34:220-223.

Table 4 - RapID STR Differential Chart (see section 14)

Organism	ARG	ESC	MNL	SBL	RAF	INU	GAL	GLU	NAG	PO ₄	TYR	HPR	LYS	PYR	HEM
β-hemolytic Streptococci:															
<i>Group A (<i>S. pyogenes</i>)</i>	99	5	17	1	1	0	9	96	72	99	90	1	99	99	97
<i>Group B (<i>S. agalactiae</i>)</i>	99	2	1	1	1	0	2	96	0	94	1	2	92	0	92
<i>Group C/G</i>	99	8	4	3	0	0	1	98	12	99	96	1	99	0	99
Enterococci:															
<i>E. avium</i>	0	96	97	93	0	0	33	1	1	0	12	56	9	99	0
<i>E. casseliflavus / mundtii</i>	48	99	98	81	94	66	95	4	96	9	2	2	90	99	0
<i>E. durans / hirae</i>	96	99	0	0	11	0	69	4	30	0	2	86	33	99	6
<i>E. faecalis</i>	97	97	95	90	1	0	1	99	95	8	17	9	90	99	7
<i>E. faecium</i>	96	97	96	0	22	0	76	1	92	1	1	36	90	99	0
<i>E. gallinarum</i>	99	99	95	0	95	94	98	5	90	8	9	12	89	99	0
<i>E. malodoratus</i>	0	99	99	94	99	0	7	0	1	1	1	14	6	99	0
<i>E. raffinosus</i>	0	96	94	91	99	0	90	0	3	6	0	26	2	99	0
Group D Streptococci:															
<i>S. bovis</i>	0	97	98	0	96	77	99	99	5	1	1	1	95	1	0
<i>S. bovis</i> var	0	98	0	0	96	0	98	99	1	1	48	1	98	0	0
<i>S. equinus</i>	0	99	88	0	6	0	0	9	85	55	0	12	81	0	0
Viridans Streptococci:															
<i>S. acidominimus</i>	98	0	0	0	0	0	92	0	0	12	15	12	98	0	0
<i>S. anginosus</i>	92	98	17	0	32	0	80	22	0	95	93	3	96	0	28
<i>S. constellatus</i>	97	89	0	4	2	0	2	88	0	93	84	1	98	0	46
<i>S. intermedius</i>	95	96	6	2	26	0	18	99	99	95	92	5	99	1	3
<i>S. mitis</i>	9	2	0	4	59	0	38	96	92	79	90	1	96	0	0
<i>S. mutans</i>	2	94	94	92	86	82	94	92	0	1	12	1	92	0	4
<i>S. salivarius / vestibularis</i>	0	91	0	0	90	40	0	96	1	10	11	1	98	0	0
<i>S. sanguis</i>	98	79	0	4	26	66	23	16	74	79	94	0	95	0	0
<i>S. sanguinis</i> ^a	17	0	0	1	96	0	86	99	39	96	89	1	96	0	0
Other:															
<i>Aerococcus</i> spp.	0	61	84	70	18	0	17	98	0	18	7	19	87	90	0
<i>Gemella</i> morbillorum	0	0	1	1	0	0	7	9	1	64	19	74	82	42	0
<i>Leuconostoc</i> citreum	0	99	29	0	2	0	2	99	0	0	0	0	29	0	0
<i>Leuconostoc</i> lactis	0	3	0	0	90	0	99	16	0	0	0	0	3	0	0
<i>Leuconostoc</i> mesenteroides group	0	88	18	0	86	0	96	80	4	0	2	0	4	0	0
<i>Listeria</i> monocytogenes	0	98	1	0	0	0	0	95	98	21	48	0	56	0	58
<i>Pediococcus</i> acidilactici	93	38	1	0	0	0	0	0	11	0	4	0	80	0	0
<i>Pediococcus</i> pentosaceus	2	96	0	0	0	0	0	0	0	0	0	0	0	18	0
<i>Streptococcus</i> pneumoniae	0	1	0	0	71	69	93	93	84	0	79	5	96	0	0
<i>Weissella</i> confusa	95	98	0	0	0	0	0	0	2	0	1	0	0	0	0

^aPreviously designated as *Streptococcus sanguis* II

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- Arduino, M.J., S.K. McAllister, S.M. Aguero, and L.A. Bland. 1994. Abstract C-138. Abstracts of the 94th General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Hinnenbusch, C.J., D.M. Nikolai, and D.A. Bruckner. 1991. Am. J. Clin. Pathol. 96:459-463.
- You, M.S. and R.R. Facklam. 1986. J. Clin. Microbiol. 24:607-611.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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- Lennette, E.H., A. Balows, W.J. Hausler, Jr., and J.P. Truant. 1980. Manual of Clinical Microbiology. 3rd ed. ASM, Washington, D.C.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- Poole, P.M. and G. Wilson. 1976. J. Clin. Microbiol. 29:740-745.
- Facklam, R.R. 1977. J. Clin. Microbiol. 5:184-201.
- Ruoff, K.L. and L.J. Kunz. 1982. J. Clin. Microbiol. 15:920-925.
- Berlitti, F., M.C. Thaller, S. Schippa, F. Pantanella, and R. Pompei. 1993. Int. J. Syst. Bacteriol. 43:63-68.
- Collins, M.D., D. Jones, J.A. Farrow, R. Klipper-Balz, and K.H. Schleifer. 1984. Int. J. Syst. Bacteriol. 34:220-223.
- Facklam, R.R. and M.D. Moody. 1970. Appl. Microbiol. 20:245-250.
- Holt J.G., N.R. Krieg, P.H. Sneath, J.T. Staley, and S.T. Williams. 1994. Bergey's Manual of Determinative Bacteriology. 9th ed. Williams and Wilkins, Baltimore, MD.
- Isenberg, H.D., E.M. Veilozzi, J. Shapiro, and L.G. Rubin. 1988. J. Clin. Microbiol. 26:479-483.
- Blazevic, D.J. and G.M. Ederer. 1975. Principles of Biochemical Tests in Diagnostic Microbiology. John Wiley & Sons, New York, NY.
- Facklam, R.R. 1976. Crit. Rev. Clin. Lab. Sci. 6:287-317.
- Gross, K.D., N.P. Houghton, and L.B. Senterfit. 1975. J. Clin. Microbiol. 1:54-60.
- Bridge, P.D. and P.H. Sneath. 1983. J. Gen. Microbiol. 129:565-597.
- Guilbault, G.G. 1970. Enzymatic Methods of Analysis. Pergamon Press, New York, NY.
- Nord, C.E., A.A. Lindberg, and A. Dahlback. 1975. Med. Microbiol. Immunol. 161:231-238.
- Facklam, R.R., L.G. Thacker, B. Fox, and L.A. Enriquez. 1982. J. Clin. Microbiol. 15:987-990.
- Norris, J.R. and D.W. Ribbons. 1976. Methods in Microbiology. Vol. 9. Academic Press, New York, NY.
- Westley, J.R., P.J. Anderson, V.A. Close, B. Halpern, and E.M. Lederberg. 1967. Appl. Microbiol. 15:822-825.
- Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller. 2007. Manual of Clinical Microbiology. 9th ed. ASM Press, Washington, D.C.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Appelbaum, P.C., M.R. Jacobs, W.M. Palko, E.E. Frauenhoffer, and A. Duffett. 1986. J. Clin. Microbiol. 23:843-846.
- Arduino, M.J., S.K. McAllister, S.M. Aguero, and L.A. Bland. 1994. Abstract C-138. Abstracts of the 94th General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Hinnenbusch, C.J., D.M. Nikolai, and D.A. Bruckner. 1991. Am. J. Clin. Pathol. 96:459-463.
- You, M.S. and R.R. Facklam. 1986. J. Clin. Microbiol. 24:607-611.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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18. bibliografie

- Lennette, E.H., A. Balows, W.J. Hausler, Jr., and J.P. Truant. 1980. Manual of Clinical Microbiology. 3rd ed. ASM, Washington, D.C.
- Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5th ed. ASM, Washington, D.C.
- Poole, P.M. and G. Wilson. 1976. J. Clin. Microbiol. 29:740-745.
- Facklam, R.R. 1977. J. Clin. Microbiol. 5:184-201.
- Ruoff, K.L. and L.J. Kunz. 1982. J. Clin. Microbiol. 15:920-925.
- Berlutt, F., M.C. Thaller, S. Schippa, F. Pantanella, and R. Pompei. 1993. Int. J. Syst. Bacteriol. 43:63-68.
- Collins, M.D., D. Jones, J.A. Farrow, R. Klipper-Balz, and K.H. Schleifer. 1984. Int. J. Syst. Bacteriol. 34:220-223.
- Facklam, R.R. and M.D. Moody. 1970. Appl. Microbiol. 20:245-250.
- Holt J.G., N.R. Krieg, P.H. Sneath, J.T. Staley, and S.T. Williams. 1994. Bergey's Manual of Determinative Bacteriology. 9th ed. Williams and Wilkins, Baltimore, MD.
- Isenberg, H.D., E.M. Veilozzi, J. Shapiro, and L.G. Rubin. 1988. J. Clin. Microbiol. 26:479-483.
- Blazevic, D.J. and G.M. Ederer. 1975. Principles of Biochemical Tests in Diagnostic Microbiology. John Wiley & Sons, New York, NY.
- Facklam, R.R. 1976. Crit. Rev. Clin. Lab. Sci. 6:287-317.
- Gross, K.D., N.P. Houghton, and L.B. Senterfit. 1975. J. Clin. Microbiol. 1:54-60.
- Bridge, P.D. and P.H. Sneath. 1983. J. Gen. Microbiol. 129:565-597.
- Guilbault, G.G. 1970. Enzymatic Methods of Analysis. Pergamon Press, New York, NY.
- Nord, C.E., A.A. Lindberg, and A. Dahlback. 1975. Med. Microbiol. Immunol. 161:231-238.
- Facklam, R.R., L.G. Thacker, B. Fox, and L.A. Eriquez. 1982. J. Clin. Microbiol. 15:987-990.
- Norris, J.R. and D.W. Ribbons. 1976. Methods in Microbiology. Vol. 9. Academic Press, New York, NY.
- Westley, J.R., P.J. Anderson, V.A. Close, B. Halpern, and E.M. Lederberg. 1967. Appl. Microbiol. 15:822-825.
- Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfaller. 2007. Manual of Clinical Microbiology. 9th ed. ASM Press, Washington, D.C.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Appelbaum, P.C., M.R. Jacobs, W.M. Palko, E.E. Frauenhoffer, and A. Duffett. 1986. J. Clin. Microbiol. 23:843-846.
- Arduino, M.J., S.K. McAllister, S.M. Aguero, and L.A. Bland. 1994. Abstract C-138. Abstracts of the 94th General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
- Hinnenbusch, C.J., D.M. Nikolai, and D.A. Bruckner. 1991. Am. J. Clin. Pathol. 96:459-463.
- You, M.S. and R.R. Facklam. 1986. J. Clin. Microbiol. 24:607-611.
- Clinical and Laboratory Standards Institute (CLSI). 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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