remel RapID™ Staph Plus System

EN

REF R8311009..... 20

1. INTENDED USE

The RapID™ Staph Plus System is a qualitative micromethod using enzyme reactions to identify clinically important staphylococcal isolates and other Gram-positive cocci grown on agar. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

A complete listing of the organisms addressed by the RapID Staph Plus System is given in the RapID Staph Plus Differential Chart.

2. SUMMARY AND EXPLANATION

The RapID Staph Plus System is comprised of (1) RapID Staph Plus Panels and (2) RapID Staph Plus Reagent. RapID Staph Plus Panels are disposable plastic trays with 18 reaction cavities, which contain dehydrated reactants. The panel allows for the simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in RapID Inoculation Fluid is used as the test inoculum, which rehydrates and initiates the test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison of results with the probability values in the Differential Chart (Table 4), or by use of the RapID ERIC™ software.

3. PRINCIPLE

The tests used in RapID Staph Plus System are based upon the microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional and single-substrate chromogenic tests and are described in Table 1.

4. REAGENTS

RapID Staph Plus Reagent (provided with kit)(15 ml/Bottle)

Reactive ingredient per liter:

p-Dimethylaminocinnamaldehyde 0.06 g

RapID Staph Plus Reagent is intended for use only in conjunction with RapID Staph Plus System panels.

RapID Inoculation Fluid (R8325106, supplied separately) (2 ml/Tube)

KCl 6.0 g

CaCl₂ 0.5 g

Demineralized Water 1000.0 ml

RapID Nitrate A Reagent (R8309003, supplied separately) (15 ml/Btl)

Sulfanilic Acid 8.0 g

Glacial Acetic Acid 280.0 ml

Demineralized Water 720.0 ml

RapID Nitrate B Reagent (R8309004, supplied separately) (15 ml/Btl)

N,N-Dimethyl-1-naphthylamine 6.0 g

Glacial Acetic Acid 280.0 ml

Demineralized Water 720.0 ml

5. WARNINGS AND PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after use. Directions should be read and followed carefully.

Non-disposable apparatus should be sterilised by any appropriate procedure after use, although the preferred method is to autoclave for 15 minutes at 121°C; disposables should be autoclaved or incinerated. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated area swabbed with a standard bacterial disinfectant or 70% alcohol. Do NOT use sodium hypochlorite. Materials used to clean spills, including gloves, should be disposed of as biohazardous waste.

Do not use reagents beyond the printed expiration dates.

Do not use if there is evidence of contamination or other signs of deterioration.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

Caution!

1. RapID Nitrate A Reagent and RapID Nitrate B Reagent may cause irritation to skin, eyes, and respiratory system.

2. RapID Staph Plus Reagent is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.

3. Refer to Safety Data Sheet for detailed information on reagent chemicals.

Composition / information on ingredients

2-Methoxyethanol 109-86-4

Acetic acid 64-19-7

Hydrochloric acid 7647-01-0

WARNING! This product contains a chemical known in the State of California to cause birth defects or other reproductive harm.

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-352-323-3500 (Call Collect)

DANGER

US ONLY

US & EU

H315	Causes skin irritation
H319	Causes serious eye irritation
H335	May cause respiratory irritation
H336	May cause drowsiness or dizziness
H360	May damage fertility. May damage the unborn child
H373	May cause damage to organs through prolonged or repeated exposure
P201	Obtain special instructions before use
P202	Do not handle until all safety precautions have been read and understood
P281	Use personal protective equipment as required
P264	Wash face, hands and any exposed skin thoroughly after handling
P280	Wear eye/face protection
P260	Do not breathe dust/fume/gas/mist/vapors/spray
P271	Use only outdoors or in a well-ventilated area
P308+313	IF exposed or concerned: Get medical attention/advice
P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing
P302+P352	IF ON SKIN: Wash with plenty of soap and water
P332+P313	IF skin irritation occurs: Get medical advice/attention
P362	Take off contaminated clothing and wash before reuse
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P337+P313	IF eye irritation persists: Get medical advice/attention
P405	Store locked up
P403+P233	Store in a well-ventilated place. Keep container tightly closed
P501	Dispose of contents/container to an approved waste disposal plant

6. STORAGE

2°C - 8°C

RapID Staph Plus System should be stored in its original container at 2-8°C until use. Allow product to equilibrate to room temperature before use. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. RapID Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

7. PRODUCT DETERIORATION

This product should not be used if (1) the expiration date has passed, (2) the plastic tray is broken or the lid is compromised, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{1,5,8}

9. MATERIALS SUPPLIED

• 20 RapID Staph Plus System Panels

• 20 report forms

• RapID Staph Plus Reagent (one plastic dropper-bottle containing reagent sufficient for 20 panels)

• 2 chipboard incubation trays

• 1 color guide

• Instructions for use (IFU)

10. CONTENTS SYMBOLS

Staph Plus Panels	Staph Plus Panels
Report Forms	RapID Report Forms
Staph Plus Reagent	Staph Plus Reagent
Incubation Trays	Incubation Trays

11. MATERIALS REQUIRED BUT NOT SUPPLIED

- Loop sterilization device
- Inoculating loop, swabs, collection containers
- Incubators
- Supplemental media
- Quality control organisms
- Gram stain reagents
- Microscope slides
- Cotton swabs
- RapID Inoculation Fluid - 2 ml (R8325106)
- RapID Nitrate A Reagent (R8309003)
- RapID Nitrate B Reagent (R8309004)
- McFarland #3 turbidity standard or equivalent (R20413),
- Pipettes
- ERIC (Electronic RapID Compendium, R8323600) (Optional).

12. PROCEDURE

Note: Do not interchange reagents among different RapID Systems.

Inoculum Preparation:

- Test organisms must be grown in pure culture, examined by Gram stain, and tested for catalase production prior to inoculum preparation. Only catalase-positive, Gram-positive cocci that morphologically resemble staphylococci should be tested using RapID Staph Plus System. Gram-positive, catalase-negative cocci in chains, resembling streptococci, and Gram-positive bacilli should not be tested.
- Test organisms should be removed from growth media commonly used for staphylococci. The following media are recommended:

Blood Agar (TSA w/ 5% Sheep Blood), Columbia Blood Agar w/ 5% Sheep Blood, or Columbia CNA Agar w/ Sheep Blood.

Notes:

- Plates used for inoculum preparation should be less

Table 1. Principles and Components of the RapID Staph Plus System

Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
1	ADH	Arginine	2%	Utilization of the amino acid substrate produces basic products which raise the pH and change the indicator.	1, 3, 4, 6, 7
2	ODC	Ornithine	2%		
3	LIP	Fatty acid ester	2%	Utilization of the fatty acid substrate produces acidic products which lower the pH and change the indicator.	1, 3, 7
4	SUC	Sucrose	2%	Utilization of the carbohydrate substrate produces acidic products which lower the pH and change the indicator.	1, 3, 4, 6, 7
5	MANO	Mannose	2%		
6	PO ₄	<i>p</i> -Nitrophenyl-phosphate	0.1%	The enzymatic hydrolysis of the colorless aryl-substituted glycoside or phosphoester releases yellow <i>o</i> - or <i>p</i> -nitrophenol which is detected by the formation of a yellow color.	1, 2, 5, 6, 9
7	αGLU	<i>p</i> -Nitrophenyl-α,D-glucoside	0.1%		
8	βGLU	<i>p</i> -Nitrophenyl-β,D-glucoside	0.1%		
9	ONPG	<i>o</i> -Nitrophenyl-β,D-galactoside	0.1%		
10	GUR	<i>p</i> -Nitrophenyl-β,D-glucuronide	0.1%		
11	NAGA	<i>p</i> -Nitrophenyl-N-acetyl-β,D-glucosaminide	0.1%		
12	URE	Urea	0.6%	Utilization of urea creates basic products that change the pH indicator.	1, 3, 4, 6, 7
13	PYR	Pyrrrolidine-β-naphthylamide	0.05%	Enzymatic hydrolysis of the arylamide substrate releases free β-naphthylamine which is detected with RapID Staph Plus Reagent.	2, 6, 9, 10
14	ARG	Arginine-β-naphthylamide	0.05%		
15	ALA	Alanine-β-naphthylamide	0.05%		
16	LEU	Leucine-β-naphthylamide	0.05%		
17	LGLY	Leucyl-glycine-β-naphthylamide	0.05%		
18	NIT	Potassium nitrate	0.015%	Utilization of nitrate results in the formation of nitrite which is detected by nitrate reagents.	1, 3, 4, 6, 7

than 72 hours old. Plates cultivated for 18-24 hours are recommended. Slow growing strains may be incubated for 48-72 hours if insufficient growth is available after 18-24 hours.

- The use of media other than those recommended may compromise test performance.
- Using a cotton swab or inoculating loop, suspend sufficient growth from the culture in RapID Inoculation Fluid (2 ml) to achieve a visual turbidity minimally equivalent to a #3 McFarland turbidity standard.


Notes:

- Suspensions less turbid than a #3 McFarland standard may result in aberrant reactions. Suspensions slightly more turbid will not compromise test performance and are recommended for stock cultures and quality control strains.
- Suspensions should be mixed thoroughly and vortexed if required.
- Suspensions should be used within 15 minutes of preparation.

- An agar plate may be inoculated to verify purity and any additional testing that may be required using a loopful of the test suspension from the Inoculation Fluid tube. Incubate the plate for 18-24 hours at 35-37°C.

Inoculation of RapID Staph Plus Panel:

- Peel back the panel lid over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
- Using a pipette, gently transfer the entire contents of the Inoculation Fluid tube into the upper right-hand corner of the panel and reseal the inoculation port by pressing the peel tab back in place.
Note: It is important that the panel receive the entire test suspension.
- After adding the test suspension and while keeping the panel on a level surface, tilt the panel back away from the test cavities at approximately 45°. See below:

- While tilted back, gently rock the panel from side to side

RapID Staph Plus Panel Test Location

Cavity #	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
Test Code	ADH	ODC	LIP	SUC	MANO	PO ₄	αGLU	βGLU	ONPG	GUR	NAGA	URE	PYR	ARG	ALA	LEU	LGLY	NIT
Reagent	None												RapID Staph Plus Reagent			RapID Nitrate A Reagent RapID Nitrate B Reagent		

Table 2. Interpretation of RapID Staph Plus System Tests*

Cavity #	Test Code	Reagent	Reaction		Comments
			Positive	Negative	
1	None	ADH	Purple, blue, bluish gray, or greenish gray	Yellow or, straw	Any purple or significant blue, or gray color should be scored positive.
2		ODC			
3	None	LIP	Yellow, yellow orange	Red, dark red orange, or orange	Any yellow or yellow orange throughout the well should be scored positive.
4		SUC			
5		MANO			
6		None			
7	None	αGLU	Any yellow or light yellow	Clear, cream or buff	Any yellow color throughout the well should be scored positive. Note, yellow color may be faint with some taxa.
8		βGLU			
9		ONPG			
10		GUR			
11		NAGA			
12	None	URE	Red or dark red orange	Yellow, yellow orange, or very pale orange	Any red or dark orange throughout the well should be scored positive.
13	RapID STAPH Plus Reagent	**PYR	Purple or dark red	Yellow, orange, pale light red or shades of pink	Only significant purple or red should be scored positive.
14		ARG	Purple, red, light red or dark pink	Yellow, orange or very pale pink	Any purple, red or dark pink should be scored positive. Yellow, pale light pinks, and orange are scored negative.
15		ALA			
16		LEU			
17	LGLY				
18	RapID Nitrate A and Nitrate B	NIT	Cherry red or dark pink	Clear, yellow, cream or very pale pink	Only a significant red or dark pink, developing immediately or within 30 seconds, should be scored positive. Pinks that develop after prolonged periods after reagent addition should be scored negative.

*NOTE: Panels are best observed by looking down through the reaction cavities against a white background.

** NOTE: PYR color development is interpreted differently from other tests requiring addition of RapID Staph Plus Reagent.

Table 3. Quality Control Chart for RapID Staph Plus Panels

Organism	ADH	ODC	LIP	SUC	MANO	PO ₄	αGLU	βGLU	ONPG	GUR	NAGA	URE	PYR	ARG	ALA	LEU	LGLY	NIT
<i>Staphylococcus haemolyticus</i> ^a ATCC® 29970	+	-	-(+)	-(+)	-	-	V	-	-	+	-	-	+	-	v	-	-	v
<i>Staphylococcus saprophyticus</i> ^a ATCC® 35552	-	-	+	+	-	+(-)	+	-	+	-	-	+	+	-	-	-(+)	-	-
<i>Enterobacter aerogenes</i> ATCC® 13048	-	+	V	+	+	+	V	+	+	-	+	-(+)	+	+	+	+	+	+
<i>Oligella urelytica</i> ATCC® 43534	-(+)	-	-	-	-	v	-	-	-	-	-	+	-	+	+	+	+	-

+ = Positive QC reaction, - = Negative QC reaction, V = equivocal or variable reaction, -(+) = Usually negative, +(-) = Usually positive

NOTE: Only reactions marked as + or - are necessary for panel quality control.

^a Key indicator strains demonstrate acceptable performance of the most labile substrate in the system and reactivity in a significant number of wells, according to Clinical and Laboratory Standards Institute recommendations for streamlined quality control.¹¹

Incubation of RapID Staph Plus Panel:

Incubate inoculated panels at 35-37°C in a non-CO₂ incubator for **at least 4 HOURS and not more than 6 HOURS**. For ease of handling, panels may be incubated in the chipboard incubation trays provided with the kit.

Scoring of RapID Staph Plus Panels:

RapID Staph Plus Panels contain 18 reaction cavities that provide 18 test scores. Reagent-requiring tests (cavities 13-18) are designated with a box drawn around them as illustrated below. Refer to Table 2 for complete information on interpreting RapID Staph Plus System tests.

- While firmly holding the RapID Staph Plus Panel on the bench top, peel off the label lid over the reaction cavities by pulling the lower right hand tab up and to the left.
- Cavities 1 (ADH) through 12 (URE) are read without addition of reagent. To the remaining cavities (13 through 18), add the following reagents as indicated:
 - Add 2 drops of RapID Staph Plus Reagent to cavities 13 (PYR) through 17 (LGLY).
 - Add 1 drop of RapID Nitrate A Reagent to cavity 18 (NIT).
 - Add 1 drop of RapID Nitrate B Reagent to cavity 18 (NIT).

Note: In most cases, immediate color development will be observed after addition of reagents. If immediate color development is not observed, allow approximately 30 seconds (but no more than 1 minute) for color development.

- After reagent addition, read and score the test cavities from left to right using the interpretation guide presented in color guide and Table 2. Panels should be read by looking down through the reaction cavities against a white background. Record the scores in the appropriate boxes on the report form.
- To identify the isolate, reference the microcode obtained in ERIC or use the Differential Chart included in this IFU.

13. RESULTS AND RANGE OF EXPECTED VALUES

The RapID Staph Plus Differential Chart (Table 4) illustrates expected results for the RapID Staph Plus System. Differential Chart results are expressed as a series of positive percentages for each system test. This information statistically supports the use of each test and provides the basis, through numerical coding of digital test results, for a probabilistic approach to the identification of the test isolate. Identifications are made using individual test scores from RapID Staph Plus panels in conjunction with other laboratory information (e.g., Gram stain, colonial morphology, coagulase test) to produce a pattern that statistically resembles known

reactivity for taxa recorded in the RapID Staph Plus System database. These patterns are compared through the use of the RapID Staph Plus Differential Chart or by derivation of a microcode and the use of ERIC.

14. QUALITY CONTROL

All lot numbers of RapID Staph Plus System have been tested using the quality control organisms listed in Table 3 and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported. Expected results for selected quality control organisms are listed in Table 3.

Notes:

- The quality control of RapID Staph Plus Reagent is accomplished by obtaining the expected reactions for tests requiring the addition of the reagent (cavities 13-17).
- Organisms which have been repeatedly transferred on agar media for prolonged periods may provide aberrant results.
- Quality control strains should be stored frozen or lyophilized. Prior to use, transfer quality control strains 2-3 times after removing from storage on agar recommended for use with the RapID Staph Plus System.
- Formulations, additives, and ingredients of culture media vary from manufacturer to manufacturer and may vary from batch to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve discrepancies.

15. LIMITATIONS

- The use of RapID Staph Plus System and the interpretation of results requires the knowledge of a competent microbiologist, familiar with laboratory procedures, who is trained in general microbiological methods and judiciously makes use of training, experience, specimen information, and other pertinent procedures before reporting the identification obtained using this system.
- The RapID Staph Plus System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.
- The RapID Staph Plus System is designed for use with the taxa listed in the RapID Staph Plus Differential Chart (Table 4). Gram-positive bacilli, Gram-negative cocci, or catalase-negative Gram-positive cocci in chains should

not be tested.

- Expected values listed for RapID Staph Plus System tests may differ from conventional test results or previously reported information.
- The accuracy of the RapID Staph Plus System is based upon the statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single test found in the RapID Staph Plus System to establish the identification of a test isolate is subject to the error inherent in that test alone.
- Results obtained using the RapID Staph Plus System are dependent upon adherence to the procedures indicated. Any change or modification in the procedure may result in aberrant results.

16. PERFORMANCE CHARACTERISTICS

The RapID Staph Plus System performance characteristics have been established by laboratory testing of clinical, reference, and stock cultures at Remel. Among 291 strains tested, 280 (96%) agreed with the reported reference identifications (with or without supplementary tests). One (0.3%) strain provided a questionable microcode that did not result in an identification. Ten strains (3.4%) did not agree with reported reference identifications.

17. BIBLIOGRAPHY

- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Guilbert, G.G. 1970. Methods of Enzymatic Analysis. Pergamon Press, New York, NY.
- Heikens, E., A. Fleer, A. Paauw, A. Florijn, and A.C. Fluit. 2005. J. Clin. Microbiol. 43:2286-2290.
- Ieven, M., J. Verhoeven, S.R. Pattyn, and H. Goossens. 1995. J. Clin. Microbiol. 33:1060-1063.
- Leber, A.L. 2016. Clinical Microbiology Procedures Handbook. 4th edition, ASM Press, Washington, D.C.
- Murray, P.R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. Manual of Clinical Microbiology. 12th Edition. ASM Press.
- Bascomb, S. and M. Manafi. 1998. Clin. Microbiol. Rev. 11:318-340.
- Norris, J.R. and D.W. Ribbons. 1976. Methods in Microbiology. Vol. 9. Academic Press, New York, NY.

11. Clinical and Laboratory Standards Institute. 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

18. PACKAGING

REF R8311009 RapID Staph Plus System 20 Tests/Kit

19. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

RapID™ is a trademark of Thermo Fisher Scientific and its subsidiaries.

ERIC™ is a trademark of Thermo Fisher Scientific and its subsidiaries.

ATCC® is a registered trademark of American Type Culture Collection.

For technical assistance please contact your local distributor

Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
 Tel: (800) 255-6730 • International: (913) 888-0939
www.oxid.com/IFU
 Europe +800 135 79 135 • US 1 855 2360 190
 CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU8311009	August 2023 Updated to meet IVDR requirements

Printed in the UK

Table 4 - RapID Staph Plus Differential Chart

Organism	ADH	ODC	LIP	SUC	MANO	PO ₄	αGLU	βGLU	ONPG	GUR	NAGA	URE	PYR	ARG	ALA	LEU	LGLY	NIT
<i>Staphylococcus</i>																		
<i>S. arlettae</i>	0	0	0	81	50	79	74	99	92	91	0	0	0	0	88	0	0	0
<i>S. aureus</i>	54	2	79	96	89	94	95	92	2	1	20	26	50	1	8	11	2	77
<i>S. auricularis</i>	44	5	42	19	0	5	39	0	4	0	0	1	99	96	11	39	4	58
<i>S. capitis ss capitis</i>	29	0	39	2	76	2	5	0	0	0	1	1	44	1	27	5	0	69
<i>S. capitis ss urealyticus</i>	13	0	17	30	84	5	90	0	0	0	0	99	77	1	22	4	0	80
<i>S. caprae</i>	99	0	93	1	80	7	22	0	8	0	0	36	99	0	1	1	0	99
<i>S. carnosus</i>	99	0	0	0	89	67	2	0	96	0	0	0	98	0	2	22	0	99
<i>S. chromogenes</i>	9	0	92	66	49	69	16	6	0	0	0	39	99	0	99	56	80	99
<i>S. coagulans^d</i>	99	50	99	17	0	99	1	0	53	0	0	99	99	0	44	80	0	99
<i>S. cohnii ss urealyticus</i>	7	0	80	1	77	72	60	19	88	90	6	94	95	0	69	2	1	5
<i>S. delphini</i>	30	0	24	47	52	99	0	0	99	0	0	90	99	0	91	90	2	9
<i>S. epidermidis</i>	47	6	77	90	18	94	88	9	29	0	2	90	0	0	2	2	0	70
<i>S. equorum</i>	0	0	1	88	90	62	2	99	27	74	2	67	99	0	87	38	0	99
<i>S. felis</i>	21	0	99	13	92	75	0	0	87	0	0	99	99	0	93	0	0	41
<i>S. gallinarum</i>	9	0	22	99	92	98	51	99	58	0	21	92	90	0	2	17	0	99
<i>S. haemolyticus</i>	95	0	71	91	1	3	97	40	11	70	50	0	99	3	74	11	5	86
<i>S. hominis ss hominis</i>	9	2	84	80	3	15	98	14	11	0	5	97	23	1	56	9	5	78
<i>S. hominis ss novo</i>	0	0	99	99	0	0	0	92	0	0	0	99	0	0	0	0	0	99
<i>S. hyicus</i>	12	0	99	85	80	55	0	27	5	89	0	0	0	9	91	48	0	91
<i>S. intermedius</i>	26	0	95	68	43	99	9	12	95	0	0	35	99	0	59	66	2	99
<i>S. kloosii</i>	0	0	98	0	0	38	99	85	83	45	0	93	93	0	9	0	0	0
<i>S. lugdunensis</i>	3	99	21	95	80	9	79	90	13	1	1	42	91	0	33	9	5	79
<i>S. muscae</i>	0	0	99	74	0	0	0	99	0	90	0	0	0	0	52	20	0	99
<i>S. pasteurii</i>	88	0	29	91	41	90	66	99	0	98	1	80	2	0	10	5	0	90
<i>S. saprophyticus</i>	1	0	92	95	2	82	88	39	92	0	9	99	74	9	44	31	2	1
<i>S. schleiferi^f</i>	99	0	99	5	50	99	4	0	0	0	9	3	92	0	47	52	0	41
<i>S. simulans</i>	95	0	12	86	24	74	91	1	96	71	15	87	95	0	11	9	0	99
<i>S. ureilyticus^g</i>	0	0	98	2	9	71	26	0	1	4	0	2	13	6	22	8	0	0
<i>S. warneri</i>	57	3	87	81	7	4	92	89	3	74	1	84	24	8	29	11	3	33
<i>S. xylosus</i>	4	0	96	93	81	80	74	93	84	83	18	94	95	6	12	15	0	88
Other Organisms																		
<i>Kocuria rosea</i>	0	0	99	82	9	0	99	18	0	0	0	0	0	44	38	55	11	83
<i>Kocuria varians</i>	0	0	0	14	40	9	99	0	13	0	0	99	26	58	95	38	12	99
<i>Kytococcus sedentarius</i>	14	0	3	0	0	51	99	0	0	0	0	0	3	89	91	80	98	4
<i>Macrococcus caseolyticus</i>	0	0	37	13	7	50	99	9	14	0	0	0	25	12	98	74	0	98
<i>Mammaliococcus lentus^b</i>	0	0	90	99	97	98	25	99	57	0	95	3	9	7	22	73	0	99
<i>Mammaliococcus sciuri^f</i>	0	0	3	98	93	95	56	99	0	16	98	2	3	0	4	21	3	98
<i>Mammaliococcus vitulinus^f</i>	0	0	0	99	18	40	6	81	0	0	0	0	0	0	0	9	0	99
<i>Micrococcus sp.</i>	4	0	9	3	1	48	92	2	0	0	0	44	92	95	95	95	90	3
<i>Rothia kristinae^e</i>	0	0	0	98	95	1	99	95	5	0	0	1	99	81	95	88	81	2
<i>Rothia mucilaginosa</i>	0	0	0	99	0	0	0	0	99	0	0	0	99	99	99	99	0	99

^a Previously designated as *S. cohnii ss cohnii*

^b Previously designated as *Staphylococcus lentus*

^c Previously designated as *Staphylococcus schleiferi ss schleiferi*

^d Previously designated as *Staphylococcus schleiferi ss coagulans*

^e Previously designated as *Staphylococcus sciuri*

^f Previously designated as *Staphylococcus vitulinus*

^g Previously designated as *Kocuria kristinae*

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17. BIBLIOGRAFIA

- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Guilbert, G.G. 1970. Methods of Enzymatic Analysis. Pergamon Press, New York, NY.
- Heikens, E., A. Fleer, A. Paauw, A. Florijn, and A.C. Fluit. 2005. J. Clin. Microbiol. 43:2286-2290.
- leven, M., J. Verhoeven, S.R. Pattyn, and H. Goossens. 1995. J. Clin. Microbiol. 33:1060-1063.
- Leber, A.L. 2016. Clinical Microbiology Procedures Handbook. 4th edition, ASM Press, Washington, D.C.
- Murray, P.R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C.
- MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
- K.C. Carroll, M.A. Pfaller, M.L. Landry, A.J. McAdam, R. Patel, S.S. Richter, D.W. Warnock. 2019. Manual of Clinical Microbiology. 12th Edition. ASM Press.
- Bascomb, S. and M. Manafi. 1998. Clin. Microbiol. Rev. 11:318-340.
- Norris, J.R. and D.W. Ribbons. 1976. Methods in Microbiology. Vol. 9. Academic Press, New York, NY.
- Clinical and Laboratory Standards Institute. 2008. Quality Control for Commercial Microbial Identification Systems; Approved Guideline. M50-A. CLSI, Wayne, PA.

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

www.thermofisher.com/microbiology

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