Key Code TSMX7826B
www.thermofisher.com

Europe + 800 135 79 135 US 1 855 236 0910
CA 1 855 805 8539 ROW +31 20 794 7071

remel EN Staphaurex Plus

1. INTENDED USE

Staphaurex™ Plus is a qualitative latex slide agglutination test for the differentiation of *Staphylococcus aureus* from other *Staphylococcus* species isolates grown on agar by the detection of clumping factor and Protein A and/or surface antigens specific to *Staphylococcus aureus*. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

2. SUMMARY AND EXPLANATION OF THE TEST

S. aureus possess a number of properties which are used to confirm identification. These include free coagulase, clumping factor (bound coagulase), thermonuclease and protein A¹. The tube coagulase test detects free coagulase and is considered as a reference test for *S. aureus*¹. This test, however, takes 4 to 24 hours and plasma may show a lot-to-lot variation². Over the past decade particle agglutination assays have been developed which give a much more rapid identification^{3,4}. These first generation assays are based on latex particles or red cells coated with either fibrinogen alone, to detect clumping factor, or fibrinogen and immunoglobulin G (IgG), to detect both clumping factor and staphylococcal protein A.

Recently it has been shown that these tests can fail to detect certain strains of *S. aureus*, particularly a proportion of methicillin/oxacillin resistant strains (MRSA)^{5,6,7}. Some of these strains may express undetectable levels of clumping factor and protein A⁸.

Two antigens, somatic type 18⁹ and capsular type 5^{10,11} have been associated with the methicillin-resistant phenotype. The incorporation of antisera to these antigens may improve the sensitivity of agglutination assays for MRSA strains. Investigations on strains that are negative in rapid assays have shown that antibodies to a single somatic or capsular antigen are insufficient to detect all strains that are negative with the first generation of particle agglutination tests. Staphaurex Plus uses latex beads coated with fibrinogen to detect the majority of clinical strains and IgG specific for a carefully selected group of strains that are negative in the first generation tests.

3. PRINCIPLE OF THE PROCEDURE

The Staphaurex Plus Test Latex consists of yellow latex particles which have been coated with fibrinogen and rabbit immunoglobulin G (IgG) specific for *S. aureus*. When a drop of the reagent is mixed on a card with *S. aureus* organisms, rapid agglutination occurs through the interaction of (i) fibrinogen and clumping factor, (ii) the Fc portion of IgG and protein A or (iii) specific IgG and cell surface antigens.

Some strains of *Staphylococcus spp.*, particularly *S. saprophyticus*,

may cause non-specific aggregation of latex particles. Therefore a Control Latex is provided to assist with the identification of non-specific reactions.

4. REAGENTS

KIT CONTENTS

| Staphaurex Plus | ZL33/R30950102 150 Tests | ZL34/R30950201 450 Tests |
|---|-----------------------------|-----------------------------|
| 1. Test Latex (Yellow cap) | 1 dropper bottle | 3 dropper bottles |
| 2. Control Latex (Grey cap) | 1 dropper bottle | 3 dropper bottles |
| 3. Disposable Reaction Cards (RT64/R30369001) | 2 packs | 6 packs |
| 4. Disposable Mixing Sticks | 3 bundles | 9 bundles |
| 5. Instructions for Use | 1 | 1 |

5. DESCRIPTION OF REAGENTS, PREPARATION FOR USE AND RECOMMENDED STORAGE CONDITIONS

See also **Warnings and Precautions**.

The latex suspensions are provided ready to use and should be stored in an upright position at 2 to 8°C, where they will retain activity at least until the date shown on the bottle label. Do not freeze. Avoid storage at room temperature (15 to 30°C). Do not stand the reagent in bright light on the bench.

TEST LATEX

Test Latex

A buffered suspension of yellow polystyrene latex particles coated with an enzymic digest of human fibrinogen (approx. 0.02% w/v) and rabbit IgG (approx. 0.02% w/v). Contains 0.05% w/v Bronidox® preservative¹².

Materials of human origin have been tested for the presence of hepatitis B surface antigen, anti-HCV and anti-HIV-1/HIV-2 and found to be negative.

CONTROL LATEX

Control Latex

A buffered suspension of yellow polystyrene latex particles with bovine serum albumin (approx. 0.2% w/v) unreactive with *S. aureus*. Contains 0.05% Bronidox® preservative¹².

Reaction Cards and Mixing Sticks should be stored at room temperature (15 to 30°C). Staphaurex Plus (ZL33/R30950102 and ZL34/R30950201) was developed using RT64/R30369001 Disposable Reaction Cards.

Do not substitute another disposable slide for the RT64/R30369001 Disposable Reaction Cards when samples are tested using Staphaurex Plus.

6. WARNINGS AND PRECAUTIONS

IVD For *in vitro* diagnostic use only. For professional use only.

Please refer to the manufacturer's safety data sheet and the product labelling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. In the event of malfunction do not use device.

HEALTH AND SAFETY INFORMATION

- CAUTION:** This kit contains human sourced components. No known test can offer complete assurance that products derived from human sources will not transmit infection. Therefore, all human sourced material should be considered potentially infectious. It is recommended that these reagents and test specimens be handled using established good laboratory working practices.
- Non-disposable apparatus should be sterilised by any appropriate procedure after use, although the preferred method is to autoclave for 15 minutes at 121°C. Disposables should be autoclaved or incinerated. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated areas swabbed with a standard bacterial disinfectant. Materials used to clean spills, including gloves, should be disposed of as biohazardous waste.
- Wear laboratory coat, disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- When used in accordance with the principles of Good Laboratory Practice, good standards of occupational hygiene and the instructions in these Instructions for Use, the reagents supplied are not considered to present a hazard to health.

ANALYTICAL PRECAUTIONS

- Do not use the reagents beyond the stated expiry date.
- Latex reagents should be brought to room temperature (15 to 30°C) before use. Latex reagents which show signs of aggregation or 'lumpiness' before use may have been frozen and should not be used.
- It is important when using dropper bottles that they are held vertically and that the drop forms at the tip of the nozzle. If the nozzle becomes wet an incorrect volume will form around the end and not at the tip; if this occurs dry the nozzle before proceeding.
- Do not touch the reaction areas on the cards.
- Do not interpret agglutination that appears after 30 seconds as a positive result. Prolonged rocking can result in false-positive reactions with some coagulase-negative isolates.
- Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.

7. SPECIMEN COLLECTION AND STORAGE

For details of specimen collection and treatment a standard text book should be consulted¹. Cultures may be tested from any of the following media:

| | |
|----------------------------------|-----------------------------------|
| Blood agar | Columbia CNA agar |
| Nutrient agar | Mueller Hinton agar with 5% blood |
| Tryptone soya agar | Baird-Parker agar |
| Tryptone soya agar with 5% blood | Mannitol-salt agar† |
| Columbia blood agar | |

†Note: specimens grown on media containing antibiotics or a high-salt-supplemented medium such as Mannitol-salt agar may give an agglutination containing stringy aggregates.

THE USE OF FRESH CULTURES GROWN OVERNIGHT IS RECOMMENDED.

8. PROCEDURE

MATERIALS PROVIDED

Sufficient materials are provided for 150 (ZL33/R30950102) or 450 (ZL34/R30950201) tests, see **Kit Contents**.

TEST PROCEDURE

Please read **Analytical Precautions** carefully before performing the test.

- Step 1** Shake vigorously and examine the latex reagents for aggregation before use. Refer to **Quality Control Section** and **Visual Inspection** for additional instructions.
- Step 2** For each test sample place one drop of **Test Latex** in one circle on a Reaction Card **1 drop** (RT64/R30369001) and one drop of **Control Latex** in a separate circle. Ensure that the dropper bottles are held vertically to dispense an accurate drop.
- Step 3** Using a mixing stick, remove sufficient growth from a pure culture or well-isolated colonies to cover the blunt end of the stick. As a guide, an amount of growth roughly equivalent to six average-sized colonies should be used.
- Step 4** Emulsify the sample of culture in the drop of **Test Latex** by rubbing with the flat end of the stick. Rub thoroughly, but not too vigorously or the surface of the card may be damaged. Some strains, particularly of species other than *S. aureus* remain difficult to emulsify and this should be noted, since lumps of unemulsified culture can make the latex appear 'rough' or 'stringy' when read. Spread the latex over approximately half the area of the circle. Discard the mixing stick for safe disposal. **Emulsify sample**
- Step 5** Using a separate stick, emulsify a similar culture sample in the **Control Latex**, as **sample** stated in Step 4. Discard the mixing stick for safe disposal. **Emulsify sample**
- Step 6** Rock the card slowly for up to 30 seconds while observing for agglutination. The card should be held at normal reading distance (25 to 35 cm) from the eyes. Do not use a magnifying lens. **Rock**
- Step 7** Discard the used Reaction Card for safe disposal.

9. RESULTS

Positive Result

Agglutination of the Test Latex accompanied by a lack of agglutination of the Control Latex indicates the presence of either coagulase, protein A or antigens commonly found on *S. aureus* in the culture under test. Most positive reactions will be almost instantaneous. False positive results can occur if the test is read after more than 30 seconds.

Negative Result

Lack of agglutination in both reagents means that the culture under test is unlikely to be *S. aureus*.

Non-interpretable Result

Visible agglutination of the Control Latex, whether stronger or weaker than the Test Latex, indicates a non-specific reaction.

QUALITY CONTROL

Quality control testing should be run with each shipment and new kit lot number received. Each laboratory should follow their state and local requirements.

Any departure from the expected results indicates that there may be a problem with the reagents, which must be resolved before further use with clinical samples.

Visual inspection

The latex suspensions should always be inspected for aggregation as they are dropped onto the Reaction Card. If there is evidence of clumping before addition of the test sample the suspension should not be used. After prolonged storage some aggregation or drying may have occurred around the top of the bottle. If this is observed the bottle should be shaken vigorously for a few seconds until resuspension is complete.

Control procedure

The performance of the Test and Control Latex reagents should be confirmed using fresh, overnight cultures of reference strains of bacteria, following the method described in **Test Procedure**. Suitable reference strains are shown below.

| SPECIES | EXPECTED RESULT | |
|--|-----------------|---------------|
| | TEST LATEX | CONTROL LATEX |
| <i>S. aureus</i> (ATCC [®] 25923 TM) | + | - |
| <i>S. epidermidis</i> (ATCC [®] 12228 TM) | - | - |

10. INTERPRETATION OF RESULTS

A positive reaction indicates the presence of one or more of clumping factor, protein A or cell surface antigens in the culture under test and a negative result indicates their absence.

11. LIMITATIONS OF THE PROCEDURE

- Specimens grown on media containing antibiotics or a high-salt-supplemented medium such as mannitol-salt agar may give an agglutination containing stringy aggregates.
- Some species of staphylococci in addition to *S. aureus* notably *S. hyicus*, *S. intermedius*, *S. lugdunensis* and *S. schleiferi*, may give positive results in coagulase tests and may also react in rapid latex procedures. If necessary these species may be identified by biochemical test procedures. *S. hyicus* and *S. intermedius* are encountered rarely in the clinical laboratory.
- Some other coagulase negative staphylococcal species, such as *S. capitis* possess plasma protein binding factors, which do not react in the Staphaurex Plus test. However, a few strains identified biochemically as *S. saprophyticus* have given weak positive reactions and further identification of urinary isolates may be required.
- Some streptococci and possibly other organisms possess immunoglobulin or other plasma protein binding factors which can react in the latex test and there are several bacteria such as *E. coli*, which are able to non-specifically agglutinate latex particles. To eliminate potential interference from these organisms a Gram stain and catalase test should be performed so that only organisms with staphylococcal morphology are tested.
- All questionable results should be checked for purity and identified by an alternative method.

12. EXPECTED RESULTS

Strong agglutination with *S. aureus* cultures, no agglutination with staphylococci which possess neither clumping factor, protein A or surface antigens characteristic of *S. aureus*.

13. SPECIFIC PERFORMANCE CHARACTERISTICS

The performance of Staphaurex Plus has been evaluated in four North American and seven European microbiological reference laboratories on a total of 1293 routine (presumed staphylococcal) clinical isolates and 820 stored cultures. The cultures were tested in parallel with the tube coagulase procedure, Gram stain and at least one alternative rapid test for the identification of *S. aureus*. The results are summarised **Tables 1 and 2**.

CLINICAL ISOLATES

Methicillin Resistant *S. aureus* (MRSA)

A total of 241 fresh *S. aureus* cultures shown to be resistant to one or more antibiotics were tested in the American and European reference laboratories. Staphaurex Plus correctly identified 240 of these isolates. The discrepant isolate was positive with a tube coagulase test and an alternative rapid latex test.

The sensitivity of Staphaurex Plus on this group of MRSA cultures is estimated to be 99.6% (240/241).

Methicillin Sensitive *S. aureus* (MSSA)

Staphaurex Plus correctly identified 700 of 703 confirmed *S. aureus* cultures from the microbiological reference laboratories. The discrepant isolates included two which also gave a negative result with the alternative rapid latex test.

The sensitivity of Staphaurex Plus on this group of MSSA cultures is estimated to be 99.6% (700/703).

Other staphylococci

A total of 349 fresh non-*S. aureus* staphylococcal isolates were also tested. Staphaurex Plus gave a negative result with 324 of these isolates which included *S. saprophyticus*, *S. epidermidis* and *S. haemolyticus*. The remaining 25 cultures which gave a positive result with Staphaurex Plus included 16 which were also positive with an alternative rapid latex test.

The specificity of Staphaurex Plus on this group of non-*S. aureus* staphylococcal cultures is estimated to be 92.8% (324/349).

Overall Performance of Staphaurex Plus in Comparison with Tube Coagulase on *S. aureus* Isolates

| | |
|----------------------|-------|
| Relative Sensitivity | 99.6% |
| Relative Specificity | 92.8% |
| Overall agreement | 97.8% |

NOTE: Staphaurex Plus gave a non-interpretable result with 0.15% (2/1295) of the fresh cultures, which has been excluded from the summary above.

STORED CULTURES

Methicillin Resistant *S. aureus* (MRSA)

A total of 336 stored *S. aureus* cultures shown to be resistant to one or more antibiotics were tested. Staphaurex Plus correctly identified 335 of these isolates. The discrepant culture was positive with a tube coagulase test and negative with an alternative rapid latex test.

The sensitivity of Staphaurex Plus on this group of MRSA cultures is estimated to be 99.7% (335/336).

Methicillin Sensitive *S. aureus* (MSSA)

Staphaurex Plus correctly identified 326 of 332 confirmed *S. aureus* cultures from the microbiological reference laboratories. The discrepant cultures included four which also gave a negative result with the alternative rapid latex test.

The sensitivity of Staphaurex Plus on this group of MSSA cultures is estimated to be 98.2% (326/332).

Other staphylococci

A total of 152 stored non-*S. aureus* staphylococcal cultures were also tested. Staphaurex Plus gave a negative result with 144 of these isolates which included *S. saprophyticus*, *S. epidermidis* and *S. haemolyticus*. The remaining eight cultures which gave a positive result with Staphaurex Plus included two which were also positive with an alternative rapid latex test.

The specificity of Staphaurex Plus on this group of non-*S. aureus* staphylococcal cultures is estimated to be 94.7% (144/152).

Overall Performance of Staphaurex Plus in Comparison with Tube Coagulase on stored *S. aureus* Cultures

| | |
|----------------------|-------|
| Relative Sensitivity | 99.0% |
| Relative Specificity | 94.7% |
| Overall agreement | 98.1% |

NOTE: Staphaurex Plus gave a non-interpretable result with 0.36% (3/823) of the stored cultures, which have been excluded from the summary above.

Table 1 Reactivity of Staphaurex Plus on Presumed Staphylococcal fresh Clinical Isolates^a

| | Staphaurex Plus result | | |
|---|------------------------|----------|-------|
| | Positive | Negative | Total |
| Methicillin Resistant <i>S. aureus</i> (MRSA) | 240 | 1 | 241 |
| Methicillin Sensitive <i>S. aureus</i> (MSSA) | 700 | 3 | 703 |
| Non- <i>S. aureus</i> isolates ^b | 25 | 324 | 349 |

^a Staphaurex Plus gave a non-interpretable result with 2 samples. These have been excluded from the table.

^b includes *S. saprophyticus*, *S. epidermidis* and *S. haemolyticus*.

Table 2 Reactivity of Staphaurex Plus on Stored Staphylococcal Cultures^a

| | Staphaurex Plus result | | |
|---|------------------------|----------|-------|
| | Positive | Negative | Total |
| Methicillin Resistant <i>S. aureus</i> (MRSA) | 335 | 1 | 336 |
| Methicillin Sensitive <i>S. aureus</i> (MSSA) | 326 | 6 | 332 |
| Non- <i>S. aureus</i> cultures ^b | 8 | 144 | 152 |

^a Staphaurex Plus gave a non-interpretable result with 3 samples. These have been excluded from the table.

^b includes *S. saprophyticus*, *S. epidermidis* and *S. haemolyticus*.

14. BIBLIOGRAPHY




- Kloos, W.E. and Lambe, D.W. (1991). *Staphylococcus. Manual of Clinical Microbiology*, 5th Ed., edited by Balows, A., Hausler, W.J., Herrmann, K.L., Isenberg, H.D. and Shadomy, H.J. American Society for Microbiology, Washington, D.C. Pages 222-237.
- Selepak, S.T. and Witebsky F.G. (1985). Inoculum size and lot-to-lot variation as significant variables in the tube coagulase test for *Staphylococcus aureus*. *J. Clin. Microbiol.*, 22, 835-837.
- Essers, L. and Radebold, K. (1980). Rapid and reliable identification of *Staphylococcus aureus* by a Latex Agglutination Test. *J. Clin. Microbiol.*, 12, 641-643.
- Orsi, A., Bartoloni, A., et al (1989). Evaluation of six different agglutination methods for rapid identification of *Staphylococcus aureus*. *Eur. J. Clin. Microbiol. Infect. Dis.*, 8, 655-656.
- Piper, J., Hadfield, T., et al (1988). Efficacies of rapid agglutination tests for identification of Methicillin-Resistant Staphylococcal strains as *Staphylococcus aureus*. *J. Clin. Microbiol.*, 26, 1907-1909.
- Ruane, P., Morgan, M.A., et al (1986). Failure of rapid agglutination methods to detect Oxacillin-Resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 24, 490-492.
- Wanger, A.R., Morris, S.L., et al (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583-2588.

- Roberts, J.I.S. and Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
- Chabbert, Y.A. and Pillet, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus aureus*. *Nature*, 213, 1137.
- Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.
- Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.
- Henkel KGaA. Manufacturers Information and Safety Data Sheet for Bronidox[®] L.

15. PACKAGING

| | | |
|-----|---------------------|-----|
| REF | ZL33/R30950102..... | 150 |
| | ZL34/R30950201..... | 450 |

16. SYMBOL LEGEND

| | |
|---|---|
|  | Catalogue Number |
|  | In Vitro Diagnostic Medical Device |
|  | Consult Instructions for Use (IFU) |
|  | Temperature Limitations (Storage temp.) |
|  | Contains sufficient for <N> tests |
|  | Not for near patient testing |
|  | Batch Code (Lot Number) |
|  | Use By (Expiration Date) |
|  | Importer |
|  | Unique Device Identifier |
|  | Authorized representative in the European Community |
|  | UK Conformity Assessed |
|  | European Conformity Assessment |
|  | Manufacturer |

Bronidox[®] is the registered trade name of Cognis UK Ltd. ATCC[®] is a registered trademark of American Type Culture Collection.

Remel Europe Ltd
Clipper Boulevard West, Crossways
Dartford, Kent, DA2 6PT, UK
www.thermofisher.com

For technical assistance please contact your local distributor

| Version | Date of modifications introduced |
|---------|---|
| X7826B | January 2024 Updated to meet IVDR requirements |

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www.oxid.com/ifu • www.thermofisher.com

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7. Wanger, A.R., Morris, S.L., et al (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583-2588.
8. Roberts, J.I.S. and Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
9. Chabbert, Y.A. and Pillet, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus*. *Nature*, 213, 1137.
10. Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.
11. Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.
12. Henkel KGaA. Manufacturers Information and Safety Data Sheet for Bronidox™ L.

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| | | |
|-----|---------------------|-----|
| REF | ZL33/R30950102..... | 150 |
| | ZL34/R30950201..... | 450 |

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Key Code TSMX7826B

www.oxoid.com/ifu • www.thermofisher.com

Evropa + 800 135 79 135
Kanada 1 855 805 8539USA 1 855 236 0910
Zbytek světa +31 20 794 7071

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- Roberts, J.I.S. and Gaston, M.A.** (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
- Chabbert, Y.A. and Pillet, J.** (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus aureus*. *Nature*, 213, 1137.
- Nelles, M.J., Niswander, C.A., et al** (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.
- Fournier, J.M., Bouvet, A., et al** (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.
- Henkel KGaA.** Informace od výrobce a bezpečnostní list pro přípravek Bronidox™ L.

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8. Roberts, J.I.S. and Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
9. Chabbert, Y.A. and Pilllet, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus*. *Nature*, 213, 1137.
10. Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.
11. Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.
12. Henkel KGaA. Producentinformation og sikkerhedsdatablad for Bronidox® L.

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⁵ Piper, J., Hadfield, T., et al (1988). Efficacies of rapid agglutination tests for identification of Methicillin-Resistant Staphylococcal strains as *Staphylococcus aureus*. *J. Clin. Microbiol.*, 26, 1907–1909.

⁶ Ruane, P., Morgan, M.A., et al (1986). Failure of rapid agglutination methods to detect Oxacillin-Resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 24, 490–492.

⁷ Wanger, A.R., Morris, S.L., et al (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583–2588.

⁸ Roberts, J.I.S. und Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837–840.

⁹ Chabbert, Y.A. und Pillet, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus aureus*. *Nature*, 213, 1137.

¹⁰ Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14–18.

¹¹ Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932–1933.

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- Wanger, A.R., Morris, S.L., et al (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583-2588.
- Roberts, J.I.S. and Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
- Chabbert, Y.A. and Pillet, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus*. *Nature*, 213, 1137.
- Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.
- Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.
- Henkel KGaA. Manufacturers Information and Safety Data Sheet for Bronidox® L.

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Remel Europe Ltd
Clipper Boulevard West, Crossways
Dartford, Kent, DA2 6PT, HB
www.thermofisher.com

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- Wanger, A.R., Morris, S.L., et al (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583-2588.
- Roberts, J.I.S. and Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
- Chabbert, Y.A. and Pillot, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus*. *Nature*, 213, 1137.
- Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.
- Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.
- Henkel KGaA. Manufacturers Information and Safety Data Sheet for Bronidox® L.

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- Wanger, A.R., Morris, S.L., et al (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583-2588.
- Roberts, J.I.S. and Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
- Chabbert, Y.A. and Pilllet, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus*. *Nature*, 213, 1137.
- Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysacchande. *Infect. Immun.*, 49, 14-18.
- Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysacchande Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.
- Henkel KGaA. Valmistajan tiedot ja Bronidox[®] L-tuotteen käyttöohjeet.

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Painettu Iossa-Britanniassa

Code d'identification TSMX7826B
www.oxid.com/ifu • www.thermofisher.com

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
- Kloos, W.E. and Lambe, D.W. (1991). *Staphylococcus. Manual of Clinical Microbiology*, 5th Ed., edited by Balows, A., Hausler, W.J., Herrmann, K.L., Isenberg, H.D. and Shadomy, H.J. *American Society for Microbiology, Washington, D.C.* Pages 222-237.
- Selepak, S.T. and Witesky F.G. (1985). Inoculum size and lot-to-lot variation as significant variables in the tube coagulase test for *Staphylococcus aureus*. *J. Clin. Microbiol.*, 22, 835-837.
- Essers, L. and Raebold, K. (1980). Rapid and reliable identification of *Staphylococcus aureus* by a Latex Agglutination Test. *J. Clin. Microbiol.*, 12, 641-643.
- Orsi, A., Bartoloni, A., et al (1989). Evaluation of six different agglutination methods for rapid identification of *Staphylococcus aureus*. *Eur. J. Clin. Microbiol. Infect. Dis.*, 8, 655-656.
- Piper, J., Hadfield, T., et al (1988). Efficacies of rapid agglutination tests for identification of Methicillin-Resistant Staphylococcal strains as *Staphylococcus aureus*. *J. Clin. Microbiol.*, 26, 1907-1909.
- Ruane, P., Morgan, M.A., et al (1986). Failure of rapid agglutination methods to detect Oxacillin-Resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 24, 490-492.

- Wanger, A.R., Morris, S.L., et al (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583-2588.
- Roberts, J.I.S. and Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
- Chabbert, Y.A. and Pillet, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus aureus*. *Nature*, 213, 1137.
- Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.
- Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.
- Henkel KGaA. Manufacturers Information and Safety Data Sheet for Bronidox® L.

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^b inclui *S. saprophyticus*, *S. epidermidis* e *S. haemolyticus*.

14. BIBLIOGRAFIA

- Kloos, W.E. and Lambe, D.W. (1991). *Staphylococcus. Manual of Clinical Microbiology*, 5th Ed., edited by Balows, A., Hausler, W.J., Herrmann, K.L., Isenberg, H.D. and Shadomy, H.J. *American Society for Microbiology, Washington, D.C.* Pages 222-237.
- Selepak, S.T. and Witebsky F.G. (1985). Inoculum size and lot-to-lot variation as significant variables in the tube coagulase test for *Staphylococcus aureus*. *J. Clin. Microbiol.*, 22, 835-837.
- Essers, L. and Radebold, K. (1980). Rapid and reliable identification of *Staphylococcus aureus* by a Latex Agglutination Test. *J. Clin. Microbiol.*, 12, 641-643.
- Orsi, A., Bartoloni, A., et al (1989). Evaluation of six different agglutination methods for rapid identification of *Staphylococcus aureus*. *Eur. J. Clin. Microbiol. Infect. Dis.*, 8, 655-656.
- Piper, J., Hadfield, T., et al (1988). Efficacies of rapid agglutination tests for identification of Methicillin-Resistant Staphylococcal strains as *Staphylococcus aureus*. *J. Clin. Microbiol.*, 26, 1907-1909.
- Ruane, P., Morgan, M.A., et al (1986). Failure of rapid agglutination methods to detect Oxacillin-Resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 24, 490-492.
- Wanger, A.R., Morris, S.L., et al (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583-2588.
- Roberts, J.I.S. and Gaston, M.A. (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
- Chabbert, Y.A. and Pillet, J. (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus*. *Nature*, 213, 1137.

¹⁰ Nelles, M.J., Niswander, C.A., et al (1985). Reactivity of type-specific monoclonal antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.

¹¹ Fournier, J.M., Bouvet, A., et al (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.

¹² Henkel KGaA. Manufacturers Information and Safety Data Sheet for Bronidox® L.

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Remel Europe Ltd
Clipper Boulevard West, Crossways
Dartford, Kent, DA2 6PT, Reino Unido
www.thermofisher.com

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^b include *S. saprophyticus*, *S. epidermidis* și *S. haemolyticus*.

14. BIBLIOGRAFIE

- Kloos, W.E. and Lambe, D.W.** (1991). *Staphylococcus. Manual of Clinical Microbiology*, 5th Ed., edited by Balows, A., Hausler, W.J., Herrmann, K.L., Isenberg, H.D. and Shadomy, H.J. *American Society for Microbiology, Washington, D.C.* Pages 222-237.
- Selepak, S.T. and Witebsky F.G.** (1985). Inoculum size and lot-to-lot variation as significant variables in the tube coagulase test for *Staphylococcus aureus*. *J. Clin. Microbiol.*, 22, 835-837.
- Essers, L. and Radebold, K.** (1980). Rapid and reliable identification of *Staphylococcus aureus* by a Latex Agglutination Test. *J. Clin. Microbiol.*, 12, 641-643.
- Orsi, A., Bartoloni, A., et al** (1989). Evaluation of six different agglutination methods for rapid identification of *Staphylococcus aureus*. *Eur. J. Clin. Microbiol. Infect. Dis.*, 8, 655-656.
- Piper, J., Hadfield, T., et al** (1988). Efficacies of rapid agglutination tests for identification of Methicillin-Resistant Staphylococcal strains as *Staphylococcus aureus*. *J. Clin. Microbiol.*, 26, 1907-1909.
- Ruane, P., Morgan, M.A., et al** (1986). Failure of rapid agglutination methods to detect Oxacillin-Resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 24, 490-492.
- Wanger, A.R., Morris, S.L., et al** (1992). Latex agglutination-negative Methicillin-Resistant *Staphylococcus aureus* recovered from Neonates: epidemiologic features and comparison of typing methods. *J. Clin. Microbiol.*, 30, 2583-2588.
- Roberts, J.I.S. and Gaston, M.A.** (1987). Protein A and Coagulase expression in epidemic and non-epidemic *Staphylococcus aureus*. *J. Clin. Pathol.*, 40, 837-840.
- Chabbert, Y.A. and Pillet, J.** (1967). Correlation between "Methicillin Resistance" and serotype in *Staphylococcus*. *Nature*, 213, 1137.
- Nelles, M.J., Niswander, C.A., et al** (1985). Reactivity of type-specific monoclonal

antibodies with *Staphylococcus aureus* clinical isolates and purified capsular polysaccharide. *Infect. Immun.*, 49, 14-18.

¹¹ **Fournier, J.M., Bouvet, A., et al** (1987). Predominance of capsular polysaccharide Type 5 among oxacillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.*, 25, 1932-1933.

¹² **Henkel KGaA.** Manufacturers Information and Safety Data Sheet for Bronidox® L.

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Remel Europe Ltd
Clipper Boulevard West, Crossways
Dartford, Kent, DA2 6PT, Regatul Unit
www.thermofisher.com

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