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- Facklam, R.R., Cooksey, R.C., et al (1979).** Evaluation of Commercial Latex Agglutination Reagents for Grouping Streptococci. J. Clin. Microbiol., 10, 641.
- Facklam, R.R. and Smith, P.B. (1976).** The Gram Positive Cocci. Human Pathology, 7, 187.
- Fuller, A.T. (1938).** The Formamide Method for the Extraction of Polysaccharides from Haemolytic Streptococci. Brit. J. Exp. Path., 19, 130.
- Harvey, C.L. and McIlmurray, M.B. (1984).** Streptococci with dual antigen specificity for Lancefield groups D and G. Eur. J. Clin. Microbiol., 3, 526.
- Lancefield, R.C. (1938).** A Micro Precipitin-Technic for Classifying Hemolytic Streptococci, and Improved Methods for Producing Antisera. Proc. Soc. Exp. Biol., N.Y., 38, 473.
- Maxted, W.R. (1948).** Preparation of Streptococcal Extracts for Lancefield Grouping. Lancet, ii, 255.
- Medrek, T.F. and Barns, E.M. (1962).** The Influence of the Growth Medium on the Demonstration of a Group D Antigen in Faecal Streptococci. J. gen. Microbiol., 28, 701.
- Nowlan, S.S. and Deibel, R.H. (1967).** Group Q Streptococci. 1. Ecology, Serology, Physiology, and Relationship to Established Enterococci. J. Bact., 94, 291.
- Parker, M.T. and Ball, L.C. (1976).** Streptococci and Aerococci Associated with Systemic Infection in Man. J. Med. Microbiol., 9, 275.
- Rantz, L.A. and Randall, E. (1955).** Use of Autoclaved Extracts of Hemolytic Streptococci for Serological Grouping. Stanford Med. Bull., 13, 290.
- Watson, B.K., Moellering, R.C., et al (1975).** Identification of Streptococci: Use of Lysozyme and Streptomyces albus Filtrate in the Preparation of Extracts for Lancefield Grouping. J. Clin. Microbiol., 1, 274.

19. Data on file

*Marque commerciale.

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Remel Europe Ltd, Clipper Boulevard West,
Crossways, Dartford, Kent DA2 6PT, UK

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