

Customer Notification

ECR1038594

Date: 16 October 2017

Dear Valued Customer,

This notification is to provide information regarding a change to our Gibco serum products.

Description of Change:

Change in the label and COA storage temperature range on our Gibco serum labels from -5 to -20°C to $\leq$ -10°C.

Purpose of Change:

This change is to align with current United States Pharmacopeia and European Pharmacopeia requirements for serum storage.

Detailed Description of Change:

The product label and associated COAs for our Gibco serum products will transition to the new storage temperature range starting in late November 2017, for product to be released late 2017. While we will begin to make these changes for product manufactured in November, please be aware that there will be a mix between existing and new packaging over several months as we exhaust our existing inventory. These changes will be applied in conjunction with the branding and intended use changes also affecting product labels and COAs notified earlier this year.

This change will impact the storage temperature range printed on COAs where applicable and product labels however it does not reflect a change to current and actual conditions for serum storage while in our facilities. Additionally this does not impact our shipping procedures which ensure product remains frozen during transit.

The above change will not impact our product quality, stability or performance and there will be no changes in the manufacturing, storage or transport processes.

Finished Products Impacted:

All Gibco Serum products.

If you have questions regarding this notification, please contact us at your convenience.

Email techsupport@thermofisher.com

In US/Canada call 800-955-6288, opt 5 ext. 40562

In Europe call 00 800 5345, option 3

For all other regions, please refer to our website for your local phone number

Best regards,



Janna Sanford
Quality Systems

Gibco sera (Catalog) SKU's impacted

For questions about custom SKU's, please contact us by email at techsupport@thermofisher.com

ITEM NUMBER	DESCRIPTION
10082139	FBS, CERT, HI, USA ORIGIN 100ML
10082147	FBS, CERT, HI, USA ORIGIN 500ML
10083129	FBS CERTIFIED - GAMMA IRRAD 1000ML
10083137	FBS CERTIFIED - GAMMA IRRAD 100ML
10083145	FBS CERTIFIED - GAMMA IRRAD 500ML
10086148	FBS QUAL USA ORIGIN G/IRRAD 500ML
10091130	FETAL BOVINE SERUM NZ ORIGIN 100ML
10091148	FETAL BOVINE SERUM NZ ORIGIN 500ML
10091155	FETAL BOVINE SERUM NZ ORIGIN 1L
10093136	FBS HI NZ ORIGIN 100ML
10093177	FBS HI NZ ORIGIN 500ML
10093188	FBS HI NZ ORIGIN
10094134	FBS GAMMA IRRAD N.ZEALAND 100ML
10094142	FBS GAMMA IRRAD NZ ORIGIN 500ML
10094159	FBS GAMMA IRRAD N.ZEALAND 1000ML
10099133	FBS QUALIFIED AUSTRALIA ORIGIN 100ML
10099141	FBS QUALIFIED AUSTRALIA ORIGIN 500ML
10099158	FBS QUALIFIED AUSTRALIA ORIGIN 1L
10100139	FBS HI AUSTRALIAN ORIGIN 100ML
10100147	FBS HI AUSTRALIAN ORIGIN 500ML
10100154	FBS HI AUSTRALIAN ORIGIN 1L
10101137	FBS GAMMA IRRAD AUSTRALIAN 100ML
10101145	FBS GAMMA IRRAD AUSTRALIAN 500ML
10101152	FBS GAMMA IRRAD AUSTRALIAN 1L
10270098	FBS SOUTH AMERICAN (CE) 100ML
10270106	FBS SOUTH AMERICAN (CE) 500ML
10371011	DONOR BOVINE SERUM WITH IRON
10371029	DONOR BOVINE SERUM WITH IRON 500ML
10371045	DONOR BOVINE SERUM WITH IRON
10437010	FBS 100ML
10437028	FBS, USDA APPD. ORIGIN 500ML
10437036	FBS 1000ML
10438018	FBS, HI 100ML
10438026	FBS, HI 500ML
10438034	FBS, HI 1000ML
10439001	FBS, ES CELL QUAL 50ML

10439016	FBS, ES CELL QUAL 100ML
10439024	FBS, ES CELL QUAL. USDA APPD. 500ML
10493106	FBS SOUTH AMERICAN HI & GI 500ML
10499044	FBS GAMMA IRRAD S.AMERICAN(CE) 500ML
10500056	FBS HEAT INACT. S.AMERICAN(CE) 100ML
10500064	FBS HEAT INACT. S.AMERICAN(CE) 500ML
10509024	FBS ULTRA LOW IGG N.AMER G/I 500ML
10603017	FBS DIALYSED GAMMA IRAD N.AM 500ML
12662002	MSC QUALIFIED FBS (USDA) 50ML
12662011	MSC QUALIFIED FBS 100ML
12662029	MSC QUALIFIED FBS 500ML
12664017	MSC QUALIFIED AUSTRALIAN FBS 100ML
12664025	MSC QUALIFIED AUSTRALIAN FBS 500ML
12763017	MSC QUALIFIED FBS(US ORIGIN) 100ML
12763025	MSC QUALIFIED FBS(US ORIGIN) 500ML
16000036	FBS, CERT, USA ORIGIN 100ML
16000044	FBS, CERT, USA ORIGIN 500ML
16000069	FBS, CERT, USA ORIGIN 1000ML
16010142	NEWBORN CALF SERUM 1L
16010159	NEWBORN CALF SERUM 500ML
16010167	NEWBORN CALF SERUM 100ML
16030033	DONOR BOVINE SERUM 100ML
16030074	DONOR BOVINE SERUM 500ML
16030108	DONOR BOVINE SERUM 1L
16030273	DONOR BOVINE SERUM
16050114	HORSE SERUM 1L
16050122	HORSE SERUM 500ML
16050130	HORSE SERUM 100ML
16070096	LAMB SERUM 500ML
16110082	CHICKEN SERUM 500ML
16120099	RABBIT SERUM 100ML
16120107	RABBIT SERUM 500ML
16125015	RABBIT SERUM HEAT INACTIVATED 100ML
16125023	RABBIT SERUM HEAT INACTIVATED 500ML
16140063	FBS, HI, QUALIFIED, USA ORIGIN 100ML
16140071	FBS, HI 500ML
16140089	FBS, HI 1000ML
16141002	ES QUALIFIED FBS 50ML
16141061	FBS, ES CELL QUAL 100ML

16141079	FBS, ES CELL QUAL, USA ORIGIN 500ML
16170060	BOVINE SERUM 1L
16170078	BOVINE SERUM 500ML
16170086	BOVINE SERUM 100ML
16175028	DONOR BOVINE SERUM HEAT INACT 500ML
16210064	GOAT SERUM 100ML
16210072	GOAT SERUM 500ML
16250078	FBS, ULTRA LOW IGG, USA ORIGIN 500ML
16250086	FBS, ULTRA LOW IGG, USA ORIGIN 100ML
26010066	NEWBORN CALF SERUM HEAT INACT 100ML
26010074	NEWBORN CALF SERUM HEAT INACT 500ML
26050070	HORSE SERUM HEAT INACTIVATED 100ML
26050088	HORSE SERUM HEAT INACTIVATED 500ML
26140079	FBS QUALIFIED, USA ORIGIN 500ML
26140087	FBS QUALIFIED, USA ORIGIN 100ML
26140095	FBS 1000ML
26170035	BOVINE SERUM HEAT INACTIVATED 100ML
26170043	BOVINE SERUM HEAT INACTIVATED 500ML
26171025	DONOR BOVINE SERUM GAMMA IRRAD 500ML
26171033	DONOR BOVINE SERUM GAMMA IRRAD 1L
26250084	PORCINE SERUM 500ML
26250095	PORCINE SERUM 100ML
26400036	FBS, DIALYZED, USA ORIGIN 100ML
26400044	FBS, DIALYZED, USA ORIGIN 500ML
30044184	FBS ES CELL QUAL NZ ORIGIN 100ML
30044333	FBS ES CELL QUAL NZ ORIGIN 500ML
30048185	FBS HI GAMMA IRRAD AUSTRALIAN 100ML
30048334	FBS HI GAMMA IRRAD AUSTRALIAN 500ML
30067185	FBS DIALYZED ORIGIN: NZ 100ML
30067334	FBS DIALYZED ORIGIN: NZ 500ML
66050022	HORSE SERUM-GAMMA RAD 500ML
1921005PG	FBS ULTRA-LOW IGG ORIGIN: NZ 100ML
1921005PJ	FBS ULTRA-LOW IGG ORIGIN: NZ 500 ML