

gibco



Custom product manufacturing

Meeting your media production needs, from
bench to full-scale production

ThermoFisher
SCIENTIFIC

Custom product manufacturing

Activity flow

Gibco™ media production capabilities provide flexibility and high-quality standards to meet your media production needs. Our pilot labs employ the same technology found

in our full-scale media manufacturing facility, allowing for seamless transfers and assurance that your quality requirements will be met at each step in the production process.

Non-GMP pilot manufacturing

Timeline: 10 business days*

Step 1

Customer confirmation of pilot manufacturing specifications

Manufacturing specifications help ensure that your media production needs are met.

Your medium is now scheduled for production.

Step 2

Medium manufacturing

A pilot is produced per your manufacturing specifications.

Your medium is ready for testing.

Step 3

Pilot testing

Standard release testing is performed and a Memorandum of Assay is prepared.

Your medium is ready to be packaged.

Step 4

Packaging

The product is packaged and labeled per your specifications.

Your medium is ready to be shipped.

Step 5

Shipping

The product is shipped.

The delivery timeline is dependent on your ship-to location.



Full-scale GMP manufacturing

Timeline: 16–20 weeks**

Step 7

Prepare for tech transfer

Our process development team prepares a technical transfer record from the pilot activity.



Step 8

Tech transfer

A cross-functional team implements the process and activates products for full-scale manufacturing.



Step 9

Full-scale manufacturing and testing

The medium is produced and tested per agreed-upon requirements.

Your full-scale medium is ready to be shipped.



Step 10

Shipping

The final product is shipped per your specifications.

Step 6



Customer evaluation and approval

Your pilot medium is tested and approved by you.

Then, full-scale manufacturing can begin.

* The timeline listed is for standard materials and testing per batch. Manufacturing timelines may be extended for formulations containing animal origin components and/or customer formulations containing new raw materials that require an additional procurement timeline.

** When leveraging suppliers qualified by Thermo Fisher Scientific.

gibco

Find out more at
[thermofisher.com/advbio](https://www.thermofisher.com/advbio)

ThermoFisher
S C I E N T I F I C

For Research Use Only. Not for use in diagnostic procedures. © 2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.
COL23243 0319