

SENSITITRE™ JUSTONE™ PLATE FORMAT

Plate Code: **JOEYE2**

	1	2	3	4	5	6	7	8	9	10	11	12
A	POS	CIP	CIP	CIP	MXF	MXF	POS	CIP	CIP	CIP	MXF	MXF
		1	2	4	0.5	1		1	2	4	0.5	1
B	ERY	ERY	ERY	ERY	NEO	NEO	ERY	ERY	ERY	ERY	NEO	NEO
	0.5	1	2	4	4	8	0.5	1	2	4	4	8
C	OXY	OXY	OXY	OXY	AMI	AMI	OXY	OXY	OXY	OXY	AMI	AMI
	0.5	1	2	4	16	32	0.5	1	2	4	16	32
D	GEN	GEN	GEN	TIC	TIC	TIC	GEN	GEN	GEN	TIC	TIC	TIC
	2	4	8	16	32	64	2	4	8	16	32	64
E	CHL	CHL	CHL	TOB	TOB	TOB	CHL	CHL	CHL	TOB	TOB	TOB
	4	8	16	4	8	16	4	8	16	4	8	16
F	OFL	OFL	OFL	OFL	POL	POL	OFL	OFL	OFL	OFL	POL	POL
	0.12	0.25	0.5	1	5	10	0.12	0.25	0.5	1	5	10
G	BAC	BAC	XNL	XNL	FAZ	FAZ	BAC	BAC	XNL	XNL	FAZ	FAZ
	2	4	2	4	8	16	2	4	2	4	8	16
H	DOX	DOX	DOX	DOX	DOX	SXT	DOX	DOX	DOX	DOX	DOX	SXT
	0.12	0.25	0.5	1	2	2/38	0.12	0.25	0.5	1	2	2/38

ANTIMICROBICS

POS	Positive Control
ERY	Erythromycin
OXY	Oxytetracycline
GEN	Gentamicin
CHL	Chloramphenicol
OFL	Ofloxacin
BAC	Bacitracin
DOX	Doxycycline
CIP	Ciprofloxacin
XNL	Ceftiofur
TIC	Ticarcillin
TOB	Tobramycin
MXF	Moxifloxacin
NEO	Neomycin
AMI	Amikacin
POL	Polymixin B
FAZ	Cefazolin
SXT	Trimethoprim / sulfamethoxazole