

# CultiMaxx shelving systems

Optimized shelving systems for Thermo Scientific™ Forma™ Steri-Cult™ CO<sub>2</sub> Incubators to accommodate large vessels







Stainless steel shelving



Retrofittable in the field



Easy to install

Description	Vessel	No. of vessels/ incubator	Advantages compared to standard shelves	Max. weight capacity (incubator)	Max. weight capacity per shelf	Cat. No.
Thermo Scientific™ CultiMaxx™ shelving system for Forma Steri-Cult CO <sub>2</sub> Incubators	 Cell culture bags	Up to 154 1 L bags*	<ul style="list-style-type: none"> <li>Up to 2.5x more culture bag capacity per unit</li> <li>Increased exposure to promote better gas exchange with rack-style design for a variety of bags</li> <li>Improved accessibility to individual bags with extending shelves</li> <li>Base set includes 8 shelves; additional shelves can be ordered (Cat. No. 332474G02)</li> </ul>	50 kg / 110 lb	3 kg / 6 lb	332474G01
CultiMaxx shelving system for Forma Steri-Cult CO <sub>2</sub> Incubators	 Thermo Scientific™ Nunc™ Cell Factory™ systems with 10 or 13 layers	Up to 9	<ul style="list-style-type: none"> <li>Cascade design offers improved stacking to allow for filter clearance on Nunc Cell Factory systems or similar vessels</li> <li>Reinforced to minimize deflection under maximum load</li> <li>Individual vessel access</li> </ul>	51 kg / 112 lb	17 kg / 37 lb	332239G01
CultiMaxx shelving system for Forma Steri-Cult CO <sub>2</sub> Incubators	 Nunc Cell Factory systems with 4 layers	Up to 15	<ul style="list-style-type: none"> <li>Cascade design offers improved stacking to allow for filter clearance on Nunc Cell Factory systems or similar vessels</li> <li>Reinforced to minimize deflection under maximum load</li> <li>Up to 25% increased capacity and individual vessel access</li> </ul>	40 kg / 88 lb	8 kg / 17 lb	332236G01
CultiMaxx shelving system for Forma Steri-Cult CO <sub>2</sub> Incubators	 G-Rex™ 500M-CS system	Up to 9	<ul style="list-style-type: none"> <li>Up to 50% more capacity</li> <li>Redesigned for greater load bearing</li> <li>Easy-to-clean reinforced design with wipeable surfaces</li> </ul>	60 kg / 132 lb	20 kg / 44 lb	332342G01

\* This is a single scenario. Max. configuration may vary based on brand of bag used and application needs.



Learn more at [thermofisher.com/cultimaxx](https://thermofisher.com/cultimaxx)

thermo scientific