

Comments for pPICZ A
3329 base pairs

5' AOX1 promoter region: bases 1-941
 5' end of AOX1 mRNA: base 824
 5' AOX1 priming site: bases 855-875
 Multiple cloning site: bases 932-1011
 c-Myc epitope tag: bases 1012-1044
 6xHis tag: bases 1057-1077
 3' AOX priming site: bases 1159-1179
 3' end of mRNA: base 1250
 AOX1 transcription termination region: bases 1078-1418
 Fragment containing TEF1 promoter: bases 1419-1830
 EM7 promoter: bases 1831-1898
 Zeocin™ resistance ORF: bases 1899-2273
 CYC1 transcription termination region: bases 2274-2591
 pUC origin: bases 2602-3275 (c)

(c) = complementary strand

* The restriction site between the *Not* I site and the Myc epitope is different in each version of pPICZ:

Apa I in pPICZ A
Xba I in pPICZ B
*Sna*B I in pPICZ C

