

# Freezer<sup>PRO</sup>

From Brooks LIFE SCIENCES

## Sample Management System



### FreezerPro®

FreezerPro® is an intuitive, fast, reliable and secure solution which enables users to know precisely where a laboratory sample is located even before opening the freezer door.

FreezerPro is a biological sample management inventory system with over 250 installations globally. It provides laboratories with a secure web based application for managing samples.

Organisations looking for a commercial off the shelf solution for freezer inventory management, with no expensive setup or implementation costs, as well as central biobanks with millions of records, have relied on FreezerPro® for the secure management of samples and sample information. The FreezerPro® Cloud uses state-of-the-art, cross-region servers and storage technologies managed directly by Brooks Life Science Systems. Customers are able to access FreezerPro Cloud via a web browser from anywhere in the world, at any time.

### Compatibility

In addition to being compatible with PC/MAC, Tablet, Smartphone and web-based operation, FreezerPro is also compatible with many barcoding and RFID hardware scanners, readers and printers.

### Compliance

- ✔ FDA 21 CFR Parts 11, 21, 58, 210, 211, 820
- ✔ cGLP/cGMP (Good Laboratory Practice and Good Manufacturing Practice)
- ✔ Section 508 Accessibility Guidelines
- ✔ HIPAA Accountability Act of 1996
- ✔ HIPAA Patient Safety and Quality Improvement Act of 2005 (Patient Safety Act)

Automated  
Storage Systems

Cryopreservation &  
Cold Chain Solutions

Informatics &  
Technical Solutions

Sample Storage,  
Lab Services & Transport

Storage Consumables  
& Instruments

Learn more – [www.brookslifesciences.com](http://www.brookslifesciences.com)

Contact us – [www.brookslifesciences.com/contact-us](http://www.brookslifesciences.com/contact-us)

© Copyright 2018 Brooks Automation, Inc. B601-17



## Features

Features	Basic	Standard	Enterprise	Elite
Truly web based - install once, access from anywhere	✓	✓	✓	✓
Complete sample management	✓	✓	✓	✓
Track sample in and out of freezers	✓	✓	✓	✓
Icon based display of box contents	✓	✓	✓	✓
Sample inventory alerts	✓	✓	✓	✓
Easy to follow workflow	✓	✓	✓	✓
Virtual freezers	✓	✓	✓	✓
Reporting and data export/import - HTML, CSV, PDF, TXT formats	✓	✓	✓	✓
Live search, instant search by word or part of the keyword	✓	✓	✓	✓
Smart search	✓	✓	✓	✓
Live Query builder and smart folders	✓	✓	✓	✓
Full audit trails	✓	✓	✓	✓
Mobile devices support	✓	✓	✓	✓
Multi-language support		✓	✓	✓
Comprehensive barcode label and RFID tag designer		✓	✓	✓
2D barcode plate scanners support (FluidX, Visionmate, Biosero, Ziath)		✓	✓	✓
Support for sample sources		✓	✓	✓
Support for sample groups		✓	✓	✓
Support for pick lists		✓	✓	✓
Report templates (visual report builder)		✓	✓	✓
Support for freezer sensors			✓	✓
Permission based access restriction to sample locations and information			✓	✓
Available shipping module			✓	✓
Available workflow module			✓	✓
Available RFID kit module			✓	✓
GLP mode for strict operations			✓	✓
Support for custom vials identifiers			✓	✓
SSL data encryption			✓	✓
HIPAA compliance				✓
User authentication via active directory domain services (LDAP)				✓
Customer postgresSQL, MYSQL, Microsoft SQL, Oracle Database support				✓
Application programming interface (API)				✓
Hyper-V support on Microsoft Windows Server				✓
Advanced workflow management				✓
Shipping module included				✓
Email tied communications and alerts				✓
Clinical data management				✓
Advanced reporting suite				✓
Configurable charting and graphing				✓

Automated  
Storage Systems

Cryopreservation &  
Cold Chain Solutions

Informatics &  
Technical Solutions

Sample Storage,  
Lab Services & Transport

Storage Consumables  
& Instruments

### Authorised distributor

#### In Australia:

For customer service, call 1300-735-292  
To fax an order, use 1800-067-639  
To email an order, [ordersau@thermofisher.com](mailto:ordersau@thermofisher.com)

#### In New Zealand:

For customer service, call 0800-933-966  
To fax an order, use 0800-329-246  
To email an order, [ordersnz@thermofisher.com](mailto:ordersnz@thermofisher.com)

**ThermoFisher**  
SCIENTIFIC