

40cc Sterile Mucus Specimen Trap for Pneumatic Tube Transport

- Eliminate the possibility of accidental tube separation specimen transport with this sterile mucus specimen trap
- Packaged with an extra cap designed to seal the trap for transport
- After collection, simply remove the standard tubed cap and replace it with the transport cap
- Each packed sterile, in easy-to-peel envelope that assures aseptic handling
- Write-on bottle

Product	Catalogue No.	Packaging
40cc Sterile Mucus Specimen Trap for Pneumatic Tube Transport	BSE405	24/Carton



40cc Sterile Mucus Specimen Traps

- See-through plastic container with gradations clearly marked
- Disposable adaptor prevents contamination from non-sterile suction tube
- Screw-on cap protects lips - helps maintain sterile pouring surface
- Each packed sterile, in easy-to-peel envelope that assures aseptic handling
- Write-on bottle

Product	Catalogue No.	Packaging
40cc Sterile Mucus Specimen Trap	BSE406	24/Carton

80cc Sterile Mucus Specimen Traps

- See-through plastic container with gradations clearly marked
- Disposable adaptor prevents contamination from non-sterile suction tube
- Screw-on cap protects lips - helps maintain sterile pouring surface
- Collapse-proof, 5" long heavy wall non-latex tubing
- Write-on bottle

Product	Catalogue No.	Packaging
80cc Sterile Mucus Specimen Trap	BSE407	25Carton

