










cerillo

co-culture duet system



features

-  observe and study individual microbial populations in real-time
-  measure the growth of one microbial population alongside other populations
-  porous barrier physically isolates populations while maintaining fluidic contact
-  multiplexed plate accommodates 18 duets
-  physically accessible wells enable population sampling during or after experiment
-  compatible with automation systems

available products

co-culture duet system



easily run comprehensive bacterial interaction studies



save valuable time and resources with this accessible, off-the-shelf solution



compare experiments and facilitate collaborations with standardized system



scale experiments more easily with automation compatibility

**with the co-culture
duet system, expand
lab capabilities while
improving accessibility,
standardization and
scalability in your bacterial
interaction studies**

All cerillo products are for research use only and are not intended for use in diagnostic procedures.

To learn more about Cerillo and our products, visit us at www.cerillo.bio or by email at info@cerillo.bio

Authorised distributor

In Australia:

For customer service, call 1300-735-292
To email an order, ordersanz@thermofisher.com
To order online: thermofisher.com

In New Zealand:

For customer service, call 0800-933-966
To email an order, ordersanz@thermofisher.com
To order online: thermofisher.com

ThermoFisher
SCIENTIFIC