

CBS™ High Security Vitrification straw*

*CE marked / FDA 510(k) premarket notification clearance # K092398 for 4-8 cell and blastocyst embryos

Exclusive High Security closed system
for oocyte and embryo vitrification!



Highest level of safety and quality in vitrification


 **CryoBio System**
Groupe I.M.V. Technologies

www.cryobiosystem-imv.com

CBS™ High Security Vitrification straw*

*CE marked / FDA 510(k) premarket notification clearance # K092398 for 4-8 cell and blastocyst embryos

1. High Security ultra-thin outer straw
Pre-sealed weighted end and flared opening
2. Capillary tube with pre-formed fully clear gutter
co-injected with its colored identification rod
3. Blue introduction device



RELIABILITY

Thermo-sealed ionomeric resin straw



SYMS III®

- SECURITY Straightforward simple use
- EFFICIENCY Safe and fast procedure
- ONE-STEP OPENING AND RETRIEVAL



- INTERNATIONAL STANDARDS COMPLIANCE
CE, FDA 510(k) cleared

- TRACEABILITY Color coding and LN2 resistant labels for easy storage and identification



- EASY STORAGE
Compatible with CBS™ storage concept



Authorised distributor

In Australia:

For customer service, call 1300-735-292
To fax an order, use 1800-067-639
To email an order, ordersau@thermofisher.com

In New Zealand:

For customer service, call 0800-933-966
To fax an order, use 0800-329-246
To email an order, ordersnz@thermofisher.com

ThermoFisher
SCIENTIFIC