



# DURAN® PURE

## DEVELOPED FOR PHARMA

**DURAN® PURE** is a range of high-quality bottles and closures developed for the needs of the pharmaceutical industry. **DURAN® PURE** meets the requirements of the international pharmacopeias and is available with the certificates required by the regulatory authorities.

### FEATURES

- **EP/USP Type I borosilicate glass**
- **Size range 25 ml - 20 l**
- **Protective cover** prevents contamination
- **Single batch delivery** CoA available
- **100% traceability**
- **100% visual inspection**
- **Change control**
- **Drug Master File**

---

### APPLICATION AREAS

- Manufacturing of pharmaceutical products
- Storage of bulk drug intermediates
- Transportation of APIs



## DURAN® PURE Bottles

Article number clear glass	Article number amber coloured	Volume [ml]	DIN thread GL	Bottles per box [pcs.]
81 801 14 01	81 806 14 06	25	25	10
81 801 24 06	81 806 24 02	100	45	10
81 801 36 08	81 806 36 04	250	45	10
81 801 44 07	81 806 44 03	500	45	10
81 801 54 03	81 806 54 08	1 000	45	10
81 801 63 05	81 806 63 01	2 000	45	10
81 801 73 01	81 806 73 06	5 000	45	1
81 801 86 06	81 806 86 02	10 000	45	1
81 801 91 05	81 806 91 01	20 000	45	1

**ALSO AVAILABLE: DURAN® PURE GLS 80®** bottles and screw caps

## DURAN® PURE Premium Screw Caps

DURAN® PURE Premium Screw Caps made from PFA (fluoropolymer) are chemically and thermally very resistant and compatible with autoclaving and other common sterilization methods. The products include silicone seals which are coated with highly inert PTFE.

### FEATURES

- Caps made from PFA with PTFE coated silicone seal
- Produced and packaged in clean rooms
- Single batch delivery CoA available
- 100% traceability
- Change control
- Drug Master File



Article number	DIN thread GL	Description	Screw Caps per PE bag [pcs.]
89 108 13 06	25	DURAN® PURE Premium Screw Cap	5
89 112 97 05	25	DURAN® PURE Premium Screw Cap, <b>supplied without cap liner.</b> Must be used with 292481456 for full functionality	5
29 248 14 56	25	DURAN® PURE Silicone Cap Liner, PTFE coated	10
89 108 86 79	45	DURAN® PURE Premium Screw Cap	5
89 886 28 02	45	DURAN® PURE Premium Lip Seal Cap	5
89 115 31 38	45	DURAN® PURE Premium Screw Cap, <b>supplied without cap liner.</b> Must be used with 292482958 for full functionality	5
29 248 29 58	45	DURAN® PURE Silicone Cap Liner, PTFE coated	10
89 108 86 78	45	DURAN® PURE Premium Pouring Ring	5

001 / EN 05-2022



Contact us at [DWK.com](https://www.dwk.com)

© 2022 DWK Life Sciences  
The trademarks used are owned by DWK Life Sciences GmbH,  
DWK Life Sciences LLC or DWK Life Sciences Ltd.

## Authorised distributor

---

### **In Australia:**

For customer service, call 1300-735-292

To email an order, [ordersanz@thermofisher.com](mailto:ordersanz@thermofisher.com)

To order online: [thermofisher.com](http://thermofisher.com)

### **In New Zealand:**

For customer service, call 0800-933-966

To email an order, [ordersanz@thermofisher.com](mailto:ordersanz@thermofisher.com)

To order online: [thermofisher.com](http://thermofisher.com)

**ThermoFisher**  
S C I E N T I F I C