



# DURAN® PURE GLS 80®

## DEVELOPED FOR PHARMA

**DURAN® PURE** is a range of high-quality bottles and closures developed for the needs of the pharmaceutical industry.

**DURAN® PURE** meets the requirements of the international pharmacopeias and is available with the certificates required by the regulatory authorities.

### FEATURES

- EP/USP Type I borosilicate glass
- Size range 500 ml - 20 l
- **GLS 80®** thread for fast closing & opening
- **Protective cover** prevents contamination
- **Single batch delivery** CoA available
- **100% traceability**
- **100% visual inspection**
- **Change control**
- **Drug Master File**

### APPLICATION AREAS

- Manufacturing of pharmaceutical products
- Storage of bulk drug intermediates
- Transportation of APIs



## DURAN® PURE GLS 80® Bottles

Article number clear glass	Article number amber coloured	Volume [ml]	Bottles per box [pcs.]
81 860 44 03	81 866 44 09	500	10
81 860 54 08	81 866 54 05	1 000	10
81 860 63 01	81 866 63 07	2 000	10
81 860 73 06	81 866 73 03	5 000	1
81 860 86 02	81 866 86 08	10 000	1
81 860 91 01	81 866 91 07	20 000	1

**ALSO AVAILABLE:** DURAN® PURE bottles and screw caps with GL 45 thread

## DURAN® PURE GLS 80® Screw Caps

DURAN® PURE GLS 80® Screws Caps made from PSU are chemically and thermally very resistant and compatible with autoclaving and other common sterilization methods. The product includes a silicone seal which is coated with highly inert PTFE.

### FEATURES

- Cap made from PSU with PTFE coated silicone seal
- GLS 80® thread for fast closing & opening
- Produced and packaged in clean rooms
- Single batch delivery CoA available
- 100% traceability
- Change control



Article number	Description	Screw Caps per PE bag [pcs.]
89 116 91 09	DURAN® PURE GLS 80® Screw Cap	1
89 116 91 82	DURAN® PURE GLS 80® Screw Cap, <b>supplied without cap liner.</b> Must be used with 292489102 for full functionality	1
29 248 91 02	DURAN® PURE Silicone Cap Liner, PTFE coated	10

001 / EN 05-2022



Contact us at [DWK.com](https://www.dwk.com)

© 2022 DWK Life Sciences  
The trademarks used are owned by DWK Life Sciences GmbH,  
DWK Life Sciences LLC or DWK Life Sciences Ltd.

## Authorised distributor

---

### **In Australia:**

For customer service, call 1300-735-292  
To email an order, [ordersanz@thermofisher.com](mailto:ordersanz@thermofisher.com)  
To order online: [thermofisher.com](http://thermofisher.com)

### **In New Zealand:**

For customer service, call 0800-933-966  
To email an order, [ordersanz@thermofisher.com](mailto:ordersanz@thermofisher.com)  
To order online: [thermofisher.com](http://thermofisher.com)

**ThermoFisher**  
S C I E N T I F I C