

HandPRO™

Cleanroom Nitrile & Vinyl Gloves

An Exceptional Glove for Every Cleanroom Need

Correct glove choice for your critical environment is imperative. HandPRO brand nitrile, latex and vinyl cleanroom gloves provide the protection you need – for your employees and your product. Manufactured to **rigorous standards**, HandPRO gloves prevent contamination and provide reliable protection.

ISO 4 | ULTRA CLEAN CLASS 10 COMPATIBLE

HandPRO 9300 NEW Cross Linking Formulation

- 100% Nitrile
- Accelerator-Free
- Sulfur-Free
- Anti-Fatigue
- Enhanced Strength
- Anti-Slip Grip
- 300mm (12")

Size	Hourglass Item No.	Fisher Scientific Cat. No.
XS	9300	17-900-900
S	9301	17-900-901
M	9302	17-900-902
L	9303	17-900-903
XL	9304	17-900-904
XXL*	9305	17-900-905

100 Pcs/Bag 10 Bags/Case
* 80 Pcs/Bag



ISO 5 | CLEAN CLASS COMPATIBLE

HandPRO 9100 NEW X-Linking Formula

- 100% Nitrile
- Accelerator-Free
- Sulfur-Free
- Enhanced Strength
- Anti-Fatigue
- Anti-Slip Grip
- 300mm (12")



HandPRO 7100 (12") HandPRO 6100 (16") Trusted Formulation

- 100% Nitrile
- Powder-Free
- Finger Textured
- White Color



HandPRO 2400 Strong & Clean Latex

- 100% Natural Rubber Latex
- Double Chlorinated
- Natural Color
- Fully Textured
- 300mm (12")



Size	Hourglass Item No.	Fisher Scientific Cat. No.
XS	9100	19-423-100
S	9101	19-423-101
M	9102	19-423-102
L	9103	19-423-103
XL	9104	19-423-104
XXL*	9105	19-423-105

100 Pcs/Bag 10 Bags/Case
*80 Pcs/Bag

Size	Hourglass Item No.	Fisher Scientific Cat. No.
XS	7100	19-423-023
	6100	19-423-050
S	7101	19-423-024
	6101	19-423-051
M	7102	19-423-025
	6102	19-423-052
L	7103	19-423-026
	6103	19-423-053
XL	7104	19-423-027
	6104	19-423-054

SERIES 7100: 100 Pcs/Bag; 10 Bags/Case
SERIES 6100: 50 Pcs/Bag; 10 Bags/Case

Size	Hourglass Item No.	Fisher Scientific Cat. No.
XS	2400	19-423-106
S	2401	19-423-107
M	2402	19-423-108
L	2403	19-423-109
XL	2404	19-423-110
XXL*	2405	19-423-111

100 Pcs/Bag 10 Bags/Case
*80 Pcs/Bag

ISO 6 | CONTROLLED ENVIRONMENT COMPATIBLE

HandPRO 1700 Strength, Quality, Comfort

- CE Bagged
- 100% Nitrile
- High Tactile Sense
- Micro-Textured
- Powder-Free
- White Color
- 240mm (9")



HandPRO 1750 Accelerator-Free Nitrile

- CE Bagged
- Enhanced Strength
- Secure Wet Grip
- 100% Nitrile
- Sulfur-Free
- Anti-Fatigue
- 240mm (9")



Size	Hourglass Item No.	Fisher Scientific Cat. No.
XS	1700	19-423-029
S	1701	19-423-030
M	1702	19-423-031
L	1703	19-423-032
XL	1704	19-423-033
100 Pcs/Bag		10 Bags/Case

Size	Hourglass Item No.	Fisher Scientific Cat. No.
XS	1750	17-500-100
S	1751	17-500-101
M	1752	17-500-102
L	1753	17-500-103
XL	1754	17-500-104
XXL*	1755	17-500-105
100 Pcs/Bag		10 Bags/Case
*80 Pcs/Bag		

CLEANROOM VINYL

HandPRO 1550 ISO 5 | Clean Class Vinyl

- 100% Clear Vinyl
- Consistent Quality
- Powder-Free
- Smooth
- 300mm (12")



HandPRO 550 ISO 6 | CE Vinyl

- 100% Clear Vinyl
- Consistent Quality
- Powder-Free
- Smooth
- 240mm (9")



Size	Hourglass Item No.	Fisher Scientific Cat. No.
S	1551	19-423-034
M	1552	19-423-035
L	1553	19-423-036
XL	1554	19-423-037
100 Pcs/Bag		10 Bags/Case

Size	Hourglass Item No.	Fisher Scientific Cat. No.
S	551	18-040-1266
M	552	18-040-1267
L	553	18-040-1268
XL	554	18-040-1269
100 Pcs/Bag		10 Bags/Case



Trademarks used are owned as indicated at www.fishersci.com/trademarks

Distributed in Australia and New Zealand by:

ThermoFisher
SCIENTIFIC

Australia: For customer service, call 1300-735-292
Visit us online at: www.thermofisher.com.au
New Zealand: For customer service, call 0800-933-966
Visit us online at: www.thermofisher.co.nz

©2015 Thermo Fisher Scientific Inc. All rights reserved.

BOXED NITRILE

AQUA1900 Premium Non-Med

- 100% Nitrile
- Consistent Quality
- Comfort & Strength
- Textured Grip
- Powder-Free
- 240mm (9")



Size	Hourglass Item No.	Fisher Scientific Cat. No.
S	1901	19-423-038
M	1902	19-423-039
L	1903	19-423-040
XL	1904	19-423-041
100 Pcs/Box		10 Box/Case

FreeStyle1100 Accelerator-Free Exam Glove

- Secure Wet Grip
- Enhanced Strength
- Sulfur-Free
- Anti-Fatigue
- Powder-Free
- 240mm (9")



Size	Hourglass Item No.	Fisher Scientific Cat. No.
XXS	1106	19-423-016
XS	1100	19-423-006
S	1101	19-423-007
M	1102	19-423-008
L	1103	19-423-009
XL	1104	19-423-010
XXL*	1105	15-423-050
100 Pcs/Box		10 Box/Case
*80 Pcs/Box		



A Thermo Fisher Scientific Brand