



## Gloves

Whether you need hand protection for a delicate lab procedure or for an ISO Class 3 or higher cleanroom, our outstanding glove line provides the consistent level of quality you demand.

**JACKSON**  
SAFETY BRAND

**KIMTECH**  
PURE\* BRAND

**KleenGuard**  
BRAND

**Kimberly-Clark**  
\*

# KIMTECH PURE\* G3 and G5 Sterile Gloves

## KIMTECH PURE\* G3 Sterile STERLING\* Nitrile Gloves

Our KIMTECH PURE\* G3 Sterile STERLING\* Nitrile Gloves are Engineered for Protection and Designed for Comfort. When double-donned, the G3 Sterile STERLING\* provides the feel and dexterity of a single glove without the risk of Type 1 latex allergic reactions.

### Product Features:

- Latex-free
- Hand-specific pairs
- Static dissipative in use
- 4.0 mil thickness
- Packaged for aseptic donning

Recommended for ISO Class 4 or higher cleanroom environments



### KIMTECH PURE\* G3 Sterile STERLING\* Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	Color	6	6.5	7	7.5	8	8.5	9	10	Pairs/Bag	Total/Case
KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves - 12"	Sterling	11821	11822	11823	11824	11825	11826	11827	11828	30/10	300

## Advantages of KIMTECH PURE\* Gloves

### Acceptable Quality Level (AQL)


- Refers to the probability of having defective gloves within a lot
- The lower the AQL number, the lower the probability of defective gloves
- KIMTECH PURE\* Gloves exceed ASTM standards for freedom from pinholes and critical defects (ASTM = AQL 2.5; KIMTECH PURE\* Gloves = AQL 1.5).
- A low AQL for pinholes is important – the lower the number of pinholes, the greater the level of barrier protection



Go to [www.kimtech.com/certificates](http://www.kimtech.com/certificates) for detailed Trend & Technical Data

### The Importance of Trend & Technical Data<sup>††</sup>

- Demonstrate consistent product quality over time
- Trend data indicating particle counts and extractable levels are available for KIMTECH PURE\* G3 (nonsterile) Latex and Nitrile Gloves



Example based on historical Critical Nitrile Data  
<sup>††</sup> Tested according to the procedure as described in IEST-RP-CC005.3: Gloves and Finger Cots Used in Cleanrooms and Other Controlled Environments

## KIMTECH PURE\* G3 Sterile White Nitrile Gloves

(formerly SAFESKIN\* Sterile Critical Nitrile Gloves)

### Product Features:

- Latex-free
- Hand-specific pairs
- Low endotoxin
- Packaged for aseptic donning


Recommended for ISO Class 4 or higher cleanroom environments

**STERILE**

Textured palm and fingertips

Available in 12" length

Beaded cuff



### KIMTECH PURE\* G3 Sterile White Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	Color	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0	Total/Case
KIMTECH PURE* G3 Sterile White Nitrile Gloves 12" Hand-Specific Pairs <i>formerly coded</i>	White	56888 HC61160	56889 HC61165	56890 HC61170	56891 HC61175	56892 HC61180	56893 HC61185	56894 HC61190	56887 HC61110	200 pairs

## KIMTECH PURE\* G3 Sterile Latex<sup>†</sup> Gloves

(formerly SAFESKIN\* Sterile Critical Latex<sup>†</sup> Gloves)

### Product Features:

- Natural rubber latex<sup>†</sup>
- Hand-specific pairs
- Low protein and endotoxin
- Packaged for aseptic donning

Recommended for ISO Class 4 or higher cleanroom environments

**STERILE**

Textured palm and fingertips

Available in 12" length

Beaded cuff



### KIMTECH PURE\* G3 Sterile Latex<sup>†</sup> Gloves

Packaged in a Class 100 Cleanroom

Description	Color	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0	Total/Case
KIMTECH PURE* G3 Sterile Latex <sup>†</sup> Gloves 12" Hand-Specific Pairs <i>formerly coded</i>	Natural	56843 HC1360S	56844 HC1365S	56845 HC1370S	56846 HC1375S	56847 HC1380S	56848 HC1385S	56849 HC1390S	56842 HC1310S	200 pairs

† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

## KIMTECH PURE\* G5 Sterile Latex<sup>†</sup> Gloves

(formerly SAFESKIN\* Sterile Controlled Latex<sup>†</sup> Gloves)

### Product Features:

- Natural rubber latex<sup>†</sup>
- Hand-specific pairs
- Packaged for aseptic donning

Recommended for ISO Class 5 or higher cleanroom environments

**STERILE**


Fully textured

Smooth grip

Ideal for double donning

Available in 12" length

Beaded cuff



### KIMTECH PURE\* G5 Sterile Latex<sup>†</sup> Gloves

Packaged in a Class 100 Cleanroom

Description	Color	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0	Total/Case
KIMTECH PURE* G5 Sterile Latex <sup>†</sup> Gloves 12" Hand-Specific Pairs <i>formerly coded</i>	Natural	56851 HC1160S	56852 HC1165S	56853 HC1170S	56854 HC1175S	56855 HC1180S	56856 HC1185S	56857 HC1190S	56850 HC1110S	200 pairs

† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

# KIMTECH PURE\* G3 Gloves

## KIMTECH PURE\* G3 EvT Nitrile Gloves

### Product Features:

- Low particles and extractables – designed specifically for use in cleanroom manufacturing environments
- Textured
- Cleaned to meet the needs of a critical cleanroom environment
- Natural rubber latex-free
- Tacky –provides good grip
- AQL 1.5 for pinholes
- Ambidextrous
- 12" length
- 5 mil thickness
- Color: White
- Value-priced



### KIMTECH PURE\* G3 EvT Nitrile Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	XX-Large	Gloves/Bag	Bags/Case	Total/Case
KIMTECH PURE* G3 EvT Nitrile Gloves 12" Ambi	White	38700	38701	38702	38703	38704	38705	250	6 (5 - XL, 2XL)	1,500 (1,250 - XL, 2XL)

## KIMTECH PURE\* G3 STERLING\* Nitrile Gloves

### Product Features:


- Comfortable fit
- Powder-free
- Nonsterile
- Excellent dexterity and tactile sensitivity
- Ambidextrous with high tack/grip
- **Superior glove strength versus competitive gloves**

*Recommended for ISO Class 4 or higher cleanroom environments*

Textured fingertips

Available in 12" length

Beaded cuff



### KIMTECH PURE\* G3 STERLING\* Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	Color	X-Small	Small	Medium	Large	Large +	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE* G3 STERLING* Nitrile Gloves 12" Ambi	Sterling	99234	99235	99236	99237	99238	99239	250/6	1,500

## KIMTECH PURE\* G3 Nitrile Gloves

(formerly SAFESKIN\* Critical Nitrile Gloves)

**Recommended for ISO Class 4 or higher cleanroom environments**

Textured fingertips

Available in 12" length

Beaded cuff

Bisque finish



### KIMTECH PURE\* G3 Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	Color	X-Small	Small	Medium	Large	Large+	X-Large	Gloves/Bags	Total Case
KIMTECH PURE* G3 White Nitrile Gloves 12" Ambi <i>formerly coded</i>	White	56880 HC61010	56881 HC61011	56882 HC61012	56883 HC61013	56884 HC69478	56886 HC61014	100/10	1,000
KIMTECH PURE* G3 NXT Nitrile Gloves 12" Ambi	White	62990	62991	62992	62993	62995	62994	100/10	1,000
KIMTECH PURE* G3 Blue Nitrile Gloves 12" Ambi <i>formerly coded</i>	Blue	56875 HC1150N	56876 HC2250N	56877 HC3350N	56878 HC4450N	—	56879 HC5550N	100/10	1,000

## KIMTECH PURE\* G3 Latex† Gloves

(formerly SAFESKIN\* Critical Latex† Gloves)

### Product Features:

- Natural rubber latex† for tactile sensitivity & comfort
- Low residual protein content

**Recommended for ISO Class 4 or higher cleanroom environments**

Textured palms and fingertips

Available in 12" length

Beaded cuff



### KIMTECH PURE\* G3 Latex† Gloves

Packaged in a Class 100 Cleanroom

Description	Color	X-Small	Small	Medium	Large	Large+	X-Large	Gloves/Bags	Total/Case		
KIMTECH PURE* G3 Latex† Gloves 12" Ambi <i>formerly coded</i>	Natural	—	56813 HC225	56814 HC335	56815 HC445	—	56816 HC555	100/10	1,000		
		6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0		
KIMTECH PURE* G3 Latex† Gloves 12" Hand-Specific Pairs <i>formerly coded</i>	Natural	56835 HC1360	56836 HC1365	56837 HC1370	56838 HC1375	56839 HC1380	56840 HC1385	56841 HC1390	56834 HC1310	50/4	200 pairs
KIMTECH PURE* G3 Latex† Gloves 12" Hand-Specific Multi-Pack <i>formerly coded</i>	Natural	56827 HC1460	56828 HC1465	56829 HC1470	56830 HC1475	56831 HC1480	56832 HC1485	56833 HC1490	56825 HC1410	50/4	200 pairs

† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

## KIMTECH PURE\* G5 STERLING\* Nitrile Gloves

**Product Features:**

- Comfortable fit
- Powder-free
- Nonsterile
- Excellent dexterity and tactile sensitivity
- Ambidextrous with high tack/grip
- **Superior glove strength versus competitive gloves**

*Recommended for ISO Class 5 or higher cleanroom environments*

Textured fingertips  
Available in 12" length  
Beaded cuff



KIMTECH PURE* G5 STERLING* Nitrile Gloves								Packaged in a Class 100 Cleanroom	
Description	Color	X-Small	Small	Medium	Large	Large+	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE* G5 STERLING* Nitrile Gloves 12" Ambi	Sterling	98184	98185	98186	98187	98188	98189	250/6	1,500

## KIMTECH PURE\* G5 Nitrile Gloves

*(formerly SAFESKIN\* Controlled Nitrile Gloves)*

*Recommended for ISO Class 5 or higher cleanroom environments*

Textured fingertips  
High tack  
Bisque finish  
Beaded cuff



KIMTECH PURE* G5 Nitrile Gloves								Packaged in a Class 100 Cleanroom	
Description	Color	X-Small	Small	Medium	Large	Large+	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE* G5 White Nitrile Gloves 10" Ambi	White	56863	56864	56865	56866	56867	56868	100/10	1,000
KIMTECH PURE* G5 Blue Nitrile Gloves 10" Ambi <i>formerly coded</i>	Blue	56858 HC150N	56859 HC250N	56860 HC350N	56861 HC450N	—	56862 HC550N	100/10	1,000

## KIMTECH PURE<sup>®</sup> G5 Latex<sup>†</sup> Gloves

(formerly SAFESKIN<sup>®</sup> Controlled Latex<sup>†</sup> Gloves)

### Product Features:

- Natural rubber latex<sup>†</sup> for tactile sensitivity & comfort

**Recommended for ISO Class 5 or higher cleanroom environments**

Textured palms and fingertips

Available in 10" & 12" lengths

Beaded cuff



### KIMTECH PURE<sup>®</sup> G5 Latex<sup>†</sup> Gloves

Packaged in a Class 100 Cleanroom

Description	Color	Small	Medium	Large	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE <sup>®</sup> G5 Latex <sup>†</sup> Gloves 10" Ambi <i>formerly coded</i>	Natural	<b>56805</b> HC220	<b>56806</b> HC330	<b>56807</b> HC440	<b>56808</b> HC550	100/10	1,000
KIMTECH PURE <sup>®</sup> G5 Latex <sup>†</sup> Gloves 12" Ambi <i>formerly coded</i>	Natural	<b>56809</b> HC2211	<b>56810</b> HC3311	<b>56811</b> HC4411	<b>56812</b> HC5511	100/10	1,000

† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

## KIMTECH PURE<sup>®</sup> G5 Co-Polymer Gloves

(formerly SAFESKIN<sup>®</sup> Critical Co-Polymer Gloves)

**Recommended for ISO Class 5 or higher cleanroom environments**

Smooth finish

Beaded cuff



### KIMTECH PURE<sup>®</sup> G5 Co-Polymer Gloves

Packaged in a Class 100 Cleanroom

Description	Color	Small	Medium	Large	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE <sup>®</sup> G5 Co-Polymer Gloves 12" Ambi	Clear	<b>61001</b>	<b>61002</b>	<b>61003</b>	<b>61004</b>	100/10	1,000

### DID YOU KNOW?

One square inch of hand surface area has an average of 10,000 microorganisms.

# KIMBERLY-CLARK\* Exam Gloves

## KIMBERLY-CLARK\* PURPLE NITRILE\* Exam Gloves

(formerly SAFESKIN\* PURPLE NITRILE\* Gloves)

### Product Features:

- 6 mil thickness
- Static dissipative in use
- Latex-free
- KIMBERLY-CLARK\* PURPLE NITRILE-XTRA\* material was tested for use with 10 chemotherapy drugs and glutaraldehyde
- PURPLE NITRILE-XTRA\* Exam Gloves meet NFPA Standards 1999:2008

Textured fingertips

Available in 9.5" & 12" lengths

Beaded cuff



### KIMBERLY-CLARK\* PURPLE NITRILE\* Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves 9.5" Ambi	Purple	55080	55081	55082	55083	55084	100/10 (XL = 90/10)	1,000 (XL = 900)
KIMBERLY-CLARK* PURPLE NITRILE-XTRA* Exam Gloves 12" Ambi	Purple	55090	50601	50602	50603	50604	50/10	500

## KIMBERLY-CLARK\* STERLING\* Nitrile Exam Gloves

### Product Features:

- 3.5 mil thickness
- Static dissipative in use
- Latex-free
- Cleared for use in chemotherapy
- STERLING\* NITRILE-XTRA\* Exam Gloves meet NFPA Standards 1999:2008

Textured fingertips

Available in 9.5" & 12" lengths

Beaded cuff



### KIMBERLY-CLARK\* STERLING\* Nitrile Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* STERLING* Nitrile Exam Gloves 9.5" Ambi (Replaces 55085-55089)	Sterling	50705	50706	50707	50708	50709	200/10 (XL = 170/10)	2,000 (XL = 1,700)
KIMBERLY-CLARK* STERLING* NITRILE-XTRA* Exam Gloves 12" Ambi	Sterling	53137	53138	53139	53140	53141	100/10	1,000

## DID YOU KNOW?

Though you can't see it happening, every minute of the day we lose about 30,000 to 40,000 dead skin cells off the surface of our skin.





## KIMBERLY-CLARK\* LAVENDER\* Nitrile Exam Gloves

### Product Features:

- 2.8 mil thickness
- Static dissipative in use
- Latex-free
- Exceptional value vs. standard nitrile

Textured fingertips

9.5" length

Beaded cuff



### KIMBERLY-CLARK\* LAVENDER\* Nitrile Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* LAVENDER* Nitrile Exam Gloves 9.5" Ambi	Lavender	52816	52817	52818	52819	52820	250/10 (XL = 230/10)	2,500 (XL = 2,300)

## KIMBERLY-CLARK\* Latex† Exam Gloves

(formerly SAFESKIN® Latex† Exam Gloves)

### Product Features:

- Natural rubber latex†
- Powder-free
- PFE – 6.7 mil thickness
- PFE-XTRA\* – 10.2 mil thickness

Fully textured

Available in  
9.5" & 12"  
lengths

Beaded cuff



### KIMBERLY-CLARK\* Latex† Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* PFE Latex† Exam Gloves 9.5" Ambi <i>formerly coded</i>	Natural	57110 110	57220 220	57330 330	57440 440	57550 550	100/10 (XL = 90/10)	1,000 (XL = 900)
KIMBERLY-CLARK* PFE-XTRA* Latex† Exam Gloves 12" Ambi	Natural	—	50501	50502	50503	50504	50/10	500 (XL = 500)


† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

For more information, visit: [www.kcproductselector.com/chemguide-gloves](http://www.kcproductselector.com/chemguide-gloves)

## KLEENGUARD\* G10 Arctic Blue Nitrile Gloves

### Product Features:

- Distinctive Arctic Blue color for easy ID
- 100% more gloves per case
- Ambidextrous



### KLEENGUARD\* G10 Arctic Blue Nitrile Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KLEENGUARD* G10 Arctic Blue Nitrile Gloves 9.5" Ambi	Arctic Blue	90095	90096	90097	90098	90099	200/10 (XL = 180/10)	2,000 (XL = 1,800)

## JACKSON SAFETY\* G40 Foam Nitrile Coated Gloves

### Product Features:

- Ecofriendly packaging – no separate insert, no printed bags, less waste
- Silicone-free
- Washable – maintains the EN388 ratings through at least five wash cycles




### JACKSON SAFETY\* G40 Foam Nitrile Coated Gloves

Description	Color	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)	Total/Case
JACKSON SAFETY* G40 Foam Nitrile Coated Gloves	Blue and Black	40225	40226	40227	40228	40229	12 pairs per bag • 5 bags per case

## JACKSON SAFETY\* G60 PURPLE NITRILE\* Level 3 Cut Resistant Gloves

### Product Features:

- Dyneema® fabric for comfortable cut protection
- Excellent dexterity
- High level of cut resistance
- Excellent flexibility




### JACKSON SAFETY\* G60 PURPLE NITRILE\* Level 3 Cut Resistant Gloves

Description	Color	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)	Total/Case
JACKSON SAFETY* G60 PURPLE NITRILE* Level 3 Cut Resistant Gloves	White and Black	97430	97431	97432	97433	97434	1 pair pack • 12 pairs per case

## JACKSON SAFETY\* G80 Nitrile Chemical Resistant Gloves

### Product Features:

- Improved more consistent cotton flock lined for added comfort
- Modified Nitrile formulation for more durability and enhanced chemical resistance
- 15 mil thickness



### JACKSON SAFETY\* G80 Nitrile Chemical Resistant Gloves

Description	Color	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)	Total/Case
JACKSON SAFETY* G80 Nitrile Chemical Resistant Gloves	Green	94445	94446	94447	94448	94449	12 pairs per bag • 5 bags per case

For more information, visit: [www.kcproductselector.com/chemguide-chemgloves](http://www.kcproductselector.com/chemguide-chemgloves)



## Masks

When it comes to facial protection, KIMBERLY-CLARK PROFESSIONAL\* has you covered. Our portfolio of masks delivers comfort, protection and performance for use in your cleanroom, laboratory/healthcare and industrial settings. You'll find a broad range of styles, features and levels of protection. One trusted source for your disposable mask and respirator needs.

**KIMTECH**  
BRAND

**Kimberly-Clark\***

## KIMTECH PURE\* M3 Pleat-Style Face Masks

### Product Features:

- Available with ties and earloops
- Low-linting polyethylene film outer layer
- Available in 7" & 9"
- BICOSOF\* Fabric inner facing provides comfort

*ISO Class 3 or higher cleanroom environments*



**DOUBLE BAGGED**

### KIMTECH PURE\* M3 Pleat-Style Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62465	KIMTECH PURE* M3 Face Mask	White	Knitted earloops	50	10	500
62466	KIMTECH PURE* M3 Face Mask	White	Soft ties	50	10	500

† Filtration efficiency studies are based on ASTM F1215-89 and MIL-M-36954C. Breathability performance is based on differential pressure testing.

**STERILE**

### KIMTECH PURE\* M3 Pleat-Style Sterile Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62494	KIMTECH PURE* M3 Sterile Face Mask	White	9" w/ties	20	10	200
62467	KIMTECH PURE* M3 Sterile Face Mask	White	7" w/ties	20	10	200
62470	KIMTECH PURE* M3 Sterile Face Mask	White	7" w/earloops	20	10	200

† Filtration efficiency studies are based on ASTM F1215-89 and MIL-M-36954C. Breathability performance is based on differential pressure testing.

## KIMTECH PURE\* M3 Pouch-Style Face Masks

### Product Features:

- Low-linting apertured polyethylene film outer layer
- BICOSOF\* Fabric inner facing provides comfort

*ISO Class 3 or higher cleanroom environments*

Large breathing chamber

Gap guard for additional protection



**DOUBLE BAGGED**

### KIMTECH PURE\* M3 Pouch-Style Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62484	KIMTECH PURE* M3 Face Mask	White	Two knitted headbands	50	6	300
62483	KIMTECH PURE* M3 Sterile Face Mask	White	Individually packaged, and two knitted headbands	20	10	200

**STERILE**

## DID YOU KNOW?

Particles exiting the nasal passages have been clocked at more than 100 miles per hour

## KIMTECH PURE\* M5 Face Masks

*ISO Class 5 or higher cleanroom environments*

Three-layer design

Polyester/Cellulose  
inner/outer facing



### KIMTECH PURE\* M5 Pleat-Style Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62692	KIMTECH PURE* M5 Face Mask	Blue	Earloops	50	10	500
62742	KIMTECH PURE* M5 Face Mask	Blue	Soft ties	50	6	300

## KIMTECH PURE\* M6 Face Masks

### Product Features:

- BICOSOF\* Fabric inner facing provides comfort

*ISO Class 6 or higher cleanroom environments*

Two-layer design



### KIMTECH PURE\* M6 Pleat-Style Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62477	KIMTECH PURE* M6 Face Mask	White	Knitted earloops	50	10	500

## KIMTECH PURE\* M7 Veils

### Product Features:

- Provides exceptional breathability
- Constructed of low-linting apertured polyethylene film

*ISO Class 7 or higher cleanroom environments*

Dual snaps for easy  
face adjustment

13" wide & 15" long



### KIMTECH PURE\* M7 Veils

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62751	KIMTECH PURE* M7 Veil	White	Snaps	50	3	150
62757	KIMTECH PURE* M7 Veil	White	Knitted headbands	50	3	150

## KIMBERLY-CLARK\* GUARDALL SHIELD\* Face Shields

**Product Features:**

- Equipped with a foam band and headband
- Can be worn with or without glasses or goggles
- Fog resistant

**KIMBERLY-CLARK\* GUARDALL SHIELD\* Face Shields**

Code	Description	Color	Eaches/Cartons	Total/Case
41204	Face Shield with Full Length Visor	Clear	40	40
41205	Face Shield 3/4 Length Visor	Clear	40	40



## KIMBERLY-CLARK\* SPLASHGUARD\* Visors

**Product Features:**

- Fog free procedure mask
- Pleat style with earloops
- Foam band
- Two levels of fluid protection

**KIMBERLY-CLARK\* SPLASHGUARD\* Visors**

Code	Description	Color	Eaches/Cartons	Total/Case
47137	Mask with Visor for Heavy Fluid Protection	Orange	25/4	100
47147	Mask with Wraparound Visor for Heavy Fluid Protection	Clear	25/4	100



## KIMBERLY-CLARK\* PFR95\* N95 Respirators and Surgical Masks

**Product Features:**

- Polyurethane headband
- Two levels of fluid protection

**KIMBERLY-CLARK\* Respirators**

Code	Description	Color	Size	Eaches/Cartons	Total/Case
62126	KIMBERLY-CLARK* PFR95* N95 Respirator & Surgical Masks for Fluid Protection	White	Regular	50/6	300
46727	KIMBERLY-CLARK* FLUIDSHIELD* PFR95* N95 Respirator & Surgical Masks for Heavy Fluid Protection	Orange	Regular	35/6	210

Large breathing chamber

Pouch-style for added comfort



## KIMBERLY-CLARK\* ASTM F2100-11\* Compliant Masks

**Product Features:**

- Kimberly-Clark is the first to offer masks that meet the ASTM F2100-11 at all three levels
- Packaging clearly identifies the level of protection for fluid resistance & particle filtration

**KIMBERLY-CLARK\* ASTM F2100-11\* Compliant Masks**

Code	Description	Color	Size	Units/Bags	Total/Case
25869	Level 1 Procedure Mask w/Earloops	Blue	Regular	50/10	500
28802	Level 1 Surgical Mask w/Ties	Blue	Regular	50/6	300
28821	Level 2 Mask w/Visor for Moderate Fluid Protection	Blue	Regular	25/4	100
28800	Level 3 Mask w/Visor for High Fluid Protection	Orange	Regular	25/4	100



For more information, visit: [www.kchealthcare.com/masks](http://www.kchealthcare.com/masks)



## Visit us at [www.kimtech.com](http://www.kimtech.com)

KIMBERLY-CLARK\*, KIMTECH\*, JACKSON SAFETY\* and the KIMTECH\* Brand offer you additional resources to help you select the right product for the right application.

KIMTECH\* Sustainability  
<http://www.kimtech.com/reducetoday>

KIMBERLY-CLARK\* STERLING\* Nitrile Exam Glove  
Space & Waste Calculator  
<http://www.kimtech.com/sterlingcalculator>

KIMTECH PURE\* A5 Cleanroom Apparel Cost Calculator  
<http://www.kimtech.com/a5calculator>

KIMTECH\* Brand Promotions  
<http://www.kimtech.com/winningworld>

KIMTECH\* Brand Sterile Bundle  
<http://www.kimtech.com/sterilebundle>

KIMTECH\* Brand Product Samples  
<http://www.kimtech.com/getsamples>

KIMTECH\* Brand Product Certificate Resource  
<http://www.kimtech.com/certificates>

KIMTECH\* Brand Interactive Scientific Catalog  
[www.kimtech.com/catalog](http://www.kimtech.com/catalog)

JACKSON SAFETY\* Interactive Product Catalog  
[www.jacksonsafety.com/catalog](http://www.jacksonsafety.com/catalog)

Scientific Product Selector Tool  
[www.kcproductselector.com/scientific](http://www.kcproductselector.com/scientific)

Latex To Nitrile Glove Conversion Site  
[www.kimtech.com/latexconversion](http://www.kimtech.com/latexconversion)

KLEENGUARD\* Chemical Resistance Guide  
[www.kcproductselector.com/chemicalguide](http://www.kcproductselector.com/chemicalguide)

## Index by Product Code Number

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
110	25	10615	6	16669	45	31697	4	40105	6	45024	9
220	25	10616	6	16670	45	33072	42	40106	6	45025	9
330	25	10621	6	16671	46	33330	32	40225	26	45026	9
440	25	10625	6	16674	45	33390	32	40226	26	45027	9
550	25	10631	6	18425	45	34015	39	40227	26	45092	9
682	10	10636	6	18425	46	34120	34	40228	26	45093	9
683	10	10712	16	18483	45	34133	34	40229	26	45094	9
684	10	10715	16	18624	45	34155	34	41025	39	45095	9
685	10	10718	16	18629	46	34256	34	41026	39	45096	9
686	10	10720	16	19639	42	34623	35	41029	39	45097	9
687	10	10724	16	19794	44	34644	35	41041	39	45232	9
688	10	10730	16	19796	44	34705	34	41043	39	45233	9
1510	40	10736	16	19799	44	34721	34	41044	39	45234	9
1770	38	10740	16	19822	44	34743	34	41048	39	45235	9
1772	38	10745	16	19823	44	34770	39	41055	39	45236	9
3086	38	10748	16	19831	44	34790	39	41204	30	45237	9
5120	38	10754	16	19876	43	34865	39	41205	30	45512	9
5123	38	10772	16	19878	43	34900	39	43744	15	45513	9
5320	38	10824	16	19879	43	34955	39	43745	15	45514	9
5322	38	10830	16	19891	43	34965	39	44302	8	45518	9
5511	35	10836	16	19892	43	35401	39	44303	8	45643	11
5514	35	10840	16	20378	43	35411	39	44304	8	45644	11
5517	35	10845	16	20379	43	35421	39	44305	8	45645	11
5800	38	10848	16	20380	43	35431	39	44306	8	45646	11
5810	38	10854	16	20381	43	35618	6	44307	8	45647	11
5812	38	10872	16	20525	46	35619	6	44312	8	45663	11
5816	38	11430	40	20526	46	35620	6	44313	8	45664	11
5820	38	11821	18	20527	46	35621	6	44314	8	45665	11
5830	38	11822	18	20539	43	35622	6	44315	8	45666	11
5841	38	11823	18	21294	44	35623	6	44316	8	45667	11
5843	38	11824	18	21296	44	35624	6	44317	8	46002	7
5930	39	11825	18	21297	44	35625	6	44322	8	46003	7
6001	37	11826	18	21298	44	36072	4	44323	8	46004	7
6006	37	11827	18	21299	44	36077	4	44324	8	46005	7
6070	32	11828	18	21310	44	36150	6	44325	8	46006	7
6179	33	12225	6	22514	42	36262	6	44326	8	46007	7
6211	37	12226	6	22515	42	36263	6	44327	8	46102	7
6411	37	12227	6	23006	44	36264	6	44332	8	46103	7
6471	37	12230	6	23007	44	36265	6	44333	8	46104	7
9812	10	12231	6	23008	44	36550	15	44334	8	46105	7
9813	10	12810	10	23009	44	36779	15	44335	8	46106	7
9814	10	12811	10	23610	15	36811	14	44336	8	46107	7
9815	10	12812	10	25627	42	36850	15	44337	8	46112	7
9816	10	12813	10	25631	42	36860	15	44442	8	46113	7
9817	10	12814	10	25638	42	36871	15	44443	8	46114	7
9818	10	12815	10	25659	42	36880	14	44444	8	46115	7
9829	10	12889	39	25676	42	36885	14	44445	8	46116	7
9830	10	12890	39	25679	42	37163	8	44446	8	46117	7
9906	40	12891	39	25685	42	37164	8	44452	8	46122	7
10019	6	12914	4	25688	42	37165	8	44453	8	46123	7
10020	12	12915	4	25714	43	37584	8	44454	8	46124	7
10021	12	12916	4	25716	43	37591	8	44455	8	46125	7
10022	12	12917	4	25796	4	38508	8	44456	8	46126	7
10023	12	12919	4	25797	4	38548	39	44480	15	46127	7
10024	12	12920	4	25869	30	38700	20	44481	15	46727	30
10029	6	12921	4	27222	14	38701	20	44492	15	47137	30
10030	12	12922	4	28615	43	38702	20	44494	15	47147	30
10031	12	14384	45	28800	30	38703	20	44495	15	49102	6
10032	12	14387	45	28802	30	38704	20	45002	9	49103	6
10033	12	14399	45	28821	30	38705	20	45003	9	49104	6
10034	12	14401	45	29094	42	39370	14	45004	9	49105	6
10039	6	14402	45	29110	42	39371	14	45005	9	49106	6
10468	6	14403	45	30707	45	39372	14	45006	9	49107	6
10606	6	14404	45	30707	46	40049	6	45007	9	49112	6
10608	6	16362	45	31683	4	40103	6	45022	9	49113	6
10610	6	16668	45	31696	4	40104	6	45023	9	49114	6



CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
49115	6	56848	19	62477	29	88805	4	3005056	45	HC1365	21
49116	6	56849	19	62483	28	88806	4	3005057	45	HC1365S	19
49117	6	56850	19	62484	28	88807	4	3005058	45	HC1370	21
50501	25	56851	19	62494	28	88808	4	3005059	46	HC1370S	19
50502	25	56852	19	62692	29	88842	4	3005066	45	HC1375	21
50503	25	56853	19	62742	29	88843	4	3009537	45	HC1375S	19
50504	25	56854	19	62751	29	88844	4	3009537	46	HC1380	21
50601	24	56855	19	62757	29	88845	4	3009654	45	HC1380S	19
50602	24	56856	19	62990	21	88846	4	3010334	45	HC1385	21
50603	24	56857	19	62991	21	88847	4	3010339	45	HC1385S	19
50604	24	56858	22	62992	21	88848	4	3010339	46	HC1390	21
50705	24	56859	22	62993	21	88849	4	3010343	46	HC1390S	19
50706	24	56860	22	62994	21	88850	4	3011373	42	HC1410	21
50707	24	56861	22	62995	21	89791	15	3011379	42	HC1460	21
50708	24	56862	22	66829	15	90018	13	3011380	42	HC1465	21
50709	24	56863	22	66857	14	90095	26	3011670	44	HC1470	21
52816	25	56864	22	68012	16	90096	26	3011671	44	HC1475	21
52817	25	56865	22	68018	16	90097	26	3011672	44	HC1480	21
52818	25	56866	22	68024	16	90098	26	3011673	44	HC1485	21
52819	25	56867	22	68030	16	90099	26	3011674	44	HC1490	21
52820	25	56868	22	68036	16	91600	40	3011681	44	HC150N	22
53137	24	56875	21	68040	16	92148	40	3011716	43	HC220	23
53138	24	56876	21	68048	16	94445	26	3011717	43	HC2211	23
53139	24	56877	21	68054	16	94446	26	3011718	43	HC225	21
53140	24	56878	21	68072	16	94447	26	3011719	43	HC2250N	21
53141	24	56879	21	68118	16	94448	26	3011720	43	HC250N	22
55080	24	56880	21	68124	16	94449	26	3013296	43	HC330	23
55081	24	56881	21	68130	16	95101	13	3013297	43	HC3311	23
55082	24	56882	21	68136	16	95111	13	3013298	43	HC335	21
55083	24	56883	21	68145	16	95121	13	3013299	43	HC3350N	21
55084	24	56884	21	68148	16	95131	13	3013710	46	HC350N	22
55090	24	56886	21	68154	16	97430	26	3013711	46	HC440	23
56805	23	56887	19	68160	16	97431	26	3013712	46	HC4411	23
56806	23	56888	19	68172	16	97432	26	3013851	43	HC445	21
56807	23	56889	19	68190	16	97433	26	3013853	43	HC4450N	21
56808	23	56890	19	68218	16	97434	26	3013854	43	HC450N	22
56809	23	56891	19	68224	16	97790	15	3013856	43	HC550	23
56810	23	56892	19	68230	16	98184	22	3016306	44	HC550N	22
56811	23	56893	19	68236	16	98185	22	3016307	44	HC5511	23
56812	23	56894	19	68240	16	98186	22	3016308	44	HC555	21
56813	21	57110	25	68245	16	98187	22	3016309	44	HC5550N	21
56814	21	57220	25	68248	16	98188	22	3016310	44	HC61010	21
56815	21	57330	25	68254	16	98189	22	3016311	44	HC61011	21
56816	21	57440	25	68272	16	99048	13	3020277	42	HC61012	21
56825	21	57550	25	68636	16	99234	20	3020278	42	HC61013	21
56827	21	58502	6	68640	16	99235	20	3023024	44	HC61014	21
56828	21	58503	6	68645	16	99236	20	3023025	44	HC61110	19
56829	21	58504	6	68648	16	99237	20	3023026	44	HC61160	19
56830	21	58505	6	68654	16	99238	20	3023028	44	HC61165	19
56831	21	58506	6	68672	16	99239	20	70128281	2	HC61170	19
56832	21	58507	6	69025	12	3000005	45	70128282	2	HC61175	19
56833	21	58512	6	69028	12	3000013	45	70148565	2	HC61180	19
56834	21	58513	6	69600	13	3000048	45	70148566	2	HC61185	19
56835	21	58514	6	69601	13	3000049	45	HC1110S	19	HC61190	19
56836	21	58515	6	69979	12	3000050	45	HC1150N	21	HC69478	21
56837	21	58516	6	69981	12	3000052	45	HC1160S	19		
56838	21	58517	6	75450	36	3000053	45	HC1165S	19		
56839	21	61001	23	75460	36	3000354	42	HC1170S	19		
56840	21	61002	23	75512	35	3000355	42	HC1175S	19		
56841	21	61003	23	76490	33	3000356	42	HC1180S	19		
56842	19	61004	23	77320	37	3000357	42	HC1185S	19		
56843	19	62126	30	88800	4	3000359	42	HC1190S	19		
56844	19	62465	28	88801	4	3003415	45	HC1310	21		
56845	19	62466	28	88802	4	3004880	42	HC1310S	19		
56846	19	62467	28	88803	4	3004882	42	HC1360	21		
56847	19	62470	28	88804	4	3004883	42	HC1360S	19		



## Choose KIMBERLY-CLARK PROFESSIONAL\* Products For All Your Needs In Clean Environments

KIMBERLY-CLARK PROFESSIONAL\* offers you convenience, reliability and technical expertise in a single source for gloves, masks, wipers and apparel for clean environments.

### Certified Quality Documentation:

For Certificates of Conformance, Certificates of Analysis and Certificates of Irradiation, visit [www.kimtech.com/certificates](http://www.kimtech.com/certificates) or e-mail [kimtech@kcc.com](mailto:kimtech@kcc.com).

Scan this QR code and view the KIMTECH\* Brand Scientific Products Catalog online. If you don't have a QR Scanner App, we suggest mobiletag. Download the mobiletag application for free on all App Stores or on your mobile browser at [m.mobiletag.com](http://m.mobiletag.com).

### WARNINGS: Improper use or failure to heed warnings may result in serious injury or death.

Protective Apparel and Hand Protection should be replaced immediately if ripped, torn, abraded or punctured.

Fabrics are not flame resistant unless specified. Keep away from sparks, flames and sources of ignition. Melting fabric can cause severe burns.

It is the responsibility of the user to assess the level of hazardous chemical exposure and then decide on the appropriate personal protective equipment needed for each circumstance.

The information provided within this literature refers to the performance of the fabric only, in the laboratory under controlled conditions.

Kimberly-Clark makes no representation that its products will provide complete protection from exposure to contaminants.

#### JACKSON SAFETY\* Eye Protection

Provides limited eye protection. Additional protective devices are required for full eye coverage and protection.

Does not protect you from all hazards. Even though the lenses in these safety glasses are impact resistant, they can break and shatter and will not protect you from explosions, all flying debris, and severe impacts.

Should not be worn to provide protection from other hazards including, but not limited to, extreme temperatures, high heat, airborne dusts or mists, liquid splashes, sports activities and high-energy optical radiation.

Do not use respirators for sandblasting or to protect against asbestos, lead-based paint, pesticides, oil aerosols<sup>1</sup>, vapors or gases. Do not wear respirators when conditions prevent a good facial seal. It is important to fit check the respirator each time it is worn.

<sup>1</sup>Except R95 or P95

### Our Guarantee

Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you† for your initial purchase, via FREE product, for up to \$1,000.

For more information on KIMBERLY-CLARK PROFESSIONAL\*, visit us online at [www.kcprofessional.com](http://www.kcprofessional.com).

† Guarantee extended to consuming end-user accounts only.

Distributed in Australian and New Zealand by:

**ThermoFisher**  
SCIENTIFIC

**Australia:** For customer service, call 1300-735-292

Visit us online at: [www.thermofisher.com.au](http://www.thermofisher.com.au)

**New Zealand:** For customer service, call 0800-933-966

Visit us online at: [www.thermofisher.co.nz](http://www.thermofisher.co.nz)

©2015 Thermo Fisher Scientific Inc. All rights reserved.

Kimberly-Clark warrants that its products (1) comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers and are warranted for the following periods from end user's date of purchase (verified by valid sales receipt) (a) five years for Balder Technology auto darkening filters; (b) two years for all other auto darkening filters; and (c) one year for powered air purifying respirators; (2) comply with all K-C labeling representations; and (3) are manufactured in compliance with all applicable federal, state, and local laws in effect at the time and place of manufacture of the products. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. K-C is not liable for any kind of special, incidental, or consequential damages. K-C's liability for breach of contract, tort or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.

®/\* Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. Marques déposées de Kimberly-Clark Worldwide, Inc. ou de ses filiales. The COLOR ORANGE, COLOR PURPLE and DAISY DESIGN are Registered Trademarks of KCWW. COLOR GREY and COLOR LAVENDER - Trademarks of KCWW. ©2012 KCWW.DUPONT, KEVLAR, TYVEK and MYLAR are registered trademarks of E.I. DuPont de Nemours and Company. DYNEEMA is a registered trademark of Royal DSM N.V. Printed in the USA. K02705 K9198-12-01 1/13

Effective January 2014