



Scientific Products Catalog



Exceptional
Workplaces*



KIMTECH PURE[®], KIMBERLY-CLARK[®]
and KLEENGUARD[®] Protective Apparel
Pages 3-16



KIMTECH PURE[®], KIMBERLY-CLARK[®],
JACKSON SAFETY[®] and KLEENGUARD[®] Gloves
Pages 17-26



KIMTECH PURE[®] and
KIMBERLY-CLARK[®] Masks
Pages 27-30



KIMTECH PURE[®], KIMTECH SCIENCE[®],
KIMTECH[®] and WYPALL[®] Wipers
Pages 31-40



JACKSON SAFETY[®] and SMITH & WESSON[®]
Eye Protection
Pages 41-46

Protective Apparel

- 4 KIMTECH PURE[®] Protective Apparel
- 6 KLEENGUARD[®] Protective Apparel
- 12 KIMBERLY-CLARK[®] Protective Apparel
- 14 KIMTECH PURE[®], KLEENGUARD[®] AND KIMBERLY-CLARK[®] Accessories
- 16 KIMGUARD[®] Sterilization Wrap

Gloves

- 18 KIMTECH PURE[®] G3 and G5 Sterile Gloves
- 20 KIMTECH PURE[®] G3 Gloves
- 22 KIMTECH PURE[®] G5 Gloves
- 24 KIMBERLY-CLARK[®] Exam Gloves
- 26 KLEENGUARD[®] and JACKSON SAFETY[®] Hand Protection

Masks, Veils, Face Shields and Visors

- 28 KIMTECH PURE[®] Masks and Veils
- 30 KIMBERLY-CLARK[®] and JACKSON SAFETY[®] Respirators

Wipers

- 32 KIMTECH PURE[®] Wipers
- 34 KIMTECH SCIENCE[®] Wipers
- 36 KIMTECH SCIENCE[®] Bench Top Protector
- 37 KIMTECH[®] Wipers for the WETTASK[®] Refillable Wet Wiping System
- 38 WYPALL[®] Wipers
- 40 KLEENEX[®] Hand Sanitizer and SCOTT[®] C-Fold Towels

Eye Protection

- 42 JACKSON SAFETY[®] Eye Protection
- 44 SMITH & WESSON[®] and JACKSON SAFETY[®] Eye and Goggle Protection

KIMTECH
BRAND

KleenGuard
BRAND

Kimberly-Clark
BRAND

Scott
BRAND

Kleenex
BRAND

WYPALL
BRAND

JACKSON
SAFETY Brand

 **Smith & Wesson**
SAFETY Brand

Exceptional Solutions

For Controlled Environments

At KIMBERLY-CLARK PROFESSIONAL*, we understand what makes your business successful – the people who make your products and processes exceptional. That's why we're dedicated to providing solutions that make their workplace environment exceptional – safer, healthier and more productive.



Did you know that reusable protective apparel may not be meeting your cleanroom standards? Be alert to this cleanroom threat!

Expert consultants will work with your employees to ensure the efficient adoption of our solutions in achieving sustainable success.



We can help you mitigate waste and cross-contamination issues in your current processes and reach Corporate Social Responsibility (CSR) goals with our innovative recycling program.



Exceptional
Workplaces*



Resolve to recycle the right way

We've resolved the recycling process with the right way to recycle. Now the garments and gloves used in your facility can be turned into a variety of useful, eco-friendly products.



Look for this symbol to see which products are eligible for recycling.

For more information, visit: www.kimtech.com/rightcycle





Protective Apparel

KIMTECH PURE* and KLEENGUARD* Protective Apparel deliver the reliable protection you're required to provide, and the peace of mind your team deserves. Both sterile and clean processed options offer an ideal balance of comfort, protection and efficiency.

For more information on apparel sizing, please visit:
www.kimtech.com/apparelsizguide

Cleanroom Classifications

ISO Class	ISO Class Limits [†]	FS209 Class Limits ^{††}
ISO 3	35	1
ISO 4	352	10
ISO 5	3,520	100
ISO 6	35,200	1,000
ISO 7	352,000	10,000
ISO 8	3,520,000	100,000
ISO 9	35,200,000	–

[†] (particles / m³ of air) for particles equal to and larger than 0.5 microns
^{††} (particles / ft³ of air) for particles equal to and larger than 0.5 microns

KIMTECH
PURE* BRAND

KleenGuard*
BRAND

Kimberly-Clark*

KIMTECH PURE* A5 Cleanroom Apparel – Clean Processed & Sterile with CLEAN-DON* Technology

Patent-Pending Highly Visible Blue Line along the inside of the garment signals the proper place to grasp while gowning, helping workers avoid touching the outside of the garment.

Cloth-Like Comfort is delivered by the SMS (Spunbond Meltblown Spunbond) material that goes into the KIMTECH PURE* A5 Coverall.

Built-In Thumb Loop keeps garment sleeve from rolling back and exposing arm or wrist.

Innovative Inside-Out Fold pattern presents the inside of the garment as the package is opened, reducing the risk of touching and contaminating the outside of the apparel.

Patent-Pending Snap Technology features built-in snaps that gather up legs and arms to lower the risk of touching the floor, then automatically release as garment is put on.



Typical Uses:

(ISO 5-8 Cleanroom Environments)

- Compounding
- Mixing
- Filling
- Cleaning



KIMTECH PURE* A5 Sterile Cleanroom Apparel**



Code	Description	Total/Case
88800	Small Coverall, Sterile	25
88801	Medium Coverall, Sterile	25
88802	Large Coverall, Sterile	25
88803	XL Coverall, Sterile	25
88804	2XL Coverall, Sterile	25
88805	3XL Coverall, Sterile	25
88806	4XL Coverall, Sterile	25
12917	5XL Coverall, Sterile	25
12914	6XL – 8XL Coverall, Sterile	25
25797	Hood with Ties, Sterile	100
88807	Hood, Sterile	100
NEW! 36072	Universal Integrated Hood and Mask, Sterile	75
NEW! 36077	Sleeve Protector, Sterile	100 Pairs
NEW! 31696	Universal Boots w/ Grasp Ties and Vinyl Edge, Sterile	100 Pairs
NEW! 31683	SM/MD Boots w/ Grasp Ties and Vinyl Edge, Sterile	100 Pairs
NEW! 31697	XL/2XL Boots w/ Grasp Ties and Vinyl Edge, Sterile	100 Pairs
88808	Universal Boots, Sterile	100 Pairs
12922	SM/MD Boots, w/ Ties, Vinyl, Sterile	100 Pairs
12920	XL/2XL Boots, w/ Ties, Vinyl, Sterile	100 Pairs

KIMTECH PURE* A5 Clean Processed Cleanroom Apparel

Code	Description	Total/Case
88842	Small Coverall	25
88843	Medium Coverall	25
88844	Large Coverall	25
88845	XL Coverall	25
88846	2XL Coverall	25
88847	3XL Coverall	25
88848	4XL Coverall	25
12916	5XL Coverall	25
12915	6XL – 8XL Coverall	25
25796	Hood With Ties	100
88849	Hood	100
88850	Universal Boots	100 Pairs
12921	SM/MD Boots	100 Pairs
12919	XL/2XL Boots	100 Pairs


** Gamma Irradiated SAL:10⁻⁶

KIMTECH PURE[®] A5 Cleanroom Apparel – Sterile with CLEAN-DON[®] Technology




KIMTECH PURE[®] Apparel SMS Fabric

Outer Layers – Cloth-like, yet tough and abrasion-resistant spunbond polypropylene



Middle Layers – Intricate web of microfibers that filter out many water-based liquids and dry particulates



KLEENGUARD* Brand Protective Apparel

KLEENGUARD* A10 Light Duty Protective Apparel

- Meets ANSI/ISEA 101-1996 sizing standards
- SMS material enhancement applies to KLEENGUARD* A10 Coveralls only

KLEENGUARD* A10 Coveralls

Zipper front and elastic wrists			Zipper front, elastic wrists, attached hood			Zipper front, elastic wrists, attached hood and boots		
Code	Size	Total/Case	Code	Size	Total/Case	Code	Size	Total/Case
10468	L	25	12231	L	25	10606	L	25
10636	XL	25	12225	XL	25	10608	XL	25
10616	2XL	25	12226	2XL	25	10631	2XL	25
10625	3XL	25	12227	3XL	25	10610	3XL	25
10621	4XL	25	12230	4XL	25	10615	4XL	25

KLEENGUARD* A10 Lab Coat

Snap front, elastic wrists, no pocket, white

Code	Size	Total/Case
40103	L	50
40104	XL	50
40105	2XL	50
40106	3XL	50



KLEENGUARD* A20 Breathable Particle Protection Apparel

- Breathable, patented MICROFORCE* Barrier SMS Fabric
- NFPA 99 Compliant Antistatic Material

KLEENGUARD* A20 Lab Coat

Snap front, chest and hip pocket, white

Code	Size	Total/Case
10019	M	25
10029	L	25
10039	XL	25
40049	2XL	25



KLEENGUARD* A20 Lab Coat

Hook & loop front, elastic wrists, no pocket

Code	Size	Total/Case
35618	S	30
35619	M	30
35620	L	30
35621	XL	30
35622	2XL	30
35623	3XL	30
35624	4XL	30
35625	5XL	30



KLEENGUARD* A20 Smock

Elastic wrists, no pockets, back opening


Code	Size	Total/Case
36150	Universal	60



KLEENGUARD* A20 Coveralls

Zipper front, elastic back, wrists & ankles


White Code	Blue Denim Code	Size	Total/Case
49102	58502	M	24
49103	58503	L	24
49104	58504	XL	24
49105	58505	2XL	24
49106	58506	3XL	20
49107	58507	4XL	20



KLEENGUARD* A20 Coveralls

Zipper front, elastic back, wrists, ankles & hood

White Code	Blue Denim Code	Size	Total/Case
49112	58512	M	24
49113	58513	L	24
49114	58514	XL	24
49115	58515	2XL	24
49116	58516	3XL	20
49117	58517	4XL	20



KLEENGUARD* A20 Frock/ Professional Jacket

Snap front, long length, elastic wrist, snap at collar, white

Code	Size	Total/Case
36262	M	30
36263	L	30
36264	XL	30
36265	2XL	30



KLEENGUARD® A30 Breathable Splash & Particle Protection Apparel

- Patented REFLEX® Design is 7½ times less likely to rip out than ANSI minimums – provides 12% more chest room and 6% longer body length
- Strong and abrasion resistant
- NFPA 99 compliant antistatic material
- Breathable, patented MICROFORCE® Barrier SMS Fabric
- Heavier & more durable than KLEENGUARD® A10 Coveralls
- 25 times more breathable with comparable particle holdout to Tyvek®

Typical Uses:

- Pharmaceutical research
- Dry particle contamination
- Aerospace
- Maintenance/General manufacturing



KLEENGUARD® A30 Coveralls

Zipper Front with 1" Flap & Elastic Back

Code	Size	Total/Case
46002	M	25
46003	L	25
46004	XL	25
46005	2XL	25
46006	3XL	21
46007	4XL	21



Zipper Front w/1" Flap, Elastic Back, Wrists & Ankles

Code	Size	Total/Case
46102	M	25
46103	L	25
46104	XL	25
46105	2XL	25
46106	3XL	21
46107	4XL	21



KLEENGUARD® A30 Coveralls

Zipper Front w/1" Flap, Elastic Back, Wrists, Ankles & Hood

Code	Size	Total/Case
46112	M	25
46113	L	25
46114	XL	25
46115	2XL	25
46116	3XL	21
46117	4XL	21




Zipper Front w/1" Flap, Elastic Back, Wrists, Ankles, Hood & Boots

Code	Size	Total/Case
46122	M	25
46123	L	25
46124	XL	25
46125	2XL	25
46126	3XL	21
46127	4XL	21



KLEENGUARD® Apparel MICROFORCE® Barrier SMS Fabric

Outer Layers – Cloth-like, yet tough and abrasion-resistant spunbond polypropylene



Middle Layers – Intricate web of microfibers that filter out many water-based liquids and dry particulates

KLEENGUARD* Brand Protective Apparel

KLEENGUARD* A40 Liquid & Particle Protection Apparel

- Better liquid barrier, better particulate barrier, lower in lint than Tyvek®
- Meets ANSI/ISEA 101-1996 sizing standards
- Passes ASTM F1670 testing for penetration of blood & bodily fluids

Typical Uses:

- Liquid handling
- Pressure washing
- Utilities
- Maintenance
- Cleanroom construction
- Medical research facilities

Breathable, microporous film laminate



KLEENGUARD* A40 Coveralls

Zipper front

Code	Color	Size	Total/Case
44302	White	M	25
44303	White	L	25
44304	White	XL	25
44305	White	2XL	25
44306	White	3XL	25
44307	White	4XL	25

KLEENGUARD* A40 Coveralls


Zipper front, elastic wrists & ankles

Code	Color	Size	Total/Case
44312	White	M	25
44313	White	L	25
44314	White	XL	25
44315	White	2XL	25
44316	White	3XL	25
44317	White	4XL	25

KLEENGUARD* A40 Coveralls with Blue Breathable Back

Zipper front, elastic wrists, ankles & hood

Code	Size	Total/Case
37591	L	25
37163	XL	25
37164	2XL	25
37165	3XL	25
38508	4XL	25
37584	5XL/6XL combo	25



KLEENGUARD* A40 Coveralls

Zipper front, elastic wrists, ankles & hood

Code	Color	Size	Total/Case
44322	White	M	25
44323	White	L	25
44324	White	XL	25
44325	White	2XL	25
44326	White	3XL	25
44327	White	4XL	25

KLEENGUARD* A40 Coveralls

Zipper front, elastic wrists, ankles, hood & boots

Code	Color	Size	Total/Case
44332	White	M	25
44333	White	L	25
44334	White	XL	25
44335	White	2XL	25
44336	White	3XL	25
44337	White	4XL	25



KLEENGUARD* A40 Lab Coat

Snap front, chest & hip pockets, white

Code	Size	Total/Case
44452	M	30
44453	L	30
44454	XL	30
44455	2XL	30
44456	3XL	30



5 Snap front, elastic wrists, no pocket, white

Code	Size	Total/Case
44442	M	30
44443	L	30
44444	XL	30
44445	2XL	30
44446	3XL	30



KLEENGUARD[®] A60 Bloodborne Pathogen & Chemical Splash Protection Apparel

- Passes ASTM F1670 & F1671[†] testing for penetration of blood, bodily fluids and bloodborne pathogens
- Patented REFLEX[®] Coverall Design is 7½ times less likely to rip out than ANSI minimums – provides 12% more chest room and 6% longer body length
- Outer layers provide protection of the film barrier against tearing and abrasion

Meets
ASTM F1671[†]
Bloodborne
Pathogens

Seamless front —————
to provide more protection in the primary exposure areas

Taped storm flap —————

Patented three-layer fabric construction features a middle layer of microporous film that allows heat and sweat vapor to escape, yet protects against many dry particulates and liquids, even under pressure



Typical Uses:

- Animal research
- Embalming/forensics
- Laboratories



Endorsed by the American
Bio-Recovery Association (ABRA)
www.americanbiorecoveryassociation.com

KLEENGUARD[®] A60 Coveralls

Zipper front

Code	Color	Size	Total/Case
45232	Denim Blue	M	24
45233	Denim Blue	L	24
45234	Denim Blue	XL	24
45235	Denim Blue	2XL	24
45236	Denim Blue	3XL	20
45237	Denim Blue	4XL	20

KLEENGUARD[®] A60 Coveralls

Zipper front, storm flap, elastic back, wrists, ankles & hood

Code	Color	Size	Total/Case
45022	Denim Blue	M	24
45023	Denim Blue	L	24
45024	Denim Blue	XL	24
45025	Denim Blue	2XL	24
45026	Denim Blue	3XL	20
45027	Denim Blue	4XL	20

KLEENGUARD[®] A60 Coveralls

Zipper front, storm flap, elastic back, wrists & ankles

Code	Color	Size	Total/Case
45002	Denim Blue	M	24
45003	Denim Blue	L	24
45004	Denim Blue	XL	24
45005	Denim Blue	2XL	24
45006	Denim Blue	3XL	20
45007	Denim Blue	4XL	20

KLEENGUARD[®] A60 Coveralls

Zipper front, storm flap, elastic back, wrists, ankles, hood & boots

Code	Color	Size	Total/Case
45092	Denim Blue	M	24
45093	Denim Blue	L	24
45094	Denim Blue	XL	24
45095	Denim Blue	2XL	24
45096	Denim Blue	3XL	20
45097	Denim Blue	4XL	20

KLEENGUARD[®] A60 Lab Coat

Snap front, chest & hip pockets, denim blue

Code	Size	Total/Case
45512	M	25
45513	L	25
45514	XL	25
45518	2XL	25



[†] Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Bloodborne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System

KLEENGUARD* A65 Flame Resistant Lab Apparel

- Comfortable, breathable and durable fabric resists tearing
- Flame resistant compliant per NFPA 701 FLAME PROPAGATION TEST (1989 edition)
- Complies with the flame resistance requirements cited by paragraph 6.1 of CAN/CGSB-155.20-2000
- NFPA 99 compliant antistatic material

KLEENGUARD* A65 Lab Coat

With hook & loop closure, no pocket

Code	Size	Total/Case
12810	M	25
12811	L	25
12812	XL	25
12813	2XL	25
12814	3XL	25
12815	4XL	25



KLEENGUARD* A70 Chemical Spray Protection Apparel

- Patented REFLEX* Design is 7½ times less likely to rip out than ANSI minimums – provides 12% more chest room and 6% longer body length
- 1.5 mil polyethylene film coated spunbond polypropylene fabric
- NFPA 99 compliant antistatic material
- Softer and quieter fabric

Taped storm flap

Liquid-resistant, extra-long zipper

A wick-away lining for enhanced comfort



Typical Uses:

- Hazardous waste remediation
- Environmental cleanup
- Petrochemical/oil refining
- Acid/caustic handling
- Biohazard cleanup
- Tank cleaning
- Pharmaceutical manufacturing



KLEENGUARD* A70 Coveralls - Bound Seams

Zipper front, storm flap, elastic wrists, ankles & hood

Code	Color	Size	Total/Case
09812	Yellow	M	12
09813	Yellow	L	12
09814	Yellow	XL	12
09815	Yellow	2XL	12
09816	Yellow	3XL	12
09817	Yellow	4XL	12
09818	Yellow	5XL	12

Zipper front, storm flap, elastic wrists, hood & boots

Code	Color	Size	Total/Case
00682	Yellow	M	12
00683	Yellow	L	12
00684	Yellow	XL	12
00685	Yellow	2XL	12
00686	Yellow	3XL	12
00687	Yellow	4XL	12
00688	Yellow	5XL	12

KLEENGUARD* A70 Smock - Back Opening

Yellow - elastic wrists, tie front

Code	Size	Total/Case
09829	44"	25
09830	52"	25



For more information, visit: www.kcproductselector.com/chemguide-apparel

KLEENGUARD[®] A80 Chemical Permeation & Jet Liquid Protection Apparel

- Made of Saranex[®] 23-P Film laminated to a durable fabric which provides a soft cloth-like feel against the skin
- Taped, sealable storm flap
- Taped seams
- Seamless front provides more protection in the primary exposure area
- Patented REFLEX[®] Design is 7½ times less likely to rip out than ANSI minimums – provides 12% more chest room and 6% longer body length



Typical Uses:

- Chemical handling/mixing/cleanup
- Hazardous material handling
- Hazardous waste cleanup
- Industrial Hazmat teams
- Petrochemicals
- Utilities

KLEENGUARD[®] A80 Apparel Fabric Properties

Physical Properties	Test Method	Results
Tensile Strength (MD) (CD)	ASTM D5034	64.4 lbs 52.8 lbs
Trapezoidal Tear (MD) (CD)	INDA IST 100.2	20.8 lbs 14.0 lbs
Hydrohead	AATCC 127-1998	703 cm
Flammability	CPSC 1610	Class 1
Blood Penetration	ASTM F1670	Pass
Bloodborne Pathogens	ASTM F1671‡	Pass

More than
25% stronger
than
Tychem[®] SL

Seams
taped with
Saranex[®]
23-P Film

‡ Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Bloodborne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

LEVEL
B/C
SUITS

KLEENGUARD[®] A80 Coveralls • Taped Seams

Zipper front, storm flap, elastic wrists, ankles & respirator-fit hood

Code No.	Size	Case Count
45643	L	12
45644	XL	12
45645	2XL	12
45646	3XL	10
45647	4XL	10



Zipper front, storm flap, elastic wrists, ankles, respirator-fit hood & boots

Code No.	Size	Case Count
45663	L	12
45664	XL	12
45665	2XL	12
45666	3XL	10
45667	4XL	10



For more information, visit: www.kcproductselector.com/chemguide-apparel

KIMBERLY-CLARK* BASIC* Plus Lab Coat

- Light fluid protection

KIMBERLY-CLARK* BASIC* Plus Lab Coat

Code	Description	Size	Total/Case
10020	White, knit collar and cuffs	S	25
10021	White, knit collar and cuffs	M	25
10022	White, knit collar and cuffs	L	25
10023	White, knit collar and cuffs	XL	25
10024	White, knit collar and cuffs	2XL	25
10030	Blue, knit collar and cuffs	S	25
10031	Blue, knit collar and cuffs	M	25
10032	Blue, knit collar and cuffs	L	25
10033	Blue, knit collar and cuffs	XL	25
10034	Blue, knit collar and cuffs	2XL	25



Knit collar and cuffs

Protective 3-layer SMS fabric

KIMBERLY-CLARK* Procedure Gowns

- Protective 3-layer SMS fabric
- Light to moderate fluid protection
- Knit cuff
- VELCRO® neck closure
- Tie waist
- Full back

KIMBERLY-CLARK* Procedure Gowns

Code	Size	Color	Gowns/Bags	Total/Case
69025	Univ	Blue	10/6	60
69028	XL	Blue	10/6	60



KIMBERLY-CLARK* Isolation Cover Gowns

- Protective 3-layer SMS fabric
- Low to moderate fluid protection
- Elastic cuff
- Tape-tab neck closure
- Tie waist

KIMBERLY-CLARK* Isolation Cover Gowns


Code	Size	Color	Gowns/Bags	Total/Case
69979	Univ	Yellow	10/10	100
69981	Univ	Blue	10/10	100



KIMBERLY-CLARK* Impervious Comfort Gowns

- Impervious film laminate
- Heavy fluid protection in front
- Knit cuff
- Perforated open back

KIMBERLY-CLARK* Impervious Comfort Gowns


Code	Image	Size	Color	Gowns/Bags	Total/Case
69600		Univ	Blue	10/10	100
69601		2XL	Blue	10/10	100



KIMBERLY-CLARK* Non-Reinforced Surgical Gowns - Sterile

- Set-in sleeves
- Adjustable neckline w/yellow net binding – hook & loop
- Book Fold

KIMBERLY-CLARK* Non-Reinforced Surgical Gowns - Sterile


Code	Image	Size	Color	Total/Case
90018		L	Blue	36
90048		XL	Blue	32



KIMBERLY-CLARK* ULTRA Surgical Gowns - Sterile

- Provides AAMI Level 3 liquid barrier protection in the front area
- Raglan sleeves
- Adjustable neckline w/yellow neck binding – hook & loop
- Knit cuff
- Generous cut
- Book fold

KIMBERLY-CLARK* ULTRA Surgical Gowns - Sterile

Code	Image	Size	Color	Total/Case
95101		Small	Blue	34
95111		Large	Blue	32
95121		XL	Blue	30
95131		2XL	Blue	28



For more information, visit: www.khealthcare.com/surgicalgowns

KIMTECH PURE* A8 Shoe Covers

- Dry particle and fluid protection
- Extra grip
- Low linting, Category 1, Helmke Drum Test
- Packed for cleanroom



KIMTECH PURE* A8 Shoe Covers

Code	Description	Size	Color	Eaches/Bags	Total/Case
NEW! 39370	A8 Shoe Covers w/ seamless bottom	Small/Med	White	100/3	300
NEW! 39371	A8 Shoe Covers w/ seamless bottom	Univ	White	100/3	300
NEW! 39372	A8 Shoe Covers w/ seamless bottom	XL/2XL	White	100/3	300

KLEENGUARD* A10 Shoe Covers

- Dry particle protection
- Packed for cleanroom



KLEENGUARD* A10 Shoe Covers

Code	Description	Size	Color	Eaches/Bags	Total/Case
27222	A10 Shoe Covers w/seamed elastic sole, toe & heel, with traction	XL/2XL	White	150/2	300
36811	A10 Shoe Covers w/seamed elastic sole, toe & heel, low skid	Univ	Blue	150/2	300

KLEENGUARD* A20 Shoe and Boot Covers

- Dry particle protection



KLEENGUARD* A20 Shoe and Boot Covers

Code	Description	Size	Color	Total/Case
36885	A20 Shoe Covers w/seamed elastic sole, toe & heel	Univ	White	300
36880	A20 Boot Covers	Univ	White	300



- Dry particle protection
- No tack
- Packed for cleanroom
- Recommended use - under boot covers



KLEENGUARD* A20 Shoe Covers

Code	Description	Size	Color	Eaches/Bags	Total/Case
66857	A20 Shoe Covers w/ seamless elastic sole, toe & heel	XL/2XL	Blue	100/3	300



Resistant to dry particulates



Repellent to liquid splashes



Traction



Durability



KLEENGUARD* A40 Shoe Covers

- Liquid splash protection
- Packed for cleanroom



KLEENGUARD* A40 Shoe Covers

Code	Description	Size	Color	Eaches/Bags	Total/Case
44492	A40 Shoe Covers w/vinyl sole, elastic toe & heel	M/L	White	100/4	400
44494	A40 Shoe Covers w/vinyl sole, elastic toe & heel	XL/2XL	White	100/4	400
36779	A40 Boot Covers w/ vinyl sole	M/L	White	100/4	400
44495	A40 Boot Covers w/ vinyl sole	XL/2XL	White	100/4	400



KLEENGUARD* Bouffant Caps

- Dry particle protection
- Latex free
- Packed for cleanroom



KLEENGUARD* A10 and A20 Bouffant Caps

Code	Description	Size	Color	Eaches/Bags	Total/Case
36850	A10 Bouffant Cap	Medium - 21"	White	100/10	1,000
36860	A10 Bouffant Cap	Large - 24"	White	100/10	1,000
66829	A20 Bouffant Cap	Large - 24"	White	100/5	500



KLEENGUARD* and KIMBERLY-CLARK* Sleeves

- Wide range of protection

KLEENGUARD* and KIMBERLY-CLARK* Sleeves

Code	Description	Size	Color	Eaches/Bags	Total/Case
23610	A10 Sleeve - Dry Particle Protection	Univ	White	200	200
44480	A40 Sleeve - Dry Particle & Liquid Protection	Univ	White	200	200
36871	A60 Sleeve - Bloodborne Pathogen & Chemical Splash Protection	Univ	Blue	200	200
89791	KIMBERLY-CLARK* - Sterile - Cleanroom Packaged	Univ	Blue	100	100



KLEENGUARD* Aprons

- Wide range of protection

KLEENGUARD* Aprons

Code	Description	Size	Color	Total/Case
43744	A10 Apron - Polyethylene Coated	28" x 36"	White	100
36550	A20 Apron	28" x 40"	White	100
43745	A20 Apron - Polyethylene Coated	28" x 46"	White	500
44481	A40 Apron - Ties in Back	28" x 40"	White	100
97790	A70 Apron - Ties in Back	44"	Yellow	100



KIMGUARD* Sterilization Wrap

Today, Kimberly-Clark has six grades of wrap fabric — from KC100 to KC600 — as well as KIMGUARD ONE-STEP*, to offer precisely the right protection for your needs. Choosing is easy; the higher the number, the stronger the wrap. Plus, every grade features our exclusive Daisy¹ bond pattern, so you'll know at a glance you've selected Kimberly-Clark KIMGUARD* Sterilization Wrap, the brand that delivers quality and protection, time after time.

- Steam and EO compatible
- Sequential (single sheet) wrap
- Also available in simultaneous (two sheets bonded together) KIMGUARD* ONE-STEP



KIMGUARD* KC100 Sterilization Wrap

KIMGUARD* KC100 Sterilization Wrap

Code	Size	Sheets per Case
10712	12" x 12" / 30cm x 30cm	1000
10715	15" x 15" / 38cm x 38cm	1000
10718	18" x 18" / 45cm x 45cm	1000
10720	20" x 20" / 50cm x 50cm	1000
10724	24" x 24" / 60cm x 60cm	500
10730	30" x 30" / 76cm x 76cm	300
10736	36" x 36" / 91cm x 91cm	300
10740	40" x 40" / 101cm x 101cm	250
10745	45" x 45" / 114cm x 114cm	250
10748	48" x 48" / 121cm x 121cm	250
10754	54" x 54" / 137cm x 137cm	100
10772	54" x 72" / 137cm x 182cm	100

KIMGUARD* KC400 Sterilization Wrap

KIMGUARD* KC400 Sterilization Wrap

Code	Size	Sheets per Case
68218	18" x 18" / 45cm x 45cm	600
68224	24" x 24" / 60cm x 60cm	250
68230	30" x 30" / 76cm x 76cm	200
68236	36" x 36" / 91cm x 91cm	150
68240	40" x 40" / 101cm x 101cm	150
68245	45" x 45" / 114cm x 114cm	100
68248	48" x 48" / 121cm x 121cm	50
68254	54" x 54" / 137cm x 137cm	50
68272	54" x 72" / 137cm x 182cm	50

KIMGUARD* KC200 Sterilization Wrap

KIMGUARD* KC200 Sterilization Wrap

Code	Size	Sheets per Case
68012	12" x 12" / 30cm x 30cm	1000
68018	18" x 18" / 45cm x 45cm	1000
68024	24" x 24" / 60cm x 60cm	500
68030	30" x 30" / 76cm x 76cm	300
68036	36" x 36" / 91cm x 91cm	300
68040	40" x 40" / 101cm x 101cm	250
68048	48" x 48" / 121cm x 121cm	100
68054	54" x 54" / 137cm x 137cm	100
68072	54" x 72" / 137cm x 182cm	100

KIMGUARD* KC500 Sterilization Wrap

KIMGUARD* KC500 Sterilization Wrap

Code	Size	Sheets per Case
68118	18" x 18" / 45cm x 45cm	600
68124	24" x 24" / 60cm x 60cm	250
68130	30" x 30" / 76cm x 76cm	250
68136	36" x 36" / 91cm x 91cm	150
68145	45" x 45" / 114cm x 114cm	100
68148	48" x 48" / 121cm x 121cm	50
68154	54" x 54" / 137cm x 137cm	50
68160	60" x 60" / 152cm x 152cm	50
68172	54" x 72" / 137cm x 182cm	50
68190	54" x 90" / 137cm x 228cm	50

KIMGUARD* KC300 Sterilization Wrap

KIMGUARD* KC300 Sterilization Wrap

Code	Size	Sheets per Case
10824	24" x 24" / 60cm x 60cm	500
10830	30" x 30" / 76cm x 76cm	300
10836	36" x 36" / 91cm x 91cm	300
10840	40" x 40" / 101cm x 101cm	250
10845	45" x 45" / 114cm x 114cm	100
10848	48" x 48" / 121cm x 121cm	100
10854	54" x 54" / 137cm x 137cm	100
10872	54" x 72" / 137cm x 182cm	100

KIMGUARD* KC600 Sterilization Wrap

KIMGUARD* KC600 Sterilization Wrap

Code	Size	Sheets per Case
68636	36" x 36" / 91cm x 91cm	144
68640	40" x 40" / 101cm x 101cm	96
68645	45" x 45" / 114cm x 114cm	96
68648	48" x 48" / 121cm x 121cm	48
68654	54" x 54" / 137cm x 137cm	48
68672	54" x 72" / 137cm x 182cm	48

For more information, visit:
www.khealthcare.com/sterilewraps



Gloves

Whether you need hand protection for a delicate lab procedure or for an ISO Class 3 or higher cleanroom, our outstanding glove line provides the consistent level of quality you demand.

JACKSON
SAFETY BRAND

KIMTECH
PURE* BRAND

KleenGuard
BRAND

Kimberly-Clark
*

KIMTECH PURE* G3 and G5 Sterile Gloves

KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves

Our KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves are Engineered for Protection and Designed for Comfort. When double-donned, the G3 Sterile STERLING* provides the feel and dexterity of a single glove without the risk of Type 1 latex allergic reactions.

Product Features:

- Latex-free
- Hand-specific pairs
- Static dissipative in use
- 4.0 mil thickness
- Packaged for aseptic donning

Recommended for ISO Class 4 or higher cleanroom environments



KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	Color	6	6.5	7	7.5	8	8.5	9	10	Pairs/Bag	Total/Case
KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves - 12"	Sterling	11821	11822	11823	11824	11825	11826	11827	11828	30/10	300

Advantages of KIMTECH PURE* Gloves

Acceptable Quality Level (AQL)


- Refers to the probability of having defective gloves within a lot
- The lower the AQL number, the lower the probability of defective gloves
- KIMTECH PURE* Gloves exceed ASTM standards for freedom from pinholes and critical defects (ASTM = AQL 2.5; KIMTECH PURE* Gloves = AQL 1.5).
- A low AQL for pinholes is important – the lower the number of pinholes, the greater the level of barrier protection



Go to www.kimtech.com/certificates for detailed Trend & Technical Data

The Importance of Trend & Technical Data^{††}

- Demonstrate consistent product quality over time
- Trend data indicating particle counts and extractable levels are available for KIMTECH PURE* G3 (nonsterile) Latex and Nitrile Gloves



Example based on historical Critical Nitrile Data

^{††} Tested according to the procedure as described in IEST-RP-CC005.3: Gloves and Finger Cots Used in Cleanrooms and Other Controlled Environments

KIMTECH PURE* G3 Sterile White Nitrile Gloves

(formerly SAFESKIN* Sterile Critical Nitrile Gloves)

Product Features:

- Latex-free
- Hand-specific pairs
- Low endotoxin
- Packaged for aseptic donning


Recommended for ISO Class 4 or higher cleanroom environments

STERILE

Textured palm and fingertips

Available in 12" length

Beaded cuff



KIMTECH PURE* G3 Sterile White Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	Color	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0	Total/Case
KIMTECH PURE* G3 Sterile White Nitrile Gloves 12" Hand-Specific Pairs <i>formerly coded</i>	White	56888 HC61160	56889 HC61165	56890 HC61170	56891 HC61175	56892 HC61180	56893 HC61185	56894 HC61190	56887 HC61110	200 pairs

KIMTECH PURE* G3 Sterile Latex[†] Gloves

(formerly SAFESKIN* Sterile Critical Latex[†] Gloves)

Product Features:

- Natural rubber latex[†]
- Hand-specific pairs
- Low protein and endotoxin
- Packaged for aseptic donning

Recommended for ISO Class 4 or higher cleanroom environments

STERILE

Textured palm and fingertips

Available in 12" length

Beaded cuff



KIMTECH PURE* G3 Sterile Latex[†] Gloves

Packaged in a Class 100 Cleanroom

Description	Color	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0	Total/Case
KIMTECH PURE* G3 Sterile Latex [†] Gloves 12" Hand-Specific Pairs <i>formerly coded</i>	Natural	56843 HC1360S	56844 HC1365S	56845 HC1370S	56846 HC1375S	56847 HC1380S	56848 HC1385S	56849 HC1390S	56842 HC1310S	200 pairs

† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

KIMTECH PURE* G5 Sterile Latex[†] Gloves

(formerly SAFESKIN* Sterile Controlled Latex[†] Gloves)

Product Features:

- Natural rubber latex[†]
- Hand-specific pairs
- Packaged for aseptic donning

Recommended for ISO Class 5 or higher cleanroom environments

STERILE


Fully textured

Smooth grip

Ideal for double donning

Available in 12" length

Beaded cuff



KIMTECH PURE* G5 Sterile Latex[†] Gloves

Packaged in a Class 100 Cleanroom

Description	Color	6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0	Total/Case
KIMTECH PURE* G5 Sterile Latex [†] Gloves 12" Hand-Specific Pairs <i>formerly coded</i>	Natural	56851 HC1160S	56852 HC1165S	56853 HC1170S	56854 HC1175S	56855 HC1180S	56856 HC1185S	56857 HC1190S	56850 HC1110S	200 pairs

† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

KIMTECH PURE* G3 Gloves

KIMTECH PURE* G3 EvT Nitrile Gloves

Product Features:

- Low particles and extractables – designed specifically for use in cleanroom manufacturing environments
- Textured
- Cleaned to meet the needs of a critical cleanroom environment
- Natural rubber latex-free
- Tacky –provides good grip
- AQL 1.5 for pinholes
- Ambidextrous
- 12" length
- 5 mil thickness
- Color: White
- Value-priced



KIMTECH PURE* G3 EvT Nitrile Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	XX-Large	Gloves/Bag	Bags/Case	Total/Case
KIMTECH PURE* G3 EvT Nitrile Gloves 12" Ambi	White	38700	38701	38702	38703	38704	38705	250	6 (5 - XL, 2XL)	1,500 (1,250 - XL, 2XL)

KIMTECH PURE* G3 STERLING* Nitrile Gloves

Product Features:


- Comfortable fit
- Powder-free
- Nonsterile
- Excellent dexterity and tactile sensitivity
- Ambidextrous with high tack/grip
- **Superior glove strength versus competitive gloves**

Recommended for ISO Class 4 or higher cleanroom environments

Textured fingertips

Available in 12" length

Beaded cuff



KIMTECH PURE* G3 STERLING* Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	Color	X-Small	Small	Medium	Large	Large +	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE* G3 STERLING* Nitrile Gloves 12" Ambi	Sterling	99234	99235	99236	99237	99238	99239	250/6	1,500

KIMTECH PURE* G3 Nitrile Gloves

(formerly SAFESKIN* Critical Nitrile Gloves)

Recommended for ISO Class 4 or higher cleanroom environments

Textured fingertips

Available in 12" length

Beaded cuff

Bisque finish



KIMTECH PURE* G3 Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	Color	X-Small	Small	Medium	Large	Large+	X-Large	Gloves/Bags	Total Case
KIMTECH PURE* G3 White Nitrile Gloves 12" Ambi <i>formerly coded</i>	White	56880 HC61010	56881 HC61011	56882 HC61012	56883 HC61013	56884 HC69478	56886 HC61014	100/10	1,000
KIMTECH PURE* G3 NXT* Nitrile Gloves 12" Ambi	White	62990	62991	62992	62993	62995	62994	100/10	1,000
KIMTECH PURE* G3 Blue Nitrile Gloves 12" Ambi <i>formerly coded</i>	Blue	56875 HC1150N	56876 HC2250N	56877 HC3350N	56878 HC4450N	—	56879 HC5550N	100/10	1,000

KIMTECH PURE* G3 Latex† Gloves

(formerly SAFESKIN* Critical Latex† Gloves)

Product Features:

- Natural rubber latex† for tactile sensitivity & comfort
- Low residual protein content

Recommended for ISO Class 4 or higher cleanroom environments

Textured palms and fingertips

Available in 12" length

Beaded cuff



KIMTECH PURE* G3 Latex† Gloves

Packaged in a Class 100 Cleanroom

Description	Color	X-Small	Small	Medium	Large	Large+	X-Large	Gloves/Bags	Total/Case		
KIMTECH PURE* G3 Latex† Gloves 12" Ambi <i>formerly coded</i>	Natural	—	56813 HC225	56814 HC335	56815 HC445	—	56816 HC555	100/10	1,000		
		6.0	6.5	7.0	7.5	8.0	8.5	9.0	10.0		
KIMTECH PURE* G3 Latex† Gloves 12" Hand-Specific Pairs <i>formerly coded</i>	Natural	56835 HC1360	56836 HC1365	56837 HC1370	56838 HC1375	56839 HC1380	56840 HC1385	56841 HC1390	56834 HC1310	50/4	200 pairs
KIMTECH PURE* G3 Latex† Gloves 12" Hand-Specific Multi-Pack <i>formerly coded</i>	Natural	56827 HC1460	56828 HC1465	56829 HC1470	56830 HC1475	56831 HC1480	56832 HC1485	56833 HC1490	56825 HC1410	50/4	200 pairs

† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

KIMTECH PURE* G5 STERLING* Nitrile Gloves

Product Features:

- Comfortable fit
- Powder-free
- Nonsterile
- Excellent dexterity and tactile sensitivity
- Ambidextrous with high tack/grip
- **Superior glove strength versus competitive gloves**

Recommended for ISO Class 5 or higher cleanroom environments

Textured fingertips
Available in 12" length
Beaded cuff



KIMTECH PURE* G5 STERLING* Nitrile Gloves								Packaged in a Class 100 Cleanroom	
Description	Color	X-Small	Small	Medium	Large	Large+	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE* G5 STERLING* Nitrile Gloves 12" Ambi	Sterling	98184	98185	98186	98187	98188	98189	250/6	1,500

KIMTECH PURE* G5 Nitrile Gloves

(formerly SAFESKIN Controlled Nitrile Gloves)*

Recommended for ISO Class 5 or higher cleanroom environments

Textured fingertips
High tack
Bisque finish
Beaded cuff



KIMTECH PURE* G5 Nitrile Gloves								Packaged in a Class 100 Cleanroom	
Description	Color	X-Small	Small	Medium	Large	Large+	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE* G5 White Nitrile Gloves 10" Ambi	White	56863	56864	56865	56866	56867	56868	100/10	1,000
KIMTECH PURE* G5 Blue Nitrile Gloves 10" Ambi <i>formerly coded</i>	Blue	56858 HC150N	56859 HC250N	56860 HC350N	56861 HC450N	—	56862 HC550N	100/10	1,000

KIMTECH PURE[®] G5 Latex[†] Gloves

(formerly SAFESKIN[®] Controlled Latex[†] Gloves)

Product Features:

- Natural rubber latex[†] for tactile sensitivity & comfort

Recommended for ISO Class 5 or higher cleanroom environments

Textured palms and fingertips

Available in 10" & 12" lengths

Beaded cuff



KIMTECH PURE[®] G5 Latex[†] Gloves

Packaged in a Class 100 Cleanroom

Description	Color	Small	Medium	Large	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE [®] G5 Latex [†] Gloves 10" Ambi <i>formerly coded</i>	Natural	56805 HC220	56806 HC330	56807 HC440	56808 HC550	100/10	1,000
KIMTECH PURE [®] G5 Latex [†] Gloves 12" Ambi <i>formerly coded</i>	Natural	56809 HC2211	56810 HC3311	56811 HC4411	56812 HC5511	100/10	1,000

† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

KIMTECH PURE[®] G5 Co-Polymer Gloves

(formerly SAFESKIN[®] Critical Co-Polymer Gloves)

Recommended for ISO Class 5 or higher cleanroom environments

Smooth finish

Beaded cuff



KIMTECH PURE[®] G5 Co-Polymer Gloves

Packaged in a Class 100 Cleanroom

Description	Color	Small	Medium	Large	X-Large	Gloves/Bags	Total/Case
KIMTECH PURE [®] G5 Co-Polymer Gloves 12" Ambi	Clear	61001	61002	61003	61004	100/10	1,000

DID YOU KNOW?

One square inch of hand surface area has an average of 10,000 microorganisms.

KIMBERLY-CLARK* Exam Gloves

KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves

(formerly SAFESKIN* PURPLE NITRILE* Gloves)

Product Features:

- 6 mil thickness
- Static dissipative in use
- Latex-free
- KIMBERLY-CLARK* PURPLE NITRILE-XTRA* material was tested for use with 10 chemotherapy drugs and glutaraldehyde
- PURPLE NITRILE-XTRA* Exam Gloves meet NFPA Standards 1999:2008

Textured fingertips

Available in 9.5" & 12" lengths

Beaded cuff



KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* PURPLE NITRILE* Exam Gloves 9.5" Ambi	Purple	55080	55081	55082	55083	55084	100/10 (XL = 90/10)	1,000 (XL = 900)
KIMBERLY-CLARK* PURPLE NITRILE-XTRA* Exam Gloves 12" Ambi	Purple	55090	50601	50602	50603	50604	50/10	500

KIMBERLY-CLARK* STERLING* Nitrile Exam Gloves

Product Features:

- 3.5 mil thickness
- Static dissipative in use
- Latex-free
- Cleared for use in chemotherapy
- STERLING* NITRILE-XTRA* Exam Gloves meet NFPA Standards 1999:2008

Textured fingertips

Available in 9.5" & 12" lengths

Beaded cuff



KIMBERLY-CLARK* STERLING* Nitrile Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* STERLING* Nitrile Exam Gloves 9.5" Ambi (Replaces 55085-55089)	Sterling	50705	50706	50707	50708	50709	200/10 (XL = 170/10)	2,000 (XL = 1,700)
KIMBERLY-CLARK* STERLING* NITRILE-XTRA* Exam Gloves 12" Ambi	Sterling	53137	53138	53139	53140	53141	100/10	1,000

DID YOU KNOW?

Though you can't see it happening, every minute of the day we lose about 30,000 to 40,000 dead skin cells off the surface of our skin.



KIMBERLY-CLARK* LAVENDER* Nitrile Exam Gloves

Product Features:

- 2.8 mil thickness
- Static dissipative in use
- Latex-free
- Exceptional value vs. standard nitrile

Textured fingertips

9.5" length

Beaded cuff

AQL
1.5

KIMBERLY-CLARK* LAVENDER* Nitrile Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* LAVENDER* Nitrile Exam Gloves 9.5" Ambi	Lavender	52816	52817	52818	52819	52820	250/10 (XL = 230/10)	2,500 (XL = 2,300)

KIMBERLY-CLARK* Latex† Exam Gloves

(formerly SAFESKIN® Latex† Exam Gloves)


Product Features:

- Natural rubber latex†
- Powder-free
- PFE – 6.7 mil thickness
- PFE-XTRA* – 10.2 mil thickness

Fully textured

Available in
9.5" & 12"
lengths

Beaded cuff

AQL
1.5

KIMBERLY-CLARK* Latex† Exam Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KIMBERLY-CLARK* PFE Latex† Exam Gloves 9.5" Ambi <i>formerly coded</i>	Natural	57110 110	57220 220	57330 330	57440 440	57550 550	100/10 (XL = 90/10)	1,000 (XL = 900)
KIMBERLY-CLARK* PFE-XTRA* Latex† Exam Gloves 12" Ambi	Natural	—	50501	50502	50503	50504	50/10	500 (XL = 500)


† CAUTION: This product contains Natural Rubber Latex which may cause allergic reactions.

For more information, visit: www.kcproductselector.com/chemguide-gloves

KLEENGUARD* G10 Arctic Blue Nitrile Gloves

Product Features:

- Distinctive Arctic Blue color for easy ID
- 100% more gloves per case
- Ambidextrous



KLEENGUARD* G10 Arctic Blue Nitrile Gloves

Description	Color	X-Small	Small	Medium	Large	X-Large	Gloves/Boxes	Total/Case
KLEENGUARD* G10 Arctic Blue Nitrile Gloves 9.5" Ambi	Arctic Blue	90095	90096	90097	90098	90099	200/10 (XL = 180/10)	2,000 (XL = 1,800)

JACKSON SAFETY* G40 Foam Nitrile Coated Gloves

Product Features:

- Ecofriendly packaging – no separate insert, no printed bags, less waste
- Silicone-free
- Washable – maintains the EN388 ratings through at least five wash cycles




JACKSON SAFETY* G40 Foam Nitrile Coated Gloves

Description	Color	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)	Total/Case
JACKSON SAFETY* G40 Foam Nitrile Coated Gloves	Blue and Black	40225	40226	40227	40228	40229	12 pairs per bag • 5 bags per case

JACKSON SAFETY* G60 PURPLE NITRILE* Level 3 Cut Resistant Gloves

Product Features:

- Dyneema® fabric for comfortable cut protection
- Excellent dexterity
- High level of cut resistance
- Excellent flexibility




JACKSON SAFETY* G60 PURPLE NITRILE* Level 3 Cut Resistant Gloves

Description	Color	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)	Total/Case
JACKSON SAFETY* G60 PURPLE NITRILE* Level 3 Cut Resistant Gloves	White and Black	97430	97431	97432	97433	97434	1 pair pack • 12 pairs per case

JACKSON SAFETY* G80 Nitrile Chemical Resistant Gloves

Product Features:

- Improved more consistent cotton flock lined for added comfort
- Modified Nitrile formulation for more durability and enhanced chemical resistance
- 15 mil thickness



JACKSON SAFETY* G80 Nitrile Chemical Resistant Gloves

Description	Color	7 (S)	8 (M)	9 (L)	10 (XL)	11 (2XL)	Total/Case
JACKSON SAFETY* G80 Nitrile Chemical Resistant Gloves	Green	94445	94446	94447	94448	94449	12 pairs per bag • 5 bags per case

For more information, visit: www.kcproductselector.com/chemguide-chemgloves



Masks

When it comes to facial protection, KIMBERLY-CLARK PROFESSIONAL* has you covered. Our portfolio of masks delivers comfort, protection and performance for use in your cleanroom, laboratory/healthcare and industrial settings. You'll find a broad range of styles, features and levels of protection. One trusted source for your disposable mask and respirator needs.

KIMTECH^{*}
BRAND

Kimberly-Clark^{*}

KIMTECH PURE* M3 Pleat-Style Face Masks

Product Features:

- Available with ties and earloops
- Low-linting polyethylene film outer layer
- Available in 7" & 9"
- BICOSOF* Fabric inner facing provides comfort

ISO Class 3 or higher cleanroom environments



KIMTECH PURE* M3 Pleat-Style Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62465	KIMTECH PURE* M3 Face Mask	White	Knitted earloops	50	10	500
62466	KIMTECH PURE* M3 Face Mask	White	Soft ties	50	10	500

† Filtration efficiency studies are based on ASTM F1215-89 and MIL-M-36954C. Breathability performance is based on differential pressure testing.

DOUBLE BAGGED

KIMTECH PURE* M3 Pleat-Style Sterile Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62494	KIMTECH PURE* M3 Sterile Face Mask	White	9" w/ties	20	10	200
62467	KIMTECH PURE* M3 Sterile Face Mask	White	7" w/ties	20	10	200
62470	KIMTECH PURE* M3 Sterile Face Mask	White	7" w/earloops	20	10	200

† Filtration efficiency studies are based on ASTM F1215-89 and MIL-M-36954C. Breathability performance is based on differential pressure testing.

STERILE

KIMTECH PURE* M3 Pouch-Style Face Masks

Product Features:

- Low-linting apertured polyethylene film outer layer
- BICOSOF* Fabric inner facing provides comfort

ISO Class 3 or higher cleanroom environments

Large breathing chamber

Gap guard for additional protection



DOUBLE BAGGED

KIMTECH PURE* M3 Pouch-Style Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62484	KIMTECH PURE* M3 Face Mask	White	Two knitted headbands	50	6	300
62483	KIMTECH PURE* M3 Sterile Face Mask	White	Individually packaged, and two knitted headbands	20	10	200

STERILE

DID YOU KNOW?

Particles exiting the nasal passages have been clocked at more than 100 miles per hour

KIMTECH PURE* M5 Face Masks

ISO Class 5 or higher cleanroom environments

Three-layer design

Polyester/Cellulose
inner/outer facing



KIMTECH PURE* M5 Pleat-Style Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62692	KIMTECH PURE* M5 Face Mask	Blue	Earloops	50	10	500
62742	KIMTECH PURE* M5 Face Mask	Blue	Soft ties	50	6	300

KIMTECH PURE* M6 Face Masks

Product Features:

- BICOSOF* Fabric inner facing provides comfort

ISO Class 6 or higher cleanroom environments

Two-layer design



KIMTECH PURE* M6 Pleat-Style Face Masks

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62477	KIMTECH PURE* M6 Face Mask	White	Knitted earloops	50	10	500

KIMTECH PURE* M7 Veils

Product Features:

- Provides exceptional breathability
- Constructed of low-linting apertured polyethylene film

ISO Class 7 or higher cleanroom environments

Dual snaps for easy
face adjustment

13" wide & 15" long



KIMTECH PURE* M7 Veils

Code	Description	Color	Style	Units/Bag	Bags/Case	Total/Case
62751	KIMTECH PURE* M7 Veil	White	Snaps	50	3	150
62757	KIMTECH PURE* M7 Veil	White	Knitted headbands	50	3	150

KIMBERLY-CLARK* GUARDALL SHIELD* Face Shields

Product Features:

- Equipped with a foam band and headband
- Can be worn with or without glasses or goggles
- Fog resistant

KIMBERLY-CLARK* GUARDALL SHIELD* Face Shields

Code	Description	Color	Eaches/Cartons	Total/Case
41204	Face Shield with Full Length Visor	Clear	40	40
41205	Face Shield 3/4 Length Visor	Clear	40	40



KIMBERLY-CLARK* SPLASHGUARD* Visors

Product Features:

- Fog free procedure mask
- Pleat style with earloops
- Foam band
- Two levels of fluid protection

KIMBERLY-CLARK* SPLASHGUARD* Visors

Code	Description	Color	Eaches/Cartons	Total/Case
47137	Mask with Visor for Heavy Fluid Protection	Orange	25/4	100
47147	Mask with Wraparound Visor for Heavy Fluid Protection	Clear	25/4	100



KIMBERLY-CLARK* PFR95* N95 Respirators and Surgical Masks

Product Features:

- Polyurethane headband
- Two levels of fluid protection

KIMBERLY-CLARK* Respirators

Code	Description	Color	Size	Eaches/Cartons	Total/Case
62126	KIMBERLY-CLARK* PFR95* N95 Respirator & Surgical Masks for Fluid Protection	White	Regular	50/6	300
46727	KIMBERLY-CLARK* FLUIDSHIELD* PFR95* N95 Respirator & Surgical Masks for Heavy Fluid Protection	Orange	Regular	35/6	210

Large breathing chamber

Pouch-style for added comfort



KIMBERLY-CLARK* ASTM F2100-11* Compliant Masks

Product Features:

- Kimberly-Clark is the first to offer masks that meet the ASTM F2100-11 at all three levels
- Packaging clearly identifies the level of protection for fluid resistance & particle filtration

KIMBERLY-CLARK* ASTM F2100-11* Compliant Masks

Code	Description	Color	Size	Units/Bags	Total/Case
25869	Level 1 Procedure Mask w/Earloops	Blue	Regular	50/10	500
28802	Level 1 Surgical Mask w/Ties	Blue	Regular	50/6	300
28821	Level 2 Mask w/Visor for Moderate Fluid Protection	Blue	Regular	25/4	100
28800	Level 3 Mask w/Visor for High Fluid Protection	Orange	Regular	25/4	100



For more information, visit: www.kchealthcare.com/masks



Wipers

Innovation and precision performance for your scientific environment.

From the very beginning, KIMTECH® Brand has offered scientific environments truly innovative wiper solutions. Over 60 years ago, we created KIMWIPES® Wipers, the first and still the Number 1 lab wiper. From sterile pre-saturated cleanroom wipers to the wiper that's the industry standard for labs and research, KIMTECH® Brand continues to provide the high-performance wiper solutions you need.

WYPALL
BRAND

KIMTECH
PURE® BRAND

KIMTECH
SCIENCE® BRAND

KIMTECH
BRAND

KIMTECH PURE* W4 Wipers



Product Features:

- 100% Polypropylene
- Thermally bonded; no glues or binders
- Acid-, base- and solvent-resistant
- Straight, clean, knife-cut edges minimize lint

Recommended for ISO Class 4 or higher cleanroom environments



KIMTECH PURE* W4 Wipers

Code	Delivery System	Color	Material	Size W" x L"	Bags/Case	Units/Bag	Total/Case
33390	 Flat/Double Bag	White	Polypropylene	9" x 9"	5	100	500
33330	 Flat/Double Bag	White	Polypropylene	12" x 12"	5	100	500

KIMTECH PURE* W4 Pre-Saturated Wipers


Product Features:

- 100% Polypropylene with 70% Isopropyl Alcohol / 30% Deionized Water
- Thermally bonded; no glues or binders
- Acid-, base- and solvent-resistant
- Straight, clean, knife-cut edges minimize lint
- Ideal one-step application of alcohol to critical surfaces

Recommended for ISO Class 4 or higher cleanroom environments



KIMTECH PURE* W4 Wipers — Pre-Saturated

Code	Delivery System	Color	Material	Size W" x L"	Bags/Case	Units/Bag	Total/Case
06070	 Resealable Pouch	White	Polypropylene	9" x 11"	10	40	400

KIMTECH PURE* W4 Sterile Pre-Saturated Wipers

Product Features:

- 100% Polypropylene with 70% Isopropyl Alcohol / 30% Deionized Water
- Thermally bonded; no glues or binders
- Acid-, base- and solvent-resistant
- Straight, clean, knife-cut edges minimize lint
- Ideal one-step application of alcohol to critical surfaces


Recommended for ISO Class 4 or higher cleanroom environments where sterility is required



STERILE

DOUBLE BAGGED
with Case Liner

KIMTECH PURE* W4 Wipers — Sterile, Pre-Saturated

Code	Delivery System	Color	Material	Size W" x L"	Bags/Case	Units/Bag	Total/Case
76490 	Resealable Pouch	White	Polypropylene	9" x 11"	10	40	400

KIMTECH PURE* W5 Wipers


Product Features:

- 50% Rayon / 50% Polyester Spunlace Blend
- Clothlike feel
- No binders or surfactants, resulting in low extractables

Recommended for ISO Class 5 or higher cleanroom environments


Excellent absorption efficiency and capacity

High liquid absorbency rate



DOUBLE BAGGED
with Case Liner

KIMTECH PURE* W5 Wipers

Code	Delivery System	Color	Material	Size W" x L"	Bags/Case	Units/Bag	Total/Case
06179 	Flat/Double Bag	White	Rayon / Polyester	9" x 9"	5	100	500

KIMTECH SCIENCE* KIMWIPES* Wipers

The industry standard for labs and research

Introduced over 60 years ago, KIMTECH SCIENCE* KIMWIPES* Wipers were the first delicate task wiper for laboratories and since then have been the market leader for wiping surfaces, parts, instruments and cleaning lenses.

Other products have tried to imitate, but have never matched the original. Don't put your delicate processes at risk!

Product Features:

- Also available in 2 and 3-ply for increased absorbency
- Gentle texture

One-at-a-time dispensing

Low lint and low extractables



KIMTECH SCIENCE* KIMWIPES* Wipers are preferred over the leading private label delicate task wipers.†

KIMTECH SCIENCE* KIMWIPES* Wipers	Leading Private Label Wipers
100% virgin fiber for purity	Mix of virgin and recycled fiber
Absorbs 18% more water	Absorbs less water
Absorbs water 24% faster	Slower rate of absorption
Preferred for softness; won't scratch delicate surfaces when used wet	Rougher feel

† According to an independent study

KIMTECH SCIENCE* KIMWIPES* Wipers

	Code		Delivery System	Color	Material	Size W" x L"	Boxes/Case	Wipers/Box	Total/Case
1 ply	34120		POP-UP* Box	White	1-ply tissue	4.4" x 8.4"	30	280	8,400
	34155		POP-UP* Box	White	1-ply tissue	4.4" x 8.4"	60	280	16,800
	34133		POP-UP* Box	White	1-ply tissue	11.8" x 11.8"	15	196	2,940
	34256		POP-UP* Box	White	1-ply tissue	14.7" x 16.6"	15	140	2,100
2 ply	34705		POP-UP* Box	White	2-ply tissue	11.8" x 11.8"	15	119	1,785
	34721		POP-UP* Box	White	2-ply tissue	14.7" x 16.6"	15	90	1,350
3 ply	34743		POP-UP* Box	White	3-ply tissue	11.8" x 11.8"	15	119	1,785

KIMTECH SCIENCE* Precision Wipes

Product Features:



- A value-conscious alternative that offers reliable wiping performance for delicate tasks
- Available in 1-ply or 2-ply

Low lint and extractables

Wipers dispense one at a time for ease of use



KIMTECH SCIENCE* Precision Wipes

	Code	Delivery System	Color	Material	Size W" x L"	Boxes/Case	Units/Box	Total/Case
	05511	POP-UP* Box	White	1-ply tissue	4.4" x 8.4"	60	280	16,800
	75512	POP-UP* Box	White	1-ply tissue	11.8" x 11.8"	15	196	2,940
	05514	POP-UP* Box	White	1-ply tissue	14.7" x 16.6"	15	140	2,100
	05517	POP-UP* Box	White	2-ply tissue	14.7" x 16.6"	15	90	1,350

KIMTECH SCIENCE* Lens Cleaning Stations

Product Features:

- Ideal compact unit for your eye safety program
- Use anywhere safety goggles are worn
- Anti-static, anti-fog; silicone-free cleaning solution removes dirt quickly and easily
- Wipers won't scratch delicate surfaces when used wet


Typical Uses:

- Labs
- Manufacturing facilities
- Welding operations
- Refineries
- Food processing
- Printing
- Medical/dental offices

Wipers dispense one at a time for ease of use



KIMTECH SCIENCE* Lens Cleaning Station

	Code	Delivery System	Color	Material	Size W" x L"	Boxes/Case	Units/Box	Total/Case
	34623	POP-UP* Box	White	1-ply tissue	4.4" x 8.4"	4	560 & 8 oz.	2,240
	34644	POP-UP* Box	White	1-ply tissue	4.4" x 8.4"	4	1,120 & 16 oz.	4,480


KIMTECH SCIENCE* Bench Top Protector

Product Features:

- Flat, stable surface helps protect experiments and delicate glassware
- Quickly absorbs spills
- Better wet durability than competitive sheet and roll products
- Slip-resistant polyethylene backing keeps mats in place
- Lies flat and remains flat even after drying

Fluid barrier resists harsh chemicals

Available in sheets or rolls



KIMTECH SCIENCE* Bench Top Protector

Code	Description	Color	Material	Size W" x L"	Wipers/Unit	Units/Carton	Product/Case	Total/Case
75460	Sheet	White/Green	Airlaid / polyethylene	18" x 19.5"	50	8	400	400
75450	Roll	White/Green	Airlaid / polyethylene	19" x 250'	-	-	500'	2 Rolls

KIMTECH® Wipers

for the WETTASK® Refillable Wet Wiping System

Choose the KIMTECH® Wiper that works best for your task & add the chemical solution you already use

Product Features:

- Reduces exposure to chemical vapors & splashes
- Eliminates potential chemical spills
- Clean, fresh wiper every time
- Enclosed system helps to avoid contamination of the wipers and cleaning solution when in use

Typical Uses:

- Surface preparation
- Surface cleaning, sanitizing, and disinfecting
- Cleaning and removing soil, oil and grease
- Solvent application

3 Easy Steps



Available in two sizes

KIMTECH® Wipers for the WETTASK® Refillable Wet Wiping System

Code	Chemical Types	Delivery System	Color	Material	Size W" x L"	Rolls/Case	Sheets/Roll	Total/Case	Dispenser	
06001	Solvents, Detergents, Water-based Cleaners		Roll	White	HYDROKNIT®	12" x 12.5"	6	60	360	 Small bucket free w/case
06006	Solvents, Detergents, Water-based Cleaners		Roll	White	HYDROKNIT®	9" x 15"	2	275	550	 Large bucket 09361
06211	Disinfectants and Sanitizers††		Roll	White	Spunlace	12" x 12.5"	6	90	540	 Small bucket free w/case
06411	Bleach, Disinfectants and Sanitizers††		Roll	White	Meltblown	12" x 12.5"	6	90	540	 Small bucket free w/case
06471	Bleach, Disinfectants and Sanitizers††		Roll	White	Meltblown	12" x 12.5"	6	90	540	No bucket included
77320	Bleach, Disinfectants and Sanitizers††		Roll	White	Meltblown	12" x 12.5"	12	35	420	 Small canister free w/case









†† Various quaternary amine disinfectants have differing shelf lives.

For information regarding the proper use of KIMTECH PREP® Wipers for the WETTASK® System, e-mail us at kimtech@kcc.com.

WYPALL* Brand Wipers

WYPALL* L10 Wipers













WYPALL* L10 Wipers

Code	Delivery System	Color	Material	Size (W x L)	Put-Up	Dispenser
05320 	POP-UP* Box	White	1-ply paper	9.0" x 10.25"	18/125	09352 
05322 	POP-UP* Box	White	1-ply paper	12" x 10.25"	18/125	09352
05120 	Banded	Blue	2-ply paper	9.3" x 10.25"	16/140	N/A
05123 	Banded	Blue	1-ply paper	9.1" x 10.25"	10/224	N/A
01770 	Banded	White	1-ply paper	9.3" x 10.25"	12/200	N/A
01772 	POP-UP* Box	White	1-ply paper	9" x 10.25"	18/110	09352 




WYPALL* L30 Wipers

WYPALL* L30 Wipers













Code	Delivery System	Color	Material	Size (W x L)	Put-Up	Dispenser
05816	POP-UP* Box	White	DRC	16.4" x 9.8"	6/120	09352
05800	POP-UP* Box	White	DRC	16.4" x 9.8"	8/100	09352 
05810 	POP-UP* Box	Blue	DRC	16.4" x 9.8"	8/100	09352
03086	POP-UP* Box	White	DRC	10" x 9.8"	10/120	73900 
05812 	¼ Fold	White	DRC	12.5" x 12"	12/90	09107 
05820 	Center Pull Roll	White	DRC	9.8" x 15.2"	2/300	09335 
05830 	Center Pull Roll Jr.	White	DRC	8" x 15"	6/150	09336 
05841 	Jumbo Roll	White	DRC	12.4" x 13.3"	1/950	80596, 80579 
05843 	Small Roll	White	DRC	11" x 10.4"	24/70	N/A



 Meets EPA Standards for minimum post consumer waste content.

WYPALL* X60 Wipers

WYPALL* X60 Wipers

Code	Delivery System	Color	Material	Size (W x L)	Put-Up	Dispenser
34015	 BRAG* Box	White	HYDROKNIT*	12.5" x 16.8"	1/180	73900
34790	 POP-UP* Box	White	HYDROKNIT*	9.1" x 16.8"	10/126	73900
38548	 POP-UP* Box	White	HYDROKNIT	9.1" x 12.5"	10/156	73900
34865	 1/4 Fold	White	HYDROKNIT*	12.5" x 12"	12/76	09107
34770	 1/4 Fold	White	HYDROKNIT*	11" x 23"	9/100	09107
34900	 Flat Sheet	White	HYDROKNIT*	12.5" x 16.8"	6/150	N/A
35401	 Small Roll	White	HYDROKNIT*	9.8" x 13.4"	12/130	N/A
35411	 Small Roll	Blue	HYDROKNIT*	9.8" x 13.4"	12/130	N/A
35421	 Small Roll	White	HYDROKNIT*	19.6" x 13.4"	6/130	N/A
35431	 Small Roll	Blue	HYDROKNIT*	19.6" x 13.4"	6/130	N/A
34955	 Jumbo Roll	White	HYDROKNIT*	12.5" x 13.4"	1/1100	80596, 80579
34965	 Jumbo Roll	Blue	HYDROKNIT*	12.5" x 13.4"	1/1100	80596, 80579

With 1-at-a-Time
Dispenser Port

WYPALL* X80 Wipers




WYPALL* X80 Wipers

Code	Delivery System	Color	Material	Size (W x L)	Put-Up	Dispenser
05930	 POP-UP* Box	Red	HYDROKNIT*	9.1" x 16.8"	5/80	73900
41048	 POP-UP* Box	White	HYDROKNIT*	9.1" x 16.8"	5/80	73900
41041	 BRAG* Box	Blue	HYDROKNIT*	12.5" x 16.8"	1/160	73900
41044	 BRAG* Box	White	HYDROKNIT*	12.5" x 16.8"	1/160	73900
41026	 1/4 fold	White	HYDROKNIT*	12.5" x 12"	4/50	09107
41029	 1/4 fold	Red	HYDROKNIT*	12.5" x 12"	4/50	09107
41025	 Jumbo Roll	White	HYDROKNIT*	12.5" x 13.4"	1/475	80596, 80579
41055	 Jumbo Roll	Red	HYDROKNIT*	12.5" x 13.4"	1/475	80596, 80579
41043	 Jumbo Roll	Blue	HYDROKNIT*	12.5" x 13.4"	1/475	80596, 80579



WYPALL* X90 Cloths

WYPALL* X90 Cloths

Code	Delivery System	Color	Material	Size (W x L)	Put-Up	Dispenser
12890	 POP-UP* Box	White	HYDROKNIT*	8.3" x 16.8"	5/68	73900
12891	 BRAG* Box	White	HYDROKNIT*	11.1" x 16.8"	1/136	73900
12889	 Jumbo Roll	White	HYDROKNIT*	11.1" x 13.4"	1/450	80596, 80579



KLEENEX® Moisturizing Instant Hand Sanitizer and Dispenser

KLEENEX® Moisturizing Instant Hand Sanitizer and Dispenser


Code	Description	Format	Fragrance/Color	Total/Case
91600	KLEENEX® Moisturizing Instant Hand Sanitizer 62% Alcohol	Liquid	Fruity Cucumber/White	2

Code	Description	Finish	Dimensions (W x H x D)	Total/Case
92148	Electronic Cassette Skin Care Dispenser (1200 mL)	Black	7.25" x 11.48" x 4.0"	1
11430	Electronic Cassette Dispenser Floor Stand	Black	17.7" x 57.4" x 1"	1



SCOTT® C-Fold Towels and Dispenser

SCOTT® C-Fold Towels and Dispenser

Code	Description	Color	Ply	Size (W x L)	Qty/Pack	Packs/Case	Qty/Case
01510	SCOTT® C-Fold Towels 	White	1	10.1" x 13.15"	200 sts.	12 pks.	2,400 sts.

Code	Description	Finish	Dimensions (W x H x D)	Total/Case
09906	WINDOWS® Universal Folded Towel Dispenser	White	13.3" x 18.9" x 5.9"	1





Eye Protection

JACKSON SAFETY® Brand Eye Protection has a legacy of providing high-quality products that fulfill customers' expectations with support from the global SMITH & WESSON® Brand of safety equipment.

JACKSON SAFETY* V10 EXPO* Safety Eyewear



Product Features:

- Lightweight, wraparound protection
- Adjustable temple lengths
- Polycarbonate lens with integral side protection and brow guard
- Lens provides 99.9% U.V. protection
- Meets ANSI Z87.1+ standards
- Custom imprinting available

JACKSON SAFETY* V10 EXPO* Safety Eyewear

Code	Old Code	Description	Pairs/Box
22514	3020277	Clear Lens, Black Frame, Uncoated	12
22515	3020278	Smoke Lens, Black Frame, Uncoated	12

JACKSON SAFETY* V10 ELEMENT* Safety Eyewear



Product Features:

- Lightweight popular styling with an economical price
- Universal nose bridge for comfort
- 100% polycarbonate (no rubber, metal or other components)
- Dielectric¹, thermoplastic construction
- Visible Light Transmittance (VLT): 55% (*applies to product code 33072 only*)
- Commonly used in sodium vapor lighting conditions to help reduce eye strain and increase comfort (*applies to product code 33072 only*)

NEW!

JACKSON SAFETY* V10 ELEMENT* Safety Eyewear

Code	Old Code	Description	Pairs/Box
25627	3004880	Clear Lens, Clear Temples, Uncoated	12
25631	3004882	Smoke Lens, Black Temples, Uncoated	12
25638	3011379	Indoor/Outdoor Lens, Black Temples, Uncoated	12
29094	3004883	Amber Lens, Black Temples, Mirror Coated	12
29110	3011380	Blue Mirror Lens, Black Temples, Mirror Coated	12
33072		Light Blue Lens, Light Blue Temples, Uncoated	12

¹Dielectric means all exposed parts are plastic or rubber.

JACKSON SAFETY* V30 NEMESIS* Safety Eyewear



Product Features:

- Sporty, flexible, lightweight design
- Soft-touch temples for added comfort
- Single-lens wraparound protection
- Every pair includes a neck cord
- Available in reader style lenses

JACKSON SAFETY* V30 NEMESIS* Safety Eyewear

Code	Old Code	Description	Pairs/Box
19639	3011373	Light Blue Lens, Blue Frame	12
25659	3000359	Amber Lens, Black Frame	12
25676 [†]	3000354	Clear Lens, Black Frame	12
25679	3000355	Clear AF Lens, Black Frame	12
25688 [†]	3000356	Smoke Mirror Lens, Black Frame	12
25685	3000357	Indoor/Outdoor Lens, Black Frame	12

JACKSON SAFETY* V30 NEMESIS* Gravity Feed Display Box

Select eyewear available in retail display box. Only part numbers with † are available in gravity feed displays

JACKSON SAFETY* V30 NEMESIS* CSA Safety Eyewear



Product Features:

- Sporty, flexible, lightweight design with soft-touch temples
- Large, single-lens wraparound protection
- Every pair includes a neck cord
- CSA Certified

JACKSON SAFETY* V30 NEMESIS* CSA Safety Eyewear

Code	Old Code	Description	Pairs/Box
20378	3013296	Clear Lens, Black Frame	12
20380[†]	3013298	Smoke Mirror Lens, Black Frame	12
20381[†]	3013299	Indoor/Outdoor Lens, Black Frame	12
20379[†]	3013297	Clear AF Lens, Black Frame	12

[†]CE Certified

JACKSON SAFETY* V40 HELLRAISER* Safety Eyewear



Product Features:

- Flex point temples for comfort and stability
- Single lens wraparound protection
- Lightweight with a sporty style
- Meets ANSI Z87.1+ standards and is CE Certified (EN 166)
- High-performance polycarbonate lens provides 99.9% U.V. protection
- Custom imprinting available



JACKSON SAFETY* V40 HELLRAISER* Gravity Feed Display Box
All standard packs of 12 are shipped in a convenient retail box.

JACKSON SAFETY* V40 HELLRAISER* Safety Eyewear

Code	Old Code	Description	Pairs/Box
20539	3013851	Clear Lens, Black Frame	12
28615[†]	3013853	Clear AF Lens, Black Frame	12
25714[†]	3013854	Smoke Lens, Black Frame	12
25716[†]	3013856	Indoor/Outdoor Lens, Black Frame	12

[†]CE Certified

JACKSON SAFETY* V60 30-06* RX Safety Eyewear



Product Features:

- Stylish, ergonomic, nylon frame
- Soft, comfortable, self-adjusting, rubber nose pad
- "Cheaters" style with diopters for vision assistance



Magnification
Inserts

JACKSON SAFETY* V60 30-06* RX Safety Eyewear

Code	Old Code	Description	Pairs/Box
19876	3011716	Clear Lens, Black Frame, +1.0	6
19878	3011717	Clear Lens, Black Frame, +1.5	6
19879	3011718	Clear Lens, Black Frame, +2.0	6
19891	3011719	Clear Lens, Black Frame, +2.5	6
19892	3011720	Clear Lens, Black Frame, +3.0	6

JACKSON SAFETY* Eyewear meets ANSI Z87.1+ 2010 Standards.

SMITH & WESSON® CALIBER® Safety Eyewear



Clear AF Lens/
Black Frame

BASE CURVE 7

Product Features:

- Ergonomic, adjustable rubber nose pads with built-in airflow holes
- Precision, cut-to-shape spherical lens
- Scratch-resistant polycarbonate lens
- Built-in hinge stops
- Contoured, soft-touch rubber temples
- Polycarbonate lens provides 99.9% U.V. protection
- Meets ANSI Z87.1+ standards and is CSA Z94.3 Compliant
- Custom imprinting available



S&W® CALIBER® Gravity Feed Display Box
All standard packs of 12 are shipped in a convenient retail box.

SMITH & WESSON® CALIBER® Safety Eyewear

Code	Old Code	Description	Pairs/Box
23006	3023024	Clear AF Lens, Black Frame	12
23007	3023025	Smoke AF Lens, Black Frame	12
23008	3023026	Indoor/Outdoor AF Lens, Black Frame	12
23010	3023028	Brown AF Lens, Brown Frame	12

SMITH & WESSON® EQUALIZER® Safety Eyewear



Clear Lens/
Gun Metal Frame

BASE CURVE 8

Product Features:

- Soft rubber nose pads for superior comfort and fit
- Sport grip temples eliminate slippage
- High-performance polycarbonate lens provides 99.9% U.V. protection
- Meets ANSI Z87.1+ standards
- Every glass comes with a SMITH & WESSON® micro-fiber carry pouch that doubles as a cleaning cloth
- Custom imprinting available



S&W® EQUALIZER® Gravity Feed Display Box
All standard packs of 12 are shipped in a convenient retail box.

SMITH & WESSON® EQUALIZER® Safety Eyewear

Code	Old Code	Description	Pairs/Box
21294	3016306	Clear Lens, Gunmetal Frame	12
21296	3016307	Clear AF Lens, Gunmetal Frame	12
21297	3016308	Smoke AF Lens, Gunmetal Frame	12
21298	3016309	Indoor/Outdoor Lens, Gunmetal Frame	12
21299	3016310	Amber Lens, Red Frame	12
21301	3016311	Blue Mirror Lens, Blue Frame	12

SMITH & WESSON® MAGNUM 3G® Safety Eyewear



Smoke Lens/
Black Frame

BASE CURVE 7

Product Features:

- Classic SMITH & WESSON® Magnum Style with built-in hinge stops to help prevent scratches
- Ergonomically dimensioned three-point fit for the ultimate in comfort; self-adjusting nose pad with rubber fingers
- Available in reduced sizes for smaller faces – Magnum Minis
- Cut lens for distortion-free optics

SMITH & WESSON® MAGNUM 3G® Safety Eyewear

Code	Old Code	Description	Pairs/Box
19794	3011670	Clear AF Lens, Black Frame	12
19796	3011671	Copper Blue Block Lens, Black Frame	12
19799	3011672	Clear Lens, Black Frame	12
19823	3011674	Smoke Lens, Black Frame	12
19831	3011681	Indoor/Outdoor Lens, Metallic Gray Frame	12

MAGNUM® MINI

19822	3011673	Mini, Clear Lens, Black Frame	12
--------------	---------	-------------------------------	----

SMITH & WESSON® Eyewear meets ANSI Z87.1+ 2010 Standards.

JACKSON SAFETY* V80 SG34 Goggle Protection



Clear Lens/
Green Frame

Product Features:

- Indirect ventilation
- Economical splash protection

JACKSON SAFETY* V80 SG34 Goggle Protection

Code	Old Code	Description	Pairs/Box
16362	3003415	Clear Lens, Green Frame	50
16674	3005066	Clear Lens, Green Frame	200

JACKSON SAFETY* V80 MONOGOGGLE* 211 Goggle Protection



Clear Lens/
Green Frame

Product Features:

- Adjustable Neoprene strap
- Indirect ventilation – splash goggle
- Soft, pliable plastic nosepiece
- Pantascope angle polycarbonate lens
- Patented vents increase airflow

JACKSON SAFETY* V80 MONOGOGGLE* 211 Goggle Protection

Code	Old Code	Description	Pairs/Box
14384	3000005	Clear Lens, Green Frame	36
14387	3000013	Clear VISICLEAR* Lens, Green Frame	36
16668	3005056	Clear VISICLEAR* Lens, Green Frame With Foam	36
16670	3005058	Clear VISICLEAR* Lens, Smoke Frame	36
16669	3005057	Clear VISICLEAR* Lens, Green Frame Respiratory	36
18425	3009537	Clear VISICLEAR* Replacement Lens	12

JACKSON SAFETY V80 REVOLUTION* OTG Goggle Protection



Clear VISICLEAR*
Lens/Blue Frame

Product Features:

- Adjustable Neoprene strap
- VISICLEAR* anti-fog coating offers anti-scratch and anti-static characteristics
- CROSS VENT* technology keeps air flowing in and splashes out
- Soft, pliable material around the eyes and nose

JACKSON SAFETY* V80 REVOLUTION* OTG Goggle Protection

Code	Old Code	Description	Pairs/Box
14399	3000048	Clear VISICLEAR* Lens, Blue Frame	30
14401	3000049	Smoke VISICLEAR* Lens, Blue Frame	30
14402	3000050	Amber VISICLEAR* Lens, Blue Frame	30
18483	3009654	Clear VISICLEAR* Lens, Black Frame	30
14404	3000053	Clear VISICLEAR* Replacement Lens	12

JACKSON SAFETY* V80 MONOGOGGLE* XTR* OTG Goggle Protection



Clear AF Lens/
Blue Frame

Product Features:

- This stylish goggle protects against splashes, while providing increased visibility due to its cylindrical lens
- A soft, pliable frame conforms to the face for comfort and is easy to clean
- Indirect ventilation
- Can be worn over most prescription eyeglasses

JACKSON SAFETY* V80 MONOGOGGLE* XTR* OTG Goggle Protection

Code	Old Code	Description	Pairs/Box
18624	3010334	Clear AF Lens, Blue Frame	6
30707	3010339	Clear AF Replacement Lens	12

JACKSON SAFETY* V80 WILDCAT* Goggle Protection



Amber AF Lens/
Black Frame

Product Features:

- Extreme heat resistance – will not melt, drip or ignite at 350° F for five minutes
- Adjustable side vents help increase airflow
- Fits naturally to the contour of your face
- Lightweight and comfortable with a pivot headband that allows for optimal placement
- CSA Z94.3 compliant

JACKSON SAFETY* V80 WILDCAT* Goggle Protection

Code	Old Code	Description	Pairs/Box
20525	3013710	Clear AF Lens, Black Frame	12
20526	3013711	Smoke AF Lens, Black Frame	12
20527	3013712	Amber AF Lens, Black Frame	12

JACKSON SAFETY* V90 MONOSHIELD* Goggle Protection



Clear VISICLEAR* Lens/
Green Frame

Product Features:

- Polycarbonate shield
- PVC gaskets snug against goggle
- Supplied with MONOGOGGLE* 211 VISICLEAR*

JACKSON SAFETY* V90 MONOSHIELD* Goggle Protection

Code	Old Code	Description	Pairs/Box
16671	3005059	Clear VISICLEAR* Lens, Green Frame	12
18425	3009537	Clear VISICLEAR* Replacement Lens	12

JACKSON SAFETY* V90 SHIELD* Goggle Protection



Clear AF Lens/
Blue Frame

Product Features:

- Features MONOGOGGLE* XTR*
- Polycarbonate shield
- Curved face shield conforms to shape of face
- Soft, pliable material around the eyes and nose
- Splash and impact protector
- Meets ANSI Z87.1+ standards
- Indirect ventilation

JACKSON SAFETY* V90 SHIELD* Goggle Protection

Code	Old Code	Description	Pairs/Box
18629	3010343	Clear AF Lens, Blue Frame	6
30707	3010339	Clear AF Replacement Lens	12



Visit us at www.kimtech.com

KIMBERLY-CLARK*, KIMTECH*, JACKSON SAFETY* and the KIMTECH* Brand offer you additional resources to help you select the right product for the right application.

KIMTECH* Sustainability
<http://www.kimtech.com/reducetoday>

KIMBERLY-CLARK* STERLING* Nitrile Exam Glove Space & Waste Calculator
<http://www.kimtech.com/sterlingcalculator>

KIMTECH PURE* A5 Cleanroom Apparel Cost Calculator
<http://www.kimtech.com/a5calculator>

KIMTECH* Brand Promotions
<http://www.kimtech.com/winningworld>

KIMTECH* Brand Sterile Bundle
<http://www.kimtech.com/sterilebundle>

KIMTECH* Brand Product Samples
<http://www.kimtech.com/getsamples>

KIMTECH* Brand Product Certificate Resource
<http://www.kimtech.com/certificates>

KIMTECH* Brand Interactive Scientific Catalog
www.kimtech.com/catalog

JACKSON SAFETY* Interactive Product Catalog
www.jacksonsafety.com/catalog

Scientific Product Selector Tool
www.kcproductselector.com/scientific

Latex To Nitrile Glove Conversion Site
www.kimtech.com/latexconversion

KLEENGUARD* Chemical Resistance Guide
www.kcproductselector.com/chemicalguide

Index by Product Code Number

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
110	25	10615	6	16669	45	31697	4	40105	6	45024	9
220	25	10616	6	16670	45	33072	42	40106	6	45025	9
330	25	10621	6	16671	46	33330	32	40225	26	45026	9
440	25	10625	6	16674	45	33390	32	40226	26	45027	9
550	25	10631	6	18425	45	34015	39	40227	26	45092	9
682	10	10636	6	18425	46	34120	34	40228	26	45093	9
683	10	10712	16	18483	45	34133	34	40229	26	45094	9
684	10	10715	16	18624	45	34155	34	41025	39	45095	9
685	10	10718	16	18629	46	34256	34	41026	39	45096	9
686	10	10720	16	19639	42	34623	35	41029	39	45097	9
687	10	10724	16	19794	44	34644	35	41041	39	45232	9
688	10	10730	16	19796	44	34705	34	41043	39	45233	9
1510	40	10736	16	19799	44	34721	34	41044	39	45234	9
1770	38	10740	16	19822	44	34743	34	41048	39	45235	9
1772	38	10745	16	19823	44	34770	39	41055	39	45236	9
3086	38	10748	16	19831	44	34790	39	41204	30	45237	9
5120	38	10754	16	19876	43	34865	39	41205	30	45512	9
5123	38	10772	16	19878	43	34900	39	43744	15	45513	9
5320	38	10824	16	19879	43	34955	39	43745	15	45514	9
5322	38	10830	16	19891	43	34965	39	44302	8	45518	9
5511	35	10836	16	19892	43	35401	39	44303	8	45643	11
5514	35	10840	16	20378	43	35411	39	44304	8	45644	11
5517	35	10845	16	20379	43	35421	39	44305	8	45645	11
5800	38	10848	16	20380	43	35431	39	44306	8	45646	11
5810	38	10854	16	20381	43	35618	6	44307	8	45647	11
5812	38	10872	16	20525	46	35619	6	44312	8	45663	11
5816	38	11430	40	20526	46	35620	6	44313	8	45664	11
5820	38	11821	18	20527	46	35621	6	44314	8	45665	11
5830	38	11822	18	20539	43	35622	6	44315	8	45666	11
5841	38	11823	18	21294	44	35623	6	44316	8	45667	11
5843	38	11824	18	21296	44	35624	6	44317	8	46002	7
5930	39	11825	18	21297	44	35625	6	44322	8	46003	7
6001	37	11826	18	21298	44	36072	4	44323	8	46004	7
6006	37	11827	18	21299	44	36077	4	44324	8	46005	7
6070	32	11828	18	21310	44	36150	6	44325	8	46006	7
6179	33	12225	6	22514	42	36262	6	44326	8	46007	7
6211	37	12226	6	22515	42	36263	6	44327	8	46102	7
6411	37	12227	6	23006	44	36264	6	44332	8	46103	7
6471	37	12230	6	23007	44	36265	6	44333	8	46104	7
9812	10	12231	6	23008	44	36550	15	44334	8	46105	7
9813	10	12810	10	23009	44	36779	15	44335	8	46106	7
9814	10	12811	10	23610	15	36811	14	44336	8	46107	7
9815	10	12812	10	25627	42	36850	15	44337	8	46112	7
9816	10	12813	10	25631	42	36860	15	44442	8	46113	7
9817	10	12814	10	25638	42	36871	15	44443	8	46114	7
9818	10	12815	10	25659	42	36880	14	44444	8	46115	7
9829	10	12889	39	25676	42	36885	14	44445	8	46116	7
9830	10	12890	39	25679	42	37163	8	44446	8	46117	7
9906	40	12891	39	25685	42	37164	8	44452	8	46122	7
10019	6	12914	4	25688	42	37165	8	44453	8	46123	7
10020	12	12915	4	25714	43	37584	8	44454	8	46124	7
10021	12	12916	4	25716	43	37591	8	44455	8	46125	7
10022	12	12917	4	25796	4	38508	8	44456	8	46126	7
10023	12	12919	4	25797	4	38548	39	44480	15	46127	7
10024	12	12920	4	25869	30	38700	20	44481	15	46727	30
10029	6	12921	4	27222	14	38701	20	44492	15	47137	30
10030	12	12922	4	28615	43	38702	20	44494	15	47147	30
10031	12	14384	45	28800	30	38703	20	44495	15	49102	6
10032	12	14387	45	28802	30	38704	20	45002	9	49103	6
10033	12	14399	45	28821	30	38705	20	45003	9	49104	6
10034	12	14401	45	29094	42	39370	14	45004	9	49105	6
10039	6	14402	45	29110	42	39371	14	45005	9	49106	6
10468	6	14403	45	30707	45	39372	14	45006	9	49107	6
10606	6	14404	45	30707	46	40049	6	45007	9	49112	6
10608	6	16362	45	31683	4	40103	6	45022	9	49113	6
10610	6	16668	45	31696	4	40104	6	45023	9	49114	6

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
49115	6	56848	19	62477	29	88805	4	3005056	45	HC1365	21
49116	6	56849	19	62483	28	88806	4	3005057	45	HC1365S	19
49117	6	56850	19	62484	28	88807	4	3005058	45	HC1370	21
50501	25	56851	19	62494	28	88808	4	3005059	46	HC1370S	19
50502	25	56852	19	62692	29	88842	4	3005066	45	HC1375	21
50503	25	56853	19	62742	29	88843	4	3009537	45	HC1375S	19
50504	25	56854	19	62751	29	88844	4	3009537	46	HC1380	21
50601	24	56855	19	62757	29	88845	4	3009654	45	HC1380S	19
50602	24	56856	19	62990	21	88846	4	3010334	45	HC1385	21
50603	24	56857	19	62991	21	88847	4	3010339	45	HC1385S	19
50604	24	56858	22	62992	21	88848	4	3010339	46	HC1390	21
50705	24	56859	22	62993	21	88849	4	3010343	46	HC1390S	19
50706	24	56860	22	62994	21	88850	4	3011373	42	HC1410	21
50707	24	56861	22	62995	21	89791	15	3011379	42	HC1460	21
50708	24	56862	22	66829	15	90018	13	3011380	42	HC1465	21
50709	24	56863	22	66857	14	90095	26	3011670	44	HC1470	21
52816	25	56864	22	68012	16	90096	26	3011671	44	HC1475	21
52817	25	56865	22	68018	16	90097	26	3011672	44	HC1480	21
52818	25	56866	22	68024	16	90098	26	3011673	44	HC1485	21
52819	25	56867	22	68030	16	90099	26	3011674	44	HC1490	21
52820	25	56868	22	68036	16	91600	40	3011681	44	HC150N	22
53137	24	56875	21	68040	16	92148	40	3011716	43	HC220	23
53138	24	56876	21	68048	16	94445	26	3011717	43	HC2211	23
53139	24	56877	21	68054	16	94446	26	3011718	43	HC225	21
53140	24	56878	21	68072	16	94447	26	3011719	43	HC2250N	21
53141	24	56879	21	68118	16	94448	26	3011720	43	HC250N	22
55080	24	56880	21	68124	16	94449	26	3013296	43	HC330	23
55081	24	56881	21	68130	16	95101	13	3013297	43	HC3311	23
55082	24	56882	21	68136	16	95111	13	3013298	43	HC335	21
55083	24	56883	21	68145	16	95121	13	3013299	43	HC3350N	21
55084	24	56884	21	68148	16	95131	13	3013710	46	HC350N	22
55090	24	56886	21	68154	16	97430	26	3013711	46	HC440	23
56805	23	56887	19	68160	16	97431	26	3013712	46	HC4411	23
56806	23	56888	19	68172	16	97432	26	3013851	43	HC445	21
56807	23	56889	19	68190	16	97433	26	3013853	43	HC4450N	21
56808	23	56890	19	68218	16	97434	26	3013854	43	HC450N	22
56809	23	56891	19	68224	16	97790	15	3013856	43	HC550	23
56810	23	56892	19	68230	16	98184	22	3016306	44	HC550N	22
56811	23	56893	19	68236	16	98185	22	3016307	44	HC5511	23
56812	23	56894	19	68240	16	98186	22	3016308	44	HC555	21
56813	21	57110	25	68245	16	98187	22	3016309	44	HC5550N	21
56814	21	57220	25	68248	16	98188	22	3016310	44	HC61010	21
56815	21	57330	25	68254	16	98189	22	3016311	44	HC61011	21
56816	21	57440	25	68272	16	99048	13	3020277	42	HC61012	21
56825	21	57550	25	68636	16	99234	20	3020278	42	HC61013	21
56827	21	58502	6	68640	16	99235	20	3023024	44	HC61014	21
56828	21	58503	6	68645	16	99236	20	3023025	44	HC61110	19
56829	21	58504	6	68648	16	99237	20	3023026	44	HC61160	19
56830	21	58505	6	68654	16	99238	20	3023028	44	HC61165	19
56831	21	58506	6	68672	16	99239	20	70128281	2	HC61170	19
56832	21	58507	6	69025	12	3000005	45	70128282	2	HC61175	19
56833	21	58512	6	69028	12	3000013	45	70148565	2	HC61180	19
56834	21	58513	6	69600	13	3000048	45	70148566	2	HC61185	19
56835	21	58514	6	69601	13	3000049	45	HC1110S	19	HC61190	19
56836	21	58515	6	69979	12	3000050	45	HC1150N	21	HC69478	21
56837	21	58516	6	69981	12	3000052	45	HC1160S	19		
56838	21	58517	6	75450	36	3000053	45	HC1165S	19		
56839	21	61001	23	75460	36	3000354	42	HC1170S	19		
56840	21	61002	23	75512	35	3000355	42	HC1175S	19		
56841	21	61003	23	76490	33	3000356	42	HC1180S	19		
56842	19	61004	23	77320	37	3000357	42	HC1185S	19		
56843	19	62126	30	88800	4	3000359	42	HC1190S	19		
56844	19	62465	28	88801	4	3003415	45	HC1310	21		
56845	19	62466	28	88802	4	3004880	42	HC1310S	19		
56846	19	62467	28	88803	4	3004882	42	HC1360	21		
56847	19	62470	28	88804	4	3004883	42	HC1360S	19		



Choose KIMBERLY-CLARK PROFESSIONAL* Products For All Your Needs In Clean Environments

KIMBERLY-CLARK PROFESSIONAL* offers you convenience, reliability and technical expertise in a single source for gloves, masks, wipers and apparel for clean environments.

Certified Quality Documentation:

For Certificates of Conformance, Certificates of Analysis and Certificates of Irradiation, visit www.kimtech.com/certificates or e-mail kimtech@kcc.com.

Scan this QR code and view the KIMTECH* Brand Scientific Products Catalog online. If you don't have a QR Scanner App, we suggest mobiletag. Download the mobiletag application for free on all App Stores or on your mobile browser at m.mobiletag.com.

WARNINGS: Improper use or failure to heed warnings may result in serious injury or death.

Protective Apparel and Hand Protection should be replaced immediately if ripped, torn, abraded or punctured.

Fabrics are not flame resistant unless specified. Keep away from sparks, flames and sources of ignition. Melting fabric can cause severe burns.

It is the responsibility of the user to assess the level of hazardous chemical exposure and then decide on the appropriate personal protective equipment needed for each circumstance.

The information provided within this literature refers to the performance of the fabric only, in the laboratory under controlled conditions.

Kimberly-Clark makes no representation that its products will provide complete protection from exposure to contaminants.

JACKSON SAFETY* Eye Protection

Provides limited eye protection. Additional protective devices are required for full eye coverage and protection.

Does not protect you from all hazards. Even though the lenses in these safety glasses are impact resistant, they can break and shatter and will not protect you from explosions, all flying debris, and severe impacts.

Should not be worn to provide protection from other hazards including, but not limited to, extreme temperatures, high heat, airborne dusts or mists, liquid splashes, sports activities and high-energy optical radiation.

Do not use respirators for sandblasting or to protect against asbestos, lead-based paint, pesticides, oil aerosols¹, vapors or gases. Do not wear respirators when conditions prevent a good facial seal. It is important to fit check the respirator each time it is worn.

¹Except R95 or P95

Our Guarantee

Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you† for your initial purchase, via FREE product, for up to \$1,000.

For more information on KIMBERLY-CLARK PROFESSIONAL*, visit us online at www.kcprofessional.com.

† Guarantee extended to consuming end-user accounts only.

Distributed in Australian and New Zealand by:

ThermoFisher
SCIENTIFIC

Australia: For customer service, call 1300-735-292

Visit us online at: www.thermofisher.com.au

New Zealand: For customer service, call 0800-933-966

Visit us online at: www.thermofisher.co.nz

©2015 Thermo Fisher Scientific Inc. All rights reserved.

Kimberly-Clark warrants that its products (1) comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers and are warranted for the following periods from end user's date of purchase (verified by valid sales receipt) (a) five years for Balder Technology auto darkening filters; (b) two years for all other auto darkening filters; and (c) one year for powered air purifying respirators; (2) comply with all K-C labeling representations; and (3) are manufactured in compliance with all applicable federal, state, and local laws in effect at the time and place of manufacture of the products. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. K-C is not liable for any kind of special, incidental, or consequential damages. K-C's liability for breach of contract, tort or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.

®/* Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. Marques déposées de Kimberly-Clark Worldwide, Inc. ou de ses filiales. The COLOR ORANGE, COLOR PURPLE and DAISY DESIGN are Registered Trademarks of KCWW. COLOR GREY and COLOR LAVENDER - Trademarks of KCWW. ©2012 KCWW.DUPONT, KEVLAR, TYVEK and MYLAR are registered trademarks of E.I. DuPont de Nemours and Company. DYNEEMA is a registered trademark of Royal DSM N.V. Printed in the USA. K02705 K9198-12-01 1/13

Effective January 2014