

View online: <http://www.labconco.com/product/clear-stoppering-chamber-2/3133>



## Overview

Drying Chambers provide a means to attach flasks to your freeze dry system or contain samples in bulk.

230V, 50/60Hz. For use with 6L, 12L and 18L Freeze Dryers.

**Catalog Number:** 7868030

## Specifications

- **Weight:** 30.0 lbs
- **Weight metric:** 13.6 kg
- **Dimensions:** 23.0" h x 13.5" diameter
- **Electrical:** 230 volts, 50/60 Hz, International
- **Region:** International
- **Applicable Products:** FreeZone 12 Liter, FreeZone 18 Liter, FreeZone 6 Liter

## For Product Assistance

### Authorised distributor

#### In Australia:

For customer service, call 1300-735-292  
To fax an order, use 1800-067-639  
To email an order, [ordersau@thermofisher.com](mailto:ordersau@thermofisher.com)

#### In New Zealand:

For customer service, call 0800-933-966  
To fax an order, use 0800-329-246  
To email an order, [ordersnz@thermofisher.com](mailto:ordersnz@thermofisher.com)

**ThermoFisher**  
SCIENTIFIC

