

Mini Stoppering Chamber

View online: <http://www.labconco.com/product/mini-stoppering-chamber/4851>



Overview

Drying Chambers provide a means to attach flasks to your freeze dry system or contain samples in bulk.

The Mini Stoppering Chamber attaches to a port on a Labconco Drying Chamber or Manifold (sold separately) to stopper small vials (up to 30ml each). Includes 17" of clear PVC tubing.

Can stopper up to 48 serum vials (2ml) or 8 serum vials (30ml) while maintaining vacuum. Spacers, tubing, and adapter for connection to the port are included. Vials and stoppers sold separately.

Catalog Number: 7569000

Specifications

- **Weight:** 15.0 lbs
- **Weight metric:** 6.8 kg
- **Dimensions:** 12" w x 12" d x 10" h
- **Dimensions metric:** 30.5cm x 30.5cm x 25.4cm
- **Applicable Products:** Drying Chambers, Manifolds

For Product Assistance

Authorised distributor

In Australia:

For customer service, call 1300-735-292
To fax an order, use 1800-067-639
To email an order, ordersau@thermofisher.com

In New Zealand:

For customer service, call 0800-933-966
To fax an order, use 0800-329-246
To email an order, ordersnz@thermofisher.com

ThermoFisher
SCIENTIFIC

