

Sani-Cloth®
DISINFECTANT WIPE

70%
WIPE



Sani-Cloth® 70

is the multi-purpose solution,
ideal for disinfection.

- + **Ready to use:** pre-dosed, low linting disposable alcohol disinfection wipes, requiring no mixing or measuring, with no smear formulation.
- + **Proven efficacy:** effective against bacteria, yeast and Enveloped viruses.
- + **Improved efficiencies:** quick drying, no residue formulation containing 70% Isopropanol alcohol allowing the disinfected area or equipment to be used shortly after application.
- + **Maximum flexibility:** can be used on non-invasive medical devices; non-porous hard surfaces.
- + **Recommended usage:** as part of a two-step process with Sani-Cloth® Detergent to clean prior to disinfection. Where surface are soiled, pre-cleaning is required.



Product Efficacy

SANI-CLOTH® 70%	TEST METHOD	ORGANISM	CONTACT TIME*	
Bacteria	EN 16615	Pseudomonas aeruginosa	60 sec (c)	2 min (d)
	EN 16615	Staphylococcus aureus	60 sec (c)	2 min (d)
	EN 16615	Enterococcus hirae	60 sec (c)	
	EN 13727	Pseudomonas aeruginosa	30 sec (c)	
	EN 13727	Staphylococcus aureus	30 sec (c)	
	EN 13727	Enterococcus hirae	30 sec (c)	
Yeasts	EN 16615	Candida albicans	60 sec (c)	2 min (d)
	EN 13624	Candida albicans	30 sec (c)	
Enveloped Virus	EN 14476	Vaccinia Virus	15 sec (c)	60 sec (d)
	EN 14476	Human Coronavirus ATCC VR-740	30 sec (c)	

*Clean Conditions (c) Dirty Conditions (d)

Point of Use Accessories



Bracket for canisters

Bracket: PDIA001GRY

Canister brackets affix securely to the wall, holding the canister in place at its location.

Product Details

Sani Cloth 70% Wipes	ITEM CODE	WIPE COUNT	WIPE SIZE	QTY PER CASE
200ct Large Canister	PDIXP00325A01	200	200mm x 200mm	6

Other formulations available:



Authorised distributor

In Australia:

For customer service, call 1300-735-292

To email an order, ordersanz@thermofisher.com

To order online: thermofisher.com

In New Zealand:

For customer service, call 0800-933-966

To email an order, ordersanz@thermofisher.com

To order online: thermofisher.com

ThermoFisher
SCIENTIFIC

Find out more about how our Layered Approach to Infection Prevention can help you.

PDI LTD, Pywell Road, Willowbrook Est. Industrial Estate, Corby - Northants, NN17 5XJ, UK

11.09.2023 V1 THERMO FISHER SCIENTIFIC

PDI
BE THE DIFFERENCE™